



MATERIAL SAFETY DATA SHEET

Product name: Pneu Neck II Portable Cervical Traction

Model Number: PNB3550

Date prepared: December 2, 2013

1. COMPANY IDENTIFICATION:

Manufacturer: Roscoe Medical Inc.
21973 Commerce Pkwy
Strongsville, OH. 44149

Department: Quality Department

Telephone Number: 800-633-4839

Fax Number: 650-697-3596

2. HEALTH HAZARD INFORMATION

This product is for home care use and specially designed for offering adjuvant therapy for people who have cervical spondylosis. It can reduce neck pain. The patient should operate the units under doctor's instruction or the product user manual.

Potential Health Effects:

Skin: Latex in product may cause an allergic reaction.

3. COMPOSITION / INFORMATION ON INGREDIENTS

Latex: 45% (May cause allergic reaction)

Nylon: 20%

Terylene: 30%

4. FIRST-AID MEASURES

Skin Contact:

If skin has an allergic reaction, please stop using the product, and seek a dermatologist's help.

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com



MATERIAL SAFETY DATA SHEET

5. FIRE-FIGHTING MEASURES

Extinguishing Measures:

Product is easily flammable. If product is to catch fire, use water or dry powder extinguisher directly.

Extinguishing Media:

Water, carbon dioxide, dry powder extinguisher.

Fire and Explosion Hazard Data:

Project Control – Please operate under medical personnel instruction or user manual.

Self Protection – Please choose the right traction pressure based on self physical condition. Make sure the user feels comfortable. Do not increase the traction pressure blindly to avoid accents.

6. ACCIDENTAL RELEASE MEASURES

Remove all sources of ignition. Wear appropriate personal protective equipment. Keep unnecessary personnel from entering the area. Use water, carbon dioxide, or dry powder extinguisher for firefighting.

7. HANDLING AND STORAGE

Storage:

Store in a dry, well-ventilated, no corrosive material location, away from where fire hazard may be acute. Keep the product away from wetness, direct sunshine, sharp objects, and fire. Storage temperature is 10°C ~ 25°C. Please eliminate the air after using.

8. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Dark blue

Odor: Natural Emulsion

Boiling Point: Unknown

Gas Tightness: 0 ~ 30kPa

pH: 7.0 ~ 8.0

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com



MATERIAL SAFETY DATA SHEET

9. STABILITY AND REACTIVITY

Chemical stability: Heat, cold temperature, or strong light can effect stability.

Auto Ignition Temperature: Unknown, but product will burn if there is an existing fire.

10. TOXICOLOGICAL INFORMATION

TMTD: 5g/kg

11. DISPOSAL CONSIDERATIONS

Dispose the product according to local environment regulation.

12. TRANSPORT INFORMATION

Avoid squeezing, heat, low temperature, direct sunshine, and rain.

13. OTHER INFORMATION

The product is designed for patients with a cervical spondylosis problem, and can remit the pain cause by neck pain. The above information is true and intended to be used as a guide. The patient should operate the product according to personal health condition or detailed environment.

To the best of our knowledge, the information contained herein is accurate. However, Roscoe Medical does not assume any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user.

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com