



**Philips France - Healthcare
Solutions Santé à Domicile**

AIR TRAVELING ATTESTATION PR1 (System One) CPAP devices

REMstar Pro C-Flex, REMstar Auto A-Flex, REMstar Auto P-Flex, REMstar Pro Check (Pro C-Flex+) and BiPAP Auto Bi-Flex

Dear Madam, Dear Sir,

The devices listed above are for the following intended use: deliver a therapy by Continuous Positive Airway Pressure for the treatment of Sleep Obstructive Apnea for adults weighting over 30 kgs and breathing spontaneously.

In accordance with what the user's manual indicates, they are compliant with the Standard RTCA/DO-160F section 21, category M; Emission of Radio Frequency Energy and therefore, they can be used in an aircraft when plugged to a source of AC or DC power.

Please note that these devices are not suitable for use on a plane if a modem or a humidifier is connected to the device.

Would you need any further information, please consult the user manual of the device or contact Philips SSD at the following number: +33 (0)2 51 89 36 00.

In Carquefou, on July 25th 2015

Sarah Bahuaud
Responsable Qualité & Affaires Réglementaires

