





Join Our Community
Join Us on Facebook | Follow Us on Twitter
View Us on YouTube

facebook

To Order: 1-800-P-SCHEIN (1-800-772-4346) 8am-9pm (et)

To Fax: 1-800-329-9109 24 Hours

www.henryschein.com/medical

Henry Schein introduces the 2012/2013 Physician Office Lab Supply Guide. This guide is intended to serve as an overview of our most popular diagnostic testing options and supplies needed at the point of care. Don't see what you are looking for, ask your Sales Consultant. From rapid test kits to large analyzers, from CLIA waived to highly complex, we have POL solutions for all your needs.

# **SIEMENS**





































If you have questions regarding diagnostic, pharmaceutical, or vaccine items used in your practice, the DxRx Solutions Hotline is here to help! The DxRx Solutions Hotline can provide CPT Code, package inserts, CLIA product status, product storage requirements, and brand-to-generic comparisons. Let us help you today!

DxRx Solutions Hotline
1-877-523-SHOT 8:30am - 5:00pm, ET
Dxrxsolutions@henryschein.com

# TABLE OF CONTENTS

OSOM	3	Specialty Testing	40-43
Respiratory	4–7	Cardiac	
Strep A	8–9	General Chemistry/Lipids	46
Diabetes Management	10–15	Coagulation	47
Kidney Function	16	Blood Collection	48–51
Uroflowmetry		Lab Equipment	52
Urinalysis		Lab Supplies	53–55
Fecal Occult Blood	24–29	MOHs Lab Solutions	56
OSOM		Cytology/Histology Supplies	57–59
Women's Health	31–33	Stains & Reagents	
Infectious Disease Testing	34–37	Plated Media	61
Drugs of Abuse Testing	38–39	Clinical Lab Solutions	62–63

**2** 12MS7336

# not everyone is this



Wouldn't it be awesome to be a part of an exclusive group of individuals? Think about wearing the Master's Green Jacket, holding the Stanley Cup, or posing with the Heisman Trophy.

**YOU** can be part of an exclusive group too. The **OSOM®** Influenza A&B Test is for specific customers who understand quality, demand affordability, and appreciate consistency. Are your rapid flu tests **OSOM®**?





For more information call 800 332 1042 or visit us at www.sekisuidiagnostics.com

© 2012 Seksul Diagnostics, U.C. All rights reserved. OSOM\* is a U.S. registered trademark of Seksul Diagnostics, U.C.



# small package

# **BIG REVOLUTION**

It will change the way you view rapid testing...



# New BD Veritor System CLIA-waived for Rapid Detection of Flu A+B



Helping all people live healthy lives


- Proven Performance vs PCR PCR is currently the highest sensitivity standard available
- Reproducibility You Can Count On Unique Adaptive Read Technology helps reduce false-positive results
- Results You Can Trust—
   Easy-to-read digital display eliminates subjective interpretation
- Convenient Workflow—
   Testing requires minimal hands-on time

Medical technology that captures the trusted answer

www.bd.com/ds/veritorsystem

D, BD Logo and Ventor are property of Secton, Dickinson and Company ©2011 BD

BD Diagnostics 7 Loveton Circle Sparks, MD 21152-0999 www.bd.com/ds



A rapid chromatographic immunoassay for use with symptomatic patients. This CLIA-

waived kit is to be used with the BD Veritor

freshly collected, processed and dispensed

nasopharyngeal swab samples that are

System Reader to test nasal and

directly onto an assay test device.



### **RSV** now available for use with Veritor Reader! Call For Details

- Proven performance vs. PCR
- Unique Adaptive Read Technology provides reproducibility you can count on • Easy-to-read digital display eliminates
- subjective interpretation and delivers the "Trusted Answer"
- Convenient workflow with minimal hands-on time

#256055, BD Veritor™ System Reader **(992-0003)**.....ea

#256045, CLIA-waived BD Veritor™ System Flu A+B Kit (992-0002) ......30 tests/pkg

CPT® Codes: Influenza A: 87804QW† / Influenza B: 87804-59†; Using the Objective and Differentiated Reader System: 87999†



# BinaxNOW® Influenza A & B Tests

Rapid test for identification and differentiation of influenza A and B, using specimen collected by nasal washes/aspirates, nasal swabs and nasopharyngeal swabs.

- CLIA waived
- . Differentiates types A and B on the same card
- Easy to use: no additional reagents needed
- On-board extraction
- 15-minute results
- Same specimen can be used with BinaxNOW® RSV test

#416-022 ( <b>153-9781</b> )	22/box
#416-110	
(348-7480)	10/box

CPT® Codes: Influenza A: 87804QW†/Influenza B: 87804-59†





# OSOM® Influenza A & B Tests

In vitro diagnostic test for the qualitative detection of influenza A and B antigens.

- Differentiates between influenza A and B
- Easy to read: 3 distinct lines
- Nasal swab specimens

- Less than 1 minute hands-on time; results in 10 minutes
- Room-temperature storage
- Moderately complex
- 2 extra test sticks included for external QC testing

#190 <b>(106-3126)</b>	.27/box
#191, Influenza A & B Control Kit	
(106-3132)	ea

CPT® Codes: Influenza A and B, 87804-59†



# **QUIDEL**

### QuickVue® Influenza Tests

An immunoassay for detecting influenza types A and B. Three simple steps with results in 10 minutes or less. Color-coded procedure card helps demonstrate each step of test. Kit contains positive and negative controls, and stores at room temperature.

- Detects both A and B types of influenza
- CLIA waived
- · Color-coded procedure card helps demonstrate each step of the test
- No instrumentation is needed
- Results in 10 minutes or less, right in the exam room
- Room-temperature storage with up to a 24-month shelf life from date of manufacture

#20183, A+B Test, Differentiates Type A or Type B (**227-5611**)......25/box #00317, Non-differential (227-9012) ......25/box #20103. Swab Pack\* with 25 Individually Pouched Sterile Foam Swabs (986-3022) ......25/pkg #20171, Transport Tubes\* with 25 Sterile Foam Swabs Inside Plastic Transport Tubes (**227-3550**)......25/box



CPT® Codes: Influenza A, 87804QW†/Influenza B, 87804-59†







### Alere™ Influenza A & B Test

This test provides a simple, rapid method to aid in the differential diagnosis of influenza A & B using nasal swab specimens.

- CLIA Waived
- 10-minute results
- Sensitivity (A/B): 93.8%/77.4%
  Specificity (A/B): 95.8%/98%

#412-000, Influenza A & B Test (512-0016)......25/box

CPT® Codes: Influenza A: 87804QW+; Influenza B: 87804-59+



# The #1 Rapid Flu Test Used in Physician Offices Nationwide\*

# Features of the Alere™ Flu Portfolio

CLIA waived • Dual reimbursement1 • Approved for nasal swabs • Differential diagnosis

# **NEW to the Portfolio!**

Alere™ Influenza A & B Test

- · Results in 10 minutes
- · Simple dipstick format
- · Approved for nasal swab specimens
- Specificity A / B: 95.8% / 98.0%
- Sensitivity A / B: 93.8% / 77.4%

Order# 5120016 (25 Tests)



Actual size

BinaxNOW\* Influenza A & B Test

- · Results in 15 minutes
- · Familiar card format
- Sample types: nasopharyngeal swab, nasal swab and nasal wash/aspirate specimens

Order# 153-9781 (22 Tests)

To learn more about rapid detection and differentiation of Influenza A & B, contact your local Henry Schein Representative.

Actual size

© 2012 Alere. All rights reserved. The Alere Logo, Alere and BinaxNOW are trademarks of the Alere group of companies. PN: 12000512-01 8/12

<sup>\*</sup>Data on file.

Actual reimbursement may vary.



### BinaxNOW® RSV Tests

- CLIA waived
- · Easy to use: no additional reagents needed
- On-board extraction
- 15-minute results
- Built-in control to ensure accuracy
- Same specimen can be used with BinaxNOW Influenza A & B Test

2	•
#430-100. Test	
(301-8899)	10/box
#430-122, Test	
(512-5059)	22/box
CPT Code: 87807QW†	





### Clearview® RSV Test

A 15-minute, four-step test for the qualitative detection of RSV. For use with nasopharyngeal swabs and nasal aspirate specimens (approved claim: in children less than 6 years and adults more than 60).

- CLIA waived
- 93.9% sensitivity; 97.7% specificity
- Room-temperature storage
- Built-in control to ensure accuracy
- Shelf life up to 12 months from date of manufacture

(512-8933) ......20/box

Contains: 20 test devices, 2 mL extraction reagent solution, 1 mL each positive & negative controls, & 40 disposable pipettes.



# BinaxNOW® S. pneumoniae Urinary Antigen Test

An in vitro rapid immunochromatographic assay for the qualitative detection of S. pneumoniae antigen in the urine of patients with pneumonia and in the cerebral spinal fluid (CSF) of patients with meningitis. In conjunction with culture and other methods, it is intended to aid in the diagnosis of both pneumococcal pneumonia and pneumococcal meningitis.

- Dependable accuracy
   No results negatively affected by previous antibiotic use
- Self-contained, swab-based device—no tubes or pipettes
   Easy to use: insert urine- or CSF-saturated swab into test card, add reagent, and read results in 15 minutes
- Both built-in and external positive and negative controls ensure test viability
- Accuracy/sensitivity: urine: 86% sensitivity/94% specificity; CSF: 7% sensitivity/99% specificity

#710-000, 22-Test Kit (512-0000).....ea

CPT® Code: 87449†



# **QUIDEL**®

# QuickVue® RSV Tests

- CLIA waived
- Easy-to-use dipstick format
- 3-step procedure
- 2-color results in 15 minutes
- Room-temperature storage
- Built-in internal controls
- Kit includes external controls

(227-0581)......20/box

CPT Code: 87807QW†





# ICON® SC Strep A Kits

Easy-to-use cassette with on-board extraction.

- CLIA waived
- Easy-to-read, 5-minute results
- 98% accuracy and sensitivity

(331-6584).......25/box

Contains: 25 test cassettes, sterile rayon swabs, vial of reagent A, vial of reagent B, positive control & negative control. CPT® Code: 87880QW†





# ICON® DS Strep A Kits

Dipstick-format test.

- CLIA waived
- Easy-to-read, 5-minute results
- 98% accuracy and sensitivity

(331-1419)......25/kit

Contains: 25 test strips, extraction reagent A, extraction reagent B, positive control, negative control, extraction tubes & throat swabs. CPT® Code: 87880QW†



# **QUIDEL**

# QuickVue In-Line® One-Step **Strep A Tests**

True one-step test with a unique in-line extraction and rapid procedure that allows testing to occur right at the point of care. Two-color results give you clear readability, as well as proven accuracy and dependable results. Kit stores at room temperature and includes positive and negative external controls.

- Built-in controls for CLIA compliance
- Easy to use and easy to read
- 2-color, rapid results in 5 minutes
- Room-temperature storage

**(227-4207)**......25/box CPT® Code: 87880QW†



# **QUIDEL**

# QuickVue® Dipstick Strep A Tests

- CLIA waived
- · 2-color endpoint: results are easy to read and interpret
- 95% overall agreement
- Results in 5 minutes or less
   Room-temperature storage
- Built-in positive and negative controls

(227-0611)......50/box

CPT® Code: 87880QW†



# **QUIDEL**

# QuickVue+® Strep A Tests

- Built-in controls for CLIA compliance
- Easy test procedure
  ± endpoint
- Built-in workstation
- 95% sensitivity; 98% specificity98% overall agreement
- Room-temperature storage
- Moderately complex

#20122

(**489-0075**) ......25/box

CPT® Code: 87880†

8 12MS7336




# OSOM® Ultra Strep A Tests

Ultrasensitive strep A test with the ease of use and rapid results your busy practice needs. Studies show the sensitivity of OSOM Ultra Strep A test and the standard single-swab culture method are not statistically different.

• CLIA waived

- Easy to use
- · Easy-to-read, 2-color results in minutes
- Room-temperature storage
- Includes positive and negative controls
  Includes 2 additional tests for external QC testing

(392-7389) ......25/box #149 (102-7617)......52/box CPT® Code: 87880QW†





### OSOM® Strep A Tests

- CLIA waived
- Easy-to-read 2-color results
- 96% sensitivity; 98% specificity
- Results in 5 minutes
- Room-temperature storage
- Includes positive and negative controls
  Includes 2 additional tests for external
- QC testing

(224-7320) ......52/box CPT® Code: 87880QW†



# ✓ HENRY SCHEIN®

# **ONESTEP+ STREP A TEST STRIPS**

- CLIA waived
- Room-temperature storage
- Canister packaging for ease of multiple testing
- Walk-away procedure

(107-4050)......25/box

CPT® Code: 87880QW†







#1146, Strep Select Agar SSA Plates, 100 mm (437-3477)......10/pkg #1160, TSA 5% Sheep Blood, 100 mm (437-5777)......10/pkg



Selective Strep Culture Plates

**(120-4787)**......10/pkg



SSA Selective Strep AGAR Plate

#221779 (603-4472) ......20/pkg



**Bacitracin Differentiation Disc** 

#2851HL, Disc

(437-0018).....ea



Taxo A Differentiation Disc

#2852, Disc

(146-6565) .....ea



# CAN REAL-TIME RESULTS LEAD TO BETTER OUTCOMES?

# The potential benefits of point-of-care A1C testing

- The ADA recommends that healthcare professionals use A1C point-of-care testing for timely decisions on therapy changes<sup>1</sup>
- Giving patients with diabetes on-the-spot feedback on their A1C number may result in as much as a 1% point reduction in their A1C over time<sup>2</sup>
- A 1% point reduction in A1C may lower the risk of diabetes complications by up to 40%<sup>3\*</sup>



# More efficient than the lab

- · Results in 5 minutes
- · Portable use in multiple exam rooms
- · Small (5µL) blood sample
- . Easy to use minimal training required
- · No maintenance

# Key features

- · Lab accuracy
- National Glycohemoglobin Standardization Program (NGSP) certified
- CLIA waived
- Easy to use minimal training required
- Reimbursable bill to CPT code 83037

# Bayer's A1CNow+® Multi-Test A1C System

Bayer's A1CNow+\* provides you with a fast and easy way to obtain accurate A1C results in your office.

This innovative technology enables you to communicate face-to-face with your patients about their diabetes control.

OFFICE SERVICE DESCRIPTION

DCD3024 A1CNow+ Mult-Test A1C System, 1 kit - 10 bx DCD3021 A1CNow+ Mult-Test A1C System, 1 kit - 20 bx



# Basic features

- Test time 5 Seconds
- Memory 480 results
- Sample size 0.6ul.

# Accuracy

- · Auto-calibration
- No manual coding
- Automatic Hematocrit correction

# Advanced features

- 2 Settings Basic (L1) and Advanced (L2 modes)
- 7-, 14-, and 30-day averages
- · 7-day HI and LO summary
- · Individual HI and LO ranges
- · Pre- and post-meal markers
- post-meal test reminder

# Bayer's CONTOUR® Blood Glucose Monitoring System

For fast, accurate results using No Coding\*\* technology, automatic underfill detection, control marking and temperature correction. PC download option available. For use with CONTOUR\* Blood Glucose Test Strips.

NOER NUMBER

DCD7189 CONTOUR' Blood Glucose Monitoring System, each

### REFERENCES

- 1. American Diabetes Association. Standards of medical care in diabetes 2011. Diabetes Care 2011; 34(1):511-561.
- Ferenczi A, Reddy K, Lorber Dt. Effect of immediate hemoglobin A1c results on treatment decisions in office practice. Endocr Pract. 2001; 7921:86-83.
- Stratton IM, Adler AI, Neil HAW, et al; on behalf of the UK Prospective Diabetes Study Group. Association of glycaemia with macrovascular and microvascular complications of type 2 diabetes (UKPDS 35); prospective ovservational study. BMJ. 2000;321(7258):405-412.

10 12MS7336

<sup>\*</sup> In non-pregnant adults with type 2 diabetes.





### **CONTOUR® Blood Glucose Meter**

Cleared for multiple-site testing, which means more choice and flexibility about where you test: finger, forearm, palm, abdomen, or thigh—it's up to you. However, be sure to talk with your health care professional first to ensure that alternate site testing is right for you.

- No coding requiredAutomatic underfill detection
- · Automatic control marking
- Automatic temperature control
- Small sample size (0.6 μL)
- Quick testing time—results in 15 seconds Saves information for up to 240 tests with time and date

#7151, Meter (147-1405)ea
#7099, Strips (147-9303)50/vial
#7109, Control–Normal, 2.5 mL (147-0874)ea
#06707466, Control–High, 2.5 mL (147-8283)ea
#06707458, Control–Low, 2.5 mL (147-5579)ea





# **BREEZE®2 Meter**

This meter makes testing simple. Its 10-test disc eliminates the struggle with small, individual test strips, and the use of No Coding<sup>™</sup> technology ensures easy accuracy every time you test.

- · Results in just 5 seconds
- Small sample size

- Large, easy-to-read display
   Ergonomic design is easy to handle
   Single-function buttons for easy handling
- Improved 420-test memory and 1-, 7-, 14-, and 30-day averages

Meter Kit (795-0077).....ea

Contains: meter, 1–10-strip cartridge, 1 lancing device, 10 lancets & 1 carrying case. BREEZE® was the first meter to receive the "Ease of Use" commendation from the Arthritis Foundation.



# BREEZE®2 Test Strip Disc Cartridges

- · Less chance of dropping or losing test strips
- Simply load and go

(**795-7966**) ......50/pkg

Contains: 10 strips per cartridge, 5 cartridges



### OneTouch® Ultra® 2 Blood **Glucose Monitor**

- Very small blood sample
- More choices for testing: 3 sites, including finger and arm
- FastDraw<sup>™</sup> design automatically draws blood into test strip
- Fastest test time available: just 5 seconds

(567-0390) .....ea



1101	mai contion s	solution.	
0 =	:0994 <b>7-6191)</b>		25/box
#02	0-244		50/box
	0-245 <b>7-0076)</b>		100/box

• Electropolished needle means a smaller

Meets OSHA guidelines: nonreusable,

maneuvering and easy training

On-demand training through Web and

· Preset gauge for no guessing, no

retractable needle

and reusing

gauge and more comfort for your patients

• Distinctive wing design prevents reloading





• 0.5 µL sample size

- Test results in As Fast As 1-2-3®
- Backlight
- Alternate site testing
- 300-test memory
   14-, 30-, and 90-day averages
   On-screen graphing
- 6 user-settable alarms
- Hypo/Hyperglycemic warning alarms

PRESTO™ PRO System Kit (Includes built-in test strip ejector for safely removing used test strips.) (476-0003).....ea Contains: meter, carrying case & owner's guide.

PRESTO™ System Kit (476-6255) .....ea

Contains: meter, carrying case, owners guide, 1 lancing device, 10 lancets, 10 test strips & 1 bottle of Control solution.

PRESTO<sup>™</sup> Test Strips (476-7212)......50/btl

PRESTO<sup>™</sup> Test Strips (476-2624)......100/btl

12MS7336

# LIFESCAN

### OneTouch® SureSoft™ Lancing Devices

Single-use lancing devices are electropolished, which results in more blood than an acid-washed needle. This means you can use a smaller gauge needle while still ensuring that a sufficient sample will be drawn for point-of-care testing.

CD-ROM programs Gentle, 28-ga Needle (567-0000) ......200/box Regular, 21-ga x 2-mm Needle (115-5267)......200/box

# Choose With Confidence. Choose OneTouch.



12





### TRUTOL® Glucose Tolerance

Screen for hyperglycemia (diabetes mellitus and gestational diabetes) and hypoglycemia (low blood sugar). Volume graduations are imprinted on the label to aid in dispensing. Pediatric dose chart lists precalculated dosage based on body

Orange, 50 g (122-0183)	ea
Fruit Punch, 100 g (332-1707)	24/case
Orange, 75 g (114-7224)	24/case
Fruit Punch, 50 g (917-1244)	24/case
Fruit Punch, 75 g (179-4800)	24/case
Orange, 100 g (122-0184)	ea
Lemon-Lime, 75 g (683-1429)	24/case
Lemon-Lime, 100 g (275-4299)	24/case
Lemon-Lime, 50 g (486-4583)	24/case

### **TRUEbalance Test Kit**

- · Patented TRUEfill design that promotes first-test success and accurate results
- High performance, low aguistion cost and no-coding technology with a small blood sample size and fast test results

H4H0181, Starter Kit <b>(890-0056)</b> ea
H3H01-81, Test Strips
( <b>890-0055</b> )50/box



**RUE**resu





# TRUE2go™ Blood Glucose Meter Kit

- Easy to use: simply twist the TRUE2go on the top of a new vial of strips
- Test results within 4 seconds
- · No coding: eliminates a step that can lead to inaccurate results
- · Convenient: the smallest meter in the world

The world's smallest meter system makes blood glucose testing easy and affordable. Just twist TRUE2go™ on to the top of a new vial of TRUEtest strips and you're ready to go.


Glucose Meter Kit (890-0020)	ea
TRUEtest <sup>™</sup> Test Strips (890-0034)	50/box

# NIZZO DIAGNOSTICS

# TRUEresul

# TRUEresult™ Blood **Glucose Meter Kit**

Proven performance and accuracy at a great value.



- · Results within 4 seconds: means virtually no wait time
- 0.5-mL sample size: less blood means less painful testing

  • 500-test memory: includes time/date and
- 7-day, 14-day, and 30-day averaging
- Affordable: no-coding convenience at an affordable price

Glucose Meter Kit (890-0019)	ea
TRUEtest <sup>™</sup> Test Strips (890-0034)	50/box



Give your patients no-coding meters with clinically proven accuracy!

# TRUEresult® and TRUE2go®

**Blood Glucose Monitoring Systems** · Results in as fast as 4 seconds

- . Tiny, 0.5 microliter blood sample

Both systems utilize TRUEtest Strips featuring Quad-Electrode Laser Accuracy.

# RUEbalance

**TRUE**balance

Blood Glucose Monitoring System

- Accurate results in 10 seconds
- . Small, 1.0 microliter blood sample Utilizes TRUEbalance" Test Strips featuring patented TRUEfill® technology

# To order, contact Henry Schein at 1-800-P-Schein.



Channey Bully of Bood Darma Markering Systems Sestation of the Philippolity, Southeath' (Stra'T, Assessor' (2005)64", and Freshpie Southern' Lite Systems." NES 1/18 C 8011 Nipo Deposito, Inc. St. Engosito, N. St. Special Note and No. 1982 (e.g. NEStature, 1982 (e.g.



# Glucose 201 Analyzer

The HemoCue Glucose Systems are based on a glucose-dehydrogenase method and consist of a small, dedicated analyzer and a unique, disposable microcuvette. They will produce lab-accurate results for diagnosing or monitoring diabetes, as well as monitoring neonatal blood glucose levels.

- More accurate and 68% higher. reimbursement than glucose meters
- FDA-compliant unit when screening and diagnosing
- Analyzer: 0–444 mg/dL (0–24.6 mmol/L)
- Only 5 µL of capillary, venous, or arterial blood is needed
- No inherent limitations related to hematocrit or known drug interferences

  No recalibration needed between
- cuvette lots
- 2-vear warranty

#120706, Analyzer (104-2290)ea
#110706, Microcurette, Glucose (439-9158)100/box
#110402, Microcuvette, Glucose (609-2005)100/case

# **Glucose Aqueous Liquid Control**

Stable for 12 months from date of manufacture and 30 days of opening. Requires refrigeration during storage and shipping. Once opened, it can be stored at room temperature or 2°–8°C.

#139027, 2 mL
<b>(989-8038)</b> ea
Contains: 1 level of ea 2-mL vial: high & low.



Offers the same important benefits as its predecessor. the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
- Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
- Safe and easy to use

Contains: 4 packs of 50.

Reimbursable     CLIA '88 waived     2-year warranty	
Analyzer (345-2344)	ea
Contains: analyzer, AC adapter, manua CD-ROM.	al & training
#3201 Microcuvettes (110-9167)	100/box
Contains: 4 packs of 25. (345-8733)	200/box



# Roche

### Accutrend® Plus Meter Kit

- CLIA waived
- Stores 100 results for each parameter (glucose and cholesterol)
- Out-of-meter dosing
- Onboard QC
- Positive strip lot ID
  Fast, reliable results
- Portable
- Room-temperature storage

Complete Kit #53467544160

(891-1038).....ea

Contains: 1 meter, 4-1.5-V AAA alkaline batteries, 1 user's manual, 1 quick-reference guide, 1 carrying case & 1 warranty card.



Cholesterol Test Strips #5213312160 2.9% low sensitivity. 173-224 mg/dL range.

(891-1036)......25/pkg

Contains: 25 strips per vial & 1 code strip.

Glucose Test Strips #11447475160 (891-0008) ......25/pkg Contains: 25 strips per vial & 1 code strip



Cholesterol Control 2 Level #5219957001

**(891-1037)**.....ea Contains: 1 of ea 4-mL cholesterol control level:



Glucose Control 2 Level #5213231160

(891-1039).....ea Contains: 1 of ea 4-mL glucose control level:

2.5% low & 2.4% high. CPT® Code: 82465QW†





### Sure Comfort® Pen Needles

- Protective tab ensures a safe and sterile needle
- Lubricated needle provides the best possible injection comfort
- Nontoxic, nonpyrogenic, and latex-free
- Compatible with a wide array of pens #24-1210, 31 ga x 3/16" (6 mm)

(120-0912)......100/box #24-1215, 31 ga x 5/16" (8 mm) (120-0913)......100/box

#24-1010, 29 ga x ½" (12 mm) (120-0911)......100/box



# CardioChek™ P•A

Perform on-the-spot testing for cholesterol and lipid levels with 1 convenient, portable, and accurate whole-blood test system.

- CLIA waived
- · Easy to use
- Single-testing and panel-testing capability
- Lipid profile results from a single finger-

stick whole blood sample	
CardioChek P•A (133-2551)	ea
Glucose Test Strips (261-5308)2	25/box
Total Cholesterol Test Strips (412-5021)25	tests
HDL Test Strips (107-4803)25	tests
HDL Cholesterol Controls, High & Lov (261-5338)	v ea
Triglyceride Test Strips (872-5107)25	tests
Lipid Panel Test Strips (261-7537)1	15/box
Multi Chemistry Control (261-8361)	



### Chemcard™ Cholesterol Test

- 3-minute whole-blood screening test
- Color interpretation from 150 mg/dL to 300 mg/dL in increments of 25
- No instrumentation required
- · Room-temperature storage, not to exceed 86°F (30°C).
- (100-7110).....24/box


# Chemcard™ Glucose Test

A simple, effective screening test that is used to identify abnormally high fasting plasma glucose (blood sugar) levels which are often a warning sign of diabetes

- 3-minute whole-blood screening test
- No instrumentation required
- Room temperature storage, not to exceed 86°F (30°C)

#47024, Glucose Test (844-1247) ......24/box


14 12MS7336





This handheld monitor provides accurate in-office A1c results,

allowing physicians to make on-the-spot treatment decisions



 In-office testing; no waiting for lab results

Results in just 5 minutes

- Hands-on procedure time is less than 1 minute
- Simple, 3-step procedure
- No calibration
- No refrigeration necessary if used within 4 months
- No capital equipment required
- Enables A1c testing in every exam room

CPT® Code: #83037QW†



# **SIEMENS**

# DCA™ Vantage Analyzer

Point-of-care immunoassay analyzer for diabetes management. Provides precise and accurate quantitative test results for HbA1c in whole blood during a patient's visit to help improve decision making. With a single urine sample, the analyzer also provides quantitative test results for low levels of albumin, creatinine, and A:C ratio in the physician's office or clinic.

#5075US, Analyzer
(110-0743).....ea

#5035C, HbA1c Reagent Kit
(147-8045)......10/pkg

#5068A, Hemoglobin Like Control Kit

(147-6067)....ea
Contains: 2 vials of ea control: high & low.



# #6011A, Microalbumin/Creatinine Reagent Kits

Detects and monitors microalbuminuria, an early indication of kidney disease. Achieve test results in just 7 minutes from a urine specimen, using the DCA<sup>™</sup> Vantage Analyzer. Existing DCA analyzers may be updated with a new program card.

(147-8917) ......10/pkg

CPT® Codes: 83036QW†/#82043†/#82570†

# CAN BE USED WITH BAYER A1C OR SIEMENS DCA

**Liquid Whole Blood hbA1C Controls** 

#NOD13122-100 4-Vial Kit, 1 mL 2 Level 1 and 2 Level 2 (111-9856).....ea Liquid Whole Blood hbA1C Controls

#NOD13111-100 2-Vial Kit, 1 mL 1 Level 1 and 1 Level 2 (111-9857).....ea

### **SIEMENS**

DCA Microalbumin/Creatinine Control Kit

#6012A

(147-8477)....ea
Contains: 2 vials of ea control (high & low) & reconstitution fluid.





# **SIEMENS**

# Multistix® 10LS Pro Reagent Strips

More reliable tool for detecting kidney disease. Test for glucose, ketone, high protein, low protein, specific gravity, creatinine, pH, blood, nitrite, and leukocytes.

- · Creatinine ratio to correct for varying urine concentrations
- Decreases false positive protein and false negative protein rates
- For over-65 patients and those with diabetes, CVD, or UTI
- Indication of a dilute requires patient to submit specimen, triggering another urine sample
- · Automatic identification of strip configuration for Clinitek® 50 users

(147-8028)......25/btl #1554 (147-2360) .......100/btl ČPT® Codes: 81000†/81001†/81002†/81003†/82570QW†



### **SIEMENS**

#2083

# **CLINITEK® Microalbumin Test Strips**

Test for semiquantitative determination of microalbumin/creatinine in urine. Semiquantitative tests for protein (2-4 mg/dl) and creatinine. Includes an À:C ratio résult.

- CLIA waived on CT status analyzer
- Dip-and-read test for albumin and creatinine in urine
- · Determines an A:C ratio for greatest accuracy
- Use a random urine sample-reduces need for a 24-hour collection
- Timeframe: 60 seconds

(147-4712).....25/btl

CPT® Codes: 82044†/82570†



# **Roche**

# Chemstrip® Micral Test Strips

Test for microalbumin. Fast, easy test to detect subclinical neuropathy.

- 1-minute test
- In office
- Accurate compared to radioimmunoassay (**891-6088**)......30/btl

CPT® Code: #82044QW†



### HemoCue® Albumin 201 System

Can be used for quantitative determination of low levels of albumin in urine for the purpose of screening, diagnosis, monitoring, and supporting clinical evidence in the treatment of microalbuminuria. The system consists of a small, specially designed analyzer and individually packaged microcuvettes containing freeze-dried reagent. The test is performed at the point-of-care, using a microcuvette filled by capillary action. To obtain a quantitative result, fill the microcuvette, place it in the analyzer, and read the result. The system is based on an immunoturbidimetric reaction using antihuman albumin antibodies.

The antigen-antibody complex creates turbidity, which is measured photometrically at 610 nm. The albumin concentration is proportional to the turbidity, and the result is shown on the display in mg/L.

Analyzer Kit with 1 Box of Microcuvettes (111-5593).....ea

Albumin Microcuvettes (111-5591) ......50/box



### OSOM® ImmunoDip® Urinary **Albumin Tests**

Rapid urinary albumin test for microalbuminuria to aid in the early detection of kidney disease.

- CLIA waived
- · Easy to use: simply place test device into urine sample
- Fast: results in 3 minutes
- Accurate: albumin, greater than 18 mg/L indicates positive result
   Specificity 99%, sensitivity 98%
- 18-month shelf life
- Room-temperature storage

#790-25 Test (901-0106) ......25/box CPT Code: 82044QW+

### #6011A, Microalbumin/Creatinine Reagent Kits

Detects and monitors microalbuminuria, an early indication of kidney disease. Achieve test results in just 7 minutes from a urine specimen, using the DCA™ Vantage Analyzer. Existing DCA analyzers may be updated with a new program card.

(147-8917) ......10/pkg CPT® Codes: 83036QW†/#82043†/#82570†

16 12MS7336



# Flo-Med 1100 Precision Uroflowmetry System by





What is the Flo-Med 1100 by Henry Schein?

It is a simple, non-invasive diagnostic screening device used to calculate the flow rate of urine over time. This aids in diagnosis of bladder and sphincter function as well as prostate conditions.

- Easy to use (hygienic & convenient)
- Wireless communication to laptop iPad, or other PDA (hardware capable)
- Compact/portable
- One-touch technology
- Easy calibration
- PDF can easily be stored into an EHR

# (112-5811)

Flo-Med Precision Uroflowmetry System, 1 Unit

# (112-5812)


Flo-Med Cup & Funnel Disposables, 24/Sets

Other accessories such as Laptops and printers are available separately.

†For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

# View Demo Video at:

http://vimeopro.com/user8610595/flo-med-1100-demo





Call Your Sales Consultant Today For More Information!

# Roche

# Urisys 1100™ Urine Analyzer

This CLIA-waived urine analyzer makes urine testing easier. With the simple press of a button, analyzer automatically evaluates urine test strips and prints a patient report. This easy-to-use system is an ideal choice for physician offices looking to standardize test results and improve office efficiency through automated urinalysis testing. Uses Chemstrip® 7, 10 MD, and 5 OB strips.

#3617556, Analyzer (**891-7039**) ......ea





# Roche

# Chemstrip® 10 MD Test Strips

Intended to be used for visual testing or for automated testing on the Chemstrip® 101 urine analyzer and Criterion II urine analyzer and Officion in different analyzer. Measures blood, protein, leukocytes, nitrate, pH, glucose, ketones, bilirubin, urobilinogen, and specific gravity.

(891-4205).....100/btl



Roche

Chemstrip® 5 OB Urine Test Strips #1893467 (891-4294)......100/btl



# **Roche**

# Chemstrip® 10 SG Test Strips

For visual or automated testing on the Chemstrip 101 urine analyzer and Criterion II urine analyzer. Measures leukocytes, nitrite, pH, protein, glucose, ketone, urobilinogen, bilirubin, blood, and specific gravity.

#417145 (891-8002)......100/btl



# **Roche**

# Chemstrip® GP Test Strips

Test for glucose and protein. #200743 (891-5552)......100/btl



# Roche

# Chemstrip® 9 Test Strips

Test for leukocytes, nitrite, pH, protein, glucose, ketone, urobilinogen, bilirubin, and blood.

(891-4564).....100/btl



# Chemstrip® 7 Test Strips

Test for leukocytes, pH, protein, glucose, ketone, nitrite, and blood.

#1008552 (891-1914) ......100/btl

18

# URISCAN 106-9405 1 KETONE STRIP (Kenned) Urine Test Strips 106-9418 1 GLUCOSE STRIP (Glucuse) DOCTORS' CHOICE SINCE 1995 257-5782 2 GP STRIP (Glucose, Promin) CLINICALLY PROVEN QUALITY 106-9420 2 GLUKETO STRIP (Glucose, Ketones) TEST FOR LESS 106-9421 3 GPH STRIP (Glucose, Promin. pH) 4 HemoGPH STRIP (Blood, Glucose, Protein, pH) 106-9425 106-9428 5 HemoKetoGPH STRIP (Blood, Kenones, Glucuse, Promin, pH) 61. STRIP (Blood, Protein, Nitrite, Glucose, Leucocytes, pH) 106-9432 7 STRIP (Blood, Bilirubin, Urobilinogen, Ketones, Protein, Glucose, pH) 106-9434 106-9435 8 STRIP (Blood, Bilirabin, Urobilinogen, Ketones, Protein, Nitrite, Glucose, pH) 10 SGL STRIP (Blood, Bilirubin, Urobilinogen, Kennes, Protein, Nitrite. 106-9436 Glucose, pH, Specific Gravity, Leacocytes)





### Uri-Check® Plus Test Kits

A semiquantitative dip-slide culture system for the screening and early detection of urinary tract infections. More than 100,000 bacteria per mL indicates infection. 1% diameter with patient label affixed to vial. Refer to Uri-Check reference chart to help identify bacterial growth present. 3-5-month shelf life. Nonreturnable. Requires Boekel's Micro Incubator #260700 (100-6303).

#3000UC, CLED Agar/EMB Agar (224-9604)......10/box



# Uri-Check® Test Kits

A semiguantitative dip-slide culture test for the detection of urinary tract infections. Kit includes 1" vials and caps, paddles slightly smaller than 1", and patient labels affixed to the outside of the vial.

#1000UC, CLED Agar/EMB Agar

(224-7373).....10/box

# **Uricult Test Kits**

A simple, self-contained urine test method for the detection of bacteria that commonly cause urinary tract infections. Uricult urine-culture paddles are capable of detecting bacteriuria

concentrations as low as 1000 CFU/mL of urine.

#1000. CLED/EMB

Contains: paddles, patient ID stickers & package insert.



# Accutest® Uriscreen™ Tests

- Rapid urinary-tract infection screen test for bacteriuria and presence of somatic cells
- · Detects true infections with high sensitivity and high NPV
- Positive screening result indicates presence of bacteria and/or somatic cells

(106-9335) ......20/box

CPT® Code: 87899QW†

# NOW, MORE THAN EVER, RELY ON HENRY SCHEIN BRAND PRODUCTS

SAVE WITHOUT COMPROMISING ON EXCELLENCE, WITH OUR OWN HIGH-QUALITY HENRY SCHEIN BRAND

The **ONE STEP** line of diagnostic tests features superior clinical performance, accuracy, and quick results



Health care professionals around the world use Henry Schein Brand Products to fulfill their practice supply needs. There are over 6,000 quality products to choose from —each bearing the Henry Schein Seal of Excellence—your guarantee of satisfaction.



Item#	Description	Size
900-4071	One Step+ H Pylori Test	20/Bx
900-4070	One Step+ Mono Test	20/Bx
900-4074	One Step+ Pregnancy hCG Urine Cassette Test	100/Bx
900-4073	One Step+ Pregnancy hCG Urine Cassette Test	25/Bx
900-4075	One Step+ Pregnancy hCG Combo Test	25/Bx
900-4072	One Step+ Pregnancy hCG Urine Strip Dip	25/Bx
107-4050	One Step+ Strep-A Te <mark>st Dip S</mark> tick	25/Bx





ONESTEP+ URINE ANALYSIS
READER & URINE STRIP COMBINATION

- 10-test urine strip with ascorbic acid monitor pad
- · Printed results in minutes eliminate human error
- Includes specific gravity
- EMR compatible
- Printing order of results can be adjusted
- CLIA waived
  One roll of paper is included (reorder 900-4163)

Analyzer **(900-7017)** .....ea

Contains: 1 OneStep+ urinalysis reader & 12 bottles of Urispec 11-Way urine strips.



# ✓HENRY SCHEIN® URISPEC™ 11-WAY URINE REAGENT TEST STRIPS

Urine Reagent Test Strips

- Unique "read all" feature at 60 seconds
- Similar to Multistix® 10 SG
   Plus Ascorbic Acid Test and all brands of 10 SG products
- Glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes
  • Reader available
- Ascorbic acid pad to ensure uncompromised results
- Room-temperature storage

(102-2285)......100/btl

CPT® Codes: 81000†/81001†/81002†/81003†



Test for glucose and protein, with an ascorbic acid test for accuracy.

- Protein and glucoseIdeal for metabolic, kidney,
- and diabetic monitoring Dip and read in 60 seconds (100-6567) ......100/btl









✓ HENRY SCHEIN® URISPEC™ PRO 2 GP URINE STRIPS Strips for Visual Read Tests for glucose and protein.

(900-4996) ......100/btl



✓ HENRY SCHEIN® URISPEC™ 4-WAY URINE REAGENT STRIPS

- Protein, nitrite, glucose, and leukocytes
  Ideal for OB/GYN
- · Infection, diabetes, and kidney function
- Dip and read in 60 seconds

(107-7157) ......100/btl







HENRY SCHEIN® URISPEC™ 5 OB URINE STRIPS

Strips for Visual Read

Tests for glucose, protein, blood, nitrites, and leukocytes.

(900-4995) ......100/btl





✓ HENRY SCHEIN®

URISPEC™ 10SG URINE STRIPS

Strips for Visual Read Tests for glucose, bilirubin, ketone, urobilinogen, specific gravity, blood, pH, protein, nitrites, and leukocytes.

(900-4994) ......100/btl



Answers for life.

# Siemens Adds Value—

# Every Step of the Way.



# Innovation

Bundled connectivity featuring RAPIDComm® Data Management and Connectivity Solution, DCA Vantage® and CUNITEK Status® Connect Analyzers, minimal preparation for tests, self-contained reagents, walk-away convenience, market leading position, and first to market innovation with connectivity

### Quality

Product exchange warranty, 24/7 call center, technical support, extended service agreements, internal electronic checks and built in error messages in analyzers, QC management and operator management for security, audit trail data



# Training

Proficiency testing, Personalized Education Plan (PEP) web site, telephone installation and training, webinars, patient education material, key opinion leaders, field support, video training, and procedure manuals

# Service

Technical support, reimbursement information, 24/7 call center, customer service specialists, government affairs team, and territory sales representatives in market

Contact your local Henry Schein representative or visit us at: www.usa.siemens.com/diagnostics

Order No. A91DX-POC-110347-UC1-4A00 | 03-2012 | All rights reserved © 2012 Siemens Healthcare Diagnostics Inc.

22 12MS7336

### **SIEMENS**

# **CLINITEK Status+ Urine Analyzer**

Provides convenient point-of-care testing with improved data integrity for more confidence in the results. Simple, fast, and accurate testing produces a result in about 60 seconds. Status analyzer can process basic urine strips, ratio tests, and hCG cassette tests. Full alphanumeric keyboard, easy-to-use touch-screen display, and automatic strip identification make this unit convenient to use. Eliminates need for visual reading, timing errors, and transcription errors. Memory to store 950 patient records. Future-ready for EMR when you are.

#1780

(**795-0055**).....ea



# Albustix® Strips

Test for proteinuria.

(147-1625) .......100/btl

# Combistix® Reagent Strips

Test for pH, glucose, and protein in urine. (147-9259) ......100/btl



# **SIEMENS**

### CLINITEK® Microalbumin Test Strips

Test for semiquantitative determination of microalbumin/creatinine in urine. Semiquantitative tests for protein (2–4 mg/dl) and creatinine. Includes an A:C ratio result.

(147-4712)......25/btl

CPT® Codes: 82044†/82570†

# Multistix® 10 SG Reagent Strips

Test for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes in urine.




(147-5913) ......100/box

CPT® Codes: 81000†/81001†/81002†/81003†

# Multistix® 10LS Pro Reagent Strips

More reliable tool for detecting kidney disease. Test for glucose, ketone, high protein, low protein, specific gravity, creatinine



- urine concentrations
- · Decreases false positive protein and false negative protein rates
- For over-65 patients and those with diabetes, CVD, or UTI
   Indication of a dilute requires patient to
- submit specimen, triggering another urine sample
- Automatic identification of strip configuration for Clinitek® 50 users

#1554

(147-2360) .......100/btl

CPT® Codes: 81000†/81001†/81002†/81003†/82570QW†

# Testing flexibility with simple, one-touch switching between different strip configurations Three sieve functions-confirmatory and

- microscopic-to flag abnormal results for auick review
- Convenient data retrieval for up to 500 patients by patient ID
- 200 control results
- Streamlined operation with immediate start-up
- identification of sample, color, and clarity values
- Flexibility in reporting options, on-screen display, printout and/or uploading to laboratory information system

# Multistix® Reagent Strips

Test for pH, protein, glucose, ketone, bilirubin, blood, and urobilinogen in urine.

(147-5635) .......100/btl



### Multistix® 5 Reagent Strips

- · Includes glucose, protein, blood nitrite, and leukocyte
- Easy, accurate, and reliable results in 60 seconds

#2309

(147-0768) ......100/btl

# Multistix® 9 SG Reagent Strips

Test for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, nitrite, and leukocytes in urine.

#2163

(147-5776) ......100/btl

# Hemastix® Reagent Strips

Test for blood in urine.

(147-1005)......50/btl

### Labstix® Reagent Strips

Test for pH, protein, glucose, ketone, and occult blood.

#2181

(147-6828)......100/btl CPT® Code: Urinalysis #81000†/81002†





# Uristix® 4 Reagent Strips

Test for glucose, protein, nitrite, and leukocytes in urine.

#2166

(147-4912) ......100/btl CPT® Codes: Urinalysis 81000†/81002†

# Uristix® Reagent Strips

Test for glucose and protein in urine.

(147-9199) ......100/btl

Advantus™ Analyzer

With Bar-Code Reader

**SIEMENS CLINITEK®** 

For laboratories or practices requiring medium to high throughput (up to 500 tests per hour). Respond to demands for higher productivity and high quality. Streamline your workflow with flexible operation. Immediate start-up. Automatic calibration, network ready, and offers the widest range of test options. Added QC features.

- Easy-to-use color touch-screen
  Ability to consolidate microscopy results onto one patient record
- QC lockout available or automatic prompting if regimen is not followed



(108-6007).....ea



# Trust Hemoccult ICT-the only

Immunochemical fecal occult blood test clinically tested over an extended duration and proven to work. Our dry-slide process has provided excellent proven results for over 40 years.

Excellent 14-day stability—developed with dry-slide technology, which increases stability of the sample—allowing for multiple-day collection and necessary transit/development time.

Higher sensitivity and specificity—clinically proven to have a relative sensitivity of 100% and apparent specificity of 99.1%.



No drug or dietary restrictions—Hemoccult ICT is specific to intact human hemoglobin—eliminating interactions from red meat, Vitamin C and other factors that lead to false positives.

The right choice—meeting the American Cancer Society and the U.S. Task Force for Colorectal Cancer (CRC) Foundation recommendation of multiple-day testing for CRC.

For information on developing a colorectal cancer screening program in your office, contact your Henry Schein representative today.

Blood Benking

Contribugation information Systems Chemistry Lab Automation

Flow Sytematry

temetry Hematology Molecular Diagnostics Renettanie Rapid Biagenstien



Beckman Coulter, the stylized logs and Hemocoult are trademarks of Beckman Coulter, Inc. and are registered in the USPTO. © Copyright 3011 Beckman Coulter, Inc.

**24** 





# ✓ HENRY SCHEIN®

# ONESTEP™ OCCULT BLOOD SINGLE SLIDES

- Detects hemoglobin through chemical reaction
- Compares to high-sensitivity products
- Room-temperature storage
- 100-2-hole slides and 100 sticks

(147-1246) ......100/box



# ✓ HENRY SCHEIN® ONESTEP™ OCCULT BLOOD PATIENT MAILERS

Effective and rapid method for detecting occult blood in stool. Aids in diagnosis of asymptomatic gastrointestinal problems such as colorectal cancer, ulcers, polyps, anemia, and diverticulosis.

- Detects hemoglobin through chemical reaction
- Compares to high-sensitivity products
- Room-temperature storage

(147-7182) ......50 kits/box

Contains: 50 triple 2-hole slides, 150 sticks, developer & 50 US postal mailers.



# Hemoccult II<sup>®</sup> Dispensapak™

100-Patient Kit

(254-5061) .....ea

Contains: 100 patient kits, 6-15 mL bottles developer & 1 product instruction.

Kit contains: 1 Hemoccult II triple-test card, 3 applicators & 1 illustrated diet/test instructions. Barrier mailing pouch may be purchased separately.



# Hemoccult® Single-Test Card

#60151A

(254-9052) .....ea

Contains: 100 Hemoccult single-test cards, 2–15-mL bottles of developer, 100 applicators & 1 product instruction manual.

Hemoccult® Tape Dispenser Kits

(254-8235) ......2/box Contains: Hemoccult tape dispenser (100 tests/dispenser), 1–15-mL bottle developer & product instructions.

Hemoccult® Barrier Mailing Pouches

Designed to meet US Postal Service regulations.

(254-6204)......100/box

# Hemoccult® Developer

Most clinically proven and widely used FOBT. Trusted for over 30 years.

15-mL Bottle (**254-2360**).....ea



# Hemoccult II<sup>®</sup> SENSA<sup>®</sup> Dispensapak<sup>™</sup> Plus

The fecal occult blood test that provides improved readability versus Hemoccult and other standard FOBT products.

- · Easier readability
- Higher clinical sensitivity
- Detects more true positives

# #64130A, 40-Patient Kit (254-4454) .....ea

Contains: 40 patient kits, 2-15-mL bottles developer & 1 product instruction manual.
Patient kit: 1 Hemoccult II® SENSA triple-test card, 3 flushable sample collection tissues, 3 applicators, 1 illustrated diet & test instructions, 1 patient pamphlet & 1 Barrier® mailing pouch.

#61200A, Hemoccult II Slides (224-9819)......102/box



# Hemoccult® SENSA® Single-Test Card

#64151A (254-3896) .....ea

Contains: 100 Hemoccult SENSA single-test cards, 2-15-mL bottles developer, 100 applicators & 1 product instruction manual.

Hemoccult® SENSA® Developer (254-8308) 15-mL....ea





# Hemoccult II<sup>®</sup> Dispensapak™ Plus

Most clinically proven and widely used FOBT. Trusted for over 30 years.

#61130A, 40-Patient Kit **(254-5748)** .....ea

Contains: 40 patient kits, 2–15-mL bottles developer & 1 product instruction manual. Kit contains: 1 Hemoccult II triple-test card, 3 flushable sample collection tissues, 3 applicators, 1 illustrated diet & test instructions, 1 patient pamphlet & 1 Barrier mailing pouch.



BECKMAN COULTER...

# Hemoccult II<sup>®</sup> SENSA<sup>®</sup> Elite Dispensapak<sup>™</sup>

- Easier diet and patient instructions
- Easier sample collection
- Easier processing (331-2723)......40/box

# Offer your patients a colorectal cancer (CRC) screening test they may actually use.



# Trust the brush and help boost patient compliance with InSure® FIT™

The InSure FIT Fecal Immunochemical Test with its unique water-based Blue Brush collection method makes it easier to screen for CRC. No fecal handling or test card smears required; instead, InSure FIT uses a gentle brushing of the stool in water. The 2 tests/one card format follows

American Cancer Society recommendations for CRC screening along with high sensitivity, no dietary or medicinal restrictions and a patient-friendly collection method for improved patient compliance,

Please visit **DiagnosisCAST.com** and click on ColonCancerWRAP for Information on new advances in colon cancer screening (CME credits available).



That's Easier to Live With

For more information, please visit www.insuretest.com.

Manufactured by Enterix Inc., Edison, NJ 06837 USA, a Quest Diagnostics Company.

ENTBA9013 06/2012







### InSure® FIT™ Fecal Immunochemical Test (FIT)

Test provides high sensitivity and ease of use with a simple, water-based collection method designed to be more patient-friendly. The patented brush sampling replaces unpleasant fecal collection or smearing. No dietary or medicinal restrictions, enhanced specificity and sensitivity for more productive colonoscopies. The InSure® FIT™ is based on the immunochemical detection of human hemoglobin (Hb) as an indicator of blood in the stool. The assay uses a monoclonal antibody to capture Hb on a test strip during test development. Since globin does not survive passage through the upper gastrointestinal (GI) tract, the presence of globin in the stool indicates bleeding in the colon or rectum. Clinical studies have demonstrated high sensitivity and specificity of the InSure® FIT™ test in colorectal cancer screening.

- 87% sensitivity for colorectal cancer\*
- 98% specificity
- 66% increase in patient compliance vs. a guaiac-based fecal occult blood test (gFOBT) was found in one study\*
- No handling/smearing of dry, solid fecal materials\*

CPT® Code: 82274QW† Screening G0328QW†



#90030, Combo Kit 25/25 **(450-0009)**.....ea

Contains: 25 test strips, 1–9-mL bottle of run buffer & 25 patient collection kits.



#90010, Patient Collection, 10-Pack (450-0003).....ea

Contains: Each envelope contains: 1–2-day sample card, 2 brushes & specimen collection instructions.



#16800, FOBT Controls (450-0010) .....ea

Contains: 2–1.5-mL bottles, 1 of ea: positive & negative.



#90035, Combo Kit 25/35 (**450-0008**).....ea

Contains: 25 test strips, 1–9-mL bottle of run buffer & 35 patient collection kits.



#90025, Patient Collection, 25-Pack **(450-0002)**.....ea

Contains: Each envelope contains: 1–2-day sample card, 2 brushes & specimen collection instructions.



#80025, Developer Kit

(450-0011) .....ea Contains: 25 test strips & 1–9-mL bottle of run buffer.

#11018, QC Pack, 50 Test Cards (450-0012) .....ea
Contains: 50-2-day sample test cards.

CPT® Code: 82274QW†



# HEMOSURE® IFOB TEST



# Colon Cancer Kills. Are You Screening Your Patients?



Colon Cancer is the #2 leading cause of Cancer related deaths in the US. Call New for more information.

Contact your Henry Schein representative or Hemosure Inc. at 888-435-6787 (Hemosure) or visit our website at www.hemosure.com

# Hemosure is the #1 choice for CRC screening

# Clinically Proven

- \* Greatly Reduces False Positives & False Negatives.
- \* No Drug or Diet Restrictions.
- Higher Sensitivity & Specificity.

# Cost Effective

- \* Higher Reimbursement (ROI).
- Actual Test Stays in Office or Laboratory.
- . One Test Required.

# Support

- \* Free Samples / Evaluation Kit.
- \* Inside & Outside Sales Support.
- www.HEMOSURE.com



The ACG supports the joint guideline recon endation that older quaice-based fecal occult blood tests be abandoned as a method for CRC screening... the ACG recommends the FIT as the preferred concer detection fest."

and by permission from Managhar Addition (Ad. The American may (2020) FED22950, maynyn (200)

# HEMOSURE\* iFOB TEST Immunochemical Fecal Occult Blood Test (FIT)



# Hemosure® iFOB Tests

- Best-value package (792-8076): 40% more collection tubes and home kits
- 5X higher reimbursement (ROI)
- · No dietary or drug restrictions
- Specific to human hemoglobin
- Higher patient compliance, due to ease of use
- Easy-to-read results, similar to pregnancy test
- User-friendly and more sanitary
- Complete with home kit
- CLIA waived

#T1-CK10, Complete with Home Kit (792-5630)10/box		#T1-TC01, Control (108-1373)ea CPT® Codes: 82274QW†/G0328W†
#T1-CK30, Complete with Home Kit (792-8076)30/box	#T1-CK50, Complete Kit (109-3145)50/box	

28 12MS7336



### Hemoccult® **ICT Tests**

- CLIA waived; easy test procedure Screening reimbursement code: G0328QW†
- #395065A, Collection Card Set (834-7610) .....ea

Contains: 100 cards, 100 applicator sticks, physician instructions & sample collection instructions

#395066A, Screening Kit (834-3249) .....ea

Contains: 1 box of 40 patient screening kits. Each patient screening kit contains: 1 dispensing envelope with patient instructions, 3 Hemoccult ICT collection cards, 3 flushable collection tissues, 3 applicator sticks & 1 mailing pouch.

#395067A, Test Kit (834-6389) .....ea

Contains: 20 test devices & 1-1-mL bottle Hemoccult ICT buffer (to be used with Hemoccult ICT collection cards & patient screening kits).

#395068A, Controls

(834-5620) ......4/box

#395261A, 2-Day Patient Screening Kit (237-0085)......50/box

Contains: 1 dispensing envelope, 2 Hemoccult ICT collection cards, 2 flushable collection tissues, 2 applicator sticks & 1 mailing pouch.







# **QUIDEL**

# QuickVue® iFOB Test

- Results in 5 to 10 minutes
- Analytical sensitivity: 50 µg hHb/mL buffer or 50 µg hHb/g fecés
- Specificity: specific to human hemoglobin
  Agreement is 98% when using quality
- control specimens
- Specimen collection tube doubles as convenient, sanitary dispenser
- Self-contained test cassette
- Bold burgundy end points for ease in reading and interpreting results

Complete Test Kit (950-2866)	<ul> <li>Hoom-temperature storage with 24-r shelf life from date of manufacture</li> <li>External controls available separately</li> </ul>	
(227-2458)		20/box
(223-4668)ea	Tray Pack (227-2458)	50/box
	(223-4668)	ea



### ✓HENRY SCHEIN®

### **ONESTEP+ IFOBT OCCULT BLOOD** MONOCLONAL ANTIBODY KITS

Innovative kit designed to avoid increased costs caused by noncompliant patients.

- Patient takes home stick and card only (instead of taking home buffer tubes)
- 30-day sample stability (compare to
- 4–8-day stability)

   Monoclonal; no dietary restrictions
- Clean, in-office sample handling
- Multiple-site accuracy

Cassettes (900-7458).....30/kit Contains: 30 cassettes & buffer tubes. 25-Test Complete Kit (900-4413) ......25/kit Contains: mailing cards kit, cassettes & buffer tubes. 20-Patient Mailer Kits (900-4416) ......20/kit

Contains: sticks, card mailer & instructions. All-inclusive Test Kit
CPT® Codes: 82274QW†/G0328QW†

# NOW, MORE THAN EVER, RELY ON HENRY SCHEIN BRAND PRODUCTS

SAVE WITHOUT COMPROMISING ON EXCELLENCE, WITH OUR OWN HIGH-QUALITY HENRY SCHEIN BRAND



# ✓ HENRY SCHEIN®

### ONESTEP+ IFOBT OCCULT BLOOD MONOCLONAL ANTIBODY KITS

Innovative kit designed to avoid increased costs caused by noncompliant patients.

- Patient takes home stick and card only
- (instead of taking home buffer tubes) 30-day sample stability (compare to 4-8-day stability)
- Monoclonal; no dietary restrictions
   Clean, in-office sample handling
- Multiple-site accuracy

Cassettes (900-7458)......30/kit Contains: 30 cassettes & buffer tubes.

25-Test Complete Kit (900-4413) ......25/kit

Contains: mailing cards kit, cassettes & buffer tubes.

20-Patient Mailer Kits (900-4416) ......20/kit

Contains: sticks, card mailer & instructions. All-inclusive Test Kit CPT® Codes: 82274QW†/G0328QW†

BRAND PRODUCTS

# your rapid tests should be





Sekisui Diagnostic's awesome family of rapid tests accommodates your patient's specific needs.

Accurate, with easy to use, easy to read results within minutes, backed by Sekisui Diagnostic's reputation as a world leader in diagnostics.

Anyway you say it, it's...OSOTT

For more information please call 800 332 1042 or visit our website at www.sekisuidiagnostics.com

OSOM\* Mono Test • OSOM\* Strep A Test • OSOM\* Influenza A & B Test • OSOM\* hCG Combo Test • OSOM\* BYBUS\* Test • OSOM\* hCG Urine Test • OSOM\* ImmunoDip\* Urinory Albumin Test • OSOM\* Tilchomonas Ropid Test OSOM\* H. Pylori Test • OSOM\* Card Pregnancy Test • OSOM\* Uttra Strep A Test • OSOM\* C. difficile Toxin A/B Test

© 2011 Sekinul Diognostics, LLC. All rights reserved.

CSCM\* is a regidered U.E. trademork of Seksul Diagnostics, LLC.

30 12MS7336





# ✓HENRY SCHEIN®

# **ONESTEP+ PREGNANCY STRIP TEST**

- Urine test
- CLIA waived
- 25 mIU/mL detection level
- Room-temperature storage
- Read time 3 minutes or less

(900-4072) ......25/box

CPT® Code: Urine 81025QW†





### ONE STEP+ PREGNANCY COMBO CASSETTE TEST

- Urine or serum test
- CLIA waived for urine only
  25 mIU/mL detection level
- Room-temperature storage

(900-4075) ......25/box

CPT® Codes: 81025QW† (Urine) / 84703† (Serum)





# HENRY SCHEIN® ONESTEP+ PREGNANCY CASSETTE TESTS

- Urine test
- CLIA waived25 mIU/mL detection level
- Room-temperature storage

25-Patient Kit

(900-4073).....ea

100-Patient Kit

**(900-4074)** .....ea CPT® Code: Urine 81025QW†



### OSOM® hCG Urine Test

Avoid the mess and hassle of manipulating urine samples while providing the flexibility and reliability your patients deserve.

- CLIA waived
- Dipstick: simple, 1-step procedure
  Fast results: as soon as 1 minute
- 25 mIU/mL sensitivity
- 2-color results
- Room-temperature storage
- More than 99% accurate in clinical studies

(236-2470) ......50/box

CPT® Code: 81025QW†



### OSOM® hCG Combo Test

Urine/Serum Pregnancy Test

- · CLIA waived for urine only
- Easy to read: black-and-white results
- Fast: one step; no reagents
  Sensitive: 10 mIU/mL serum,
- 20 mIU/mL urine
- Room-temperature storage; 24-month shelf life
- Easy to run

#124, Combo Test

benefit of easy-to-read, two-color results. Add 3 drops of urine and read result. Results in as soon as 1 minute
 25 mIU/mL sensitivity

OSOM® Card hCG

Pregnancy Test
Traditional physician-office

hCG test with the added

- More than 99% accurate in clinical studies
- Easy to run: no-drop counting pipette
   CLIA waived

#102. Test

(**327-0743**)......25/box

CPT® Code: 81025QW†



### Alere<sup>™</sup> hCG Combo Cassette Test (20/10 mIU/mL)

The hCG Combo Cassette test provides superior performance for the detection of human chorionic gonadotropin (hCG), a hormone produced during pregnancy.

- Cut-off level: 20 mIU/mL in urine; 10 mIU/mL in serum
- Sensitivity: > 99%; Specificity: > 99%
  Single-step test provides positive results in as little as 3 minutes for urine and 5
- minutes for serum
- CLIA-waived (for urine only),
  CPT® Code: Urine hCG 81025QW†;
  Serum hCG 84703† (Combo tests are classified as moderate when testing serum samples.)

#92215, hCG Combo Cassette Test (512-0014)......40/box



# **Alere**

# Alere™ hCG Cassette Test (25mIU/mL)

A one-step test for the early detection of human chorionic gonadotropin (hCG), a hormone produced during pregnancy. The test uses a patented, innovative technology to provide highly sensitive and specific results for pregnancy.

- Cut-off level: 25 mlU/mL
- Sensitivity: > 99%; Specificity: > 99%
  Single-step test provides positive results in as little as 3 minutes
- CLIA-waived, easy-to-use lateral flow urine test
- Shelf life: 24 months from date of manufacture
- CPT® Code: 80125QW†

#92217, hCG Cassette Test

(512-0013)......40/box



# Alere™ hCG Dipstick Test

The Alere™ hCG Dipstick test provides qualitative determination of human chorionic gonadotropin (hCG) in urine for early detection of pregnancy. A simple procedure provides easy-to-read results in 3 minutes at a detection level of 25 mIU/mL.

- Sensitivity: > 99%; Specificity: > 99%Cut-off level: 25 mlU/mL
- · Single-step test provides positive results in as little as 3 minutes
- CLIA-waived, easy-to-use lateral flow test
  Shelf life: 24 months from date of manufacture
- CPT® code: 81025QW†

#92211, hCG Dipstick Test (512-0012).....50/box

To Order: 1.800.P.SCHEIN (1.800.772.4346) 8am-9pm, et • To Fax: 1.800.329.9109 24 Hrs 12MS7336 www.henryschein.com/medical



BECKMAN COUITER.

# ICON® 25 Hcg Pregnancy Test

A one-step, rapid, combination serum/urine test designed for the early detection of pregnancy.

- 25 mIU/mL sensitivity, >99% specificity
- · 3-minute results for urine, 5-minute results
- Requires just 3 drops of specimen
  CLIA waived for urine

#43025A (254-2726) ......25/box



# 

# QuickVue® One-Step hCG Urine Test

- Results in 3 minutes or less
- 25 mIU/mL sensitivity
- 99% sensitivity, specificity, and accuracy
- · Easy to use
- Built-in controls for CLIA-waived compliance
- CLIA waived (urine only)

(227-6481).....25/box CPT® Code: 81025QW†



# OSOM® BVBLUE® Test

For Gardnerella vaginalis

The only CLIA-waived rapid test for detecting bacterial vaginosis, including gardnerella, bacteroides, prevotella, and mobilincus.

- Instant color change provides clear, easy-to-read results
- Sensitivity: 92.8% vs. Gram stain
- More sensitive and specific than wet mount
- One-minute hands-on time
- CLIA waived

#183, Test <b>(414-4170)</b>	25/box
#184, Control	
(327-4117)	ea
CPT® Code: 87905QW†	



# ICON® 20 hCG

A rapid test for detecting hCG qualitatively. #395097A

(331-0431).....25/box



# **QUIDEL**

# QuickVue+® One-Step hCG Combo Tests

- 20 mIU sensitivity in urine, 10 mIU sensitivity in serúm
- Results in 5 minutes for serum, 3 minutes for urine (CLIA waived urine only)

- 99% sensitivity, 99% specificity
   Room-temperature storage
   Built-in controls for CLIA compliance
- One-step ± endpoint

(371-6739).....30/box (371-1811) ......90/box CPT Codes: Urine, 81025†/Serum, 84703† (Combo tests

are classified as moderate when testing serum samples.)



# OSOM® Trichomonas Rapid Test

The only CLIA-waived rapid test for Trichomonas infection. Detects antigen directly from vaginal swab.

- More sensitive than wet mountEasy-to-run dipstick
- Easy-to-read 2-color results
- · Results in 10 minutes or less
- Included in CDC and ACOG guidelines
- Includes QC controls

#181, Test (406-7631) ......25/box #182 Control (611-1157) .....ea CPT® Code: 87808QW†



### Impact RPR® Tests

For Syphilis

Rapid detection of syphilis in serum or plasma.

- Simple and accurate
- 100% agreement with reference FTA testing study
- No incubation
- · Read in 8 minutes after rotation on slide rotator

Test <b>(512-1386)</b>	50/box
Dry Control Card (512-4000)	10/box
Liquid Control (512-3595)	3/box



# **QUIDEL**

### QuickVue® One-Step hCG Combo Tests and Control Sets

- · Results in 5 minutes (serum), 3 minutes (urine)
- 25 mIU/mL sensitivity (urine), 10 mIU 99% sensitivity (serum)
  99% sensitivity, 99% specificity
  Built-in controls for CLIA compliance

- CLIA waived (urine only)

#20110 Combo Test (227-7244)	50/box
#00272, Control Set (Urine) (227-4849)	ea

#00281, Control Set (Serum) (**227-5659**) .....ea CPT® Codes: Urine hCG, 81025†/Serum hCG, 84703†

(Combo tests are classified as moderate when testing serum samples.)



# VS-SENSE<sup>™</sup> Test

A rapid acidity test for the evaluation of rapid activity test to the evaluation of vaginitis (bacterial vaginosis and trichomonas). This 3-second, user-friendly test accurately detects mixed vaginal infections, which allows for immediate treatment, resulting in healthcare provider and patient satisfaction.

- Easy to use and operate
- Immediate results
- No need for a color scale
- No need to send samples to a lab and wait for results
- Sensitivity: 86.32%Specificity: 93.43%
- 3-5 second timeframe
- Room-temperature storage

#FPHSCSVS25 (411-0001).....25 swabs/box CPT Code: 83986QW†

32 12MS7336

# WOMEN'S HEALTH HCG AND SPECIALTY



# pHizatest™ Nitrazine Indicator Paper

Compare to NitraTest™

- Economical
- · Instant, easy-to-read results
- Reliable: manufacturer has made nitrazine-indicator paper for more than 50 years

(122-0181) .....ea

CPT® Code 83986†



### Affirm™ VPIII Kit

BD

- Accuracy up to 30% higher than slide microscopy
  Less than 2 minutes hands-on
- Less than 2 minutes hands-on time required
- Easy-to-read results
- Built-in procedural controls for each sample

(987-7212).....24/box

# **CYTOLOGY BRUSHES**

# Histobrush®

Cytology Brush (100-8560) ......100/bag



# Cytology Brushs

#60317, 2.3 x 240X 2.0 (889-0701).....10/case

# **Cytology Wood Scraper**

Non-sterile

#2173

(113-5208)......5000/case

# MICROSCOPE SLIDES & COVERS



# ✓ HENRY SCHEIN® GLASS MICROSCOPE SLIDES

22 mm x 22 mm; thickness: 0.17 mm–0.25 mm. 1-oz package.

Plain Slide (102-7277)...

(102-7277)......72/box

Frosted End Slide

(102-5605)......72/box Cover Glass #2, 1-oz Package

(102-6803) .....ea

# **CYTOLOGY FIXATIVE**



# Health Link

# Histology/Cytology Solutions and Stains

#400631, Cytology Fixative, 4-oz Pump Bottle (437-7476).....ea

#400600, Eosin Y Stain 1% (Alcoholic), 16-oz Bottle (113-4902).....ea

#400604, Formalin 10% Buffered, 1-gal Bottle (108-2287).....ea

#400605, Formalin 10% Buffered, 5-gal Bottle (108-2289).....ea

Spray-Cyte® Cytological Fixative

1½ oz

(168-9697) .....ea

Aerocell™

3½-oz Cytology Fixative Spray (100-0370).....ea

# BD Affirm™ VPIII

Microbial Identification System



# The Problem:

Mixed vaginal infections.

One of every four vaginal infections are caused by two or more organisms: Candida species, Gardnerella vaginalis or Trichomonas vaginalis.

Failure to diagnose these common causes of vaginitis can carry significant risk for your patients. Trichomoniasis and bacterial vaginosis (BV) are associated with increased risk of HIV infection, delivering preterm, low birth weight babies and, if untreated, PID.<sup>12</sup>

Sichermouse, CDC, Ox. STD. Fact Shares Sept 2010. Becterial visplossis, CDC, Div. STD. Eact Shares Sept 2000. Data from package maint.

80, 50 Jugo and Afforts are trademarks of Section, Ordering

# The Solution:

1 molecular test that detects & differentiates 3 organisms.

The BD Affirm" VPIII system quickly identifies and differentiates three pathogens from one single test. The Affirm VPIII technology provides RNA accuracy to enable targeted therapies.

- Results in < 45 minutes</li>
- Enables proper therapy, the first visit
- Unaffected by difficult specimens: menses, sperm or creams from self medication.<sup>3</sup>

Learn more about Affirm VPIII testing in your laboratory, call 800.553.3470.

BD Diagnostics 800 638.8663 www.bd.com/ds



# STDS—HERPES



### biokitHSV-2 Rapid Test

For Herpes Simplex Virus 2

A type-specific glycoprotein G (gG-2) immunoassay kit—a leading aid in the diagnosis of genital herpes. Two-reagent

- Results in less than 10 minutes
- Finger-stick blood or serum specimens
- Built-in procedural codes

#3000-28001

(630-3080) ......20/box

CPT Code: 86696†

# **SYPHILIS**

# **RPR Card Tests** for Syphilis

RPR (rapid plasma reagin) macroscopic nontreponemal

flocculation test kits for the qualitative and semiquantitative detection of reagin antibodies in human serum. The macroscopic RPR carbon antigen enhances the visual discrimination between reactive and nonreactive results. Kits have a shelf life of up to 24 months from date of manufacture.

- Easy-to-interpret results available in 8 minutes
- Includes controls

Syphilis Test Cards	
(653-4966)	100/box
(117-8951)	
(115-5717)	5000/box

Contains: Each kit contains RPR carbon antigen, reactive control, weak reactive control, nonreactive control, a dropping bottle, a 20-g dispensing needle, 10- or 30-well test cards & disposable stirrer pipettes.

# **Alere**

# Impact RPR® Tests

For Syphilis

Rapid detection of syphilis in serum or plasma.

Test (512-1386)	.50/box
Dry Control Card (512-4000)	.10/box
Liquid Control (512-3595)	

# C. DIFFICILE



### OSOM® C. Difficile Toxin A/B Test

- Detects both toxins: A & B
- · Excellent sensitivity and specificity
- Just 4 simple steps
   Includes 2 additional test sticks for external QC testing
   Less than 1 minute hands-on time; results
- in 20 minutes or less

Difficile Toxin A/B Test (309-0085) ......27/box CPT Code: 87324†

# **CHLAMYDIA**



# **QUIDEL**

# QuickVue® Chlamvdia Test

- Results in 12 minutes or less
- 92% sensitivity and 98.6% specificity
  Simple extraction procedure requires only 2 reagents
- Everything needed is in the kit
- Room-temperature storage
- Built-in controls, external controls are included

#B006

(227-0527).....25/box

ČPT® Code: 87810†

# HIV



# Uni-Gold™ Recombigen® HIV-1 Test

Single-use rapid immunoassay for qualitative detection of antibodies to HIV-1 in serum, plasma, and whole blood (venipuncture and finger stick). Intended for ùse in point-of-care settings as an aid in diagnosis of infection with HIV-1.

CLIA Waived

**Alere** 

- Results in 10 minutes or less
- 100% sensitivity

#1206506, Test (104-3399)......20/box #1206530, Test Control (174-3798) .....ea CPT® Code: 86703QW†



# Clearview® Complete HIV 1/2 Test

A unique 3-in-1, self-contained barrel format that handles the collection, processing, and analysis of whole blood, serum, or plasma for detecting HIV-1 and HIV-2 antibodies.

- CLIA waived for fingerstick and venipuncture whole blood
- Includes a bandage and lancet with each device
- 15 minutes to results with 99.7% sensitivity and 99.9% specificity
- 24-month shelf life from date of manufacture

#92111 CPT® Code: 86703QW+

(512-1454)......25/box



### Clearview® HIV 1/2 Stat-Pak®

A simple cassette-based format for the detection of antibodies to HIV-1 and HIV-2 in human whole blood, serum, and plasma.

- CLIA-waived for fingerstick and venipuncture whole blood
- 15 minutes to results with 99.7% sensitivity and 99.9% specificity
- 24-month shelf life from date of manufacture

#92110 (512-9386) ......20/box CPT Code: 86703QW†

34

# **HEPATITIS C**



# OraQuick™ HCV Rapid Antibody Test Kit

- The first and only FDA-approved CLIA-waived HCV rapid antibody test
- Lab accuracy (greater than 98% accurate)
  Tests for multiple HCV genotypes
- Easy to use, simple platform
- Rapid results in 20 minutes
- Fingerstick and venipuncture whole blood collection

#1001-0181	
(810-0012)	25/kit
#1001-0180	
(810-0013)	100/kit
#1001-0182, Oraquick HCV Rapid C	
(810-0011)	ea
Contains: 1 bottle of ea: positive & negat	ive.

# HIV



- Rapid human immunodeficiency virus (HIV) types 1 and 2 (antibody test is a single-use, qualitative immunoassay)
- Detects antibodies to HIV
- Only rapid HIV test for use with specimens for oral fluids, venipuncture whole blood, fingerstick whole blood, and plasma
- CLIA waived for specimens for oral fluids, venipuncture whole blood, and fingerstick; whole blood only
- Simple, 2-step, oral-fluid procedure
  Reliable results in just 20 minutes
- Greater than 99% accuracy

#1001-0079, Test (810-3221)	25/box
#1001-0078, Test (810-6864)	100/box
Control (810-3390)	ea







# SPECIAL HENRY SCHEIN BRAND PRODUCT OFFERS FOR PRIVILEGES MEMBERS ON FACEBOOK





www.facebook.com/hsmprivileges
1.866.MED.VIPS

Not a member? Join today. Visit us on Facebook for details.



### **MONONUCLEOSIS**





✓HENRY SCHEIN®

### **ONESTEP+ MONO CASSETTE TEST**

A rapid chromatographic-immunoassay for the qualitative detection of infectious mononucleosis in whole blood, serum, or plasma.

- CLIA waived
- Finger stick or collection tube
- Room-temperature storage
- Includes controls

Complete Kit

(900-4070) ......20/box

Contains: 20 pouched test devices with disposable sample droppers; disposable, heparinized capillary tubes with dispensing bulb; 1 of each control (positive & negative); 1 sample buffer; 1 procedure card & 1 package insert.

CPT® Code: 86308QW†





NO AGE RESTRICTIONS!

### ✓ HENRY SCHEIN®

### ONE STEP PRO+ MONONUCLEOSIS TEST KIT

- · CLIA waived for whole blood
- · No age restrictions
- Measures heterophile antibodies in whole blood, serum and plasma
- · Uses best technology of bovine derivatives, which eliminate cross reactivity
- Room temperature or refrigerated storage
- Two internal controls exist to avoid additional cost and ensure good laboratory practice

  Only 25 µL blood needed

Mono Test Kit

(112-5547) ......20/box

Contains: 20 test devices, 1 developer solution, 1 negative control, 1 positive control, 40 sample transfer pipettes (20 for serum/plasma & 20 for whole blood), 1 procedure card & 1 package insert.

CPT® Codes: 86308QW† (Whole Blood) / 86308† (Serum/Plasma)



### **QUIDEL**®

### QuickVue+® Infectious Mononucleosis Test Kit

- IgM-specific test
- Results in 5 minutes
- Easy to read +/- test result
- 99% sensitivity and 99% specificity
  CLIA waived for whole blood
- Flexible specimen-collection options

(227-0024).....20/box

Contains: 20 test cassettes, 1 developer, 1 of ea control (negative & positive), 20 sample pipettes & 20 capillary tubes. CPT® Codes: Whole blood: 86308QW†; Serum: 86308†





### **ICON®** Mono Test Kit

Approved for testing on all ages. Combines fingerstick whole-blood convenience with EIA membrane technology for fast, accurate results. Uses direct solid-phase immunoassay technology for early qualitative detection of IM heterophile antibodies in whole blood, serum, or plasma. Serum/plasma is moderately complex.

- Whole-blood-test waived
- 99.9% relative sensitivity and 98.8% relative specificity of serum/plasma
- 98.3% relative specificity of wholeblood assay

(111-7123) ......30/box

Contains: 30 test cassettes, 1 developer solution, 1 of each control (negative & positive), 1 bulb & capillary tubes

CPT® Code: 86308†



### Acceava® Mono Cassette

The high-performance mononucleosis diagnostic test for whole-blood, serum, or plasma samples.

- CLIA waived for whole-blood procedure only
- FDA cleared with no age limits
- Positive results in 8 minutes
- Superior sensitivity at 99.9% for earlier IM detection
- Specificity of 98.8%
- External positive and negative controls provided with your kit-no additional purchase is necessary

#92407

(512-0004) ......25/pkg





Simple in-office test for the detection of mononucleosis.

- Easy to use
- Two-color results in 5 minutes
- Includes 2 additional tests for external QC testing
- 100% sensitivity, 96% specificity
- · Clinical studies included children 4 and up
- · CLIA waived for whole blood

(**327-3297**)......27/box

CPT® Code: 86308QW†



### Alere Toxicology

### iScreen™ DX Drug Tests

- CLIA waived
- Simple dip-and-read procedure
- Results in minutes
- Laboratory accurate

....25/box

Single Dip Devices	
Specify:	
#IS AMP-DX,	
Amphetamine (AMP)	(395-0031)
#IS1 BAR-DX,	(005 0000)
Barbiturate (BAR)	(395-0036)
#IS1 BZO-DX, Benzodiazepine (BZO)	(205 0022)
#IS1 COC-DX.	(393-0033)
Cocaine (COC)	(395-0035)
#IS1 MET-DX,	(000-0000)
Methamphetamine	
(mAMP)	(395-0037)
#IS1 MTD-DX	` '
Methadone (MTD)	(395-0038)
#IS1 OPI-DX,	
Opiate (OPI)	(395-0039)
Medicare/Medicaid CPT® Cod	e: G0431QW†

#IS1 OXY-DX. Oxycodone (OXY)..... (459-0011) #IS1 MDM-DX, Ecstasy (395-0043) (MDMA) #DPC-114-DX, Phencyclidene (PCP)..... #IS1 TCA-DX, Tricyclic Antidepressants (TCA).... (395-0040) (395-0041) #IS1 THC-DX. Marijuana (THC) .... (395-0042) #IS5M-DX, COC, THC, OPI, AMP, mAMP...... (395-0044) #IS10-22015-DX, COC THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD, ....(395-0045)

Private Pay CPT® Code: 80101†

AMP: Amphetamine BZO: Benzodiazepines THC: Marijuana mAMP: Methamphetamine OPI: Opiate 2000 PCP: Phencyclidine

ID CODE TRANSLATIONS: IDE TRANSLATIONS:
BAR: Barbiturates
COC: Cocaine
MTD: Methadone
MDMA (Ecstasy): Methylenedioxymethamphetamine

OXY: Oxycodone


TCA: Tricyclic Antidepressants



### Alere Toxicology

### iCup™ DX Drug Screen Cups

- CLIA waived
- · Optional sample validity tests
- Simple process
- All inclusive
- · No urine exposure
- Tester activates test when ready



Tests up to 8 drugs: COC, THC, OPI, AMP, mAMP, BZO, OXY, and MTD. (852-0015).. ......25/box

I-DXA-167-01, 6 Panels

Tests for up to 6 drugs: COC, THC, OPI, mAMP, BZO, OXY; plus SVT adulteration strips. 

### DON'T SEE WHAT YOUR LOOKING FOR? Call us today for more options!

# Substance abuse testing with more substance.



products 800.340.4029 | aleretoxicology.com

02012 Alere. All rights reserved. Alere is a tredemark of the Alere group of compenses. 810104 REV1 7/12

### Laboratory confirmation has never been easier.

«Lab" confirmation kits bring the laboratory right to you... anytime, anywhere. ILab" is a prepaid self contained laboratory confirmation kit for your presumptive positive urine. or oral fluid point of care test devices. (Lab." kits contain everything you need to ship your specimen to our laboratory. and include detailed instructions on how to log onto our website to track your test and print laboratory reports.

Smart. Whether your needs include servicing Occupational Health, Corrections, Propostory, Physician Office Late, Clinical and Therapeutic Montturing, or CLIA Walved or Moderate Complexity testing we have a solution for you.

Connected, the other but become services including presumption positive confirmations and a broader range of drugs and configurations to service the needs of all of your accounts.

Empowering, Our new Part/Pysocribet Soutions. Group and our Employer Services teams will work with you to other a full service solution to your expounds.



Alere Toxicology

### iCassette™ Dx Test

CLIA waived, easy to use, and screens for 11 drugs on one device.

#I-DCB-1115-011: COC, THC, OPI, AMP, PCP, BZO, BAR, MDMA, PPX, OXY, TCA (459-0007) ......25/box



Alere \_ Toxicology

### iCup® AD

All-Inclusive **Drug-Screening Test** 

- Zero exposure
- Built-in validity test
- 3–13 drug configurations\*
   Results can be photocopied
- Collect and read
- 24-month shelf life

#I-DUA-157-013, 5 Panel (CLIA Waived): COC, THC, OPI, AMP, and mAMP (195-4184)......25/box

#I-DUA-157-023, 5 Panel (CLIA Waived): COC, THC, OPI, AMP, and PCP

(673-2588) ......25/box #I-DUA-157-034, 5 Panel (CLIA Waived): COC, THC, OPI, mAMP, and PCP

(204-9814) ......25/box

#I-DUA-167-012, 6-Panel (CLIA Waived): COC, THC, OPI, AMP, mAMP, PCP (OX, SG, PH)

(604-7353) ......25/box

#I-DUA-167-022, 6 Panel (Moderate): COC, THC, OPI, AMP, mAMP, and BZO (104-2333)......25/box

#I-DUA-167-291, 6 Panel (Moderate): COC, THC, OPI, AMP300, OX4, MDMA (OX, SG, PH)

(874-7298)......25/box #I-DUD-187-013, 8 Panel (Moderate):

COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR (OX, SG, PH, NI, GL, CR)

(303-6739) .....

#I-DUE-187-071, 8 Panel (Moderate): COC, THC, MOP, AMP, mAMP, PCP, BZO, OXY (OX, CR, PH)

(113-5160) ......25/box #I-DUE-1107-141, 10 Panel (Moderate):

COC, THC, OPI, AMP, mAMP, BZO, BAR, OXY, MTD, PPX (OX, CR, PH)

(113-2737)......25/box #I-DUE-1127-012, 12 Panel (Moderate):

COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD, TCA, OXY, and PPX (110-5436)......25/box

Not all configurations are CLIA waived.



Alere Toxicology

### iScreen™ One-Step Drug Card

Efficiently detects the presence of drugs in minutes.

- On-site testing
- Saves you time and money
   Simple dip-and-read procedure
   Room-temperature storage
- Laboratory accurate
- Results in minutes

#IS4A, 4-Test Card: COC, THC, OPI, and AMP (Moderate) (108-1504)......25/box #IS5A, 5-Test Panel: COC, THC, OPI, AMP,

and PCP (Moderate) (**827-2314**)......25/box

#I55MP, 5-Test Card; COC, THC, OPI,

mAMP, and PCP (Moderate) (**395-8161**) ......25/box

#DOA-164-551, 6-Test Panel; COC, THC, OPI, AMP, mAMP, and OXY (Moderate) (110-8078)......25/box

#I-DOA-364, 6-Test Card: COC, THC, OPI, AMP, BZO, and MTD (Moderate)

(680-1961) ......25/box #IS9, 9-Test Panel: COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, and MTD

(Moderate)

(172-1463)......25/box #IS-10, 10-Test Panel: COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD and

TCA (Moderate) (392-1003) ......25/box

#DOA-1124-011, 12 Panel (Moderate): COC, THC, OPI, AMP, mAMP, PCP, BZO, BAR, MTD, MDMA, OXY, PPX (106-5613)......25/box

ID CODE TRANSLATIONS:
AMP: Amphetamine; BAR: Barbiturates; BZO:
Benzodiazepines; COC: Cocaine; THC: Marijuana; MTD:
Methadone; mAMP: Methamphetamine; MDMA
(Ecstasy): Methylenedioxy-methamphetamine; OPI:
Opiate 2000; OXY: Oxycodone; PCP: Phencyclidine; and TCA: Tricyclic Antidepressants



### **Buprenorphine (BUP)** Drug-Screen Dip

Single-Dip Device (Moderate) Buprenorphine is used in the treatment of opioid addiction. Therapeutically, buprenorphine is used as a substitution treatment for opioid addicts. It is as effective as methadone but demonstrates a lower level of physical dependence.

#DBU-114

(106-8526)......25/box



### AlcoMate Premium

- CLIA status: Moderate
- 4-digit display (clear LCD)DOT approved
- · Low-battery indicator/auto power-off
- Internal mouthpiece storage slot (x2)
- Cumulative test counter
- · No calibration service required

#ALF7000, Kit (114-2179).....ea Contains: 1 AlcoMate, 50 mouthpieces, 1 extra

sensor, 1 hard carry case & 1 car charger. #ATM-Mouthpiece

(114-3049).....50/box

Q.E.D.® A150 Alcohol Saliva Screen

- DOT approved
- CLIA waived
- · Quantitative results
- Saliva test
- Fast, reliable results
- Requires no calibration or controls (102-5216)......10/box

CPT Code: 82055†



Alere Toxicology

### iCassette™

One-Step **Drug-Screening Test** 

 Minimal sample volume required

• 1–11 drug configurations

- Results in about 5 minutes
- 24-month shelf life
- Simple procedure

#I-DOA-2155, 5 Panel, CLIA Waived (447-2471)......25/box

#I-DOA-3155, 5 Panel, CLIA Waived (**724-9019**)......25/box

#1-DOA-1165, 6 Panel, CLIA Waived (922-8619) ......25/box



**Alere** Toxicology

**NEW!** 

### iScreen® OFD-Cotinine

The iScreen® OFD-Cotinine device is a rapid, oral fluid screening test that can be performed any time without the use of an instrument. The test detects elevated levels of cotinine (nicotine metabolite) in oral fluid for accurate determination of smoking status.

- Accurately screen for cotinine (primary)
- metabolite of nicotine)

   Detects up to 1-2 days after nicotine use · Qualitative detection of cotinine in oral
- fluids at a cutoff of 30 ng/mL Simple collection and testing procedure
- · Results in as little as 10 minutes
- Integrated internal procedural control #I-DCT-B702, Nicotine Test



### ✓HENRY SCHEIN®

### **ONESTEP+ H. PYLORI CASSETTE TEST**

- · Whole blood, serum, or plasma
- CLIA waived for whole blood
- Simplified
- Results in 10 minutes or less
- Allows for rapid treatment decisions
- Controls are included

Complete Kit

(900-4071) ......20/box

Contains: 20 pouched test devices with Contains: 20 pouched test devices with disposable sample droppers; disposable, heparinized capillary tubes & dispensing bulb; 1 of each control: (positive & negative); 1 sample buffer; 1 procedure card & 1 package insert.

CPT® Code: H. pylori: 86318QW†



### **W**Kimberly-Clark

### H. pylori Rapid Test MicroCOUNT Lite Meter and Accessories

- Up to 96% sensitivity and 96% specificity
- Active infection determination versus the presence of antibodies, as measured by serology
- Noninvasive, easy-to-swallow gelatin capsules
- Onsite test administration and analysis that can be completed in less than 20 minutes
- Testing not regulated by CLIA
- Reimbursement by Medicare and many

other carriers
#60495, Lite Liquid, Single-Sample Liquid Scintillation Counter <b>(108-1365)</b> ea
#60496, Lite Printer (108-1366)ea
#60442, PYtest Capsules <b>(370-9222)</b> 10/pkg
#60443, PYtest Capsules (106-0064)100/pkg
#60439, PYtest Startup Package

(106-6419) .....ea



### **QUIDEL**

### Quickvue® H. pylori GII™ Tests

For H. pylori; Serum, Plasma, and Whole Blood

- Simple one-step procedure
- Results in 5 minutes or less
- 94% sensitivity, 90% specificity92% overall agreement with EIA
- Room-temperature storage
- Built-in controls for CLIA compliance Positive and negative controls included
- CLIA-waived for whole blood

(227-2170)......30/box (227-5158) ......10/box

CPT® Codes: 86318† (Moderately Complex)/ 86318QW† (Whole Blood Only)



### OSOM® H. pylori Test

For use with whole blood, serum, or plasma.

- · Easy to use: single-reagent procedure fits easily into any practice or lab
  • Fast: results in 10 minutes or less
- Easy to read: 1 pink line indicates a negative result; 2 pink lines indicates a positive
- Small sample size: test only requires 10 μL of serum or plasma; 25 μL of whole
- · Accurate: 95.9% sensitivity

Room-temperature storage	
#175, Test	_
(309-0040)	.25/box
Control Kit	
(309-0074)	ea

Contains: 2-2-mL control vials (1 of ea: positive

& negative). CPT® Codes: 86677QW† (Serum/Plasma); 86318QW† (Whole Blood)





### Clearview® H. pylori Rapid Test

For the visual qualitative detection of antibodies specific to Helicobacter pylori in human serum, plasma, or whole blood. Simple, two-step procedure.

- · CLIA waived for whole blood only
- 88% sensitivity; 88% specificity
  Results in 10 minutes or less
- · Room-temperature storage for up to 18 months
- · Built-in procedure controls to ensure accuracy

#92405

(512-4160)......30/box

CPT® Codes: 86318QW†/86318†

### **THYROID**



Elevated TSH is the first clinical indicator of hypothyroidism. This test allows you to know in minutes if your patient's TSH is greater than 5 mLU/L. Requires 1 drop of blood from a fingerstick or venous draw.

ThyroChek® Test Kit (**818-0440**) .....ea Contains: 20 tests, pipettes, buffer, & instructions. #120C, Control Set (818-0001).....2 vials CPT® Code: 84443QW†

### DVT





### Clearview® Simplify D-dimer Test

10-Minute Rapid Test

Rapid lateral flow test for qualitative detection of D-dimer to aid in the diagnosis of deep vein thromboses (DVT), pulmonary embolism (PE), and disseminated intravascular coagulation (DIC).

- Fast: 10-minute results
- Easy: 2 easy steps
- Flexible: uses whole blood or plasma
- · Reliable: built-in control

(891-9725)......10/box

CPT® Code: 85380†

### **Microhematocrit Capillary Tubes**

Red Tip, Heparinized (130-3223)	200/vial
Blue Tip, Plain (227-5359)	200/vial

### **CONJUNCTIVITIS**



NEW AND IMPROVED!

### RPS AdenoPlus™ for Pink Eve

CLIA-waived, 10-minute, in-office test used ocia-waived, 10-minute, in-ollice test used for conjunctivitis (pink eye). Working from an easy-to-collect eye-fluid sample, AdenoPlus™ detects all known serotypes of adenovirus. Knowing the underlying infection source helps clinicians determine the most appropriate treatment.

- Reimbursed
- Proven accuracy
- Reduces antibiotic resistance

#RPS-AD. Test (770-0028) ......10/box

CPT® Code: 87809QW†

### **BLADDER CANCER**



### Status BTA Test

A Rapid Test for the Qualitative Detection of Bladder-Tumor-Associated Antigen in Human Urine

- CLIA-waived
- Uses monoclonal antibodies to detect the presence of bladder-tumor-associated
- Accurate test results with superior sensitivity and specificity
- Detects early stage recurrence cancer that cytology may miss

  Adjunct test to cystoscopy

  Simple 1-step procedure

  Results available in 5 minutes

(641-0001)......25/box

### **SEDI-RATE**



### Dispette 2 Sedimentation Rate Kit

With improved safety features plus greater efficiency. Modified Westergren sedimentation rate. Includes 100 ESR disposable plastic pipette tubes and 100 prefilled, pierceable cap reservoirs. (668-9703) .....ea

### DTM



### **DTM Test Agar Dermatophyte Test**

Medium Vial

Selective isolation of fungi from cutaneous (skin) samples.

(**224-9011**)......20/box

Item stored under refrigeration. May be shipped separately.





### Sedi-Rate™ Westergren ESR System and Accessories

The erythrocyte sedimentation rate (ESR) system is a nonspecific test that alerts physicians to the presence of disease or abnormality in the body. Sedi-Rate's pipette features a unique self-sealing aerosol barrier plug designed to seal upon contact with blood, eliminating exposure to blood aerosols. Its polypropylene vial is prefilled with 0.25 mL of 3.8% trisodium-citrate diluent and features a self-sealing stopper. Test results are ready in 60 minutes and can be read directly from the imprinted graduation marks on the pipette.

#3467, Starter Kit **(217-0930)**....ea

Contains: 100 Sedi-Rate pipettes & vials, 1–10-place rack with pipette support clamps, 1 leveling platform, 100 piercing funnels & 100 ESR transfer pipettes.

#3469, System (368-5132).....100/box Contains: 100 tests (100 pipettes & 100 vials). #3479, Pipettes (**600-7387**) ......100/box


**Individual Kit Components** #3459, Acrylic Leveling Platform (**691-4303**) .....ea

#3452, Rack with 10 Pipette Support Clamps, Holds 10 Tests (**496-3393**).....ea

# 50/0 Off! YOUR NEXT REPAIR







# Low rates do not always mean big savings with credit card processing.

Through our Complimentary Rate Analysis,
Henry Schein Financial Services can help you
understand what your current rate and overall costs are...

Sign up today and receive a \$100 Activation Bonus!

### Plus, benefit from:

- \$500 Rate Guarantee!\*
- Exclusive savings
- All major credit cards accepted
- Debit cards and e-checks accepted
- Payments through website and iPhone accepted
- Only \$1 for monthly terminal rental
- No application or set-up fees
- Detailed online reporting
- Fraud-monitoring services

1-800-443-2756

or e-mail: hsfs@henryschein.com www.henryschein.com/hsfs

42



### HITACHI Inspire the Next

### **CLA-1™ Luminometer**

This compact luminometer automatically prints a complete report that lists the tested allergens and the severity of the patient's response to each. Easy to use, the system is the simplest, quickest, and most convenient allergy testing method available today. It can analyze up to 180 allergens in 8 minutes.

(680-1614) .....ea

### CLA® Allergy Tests

20/box
20/box
20/box
20/box
IgE
20/box

(406-3895)20/box Contains: 20 CLA® Pettes & reagents.
Southern, Combines IgE Foods and Inhalants (259-4912)
Southwestern Inhalant (406-5585)20/box Contains: 20 CLA® Pettes & reagents.
Western Inhalant (168-6399)20/box Contains: 20 CLA® Pettes & reagents.

Southeastern Inhalant



Hb 201+ Handheld

Hemoglobin Analyzer
Offers the same important benefits as its predecessor, the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
   Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
  • Safe and easy to use
- Reimbursable
- CLIA '88 waived
- 2-vear warranty

* 2-year warranty	
Analyzer <b>(345-2344)</b>	ea
Contains: analyzer, AC adapter CD-ROM.	r, manual & training
#3201 Microcuvettes (110-9167)	100/box
Contains: 4 packs of 25. <b>(345-8733)</b>	200/box
Contains: 4 packs of 50.	



Contains: 20 CLA® Pettes & reagents.



### LeadCare® II Analyzer

Accurate, On-the-Spot Lead Testing

- Perform blood-lead measurements with only 2 drops of blood
- In-office immediate results in 3 minutes
- · Accuracy and sensitivity of a reference lab test (from capillary or venous blood sample)
- Two-drop sample size requires only a finger stick; allows combination of LeadCare II with other tests
- · CLIA-waived, reimbursable test

Analyzer	
(106-4441)	ea
Test Kit	
(106-5736)	48/box
CPT® Code: 83655†	



# **Empowering** Cardiovascular Alere Treatment Decisions

The right diagnostic tools provide the timely, actionable data you need to make better in-office decisions that can improve your patient's quality of life.





### Complete Lipid Profile with the Cholestech LDX® System

Get your patients to goal, and experience a smarter, faster workflow that improves office efficiencies. CLIA waived Lipid Profile, Glucose, and ALT+AST results in 5 minutes.



### Track Coumadin® Patients with the INRatio®2 System

Empower face-to-face dosage change decisions that can help your patients stay in therapeutic range. CLIA walved, fingerstick PT/INR results in about 1 minute.



### Manage Heart Failure with the Triage® Meter

Diagnose, assess severity, and risk stratify your patients at risk for heart failure in approximately 15 minutes, right in your office. BNP test is CLIA waived for whole blood.

# Be empowered. Connect to rapid results. Make smarter decisions.

Contact your sales representative today!

alere.com

635112 Asse, All rights received. Cholestech LDK, RiPlatia, Trage, and the Aline Logo, are trademarks of the Alive group of companies. Counselin is a treasment of Bristo-Myers Source Pharms Company, 1000017vA



### BISSITE

### Triage® Profiler S.O.B. (Shortness of Breath) Panel

Rapid fluorescence immunoassay of the quantitative determination of creatine kinase MB (CK-MB), myoglobin, troponin I, B-type natriuretic peptide (BNP) and Ddimer in whole blood or plasma. Test devices include 2 levels of built-in controls. Results in approximately 15 minutes, read by Triage® MeterPro™ meter.

Measurement of 5 markers:

- Troponin I
- CK-MB
- MyoglobinBNP
- D-dimer

#97300. Panel (241-0108) ......25/box



### BIGSITE

### Triage® Cardiac Panel

Measures 3 key cardiac markers: myoglobin, CK-MB, and troponin I.

- · Easy to use, with results in approximately 15 minutes
- Results read by Triage<sup>®</sup> MeterPro<sup>™</sup> meter
- Whole blood or plasma

#97000HS

(131-1964).....25/box

### BI@SITE

### Triage® Control 5 Calibration Verification

2 x 5 Levels

#88755

(111-9927) ......5/box

### BIGSITE

### Triage® Total Control 5, Control 2

For Cardio Profiler, S.O.B. and Cardiac Panel 5 x 0.25 mL

#88754

(111-9843) ......5/box

### BISSITE

### Triage® Total Control 5, Control 1

For Cardio Profiler, S.O.B. and Cardiac Panel

5 x 0.25 mL

#88753 (111-9841) ......5/box

# BIGSITE

### Triage® Cardio Profiler Panel

Measurement of 4 cardiovascular markers: troponin, CR-MB, myoglobin, BNP. Results in approximately 15 minutes. Results read by Triage MeterPro™ meter.

#97100. Panel

(706-7343) ......25/box

### BISITE

### Triage® MeterPro™ Meter

High-performance, comprehensive cardiovascular and drugs-of-abuse testing in the laboratory or at the point of care (POC). Designed to provide diagnostic results quickly and easily. Quantitative diagnostic panel results are displayed and printed in easy-to-read numbers and highlighted when the result is above the cutoff. Qualitative diagnostic panel results are displayed as positive (POS) or negative (NEG), eliminating the need for interpretation of lines that often hamper qualitative diagnostics.

- 16 immunoassays availableResults in about 15 minutes
- · Whole blood, plasma, or urine
- Customizable test panels
- Onsite training and implementation support
- Results printed, displayed, and interfaceable
- Full connectivity support



#55070MP ( <b>496-3405</b> )ea	a
#52031, Paper for Triage Meter and MeterPro, 1 Roll	
(111-6321)ea	a



### BISSITE

### Triage® D-dimer Test

Rapid, quantitative immunoassay to aid in the assessment and evaluation of patients suspected of having disseminated intravascular coagulation (DIC) and thromboembolic events, including pulmonary embolism (PE) and deep-vein thrombosis (DVT).

Performed on Triage® meters
Provides rapid, quantitative results in approximately 15 minutes

Utilizes highly sensitive fluorescence immunoassay technology

Utilizes the 3B6 D-dimer antibody, which demonstrates high clinical specificity for PE and DVT

#98100, Test (131-0268)	.25/box
#98113, Control Set, Control 1 (5 x 1 (176-9303)	
#98114, Control Set, Control 2 (5 x 1 (241-0075)	



### BISSITE

### Triage® BNP Test

Rapid, objective point-of-care immunoassay is used:

- As an aid in the diagnosis and severity of CHF
- · For the risk stratification of patients with ACS and HF

#98000XR, BNP Test <b>(131-7507)</b> 2	25/box
#98013XR, BNP Control Set, Level 1 (5 Level 1) (219-6429)	ea
#98014XR, BNP Control Set, Level 2 (5 Level 2) (219-6475)	
#98015XR, BNP Calibration Verification	n

Control Kit (1 per Level, 5 Levels) (131-9669) .....ea #98200. BNP Reagent for Beckman Coulter Immunoanalyzers

(104-4059)......100/box #98201, BNP QC Controls for Beckman

Coulter Immunoanalyzers **(478-8086)** .....ea #98202, BNP Calibrator for Beckman

Coulter Immunoanalyzers (640-1104).....ea

### BISSITE

### Triage® TOX Drug-Screen

Simultaneously detects 11 drugs in one rapid, easy-to-use test that is read on the Triage® meter.

#94400, Test	
(131-0824)	25/pkg
#94413, 0.25-mL Control Level 1	
(131-7103)	5/pkg
#94414, 0.25-mL Control Level 2 (131-0969)	5/pkg
(101 0000)	o, prig



### Cardiac Reader Analyzer/D-Dimer Systém

Cardiac Reader Analyzer



(100-3004)	<del>e</del> a
Cardiac T Quantitative Assay	
(106-5803)	10/box
Cardiac M Assay	
(106-5802)	20/hox





### Piccolo® Point-of-Care **Chemistry Analyzer Accessories**

The first broad-menu clinical chemistry analyzer designed for near-patient testing. Basic Panel contains BUN, Ca, Cl-, CRE, GLU, K+, Na+, and tCO<sub>2</sub>. Comprehensive Panel contains ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO<sub>2</sub>, and TP. Lipid Panel contains CHOL, CHOL/HDL (calculated), HDL, LDL (calculated), TRIG, and VLDL (calculated).

- General chemistry, liver, electrolyte, BMP. CMP, lipid, renal, and hepatic panels available
- · No daily controls, calibrators, or maintenance
- Intelligent Quality Control (iQCT) verifies correct performance
- No special training required
   O.1cc whole blood, serum, or plasma
- Single-use reagent discs
- Compact and portable
- Built-in QC and calibration
- Lipid panel is CLIA waived
- Whole blood, serum, or plasma · Less than 2 minutes hands-on
- Results in under 13 minutes

#100-0000, Accessory Kit (106-6368) .....ea

Contains: 1-0.1cc Minipette, 1-96-pack of Contains: 1–0.1cc Minipette, 1–96-pack of disposable Minipette tips, 1 operator's manual, training CD, 1 power cord, 1 power supply, 10 lithium-heparin tubes, 1 screwdriver & 2–12-packs of results cards.

#400-0024, Basic Metabolic Panel (137-3987)	10/box
#400-0028, Comprehensive Metabolic Panel (174-6978)	10/box
#400-0025, Lipid Panel	

(146-0189) ......10/box

### **CHOLESTECH ₹**

### Cholestech LDX® System

Complete Lipid Profile, Glucose

- CLIA waived\*\*
- Procedure saves time and enhances productivity
- Small, tough, and portable
- · Eliminates anxiety of venipuncture; requires only a fingerstick
- On-the-spot adjustment of therapy and improved patient compliance
- Accuracy comparable to laboratory tests
   Saves time: fewer follow-up visits required
   Enhances compliance through teachable moments between healthcare providers and patients

### #10-959, Analyzer

### (**738-2917**).....ea

Contains: 1 of ea: analyzer, printer, power supply, starter pack, optics check cassette, printer cable, user's manual, procedure manual & training video. Starter pack contains 200 lancets, 200 capillary tubes, 200 capillary plungers, 1
Minipette pipette, 50 pipette tips, 1 accessory tray & 1 fingerstick tip pad.

#11-785, Cholestech Thermal Printer Labels

White labels for use with Cholestech Thermal Lahel Printer

mema Laberi milei.	
(803-6378)	.400/roll

#10-989 Test for Lipid Profile
(116-9995)10/box

#10-986. Test for TC (420-7892)......10/box

#10-988. Test for TC and Glucose (956-0149).....10/box #10-987, Test for TC and HDL Panel


(600-5779) ......10/box #10-990, Test for TC, HDL, and

Glucose Panel (600-5781) ......10/box

#10-991, Test for Lipid Profile and Glucose (694-9360)......10/box #12-788, Test for ALT/AST Alanine Aminotransferace







(116-4232).....ea Contains: 1 vial each levels 1-2. #12-807, High-Sensitivity C-Reactive Protein (hs-CRP)

(**392-1076**)......10/box Contains: 1-2-mL vial of ea level (1-4). #11-255. Calibration Verification Material

10/hox

(Lipids/Glucose), 2 mL (116-7301).....ea Contains: 1 vial each levels 1–4.

#11-010, Disposable Pipette Tips For use with 35µL, 40µL, and 50µL MiniPet Pipettes. (**267-4596**)......50/box

#10-311, Capillary Plungers For use with all types of Cholestech capillary tubes.

(960-1587) ......50/box

CPT® Codes for lipid profile: #80061QW†; for lipid profile/glucose: 80061QW†/839471QW†; for ALT/AST: 84460QW†

### **CHOLESTECH**

### **Printer Cable**

Connects the Cholestech LDX to the Thermal Label Printer.

#11-787

(803-3213) .....ea



Abbett Disbetes Care

### i-STAT® 1 Handheld System

Includes these important features:

- · Easy to use: Cartridges are self-contained and the system is simple to learn and operate
- Lab-áccurate results: Advanced biosensor technology and real-time quality assurance
- provide extremely accurate testing Small sample size: Just a few drops of blood per cartridge
- Ready when you need it: Test results in as little as 2 minutes

i-STAT Analyzer (111-3324)	.ea
Cartridges Chem8+ (waived) (111-3408)25/k	хос
Creatine (waived) (111-3394)25/k	оох
CG8+ (waived) (115-6662)25/b	хос
PT <b>(111-3385)</b> 24/t	хос
BNP (111-3158)25/k	хос
Accessories i-STAT Serial Recharger (111-3403)	.ea
i-STAT Serial Downloader (111-3401)	.ea
i-STAT Printer Paper	

(111-3398) ......6/box



### CoaguChek XS® System PT/INŘ

Monitoring Solution The next generation in point-of-care prothrombotin (PT) and international normalized ratio (INR) testing, this new system

uses exclusive smart technology to help ensure accurate results and offers the flexibility and control that healthcare professionals need to confidently manage a wide range of anticoagulated patients. It is the only system to perform onboard quality control and determine patient results in a single test chamber.

- CLIA waived
- Neutralizes therapeutic levels of heparin and LMWH
- INR corrected for hematocrit within specified range; manage a broad array of patients with a single device



- Small sample size (10 μL) for easy dosing
   Strip allows dosing from either top or side for simple blood application
   Small, battery-powered, handheld meter for
- portability and efficiency
- Onboard controls: no external QC necessary
- Auto-on with strip insertion and auto-off capabilities
- 18-month strip shelf life; no refrigeration needed

3-Box Starter Kit (108-1861)ea
Contains: XS Analyzer, 3 boxes of CoaguChek XS PT strips (144 total) & CoaguChek Lancets.

5-Box Starter Kit (108-7757) .....ea

Contains: XS Analyzer, 5 boxes of CoaguChek XS PT test strips (240 total) & CoaguChek Lancets

CoaguChek XS® Plus System	WAIVED
Battery	WAIVED
Battery (109-9303)	ea
10-Box Starter Kit	
(109-3345)	ea

NOW CLIA

Contains: CoaguChek® XS Plus Analyzer, 10 boxes (480 total) of CoaguChek XS PT test strips & CoaguChek lancets.



Meter #PRO3-25, Disposable Cuvette,

(487-7457) ......25/box

5 Channels

CPT® Code: 85610QW†



### CoaguChek XS® Test Strips

For XS and XS Plus Meters (891-7593)......48/box

Contains: 24 strips per bottle, 2 bottles per box.



### VerifyNow™ Aspirin Assay

A qualitative assay to aid in the detection of platelet dysfunction due to aspirin ingestion incitrated whole blood for the point-of-care or laboratory setting.

- Simple: 1-step, closed-tube sampling
  Rapid: results in 5 minutes
  Accurate: >99% specificity; 91.4% sensitivity
- · Cost-effective: eligible for Medicare and private-payer reimbursement with CPT<sup>®</sup>
  Code 85576QW†; platelet aggregometry-arachidonic acid (agonist) (113-2708)......ea

### VerifyNow® Aspirin ASA

**CLIA Waived** (369-4955) ......25/box CPT Code: 85576QW+ (-91 Modifier)

### VerifyNow® Plavix® P2Y12 Assay Moderately Complex (261-0249)......25/box

CPT Code: 85576† (-91 Modifier)



### **ProTime® Monitor**

- Accepts and stores patient and/or operator ID for reduced transcription errors

  • Automatic: sends results quickly and easily
- to printer, computer, or LIS

  Only system that offers flexibility of both
- onboard and optional external quality controls
- Thromboplastin reagent of ISI=1.0 (AHA, CAP, and WHO recommended)
   Off-meter sampling provides safety against
- disease transmission
- Bar-coded cuvette: no coding necessary
- Longer open-pouch stability: up to 16 hours





### HernoSonse...

### INRatio2® Monitor

Allows you to obtain a prothrombin time (PT) and international normalized ratio (INR) result using 1 drop of blood from a fingerstick.

- CLIA waived
- Easy to use
- Single-button operation
- 1-step auto turn-on with strip insertion
  Larger LCD screen and easy-to-read digits
  Small, 1-drop (15 μL), sample size
- · Wider test-strip area for easier sample application
- Test results in under 1 minute
- Test strips do not require refrigeration
- 2-channel, on-board quality control runs on every test
- No need for external liquid quality control or extra strips
- Test strips are individually wrapped and stored at room temperature

### #INR2QS1 Quick Start #1 (111-3840) .....ea

Contains: analyzer, 1 box of 48 INRatio test strips, 50 Unistik® single-use lancets, 100 alcohol pads & 250 gauze brick pads.

### #INR2QS3 Quick Start #3 (111-3844).....ea

Contains: analyzer, 3 boxes of 48 INRatio test strips (144 total, 150 Unistik® single-use lancets, 200 alcohol pads & 250 gauze brick pads.

### #INR2QS5 QuickStart #5 (111-6578).....ea

Contains: analyzer, 5 boxes of 48 INRatio test strips (240 total), 250 Unistik® single-use lancets, 300 alcohol pads, 5 boxes of capillary tubes &

### 500 gauze brick pads. #INR2-6PK 6-Box Promo (112-3230).....ea

Contains: analyzer, 6 boxes of 48 INRatio test strips (288 total), 10 single-use lancets & 10 capillary tubes. CPT® Code: #85610QW†

# **BD Diagnostics** – Preanalytical Systems

# **Technology Built On Trust**

### The BD Vacutainer® Push Button Blood Collection Set

- The only blood collection system with push-button technology
- Split-second activation for improved safety
- In-vein activation assures retraction of the needle directly from the vein
- Intuitive design reduces training



HS Reference No.	BD Reference No.	Product Description
987-9975	367326	PBBCS 12" w/o Luer 21 x.75
987-2637	367324	PBBCS 12" w/o Luer 23 x .75
987-3303	367323	PBBCS 12" w/o Luer 25 x .75
987-0431	367338	PBBCS 7" with Luer 21 x .75
987-6287	367336	PBBCS 7" with Luer 23 x .75
987-8340	367335	PBBCS 7" with Luer 25 x .75
987-3358	367344	PBBCS 12" with Luer 21 x .75
919-4247	367342	PBBCS 12" with Luer 23 x .75
987-9168	367341	PBBCS 12" with Luer 25 x .75

### BD Microtainer® Contact-Activated Lancets

- Designed with a Positive Patient Experience in Mind
- Activates on Contact
- Improves Puncture Accuracy
- Simplifies the Procedure



The value of any rebates, discounts, incentives, or product provided to customer may constitute a "discount or other reduction in price" under Section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. Sec. 1320a-76(b)(3)(A)). Customer shall satisfy any and all requirements imposed on buyers relating to discounts or reductions in price, including, when required by law, to disclose all discounts or other reductions in price received from BD and to accurately report under any state or federal health care program the net cost

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2009 BD DIST072

HS Reference No.	BD Reference No.	Dimensions (Width x Depth)	Blood Volume
106-5174	366592	30 G x 1.5 mm, Purple	Low Flow
106-5210	366593	21 G x 1.8 mm, Pink	Med. Flow
106-5209	366594	1.5 mm x 2.0 mm, Blue	High Flow


View the product demonstration video!

Please visit: www.bd.com/vacutainer



**BD** Diagnostics

Preanalytical Systems 1 Becton Drive Franklin Lakes, NJ 07417 www.bd.com/vacutainer





Finger activation

- · Needle retracts automatically after use
- Cannot be reloaded after use
- · Safe for patient and user

28 ga, Fine (900-4433)	200/box
23 ga, Medium (900-4434)	200/box
17 ga, Coarse <b>(900-4435)</b>	200/box



Tenderfoot® Lancets Infant heel blood sample.

(107-9563)......50/box



Unistik® 2 Safety Lancets

Normal, for Standard Tests, Yellow, 2.4 mm (668-3836) ......50/box Normal, for Standard Tests, Yellow, 2.4 mm (668-5603) ......200/box Extra. for Larger Samples, Orange, 3 mm (159-3956).....100/box Extra, for Larger Samples, Orange, 3 mm (668-6699) ......200/box



### OneTouch® SureSoft™ Lancing Devices

- Electropolished needle means a smaller gauge and more comfort for your patients
- Meets OSHA guidelines: nonreusable, retractable needle
- Distinctive wing design prevents reloading and reusing
- Preset gauge for no guessing, no maneuvering and easy training

Gentle 28-ga Needle (567-0000) ......200/box

SureSoft Regular, 21 ga x 2 mm

(115-5267)......200/box





environment. Won't cause an allergic reaction in latex-sensitive patients or employees. Convenient packaging allows for easy, one-at-a-time dispensing that facilitates a single-use policy that will reduce danger of cross contamination between patients and healthcare workers.

#367203

(270-3884) ......25/box



Blood-Collection Needles

Thin Wall, 11/4"

Safety-engineered, multisample needle offers a simple, effective way to collect blood while it reduces the possibility of needlestick injuries. It features a patented shielding that allows for one-handed activation to cover needle immediately upon withdrawal from vein and confirm proper activation with an audible click. New needles with preattached holders are ready to use with no assembly required.

#368608, 22 ga ( <b>987-2645</b> )	48/box
#368607, 21 ga ( <b>987-6956</b> )	48/box
#368651, with Preattached H (987-6134)	
#368650, with Preattached H (987-5021)	

### **BD** BD™ Microtainer® Contact-Activated Lancets



Designed with a positive patient experience in mind. The contact activation method facilitates a consistent puncture depth and minimizes the likelihood of having to repeat the puncture. It covers only a small area at the contact point, resulting in improved visibility of the puncture site for the clinician and greater accuracy of lancet positioning when performing the puncture. Its innovative ergonomic design allows for a more comfortable grip. The lancet automatically retracts into the device, which prevents the lancet from being reused. In addition, the lot number is laser-etched on each lancet for easier tracking.

#366592, Low Blood Flow (Single Drop), Purple, 30-ga x 1.5-mm-D Needle (106-5174)200/box
#366593, Medium Blood Flow, Pink, 21-ga x 1.8-mm-D Needle (106-5210)200/box
#366594, High Blood Flow, Blue, 1.5mm-W x 2-mm-D Blade (106-5209)200/box





shield that allows for one-handed activation immediately after use to help minimize the risk of needlestick injuries and to offer clear visibility of blood flashback. #367294, 7" Tube with Luer, 25 ga (987-5250)......50/box #367292, 7" Tube with Luer, 23 ga x 3/4" (987-0529) .....ea #367287, 7" Tube with Luer, 21 ga (958-5669) ......200/case #367285, 12" Tube with Luer, 25 ga (840-3447) .....ea #367283, 12" Tube with Luer, 23 ga (987-1301).....50/box #367281, 12" Tube with Luer, 21 ga (987-8393).....50/box #367298, 12" Tube without Luer, 25 ga (**623-4670**) ......50/box #367297, 12" Tube without Luer, 23 ga (987-4692).....50/box #367296, 12" Tube without Luer, 21 ga (987-9967).....50/box #368653, 12" Tube with Preattached Holder, 23 ga (987-8940)......25/box #368652, 12" Tube with Preattached Holder, 21 ga (111-7212)......200/case



### BD Vacutainer® **Push-Button Blood** Collection Sets

Set offers split-second

protection for that single moment that could potentially change your life. Push-button safety mechanism instantly helps protect you against needlestick injury. BD offers the only blood-collection system with pushbutton technology. Set is designed to be simple, easy to use, ergonomic and, most

important, helpful in protecting you from accidental needlestick injuries.  Collection Sets50/box		
Specify: #367323, without Luer Adapter, 12" Tubing, 25 ga x ¾"	(987-3303)	
#367324, Light Blue without Luer, 12" Tubing, 23 ga x ¾"	(987-2637)	
#367326, Green without Luer, 12" Tubing, 21 ga x ¾"	(987-9975)	
#367341, Royal Blue with Luer, 12" Tubing, 25 ga x ¾"	(987-9168)	
#367336, Light Blue with Luer, 7" Tubing, 23 ga x ¾" #367338, Green with Luer,	(987-6287)	
7" Tubing, 21 ga x ¾"	(987-0431)	
#367342, with Luer, 12" Tubin (919-4247)	g, 23 ga x ¾" 50/box	
#367344, with Luer Adapter, 1 21 ga x ¾"	2" Tubing,	
(987-3358)	50/box	
#367344, with Luer, Green, 12 21 ga x ¾"	" Tubing,	
(405-4902)	200/case	
#367341, with Luer, Royal Blu 12" Tubing, 25 ga x ¾"	e,	
(679-6358)	200/case	





### BD™ Vacutainer® Plus Plastic **Blood-Collection Tubes**

Offer a safer method of blood collection. Reduce risk of tube breakage and specimen spillage. Safe, simple disposal in accordance with EPA guidelines. #367981, 3.5-mL SSTT Tube with Clot

Activator/Gel for Serum Separation, Red/Gray Conventional Top, Paper Label, Silicone-coated Interior, 13 mm x 75 mm (987-1419) ......100/box

#367986, 5-mL SST<sup>™</sup> Tube with Clot Activator/Gel for Serum Separation, Gold Hemogard Top, Paper Label, Silicone-coated Interior, 13 mm x 100 mm (987-0936).....100/box

#367988, 8.5-mL SST™ Tube with Clot Activator/Gel for Serum Separation, Red/Gray Conventional Top, Paper Label, Silicone-coated Interior, 16 mm x 100 mm (987-2786).....100/box

#367856, 3-mL EDTA Tube/K2EDTA, Lavender Hemogard<sup>™</sup> Top, Paper Label, 13 mm x 75 mm (127-4412) ......100/box

#367861, 4-mL EDTA Tube/K₂EDTA, Lavender Hemogard™ Top, Paper Label, 13 mm x 75 mm (987-7504)......100/box

#364300, 2-mL Glass Tube/K2EDTA, Lavender Conventional Top, Paper Label, 10.25 mm x 64 mm (908-1970)......1000/case

49

# Making Blood Collection Safe, Precise, Simple



### Capiject® Safety Lancets

High quality, single-use lancets offer the form, fit and function you've come to rely on.

- · Never see or handle the sharp
- · One-step activation
- · No special staging required
- · Color-coded sizes

### Surshield® Safety Blood Collection Sets

Help shield your staff against accidental needlesticks with Surshield.

- · Simple, one-handed activation
- · Audible click for confirmation that needle is securely locked
- . Available in 21G and 23G in 7" and 12" tubing, and 25G with 12" tubing,

DVLECT is a trainment of Tenuno Medical Corporation, Semented, New Jersey SURTHELLS 😚 and Tainmin, are registered trainments of Tenuno, Corporation, 25 2015 Senson Markey Committee, 1955. All rights emplaint 2020



TERUMO MEDICAL PRODUCTS



### TERUMO

### Capiject® Safety Lancets

1.5 mm x 2.0 mm, Blue (110-8852)	200/box
1.5mm x 1.5mm, Green (110-8853)	200/box
1.5 mm x 1.0 mm, Pink (110-8854)	200/box
28 ga, 1.25 mm, Purple (110-8541)	200/box
23 ga, 2.25 mm, Orange (110-8539)	200/box

# TERUMO\*

### Surshield™ Safety Wing **Blood-Collection Sets**

One-handed, easy-to-activate safety device. Requires no change in needle

insertion technique. Built-in safety design is easy to use. Latex-free.

Long, 21 ga (315-7504)	50/box
Long, 23 ga (315-0039)	50/box
<b>12" Tubing</b> Long, 25 ga <b>(315-0434)</b>	50/box
<b>7" Tubing</b> Short, 21 ga <b>(315-5499)</b>	500/case
Short, 23 ga (315-6596)	50/box



### TERUMO

### Capiject®

- Extra-wide, built-in collection lip for enhanced specimen flow
- Single-piece, tapered design for easy collection and handling
- Color-coded fill lines and caps for volume accuracy and easy identification
   Lot number and expiration date printed on each tube for populidatification.
- each tube for easy identification
- Silicone-coated anticoagulant tubes designed to reduce clotting
- Olefin-based gel for clean separations
- Tube nests in cap, allowing tube to stand alone
- · Certified for lead testing
- Latex-free

Activator (Silica Particles and Inert Gel) (315-9272)250/box
#T-MQK, Lavender, K2-EDTA (0.78 mg) (315-0280)250/box
#T-M, Brick Red, No Additive/Clot Tube (335-9969)1000/case

#T-MLH, Green, Lithium Heparin (12.5 IU) (315-0873)250/box
#T-MLHG, Mint Green, Lithium Heparin with Gel (12.5 IU) (315-1116)250/box
#T-MPS, Gray, Sodium Fluoride (0.72 mg), Potassium Oxalate (0.58 mg) (476-8940)1000/case
#T-MGA, Gold-Amber Tube with Gel Barrier and Clot Activator (Silica Particles and Inert Gel) (609-2050)1000/case
#T-MLHGA Mint Green-Amber Tube with

Lithium Heparin and Gel Barrier (12.5 IU) (315-7531).....1000/case

# Unistik 3 Safety Lancets



Unistik 3 safety lancets are pre-set, easy-to-use, OSHA compliant and feature Comfort Zone Technology for a virtually pain-free sampling.



	Unistik 3 Comfort	Unistik 3 Normal	Unistik 3 Extra	Unistik 3 Dual
Item Number:	668-6373 / 100ct.	668-6939 / 100ct.	668-4166 / 100ct.	668-0329 / 100ct.
	113-2969 / 200ct.	106-1400 / 200ct.	644-5823 / 200ct.	668-4594 / 200ct.
Specification:	Gauge: 28g,	Gauge: 23g,	Gauge: 21g,	Gauge: 18g ,
	Depth 1.8mm	Depth 1.8mm	Depth 2.0mm	Depth 1.8mm
Average Blood Volume:	5μl - 30μl	30μl - 75μl	75μΙ - 125μΙ	200μΙ+

For additional information contact Owen Mumford Customer Service at 800-421-6936 or info@owenmumfordinc.com






### UNICO

### Tube-Cubes®

Available in 3 configurations to accommodate practically all of your needs. Each cube measures  $2\%^{\circ}\text{H} \times 2\%^{\circ}\text{W} \times 2\%^{\circ}\text{D}$ . The 9-place version holds 8-mm to 16-mm/10-mL to 15-mL tubes. The 4-place version holds 26-mm to 30-mm/50-mL tubes. The single place version holds 66-mm cups, alcohol dispensers, small flasks, and more.

- Invert and pour: tubes will stay in place
- Great for phlebotomy trays
- Organize and store samples in refrigerator or freezer

#54900, 9 Place, for 8-mm-16-mm Tubes (112-4212)4/pkg
#54400, 4 Place, for 26mm–30-mm Tubes (112-4210)4/pkg
#54100, Single Place, for Up to 66-mm Tubes
(112-4214)4/pkg



### UNICO

### **Tube Rack**

Starting with a single 50-place T-Rack, you can attach as many racks as space allows. These patented racks attach together quickly and easily with their unique built-in T-connector, which allows the flexibility to attach as many racks as needed.

- Each rack holds up to 50 tubes of varying sizes in two 25-hole sections
- Top shelf has cross-referencing alphanumeric designation in order to track exact location of specific tubes in rack
   Attach additional T-Racks together; make
- Attach additional T-Racks together; make a single large rack as big as needed
- Made of durable, autoclavable polypropylene; will not break or become brittle in low-temperature refrigerators
- Easy to clean and disinfect
- Arrives disassembled and can be snapped together in seconds
- Dimensions: 8¾" x 3¾" x 5¼"

#54500, 50 Place (Expandable) (111-7761) 13-16mm .....ea



### Tube-Cube® Carrier

Durable acrylic, and designed to hold any combination of Tube-Cubes. Available in 2-, 3-, and 4-place configurations. Provides a safe, convenient way to carry patient samples to the lab.

#54040 (112-4213).....ea





### SmoothRack Test Tube Racks

Features include:

- 72-tube capacity9.7" x 5.0" x 2.6"
- Autoclavable polypropylene
- Accommodates test tubes 10 mm to 16 mm in diameter

TO THIS IS GIASTICIO	
Blue	
(720-0000)	ea
White	
(720-0001)	ea

# unico

### PowerSpin" **FX Centrifuge**

Features a 6-place angled rotor designed for 2–10-mL tubes (maximum tube size up to 16 mm x 110 mm) and a fixed speed of 3400 rpm, generating 1350 RCF (g) for superior specimen separation. Set the 0–30-minute timer to desired run length.

- Maintenance-free, brushless motors run quietly at a fixed speed of 3400 rpm. emitting only 58dB, all with a low temperature rise
- Suction-cup rubber feet hold unit firmly in place while virtually eliminating vibration
- Lid features a gasket to reduce aerosol emissions, a window for strobe speed checks, and a safety interlock switch that cuts power to motor if lid is raised

  • Accessories available on request: 8-place
- rotor and pediatric tube adapters
- 1-year manufacturer's warranty
  Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb 115 V/60 Hz

#C806, Fixed Speed, 6 Place (715-5981) .....ea



# UNICO

### PowerSpin™ LX Centrifuge

Spin blood, urine, fecal materials, or semen spiri blood, uffie, fecal friaterials, of serifer at the best speed for each type of specimen with this powerful, variable-speed, 6-place centrifuge. Maintenance-free, brushless motor runs quietly at a fixed speed up to 4000 rpm, emitting only 58dB, all with low-temperature rise. Suction-cup rubber feet hold unit firmly in place and virtually eliminate vibration. Lid features a gasket to reduced aerosol emission, a window for strobe speed check, and a safety interlock switch that cuts power to the motor if lid is raised. Angled rotor accepts 2-10-mL tubes (maximum tube size 16 mm x 110 mm). Accessories available upon request: 8-place rotor and pediatric tube adapters.

- Speed range: 300-4000 rpm, generating up to 1850 RCF (g)
- 30-minute timer
- 1-year warranty
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb 115 V/60 Hz

#C-856, Variable Speed, 6 Place (715-2656) .....ea



### UNICO

### Micro-Hematocrit Centrifuge

This powerful centrifuge, with all-metal housing, spins up to 24 microhematocrit tubes at 12,000 rpm. Tubes are held in place with a screw-down metal lid and replaceable rubber gasket/ring.
The 15-minute mechanical timer has a bell tone to signal a completed run, while the electronic brake can stop the unit on demand. Metal lid features a safety-lid latch that cuts engine power if lid is lifted. Unit comes complete with a metal hematocrit reader tray. 1-year warranty.

#C-MH30

(715-0877) .....ea



# UNICO

### Vortex Mini Mixer

This economical vortex mixer is rugged and sturdy for the most demanding applications. Dual-mode switch

lets you vortex continuously or intermittently using the touch mode: simply press your tube to the soft rubber mixing cup. Features speed control to adjust the vortexing velocity to your requirements.

- Powerful and quiet
- Variable speed up to 3000 rpm
   Switch from Auto to Touch mode
   Power On LED indicator
- 1-year limited warranty

#L-VM1000

(110-0228).....ea



### Rock-IT™ Jr. Blood-Tube Rocker

A versatile rocker for hematology (blood) tubes prior to counting; also excellent for vials and gentle mixing of injectables. Nearly silent while gently completing 25 mixing cycles per minuté. Platform is designed for 8 tubes and has a reversible pad that can be autoclaved. 1-year warranty.

Dimensions: 8.6" x 4.7" x 5"

#LTTR100. 8 Place

(159-0730) .....ea



### UNICO

### Rock-IT™ Blood Tube Rocker

A versatile rocker for hematology (blood) tubes prior to counting; also excellent for vials and gentle mixing of injectables. Nearly silent while gently completing 25 mixing cycles per minute. Platform is designed for 15 tubes and has a reversible pad that can be autoclaved. 1-year

Dimensions: 13.7" x 4.7" x 5"

#L-TTR200, 15 Place

(715-2650) .....ea

unico

### IP730 Series Microscopes

Designed to provide crisp, clear imaging on a versatile and ergonomic platform.
Available with 30° inclined binocular or trinocular heads. Features Seidentopf optical design with 55mm–75mm interpupillary distance adjustment Both heads rotate 360° and lock in place with a thumbscrew. Eyepieces are standard with 10X 22mm extra-wide field eyepieces for viewing comfort and clarity.

#IP730. Advanced Binocular Microscope. Achromat

Objectives: 4X, 10X, 40X, 100X Achromat (106-9122) .....ea

#IP730PL, Advanced Binocular Microscope, Plan Objectives: 4X, 10X, 40X, 100X Plan

Illumination: 6 V/20 W Halogen

(**109-9806**) .....ea



### #132000 Incubator

An excellent incubator for all in-office culture needs. Large capacity (0.8 cu ft) for busiest physicians' offices.

- Interior dimensions: 111/2"W x 101/2"H x 12"D; exterior dimensions: 13"W x 161/4"H x 14"D
- See-through acrylic plastic door
  Magnetic door latch
- · Pilot light and on/off switch
- · Adjustable nontip shelves
- Maximum temperature: 60°C (100-8790) .....ea





### ✓ HENRY SCHEIN®

### **LENS PAPER**

For fine lenses and optical instruments.

- · Special high wet strength
- Soft and dust-free
- · Will not scratch glass
- Ideal if liquid cleaner is used

4" x 6"

(100-8362) ......50/pkg

NOW, MORE THAN EVER, RELY ON

# **HENRY SCHEIN BRAND PRODUCTS**

SAVE WITHOUT COMPROMISING ON EXCELLENCE, WITH OUR OWN HIGH-QUALITY HENRY SCHEIN BRAND





### **✓** HENRY SCHEIN® **BRAND PRODUCTS**

### **GLASS MICROSCOPE SLIDES**

22 mm x 22 mm; thickness: 0.17 mm-0.25 mm

Plain Slide

(102-7277).....72/box

Frosted End Slide	
(102-5605)	.72/box

Cover Glass #2, 1-oz Package (102-6803) .....ea



### **Research-Grade Centrifuge Tubes**

- Acrylic spin rate is 5,000 rpm
- Polypropylene spin rate is 10,000 rpm
- All caps are made from high-density polyethylene (HDPE)

#6267. Acrylic. Non-sterile. Bulk (102-6540)......500/case

#6264. Polypropylene. Non-sterile. Bulk (280-5289) ......500/case

#6270, Acrylic, Sterile

(185-8503)......500/case Contains: 25 per bag.

#6273, Acrylic, Sterile (424-7963) ......500/case

Contains: 50 per rack. #6269, Polypropylene, Sterile

(424-6254) ......500/case

Contains: 25 per bag.

#6271, Polypropylene, Sterile

(243-5517)......500/case

Contains: 50 per rack.





### **Transfer Pipettes**

- Packaged in convenient dispenser boxes
- · Available in sterile and non-sterile versions
- Nontoxic and inert

#137135, 3 mL, Sterile, Graduated to 1 mL (119-0397) ......100/box

#137038, 5 mL, Sterile, Graduated to 1 mL (606-9661) ......400/case

#135038, 7 mL, Sterile, Graduated to 3 mL (990-1310)......400/case

#137035, 3 mL, Non-sterile, Graduated to 1 mL, Dispenser Box

(185-8544) ......500/box #136030, 4.5 mL, Non-sterile, Narrow Stem, Nongraduated, Dispenser Box (904-6934) ......500/box

#137030, 5 mL, Non-sterile, Graduated to 1 mL, Dispenser Box

(120-9438) ......500/box

#135030, 7 mL, Non-sterile, Graduated to 3 mL. Dispenser Box

(120-8751) ......500/box





### **Research-Grade Centrifuge Tubes**

50 mL. 30-mm x 116-mm Conical Centrifuge Tubes with Screw Caps

- Acrylic spin rate is 5000 rpm

Advince spiritate is 3000 rpm
 Polypropylene spiritate is 10,000 rpm
 All caps are made from high-density polyethylene (HDPE)

#6244, Polypropylene, Non-sterile, Bulk (424-5618).....500/case

#6242. Polypropylene, Non-sterile (424-9912).....500/case

Contains: 25 per rack.

#6240, Polypropylene, Non-sterile (424-8018)......500/case

Contains: 25 per bag.

#6243, Polypropylene, Sterile

(424-8999) ......500/case

Contains: 25 per rack.

#6241, Polypropylene, Sterile

(424-5191).....500/case

Contains: 25 per bag.





### **Urine Centrifuge Tubes**

PS, Graduated, 12 mL

• Tube dimensions: 20 mm x 115 mm

 Polystyrene spin rate: 1500 rpm #112017, Ames-Type

(120-9178).....500/pkg #112015, Kova Type

(424-0502) ......2000/case #112030, Sediment Bulb, Whale-Type

(424-0054) ......1500/case #113136, Blue Snap Caps for 12-mL Flared-Top Centrifuge Tubes

(**120-6724**) ......500/pkg #113137, Yellow Snap Caps with Sanitary

Grip for 12-mL Flared-Top Centrifuge Tubes

(120-8588) ......500/pkg





### Transport Tubes, Polypropylene

- Tubes are self-standing with conical
- Available with caps attached or packaged separately

#6101, 5 mL, Non-sterile, with Separate Screw Cap

(110-6203) ......1000/case

#6101C, 5 mL, with Attached Screw Cap, Non-sterile

(356-3969)......1000/case

#6101S, 5 mL, Sterile, with Attached Screw Cap (424-6233) ......500/case

#6102, 10 mL, Non-sterile, with Screw Cap (120-1971)......1000/case

#6102C, 10 mL, Non-sterile, with Attached Screw Cap

(120-6506). .....1000/case

#6102S, 10 mL, Sterile, with Attached Screw Cap

(945-1766)......500/case #6255, 50 mL, Non-sterile, with Screw Cap, Bulk

(424-6071) ......500/case #6256, 50 mL, Sterile, with Attached Screw

(999-6888) ......500/case



### Polvester-Tipped Swabs

Sterile

#25-800 D50, Aluminum Shaft, 5.5" x 0.035"

(120-0208) ......50/box

#25-800-1PD, Plastic Shaft, 6"

(120-0366) ......50/box





### **Specimen Transport Bag**

These specimen bags are available with a zipper closure or tamper-evident glue seal. The bags are printed with a biohazard symbol and feature an attached document pouch that keeps important paperwork safe.

#4918, 6" x 9" Bag with Document Pouch (106-5227)......1000/case



### **Universal Viral Transport**

#220222, Combo Kit Complete with 2 Swabs, 1 Plastic Shaft and 1 Aluminum Shaft, and 3-mL Tube (113-3040)......50/box

#220240, Swab Only

1 Aluminum and 1 Plastic Shaft (113-2874) ......100/box

54





### **Borosilicate Glass Culture Tubes**

Globe Scientific introduces a new line of high-quality borosilicate glass culture tubes. These tubes offer exceptional quality at great prices. They are neatly packaged in heavy-duty boxes of 250 and are shrink wrapped to keep them particle-free. The boxes are compact and designed to fit into most test-tube dispensers.

3 mL, 10 mm x 75 mm <b>(424-0058)</b>	250/box
5 mL, 12 mm x 75 mm <b>(424-0059)</b>	250/box
7 mL, 13 mm x 100 mm (116-7973)	1000/case
10 mL, 16 mm x 100 mm (424-0060)	250/box
12 mL, 16 mm x 125 mm (424-0062)	1000/case
15 mL, 16 mm x 150 mm (424-0061)	1000/case

# Urine Collection Containers and Cups The 4-oz (120-mL) containers are available in both sterile and non-sterile versions. Non-sterile containers have a frosted writing area for patient information. Sterile containers are individually wrapped and have a trilingual patient ID

containers are individually wrapped and have a trilingual patient ID label. They also have easy-to-turn, leak-resistant screw caps and molded graduation marks. The 6.5-oz (190-mL) cups come in both polypropylene and polystyrene versions and are available with or without paper lids.

#5912, Trilingual Patient ID La Individually Wrapped	abel, Sterile,
(606-6028)	100/case
#5913, 4-oz Container, Sterile Individually Wrapped	,
(868-0112)	100/case
#5915, 4-oz Container, Non-s Writing Space, Bulk	sterile, Etched
(642-6596)	500/case
#5916, 4-oz Container, Non-s Affixed Thermometer Strip	sterile, with
(413-4578)	500/case
#5920, 6.5-oz Specimen Colle Cup, PP	ection

(868-8877).....25/pack



### Microscope Slides, Glass

25 mm x 75 mm, 90° Ground Edges, Safety Corners

Premium-quality microscope slides are precleaned and ready to use. They feature finely-ground edges and slightly clipped corners for comfortable handling. Each box of slides contains a moisture-absorbent sheet to prevent sticking.

sheet to prevent sticking.	
<b>Plain (Colorless)</b> #1321, Unfrosted (868-0622)	.72/box
#1324, Frosted (803-6390)	.72/box
Frosted, Colors-Color on 1 End/1 S	Side
#1324W, White (868-6492)	.72/box
#1324N, Orange ( <b>424-0005</b> )144	40/case
#1324Y, Yellow ( <b>424-0007</b> )144	40/case
#1324P, Pink ( <b>424-0006</b> )144	40/case
#1324B, Blue ( <b>424-0003</b> )144	40/case
#1328, 2-sides, Frosted (868-0123)	.75/box

# Producers of exceptional quality laboratory supplies.

### **EXTENSIVE SELECTION - ONE TRUSTED SOURCE**

Microhematocrit Tubes

Pathology Containers

Serological Pipettes

Specimen Containers

Test Tube Racks


Pipettors & Pipette Tips

Microscope Slides & Accessories

Analyzer Consumables
Beakers, Bottles, Cylinders
Biohazard Specimen Bags
Blood Collection Products
Centrifuge Tubes
Collection Cups
Cytology Products
Drug Test Containers
ESR Systems

Drug Test Containers Test Tubes, Vials & Caps
ESR Systems Transfer Pipets
Formalin Filled Containers Transport Tubes
Funnels Urinalysis Products
Histology Cassettes Weighing Dishes

for more information or to request FREE sampless contact your Henry Schein representative





## MOHS LABORATORY SUPPLIES

MEDICAL | DERMATOLOGY

### **EMBEDDING MEDIA:**

Crvo O.C.T. (102-3734)4 oz. (116-1233) Paraplast® Plus 8x1kg /Case (106-4135)Tissue-Tek® O.C.T. Compound 4 07.

(Frozen Tissue Matrix)

### TISSUE FREEZING PREPARATIONS:

(104-3248)	Cytocool® II	11 oz.
(406-1248)	Flash Freeze Spray	11 oz.
(106-3127)	Neg-50 Clear	4oz.
(758-1368)	Neg-50 Blue	4 oz. x 2/Box
(106-7137)	Tissue Freezing Media Clear	4 oz. x 12/Cs.

### MOUNTING MEDIA:

Cvtoseal® 60 Low Viscosity (334-7213) 4 oz. (Toluene-Based For Use with Cover Glass)

(122-6901)Permount® 100 mL

(Toluene-Based For Hand Coverslippers)

(106-7123)SHUR/Mount™ 12 x 4 oz./Case

(Xylene-Based For Automated Coverslippers)

### **FIXATIVES:**

(108-3646)	Acetic Acid 10%	16 oz.
(120-1795)	Glacial Acetic Acid	16 oz.
(Also Known A	s Ethanoic Acid)	

(120-0443)Acetone 16 oz. (120-2827)Gallon Acetone (120-8724)Ethanol Alcohol 95% 200 Proof Gallon

(Also Known As Ethanol, Denatured Alcohol, Algarin Denatured, Ethyl Hydrate Denatured, Grain Alcohol Denatured)

(364-6848) Ethanol 100% 16 oz. (600-4289)Formalin Aceto-Alcohol Gallon (FAA Solution or Carnoy's Solution)

Formalin, 10% Buffered 10 Liter

(NBF or Neutral Buffered Formalin)

(110-5855)Pen Fix 4 Gallon/Case (106-3078)SafeFix II Gallon

(Formalin Substitute)

(234-6266)Zinfix® Zinc Formalin Fixative Gallon

### **DEHYDRATING REAGENTS:**

(113-7936)Flex 100% Alcohol 4 Gallons/Case (Premixed Blend of Methyl and Isopropyl Alcohol For H&E Staining I) (120-9026) Ammonium Hydroxide 500 mL (For H&E Staining)

(120-7599) Deionized Water (Di) 20 Liters/Box

### **BLUING AGENTS:**

(608-5631) Substitute for Scott's Tap Water Gallon (106-2876)Gallon Bluing Reagent (For H&E Staining)

### **HISTOLOGY STAINS:**

(106-3048)	Clarifier 2 Reagent	Gallon
(799-1780)	Clarifier 2 Reagent	Gallon
(For Improved Con	trast & Increased Cell Transparency)	

(831-6922)Diff-Quik Stain Set 3 x 500 mL (Also known as Romanowski Stain)

Eosin-Y Stain Solution 1% (113-4902)16 oz. (906-6070) **Eosin-Y Working Solution** Gallon (Also Eosin Y & Alcohol)

Eosin-Y Stain Solution 1% 32 oz. (652-8573) (Also Eosin Y & Alcohol)

(424-6155) Harris Hematoxylin Stain 500 ml (Mercury-Free)

(398-4749) Gills I Hematoxylin Gallon Gills II Hematoxvlin 2 Pint 4 x 16 oz./Cs. (106-2893)(487-7361)Gills III Hematoxylin 33.8 oz. (230-4808)Gills III Hematoxylin Gallon Tissue Marking Dye 6-Color Kit (106-7191)(925-7517)Toluidine Blue 0.1% 8 oz. Toluidine Blue 1% (925-4103) 16 oz.

(Also Called Diff-Quik or Romanowski Stain)

### **DECALCIFYING SOLUTION:**

(120-1248)**HCI-EDTA** Rapid Decal Gallon

### **CLEARING AGENTS:**

(111-8126)Clear-Rite 3™ Solvent 4 Gallons/Case (Can Be Used In All Tissue Processors and Automatic Stainers)

(357-0333) d-Limonine. Purified 4 Gallons/Case (Citrus-Based Xylene Substitute)

(120 - 3130)Safe-Clear Gallon (114-7751)Histo-Clear™ Gallon (120-4629)**Xylene** Gallon

### TRANSPORT & LAB SUPPLIES:

(369-4647)Petri Dish 4" (100mm) 25/Pkg. Petri Dish 4" (100mm) Glass (122-7441)12/Pkg. (115-1911) Prefilled Formalin Jar 96 x 20 mL/Case (115-1912)Prefilled Formalin Jar 96 x 40 mL/Case

We stock over 100,000 products including a full line of Histology and Cytology Containers Including: Tissue and Biopsy Cassettes, Slide Storage/Staining/Mailing; Sample containers, etc. If you don't see what you need for your practice, call your Henry Schein Sales Consultant or visit: www.henryschein.com/mohs



## Simport

### Histosette® I Tissue Processing/Embedding Cassette

Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. Cassettes can be opened or closed as often as necessary and always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassette. Anterior writing area at a 45° angle. The more acute angle makes cassettes more suitable to be used with some models of cassette-labeling instruments.

Casselles	1500/ca	se	
Specify: #M498-10, Lilac	(115-1678)	#M498-5, Yellow	(115-1684)
#M498-11, Orange	(115-1679)	#M498-6, Blue	
#M498-12, Aqua	(115-1680)	#M498-7, Peach	
#M498-2, White	(115-1681)	#M498-8, Tan	(115-1686)
#M498-3, Pink	(115-1682)	#M498-9, Gray	(115-1687)
#M498-4, Green	(115-1683)		. ,



## Simport

### Histosette® I Tissue Processing/Embedding Cassettes

Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. Just snap apart and conveniently lock the lid into the base of the cassette. They can be opened or closed as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassette for your convenience. Anterior writing area at a 30° angle. Not suitable for automated printers.

Cassettes	1500/c	case	
Specify:			
#M490-10, Lilac	(641-3848)	#M490-5, Yellow	(340-4586)
#M490-11, Orange	(630-5451)	#M490-6, Blue	
#M490-12, Aqua	(642-0225)	#M490-7, Peach	(115-1636)
#M490-2, White	(120-5194)	#M490-8, Tan	(115-1637)
#M490-3, Pink	······(600 <b>-</b> 9085)	#M490-9, Gray	(115-1639)
#M490-4, Green	(600-2992)	-	(**************************************

# Simport



### Histosette® II Tissue Processing/Embedding Cassettes

Specify:	
#M492-10, Lilac(	(115-1649)
#M492-11, Orange	(115-1650)
#M492-12, Aqua	(115-1652)
#M492-2, White	(115-1653)
#M492-3, Pink	(115-1654)
#M492-4, Green	(115-1655)
#M492-5, Yellow	(115-1656)
#M492-6, Blue	114-8888
#M492-7, Peach	(115-1657)
#M492-8, Tan	(115-1659)
#M492-9, Gray	(115-1660)





### Histosette® II Biopsy Processing/Embedding Cassette

Similar to Histosette® II Tissue Processing Cassette but specially designed to hold biopsy specimens during the embedding process, as well as in a storage cabinet. Anterior writing area is at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments

Cassettes	1500/case
Specify:	
#M493-10, Lilac	(115-1661)
#M493-11, Orange	(115-1662)
#M493-12, Aqua	
#M493-2, White	(115-1666)
#M493-3, Pink	
#M493-4, Green	
#M493-5, Yellow	
#M493-6, Blue	(115-1672)
#M493-7, Peach	
#M493-8, Tan	
#M493-9, Gray	
mivi ioo o, aray	(113-10/3)



# Simport

### Histosette® II Tissue Processing/Embedding Cassette

Most convenient for cassette-labeling instruments since covers are already removed from cassettes and are packaged separately in the case. Disposable plastic cassettes hold tissue specimens during processing and embedding, as well as in storage. Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

Cassettes	1000/case
Specify:	
#M485-10, Lilac	(115-1613)
#M485-11, Orange	
#M485-12, Aqua	(115-1615)
#M485-2, White	(115-1616)
#M485-3, Pink	
#M485-4, Green	
#M485-5, Yellow	(115-1619)
#M485-6, Blue	
#M485-7, Peach	
#M485-8, Tan	(115-1021)
#M485-9, Gray	
#IVI-100-3, GIAY	(115-1623)



# Simport

### Histosette<sup>®</sup> II Biopsy Processing/Embedding Cassette

	• • • • • • • • • • • • • • • • • • • •
Specify:	
#M486-10, Lilac	(115-1624)
#M486-11, Orange	(115-1625)
#M486-12, Aqua	(115-1627)
#M486-2, White	115-1628)
#M486-3, Pink	(115-1629)
#M486-4, Green	(115-1630)
#M486-5, Yellow	(115-1631)
#M486-6, Blue	(115-1632)
#M486-7, Peach	(115-1633)
#M486-8, Tan	(115-1634)
#M486-9, Gray	(115-1635)
	(113-1003)

# Simport

### Unisette™ Tissue Processing/Embedding Cassette

Specify:	
#M505-10, Lilac	(115-1709)
#M505-11, Orange	(115-1710)
#M505-12, Aqua	(115-1711)
#M505-2, White	(115-1712)
#M505-3, Pink	(115-1713)
#M505-4, Green	(115-1714)
#M505-5, Yellow	(115-1717)
#M505-6, Blue	
#M505-7, Peach	(115-1719)
#M505-8, Tan	(115-1720)
#M505-9, Gray	

# Simport

### Unisette™ Biopsy Processing/Embedding Cassette

Molded from a special high-density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap-latch and hinge-lock design of the Unisette™ prevents early separation of base and lid and allows for one-hand operation. Lids can be opened and closed as often as necessary and always relock securely without danger of specimen loss. Large labeling areas are located on two sides of the cassettes. Anterior writing area is at a 35° angle. Specially designed to hold biopsy specimens during the embedding process. Cassettes ......1500/case

Specify:	
#M506-10, Lilac	(115-1722)
#M506-11, Orange	(115-1723)
#M506-12, Aqua	(115-1724)
#M506-2, White	(115-1725)
#M506-3, Pink	(115-1726)
#M506-4, Green	(115-1727)
#M506-5, Yellow	(115-1728)
#M506-6, Blue	(115-1729)
#M506-7, Peach	(115-1730)
#M506-8, Tan	(115-1731)
#M506-9, Gray	(115-1732)
, ,	(



# Simport

### Micromesh™ Biopsy Processing/Embedding Cassettes

Specify:	
#M507-10, Lilac	(115-1733
#M507-11, Orange	
#M507-12, Aqua	(115-1735
#M507-2, White	(115-1736
#M507-3, Pink	
#M507-4, Green	(115-1738
#M507-5, Yellow	115-1739
#M507-6, Blue	109-7544
#M507-7, Peach	(115-1740
#M507-8, Tan	(115-1741
#M507-9, Gray	



# Simport

### **Embedding Cassette**

Specify:	
#M480-10, Lilac	(487-7323)
#M480-11, Orange	(114-9995)
#M480-12, Aqua	(114-9996)
#M480-2, White	(281-4837)
#M480-3, Pink	(769-5292)
#M480-4, Green	(120-0469)
#M480-5, Yellow	(114-9997)
#M480-6, Blue	(630-2264)
#M480-7, Peach	(114-9998)
#M480-8, Tan	(114-9999)
#M480-9, Gray	(115-1611)

**58** 12MS7336



# Simport

### **Biopsy Foam Pads**

Used to hold biopsies in place and prevent them from being lost during processing. Specially formulated foam is always verified for consistency throughout to achieve optimum solvent flow. Biopsy samples are sandwiched between two foam pads and placed either in tissue capsules or cassettes with metal or plastic lids. Resists temperatures from -40°C to 121°C. For use with cassettes.

30.2 mm x 25.4 mm x 2 mm (120-6138)......1000/case




# Simport

### **Disposable Base Molds**

Disposable base molds offer ease and convenience of operation. They are inexpensive enough to be discarded after use, yet strong enough to be reused. Specially formulated plastic material provides excellent thermal exchange. Smooth interior finish and rounded corners facilitate specimen removal. Available in the same variety of sizes as metal molds and can be used with the same styles or types of cassettes and embedding rings.

Base Molds ......1000/case

Specify: #M475-1 7 mm x 7 mm x 5 mm (114-9988)#M475-2 15 mm x 15 mm x 5 mm (114-9990)#M475-3 24 mm x 24 mm x 5 mm (114-9991)#M475-4 30 mm x 24 mm x 5 mm (114-9992)#M475-5 37 mm x 24 mm x 5 mm (114-9993)




# Simport

### **Modular Storage Drawer**

This drawer provides permanent storage and identification of up to 165 embedding rings or 250 cassettes per drawer. It is stackable to any convenient height, thanks to interlocking ridges on top and bottom. Made of high-impact resistant plastic. Identification labels included. 15%"W x 9%"D x 27%"H

#495-6, Storage Drawer (115-1677)......6/case



# Simport

### SlideFolder™ Microscope Slide Holder

The SlideFolder™ will hold up to 20 standard microscope slides 75 mm x 25 mm (3" x 1") and is made of 2 parts: a base holding the slides horizontally, offering numbered spaces for easy identification, and transparent doors that can either cover the slides or be swung behind the SlideFolder™ for space-saving purposes. The 2 transparent hinged covers offer a full view of each slide without having to remove it from its position in the SlideFolder™, and allow easy reading of ID labels with or without an optical bar-code reader. All units are stackable and take minimal space on laboratory tables or shelves.

Assorted, 2 of Each Color ...(115-1875)
#M750-20B, Blue ...............(115-1876)
#M750-20G, Green ............(115-1877)
#M750-20GY, Gray .............(115-1878)
#M750-20W, White .............(115-1879)
#M750-20Y, Yellow .................(115-1880)



# Simport

# StainTray<sup>™</sup> Slide Staining System The StainTray<sup>™</sup> has a black base made of

rough ABS plastic that can withstand a wide range of chemicals (avoid chlorinated hydrocarbons). It will accept up to 20 slides on four plastic rails covered with a polymer strip to hold slides perfectly even if tray is held at an angle. When humidity is needed, wells between rails will hold up to one mL of water securely without splashing. Middle wells will hold up to 2 mL each. Rails are raised not only to avoid water touching the slides but to make them more easily retrievable. The base will also hold excess stain solution dripping from the slides. Four rubber feet ensure greater base stability. Units are stackable for space-saving purposes.

Slide-Staining System.....ea
Specify:
#M920-1, Base with Clear
Lid .............(115-1904)
#M920-2, Base with Black



# Simport

### EasyDip™ Slide Staining Rack

Rack will hold up to 12 microscope slides with dimensions such as 75 mm x 25 mm, (3" x 1") and even 76 mm x 26 mm and with a thickness of 1.0 mm to 1.2 mm. The slides fit into individual slots for free passage and rapid drainage of staining fluids. Since they are placed vertically in the rack and not horizontally, their writing area will not be stained by the fluid, allowing their removal without the use of forceps. The lid completely covers the EasyDip" Slide Staining Jar to minimize spill and evaporation. A handle is permanently attached to the rack for easy insertion and removal of slides without the user's fingers touching the solution.

#M905-12DGY, Dark Gray (115-1897)......6/case



# Simport

### EasyDip™ Slide Staining System

The system has two components: a square staining jar and a 12-position vertical slide rack. Made of resistant acetal plastic, the staining jar will not break like most glass jars do. It will resist most staining agents, including alcohol and xylene (but not phenol, iodine or ferric chloride). The wide stable base offers greater stability, while the inside is recessed, allowing for a smaller reagent volume of only 80 mL. Easy to clean and no metals to corrode. Ideal for special stains, frozen sections and special processes.





### **Gram Stain Kits and Reagents**

Used to stain for microorganisms from cultures or specimens by the differential Gram method.

#212539, Gram Stain Kit (106-2056)......4/pkg Contains: 1–250-mL bottle ea of Gram crystal violet, Gram iodine (stabilized), Gram decolorizer,

& Gram safranin. #212524, Gram Stain Kit (645-2264)......4/pkg

Contains: 1-250-mL bottle of ea: Gram crystal violet, Gram iodine diluent with iodine concentrate, Gram isopropyl decolorizer & Gram safranin.

#212525, Gram Crystal Violet (116-6788)4/pk	g
#212527, Gram Decolorizer (224-5022)4/pkg	g
#212531, Gram Safranin (343-3274)4/pk	g
#212542, Gram Iodine (Stabilized) (613-3242)4/pkg	g
#212529, Gram Iodine (Unstabilized) (264-9099)4/pkg	g
#212536, Acridine Orange Stain (610-4654)ea	a



### Sedi-Stain™ Concentrated Reagent Stain

#421570, 12.5-mL Bottle with Dropper

Contains: Crystal violet 0.10%, safranine O 0.25%, ammonium oxalate 0.03%, ethyl alcohol (SD-3A) 10%, & water & stabilizers 89.62%.



A fast and effective biological and chemical stain remover that removes most medical and histological stains from hands, plastic, glassware, lab coats, floors, and laboratory equipment and surfaces. Phosphate-free,

and surfaces. Phospirate-free, biodegradable, and pH neutral (7–8). Easy to apply—simply rub in gently and stain is gone in a matter of seconds. Rinse with water.

### Erado-Sol removes:

- · Most medical and histological stains
- Wright stain Fresh blood
- lodine

Lugol's stain

#207-32, 32 oz

(154-6143) .....ea

### Isopropyl Alcohol

99% Purified

(120-8443) .....ea

### Alcohol-Anhydrous


**ACS Reagent** 

Suitable for histology. Contains 90.5% ethyl alcohol, 4.5% methyl alcohol, and 5% isopropyl alcohol.

#1200-1GL, 1-gal Bottle

(120-8724) .....ea





- Used as a styptic in its undiluted form
- Soluble in water or alcohol

8-mL Unidose Vial (812-0005) ......12/box



### Lugol's (Rx) Solution

- · Used for staining cervical and vaginal mucosa during colposcopic evaluations of abnormal Pap smears
- Less painful or irritating than alcoholic solutions when applied to cuts
- · Potassium iodide added to increase the solubility of the iodine

#0111-8, 8-mL Unidose Vials (812-0004) ......12/box



### **Reagents Solutions**

#400460, Acetone, 16 oz (309-2688) .....ea #400530, Potassium Hydroxide 10%, 4 oz (958-8858) .....ea

### Formaldehyde 37%, USP

(120-8716) .....ea



### **⇒**HealthLink

### **Hematology Solutions and Stains**

	QuickLink I Wright's Stain,	
	QuickLink I Wright's Stain, 3	
3 oz	QuickLink I Wright-Giemsa S	,
(103-0691	)	ea
	QuickLink III Fixative, 1 gal	ea
	QuickLink III Set, 16 oz (3)	ea
Compares	Eosinophil Stain, 16 oz s to: Hansel's <sup>®</sup> Stain <b>)</b>	ea
•	•	
	Wright Stain Kit, 8 oz )	3/box

### **→**Health**Link**

Microbiology Solutions and Stains
#400308, Chlorazol Black E, 2 oz (792-4876)ea
#400331, Decolorizer, 16 oz (108-9429)ea
#400328, Decolorizer Gram Stain, 1 gal (108-9430)ea
#400310, Gram Stain Set, 2 oz (922-5498)4/set
#400340, Gram Stain Set, Stabilized IOD, 8 oz (107-4395)4/set
#400344, lodine Stabilized, 8 oz (108-5887)ea
#400351, Lugol's Solution, 2-oz Dropper (106-0788)ea
#400350, Lugol's lodine, 16 oz (108-2277)ea
#400355, Lugol's Iodine Concentrate, 2-oz Dropper (108-2283)ea
#400353, Lugol's Iodine Concentrate, 16 oz (108-2293)ea

#400334, Safranin Replacement, 8 oz

(110-0522).....ea



Chocolate Agar Media BBL (101-3782) ......100/box

**⇔** BD

MacConkey Agar Plated (168-8219).....100/case

BD

TSA Plates with 5% Sheep's Blood (987-3908) ......100/case

**⇔**BD

Columbia CNA Agar with 5% Blood #221353 (101-3786) ......100/box

**BD** 

Culture Blood Agar 5% Media #221239

(**224-2775**) ......20/box

**BD** 

MacConkey Agar II

#221172 (**224-5653**)......20/box

BD

Mueller-Hinton Agar Plates, 150 mm

#221800 (**224-0871**) ......24/box

BD

Selective Strep Culture Plates

#221935 (**120-4787**)......10/pkg

**⊕** BD

Trypticase–Soy Agar with Blood

#221291 **(423-9799)**......100/box

BD


Trypticase Soy Broth Medium

#221815, 2 mL (101-3799)......100/pkg





**Plated Media** 

Taxo™ Differentiation Discs

#231552, A Discs
For differentiating group A streptococci.
(224-1056)......50/box
#231042, A Discs
For differentiating group A streptococci.

(693-1886) ......50/box

**BD**Sensi Disc

BD

Sulfamethoxazole with Trimethoprim

#231536, 23.75/1.25 mcg (100-6160).....50/box



The Henry Schein LabTeam is a highly trained group of Clinical Laboratory Specialists offering expert lab consultation and analysis to assist your practice in advancing patient care, maximizing revenues, and cost containment.

A broad portfolio of the leading laboratory systems available in the market, offering you the best solution for your individual practice needs.

As a distributor of a wide selection of the leading Laboratory Systems, a Henry Schein Clinical Laboratory Specialist offers you the best solution for your practice from the industry instead of just one brand.

Financial analysis for capital purchases of laboratory equipment providing select recommendations and advantageous acquisition programs.

A Henry Schein Clinical Lab Specialist will evaluate your current testing volumes to provide you a financial analysis to maximize your Laboratory solutions. In addition, Henry Schein Financial Services offers financial solutions to help you acquire the system your practice needs.

Testing and volume analysis to maximize your lab efficiency.

Our Clinical Lab Specialists can provide laboratory design and planning assistance to help you maximize your lab testing and revenues.  Professional consultation for CLIA certification.

The Henry Schein Lab Team can provide access to Laboratory consultants to help your in office achieve CLIA certification and OSHA compliance.

Knowledgeable, ongoing customer support from our dedicated Laboratory support group.

Our dedicated and experienced laboratory customer service will ensure that your ongoing ordering of consumables and reagents needs are efficiently satisfied.

Laboratory custom solutions are available based on customer needs.

Lot sequestering for reagents, calibrators, and controls are available for selected systems to keep your laboratory working smoothly.

Contact Henry Schein today for a comprehensive laboratory solution analysis with a Clinical Laboratory Specialist!

12MS7336

Please contact our Henry Schein Clinical Lab Specialists: Phone: 1-800-428-9200 • E-mail: labquestions@henryschein.com






### **Laboratory Connectivity**

Henry Schein's Clinical Laboratory Specialists will help you connect all of your Laboratory Equipment to a Laboratory Information System (LIS) to help you:

- Consolidate lab reports
- ❖ Archive and trend patients
- Store QC
- Seamlessly integrate to a third-party reference lab
- Access the results via the Web

The LIS can be connected to the EHR to simplify patient management, capture billings, and update patients' records.



Please contact our Henry Schein Clinical Lab Specialists: Phone: 1-800-428-9200 • E-mail: labquestions@henryschein.com

Call:
1-800-P-SCHEIN
(1-800-772-4346)


8am-9pm (et)

Fax: **1-800-329-9109**24 Hours

www.henryschein.com/medical

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products, excluding flut vaccine and Synagis." Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2012 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. Ffor informational purposes only Source: ReimbursementCodes com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance

12MS7336








# **Henry Schein DxRx Solutions**

### We're there when you need us with our DxRx Solutions Hotline


Do you have questions regarding the diagnostic, pharmaceutical and vaccine items that you purchase from Henry Schein? Need help finding a CPT code, package insert, or looking for brand to generic comparisons?

The Henry Schein DxRx Solutions Hotline provides information on the following topics for diagnostics, pharmaceuticals and vaccines:

- Indications and usage as provided on a package insert
- Package inserts
- Brand to generic comparisons
- CPT codes and J codes†



\*For informational purposes only; does not constitute legal advice. Source: www.immunize.org



• Storage Requirements

- Latest updates on supply issues, and much more







**1.877.523.SHOT** 8:30am–5:00pm, est 1.877.523.7468 dxrxsolutions@henryschein.com



