

Masimo® Product Catalog

2011



Masimo Corporation
40 Parker
Irvine, CA 92618
Tel: (949) 297-7000
Fax: (949) 297-7001
www.masimo.com

© 2011 by Masimo Corporation.

This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent. Product availability is subject to change without notice.

Masimo, Rainbow, SpCO, SpHb, SpMet, Cleanshield, DCI, LNCS, LNOP, Rad, Radical, , SatShare, SET, Signal Extraction Technology, Signal I.Q., Rad-5, Rad-8, Rad-9, RadLink are registered trademarks of Masimo Corporation.

M-LNCS, Blue, Radical-7, Rad-57, Pulse CO-Oximeter, Red, TF-I, RadNet, TrendCom and PPO+ are trademarks of Masimo Corporation.

The devices advertised in the catalog may not have been licensed in accordance with Canadian law.

Monitoring Devices

Masimo SET® Pulse Oximetry Devices				
Masimo PN	Description	UM	Product Image	Notes
9021	Radical® Docking Station, RDS-1 <i>Comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare cable connector.</i>	EA		
1311	Radical® Docking Station, RDS-1B <i>With Extended Battery comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare cable connector.</i>	EA		
1392	Radical® Docking Station, RDS-2 <i>Comes with choice of Power cord (International Only), and no communication capabilities.</i>	EA		
9023	Radical® Docking Station, RDS-3 limited communication capabilities <i>Comes with choice of power cord and serial output connector and Analog output / nurse call connector.</i>	EA		
9190	Rad-8®, Horizontal <i>Comes with either an LNOP or LNCS Patient Cable and a power cord.</i>	EA		
9191	Rad-8®, Vertical <i>Comes with either an LNOP or LNCS Patient Cable and a power cord.</i>	EA		

Masimo SET® Pulse Oximetry Devices

Masimo PN	Description	UM	Product Image	Notes
9230	Radical® (20-pin) Handheld, Color Screen <i>Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or Rainbow Patient Cable.</i>	EA		Does not include Rainbow Activation License. Device cannot be upgraded with Rainbow parameters.
9196	Rad-5® Pulse Oximeter <i>Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>	EA	 Shown in optional protective boot and with DCSC Sensor	Please refer to Handheld Device accessories for available colors.
9198	Rad-5® Pulse Oximeter with Sample Pack <i>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9112	Rad-5® Pulse Oximeter (International Only) <i>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9197	Rad-5v™ Pulse Oximeter <i>Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>	EA	 Shown in optional protective boot and with DCSC Sensor	Please refer to Handheld Device accessories for available colors.
9199	Rad-5v™ Pulse Oximeter with Sample Pack <i>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9113	Rad-5v™ Pulse Oximeter (International Only) <i>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			

Masimo SET® Pulse Oximetry Devices

Masimo PN	Description	UM	Product Image	Notes
9158	Bitmos 901-M Pulse Oximeter with Masimo SET® <i>Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters. Provides continuous monitoring with alarms, internal memory and battery or mains power operation. Includes Instructions for Use and either a LNOP DCSC Sensor or a LNCS Patient Cable with either an Adult LNCS DCI Sensor, LNCS DCIP Sensor, LNCS YI Sensor, LNCS TC-I Sensor or LNCS Sensor Sample pack.</i>	EA		European Distribution Only.
9160	Bitmos 901-ML Pulse Oximeter with Masimo SET® <i>Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters. Provides continuous monitoring with alarms, internal memory and battery (Lithium Ion-rechargeable) or mains power operation. Includes Instructions for Use and either a LNOP DCSC Sensor or a LNCS Patient Cable with either an Adult LNCS DCI Sensor, LNCS DCIP Sensor, LNCS YI Sensor, LNCS TC-I Sensor or LNCS Sensor Sample pack.</i>			

Masimo Rainbow SET® Devices				
Masimo PN	Description	UM	Product Image	Notes
9086	Radical-7™ Handheld, Color Screen Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNCP, Red LNCS or Rainbow Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Desat Index 3D alarm and PI Delta 3D alarm.	EA		Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers. The Radical-7 is compatible with all Radical Docking Stations.
9153	Radical-7™ Handheld, Blue Screen Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNCP, Red LNCS or Rainbow Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Desat Index 3D alarm and PI Delta 3D alarm.	EA		
9226	Rad-87™, Horizontal Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters. Includes either Red LNCS, Red LNCP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRA). Comes with serial output connector and Analog output / Nurse call connector.	EA		Rainbow Acoustic Monitoring capable. Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers.
9227	Rad-87™, Horizontal with Wireless Radio Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters and integrated wireless radio. Includes either Red LNCS, Red LNCP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRA). Comes with serial output connector and Analog output / Nurse call connector.			

Masimo Rainbow SET® Devices

Masimo PN	Description	UM	Product Image	Notes
9228	Rad-87™, Vertical Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters. Includes either a Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa). Comes with serial output connector and Analog output / Nurse call connector.	EA		Rainbow Acoustic Monitoring capable. Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers.
9229	Rad-87™, Vertical with Wireless Radio Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters and integrated wireless radio. Includes either a Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa). Comes with serial output connector and Analog output / Nurse call connector.			
9221	Pronto-7™ Spot Check Pulse CO-Oximeter Includes SpO2, Pulse Rate, Perfusion Index (PI) and SpHb parameters. Rainbow 4D DC Sensor not included.	EA		Rainbow 4D DC Reusable Sensor must be ordered separately.

Masimo Rainbow SET® Devices

Masimo PN	Description	UM	Product Image	Notes
9216 - U 9216	Rad-57™ Configurable Pulse CO-Oximeter <i>Includes SpO2, Pulse Rate and Perfusion Index (PI) parameters.</i> <i>Rainbow Sensors not included.</i>	EA	 <i>Shown with Rainbow Adult Reusable Sensor and in optional protective boot</i>	<i>Please refer to Handheld Device accessories for available colors.</i> <i>9216-U is for SpO2 only Devices.</i> <i>Hb can be included on 9216.</i>

Masimo Rainbow SET® Device Options

Masimo PN	Description	UM	Product Image	Notes
9081	Field Upgrader Used to load optional Rainbow parameters and features onto compatible monitors. Current available options include: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Desat Index 3D alarm and PI Delta 3D alarm. Each optional parameter and feature sold separately.	EA		Maximum quantity of 25 of any individual parameter or feature per Field Upgrader. Includes: 1 Field Upgrader, 1 to 5 Operator's Manuals, parameter upgrade indicator sticker(s), and postage paid envelope to return Field Upgrader to Masimo.
2630	Total Hemoglobin (SpHb) parameter Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	
2296	Carboxyhemoglobin (SpCO) parameter Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	
2297	Methemoglobin (SpMet) parameter Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	
2313	Pleth Variability Index (PVI) Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	
2631	Carboxyhemoglobin (SpCO), Fixed Term License Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	
2632	Methemoglobin (SpMet), Fixed Term License Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.	EA	No Image	

Masimo Rainbow SET® Device Options

Masimo PN	Description	UM	Product Image	Notes
2633	Total Hemoglobin (SpHb), Fixed Term License <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2636	PI Delta 3D Alarm, Fixed Term License <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2637	Desat Index 3D Alarm, Fixed Term License <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2762	Full Rainbow® License (SpHb, PVI, SpCO, SpMet) <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2946	Rainbow® Activation License <i>Allows Radical (20-pin) part number 9230 to be upgraded to Rainbow parameters.</i>	EA	No Image	
2940	Acoustic Respiration (RRa) <i>Can be ordered on a Rad87 (part number 9226, 9227, 9228, 9229) or on a Field Upgrader (part number 2926) to upgrade a compatible monitor in the field</i>	EA	No Image	
2340	Demonstration Tool <i>Used to exercise monitors alarm default setting and demonstrate optional parameters (SpHb, SpCO, SpMet and PVI) when optional parameters are installed on a compatible monitor.</i>	EA		

Philips® Compatible Devices

Masimo PN	Description	UM	Product Image	Notes
2267	Masimo SET® IntelliVue® Module <i>Note: Patient Cable must be ordered separately.</i>	EA		
1881	Masimo SET® VueLink® Module <i>Comes with VueLink cable.</i>	EA		
1841	VueLink® Cable <i>For use with Masimo SET VueLink Module, 3 ft.</i>	EA		

Philips®, IntelliVue® and VueLink® are registered trademarks of Philips Corporation.

Patient SafetyNet™ System				
Masimo PN	Description	UM	Product Image	Notes
9095	Masimo Patient SafetyNet™ System <i>Configurable for customer needs.</i>	System		
2146	Masimo Patient SafetyNet™ Compatible Radical Docking Station, RDS-1 <i>Comes with power cord and communication capabilities.</i> <i>Note: SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA		
2147	Masimo Patient SafetyNet™ Compatible Radical Docking Station, RDS-3 <i>Comes with power cord and communication capabilities.</i> <i>Note: SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA		
14422	3Com Workgroup Bridge	EA		(1) Required per Pulse Oximeter.
15007	Colubris Workgroup Bridge	EA		(1) Required per Pulse Oximeter; alternative for 3Com Workgroup Bridge.
14423	3Com Ethernet Workgroup Bridge Antenna <i>Optional accessory to 14422.</i>	EA	No Image	
15275	NET232-DTE (male) RS232 Serial to Ethernet Intelligent Cable Adapter	EA		
14530	Serial Null Modem DB-9 male to DB-9 male Cable	EA	No Image	


Patient SafetyNet™ System

Masimo PN	Description	UM	Product Image	Notes
32054	Bracket A Main <i>Includes bracket and hardware for mounting the 3Com Radio to the Radical-7 Pulse Oximeter.</i>	EA		
32055	Patient SafetyNet™ Bracket B Main <i>Includes bracket and hardware for mounting the 3Com Radio to the Rad-8 Pulse Oximeter.</i>	EA		
14256	Cisco 1242 AP Wireless Access Point <i>Includes Cisco Access Point, power cables and basic mount. Must also include (1) 14260 Antenna for each access point.</i>	EA		
14260	Cisco ANT-5353 Antenna <i>2 dB Wireless Antenna for 14256 Cisco Wireless Access Point (802.11 b/g).</i>	EA		(1) 14260 or 14261 Antenna required per 14256 Cisco Wireless Access Point.
14261	Cisco ANT-2012 Antenna <i>6 dB Wireless Antenna for 14256 (802.11 b/g).</i>	EA		May be required based on Wireless Site Survey.
14572	Cisco ANT-5135D-R Antenna <i>(802.11 a).</i>	EA		
14573	Cisco ANT-5145V-R Antenna <i>Use with 14256 (802.11a Networks).</i>	EA		
14574	Cisco ANT-5195P-R Antenna <i>Use with 14256 (802.11a Networks).</i>	EA		
14858	SMC Networks SMC6128PL2 TigerSwitch 10/100 Ethernet Switch	EA		
33612	Masimo Patient SafetyNet™ Server	EA		

Patient SafetyNet™ System

Masimo PN	Description	UM	Product Image	Notes
14804	HP R1500 G2 UPS	EA		
14857	XRackPro2 Rackmount 4U Cabinet, Noise Reduction Enclosure	EA		
14868	Terminal PC	EA		
16170	HP Touchsmart Terminal PC	EA		
14424	Paging Transmitter <i>Frequency Range (450 - 470 MHz).</i>	EA		(1) required per Floor in Hospital
14253	Pager 456 - 460 MHz <i>Pager tunable in the 456 - 460 MHz Range.</i>	EA		
14254	Pager 460 - 464 MHz <i>Pager tunable in the 460 - 464 MHz Range.</i>			
14255	Pager 464 - 468 MHz <i>Pager tunable in the 464 - 468 MHz Range.</i>			




Sensors & Cables

Rainbow® Acoustic Monitoring Sensor and Cables (for use with Masimo Rainbow Acoustic Monitoring Devices with RRA parameters)				
Masimo PN	Description	UM	Product Image	Notes
2759	Acoustic Respiration Sensor Adult Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.	Box		
2757	Dual-Channel Rainbow® Acoustic Monitoring Cable Channel 1 - Accepts any Masimo 20-pin cable or direct-connect sensor Channel 2 - Accepts the RAM Patient Cable only (PN 2758).	EA		
2758	Acoustic Respiration Patient Cable Rainbow 14-pin Patient Cable, 10 ft. For use with Acoustic Respiration sensors.	EA		

Rainbow® Adhesive Sensors (for use with Masimo Rainbow SET Devices with SpHb, SpO₂, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2416	Rainbow® R1 25 Adult Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex.	Box	 A white adhesive sensor with a long, thin white cable. The cable has the word "rainbow" printed on it. The sensor head is white with a red Masimo logo and a black square. The adhesive pad is white with a red border.	
2417	Rainbow® R1 20 Pediatric Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box	 A white adhesive sensor with a long, thin white cable. The cable has the word "rainbow" printed on it. The sensor head is white with a red Masimo logo and a black square. The adhesive pad is white with a red border.	
2414	Rainbow® R1 25L Adult Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex. Use replacement tape 2623.	Box	 A white adhesive sensor with a long, thin white cable. The cable has the word "rainbow" printed on it. The sensor head is white with a red Masimo logo and a black square. The adhesive pad is white with a red border.	
2415	Rainbow® R1 20L Infant Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement tape 2624.	Box	 A white adhesive sensor with a long, thin white cable. The cable has the word "rainbow" printed on it. The sensor head is white with a red Masimo logo and a black square. The adhesive pad is white with a red border.	




Rainbow® Starter Packs (for use with Masimo Rainbow SET Devices with SpHb, SpO₂, and SpMet parameters) Starter Packs ordered with 90 day license pricing have a 3 pack minimum order requirement

Masimo PN	Description	UM	Product Image	Notes
2825	Rainbow® R1 25 Adult Starter Pack 4 Rainbow R1 25 SpHb, SpO ₂ , SpMet Sensors	EA		For purchase only, not to be used for Evals or Demos.
2826	Rainbow® R1 20 Pediatric Starter Pack 4 Rainbow R1 20 SpHb, SpO ₂ , SpMet Sensors	EA		
2823	Rainbow® R1 25L Adult Starter Pack 4 Rainbow R1 25L SpHb, SpO ₂ , SpMet Sensors	EA		For purchase only, not to be used for Evals or Demos.
2824	Rainbow® R1 20L Infant Starter Pack 4 Rainbow R1 20L SpHb, SpO ₂ , SpMet Sensors	EA		





Rainbow® Adhesive Sensors (for use with Masimo Rainbow SET Devices with SpO₂, SpCO, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2221	Rainbow® R25 Adult Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex.	Box	 A white, Y-shaped adhesive sensor with a central cable labeled 'rainbow'. The cable has a Masimo logo at the top and a small black square at the bottom. The sensor is designed for adult use.	
2222	Rainbow® R20 Pediatric Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box	 A white, Y-shaped adhesive sensor with a central cable labeled 'rainbow'. The cable has a Masimo logo at the top and a small black square at the bottom. The sensor is designed for pediatric use.	
2219	Rainbow® R25-L Adult / Neonatal Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight <3 kg or >30 kg This product does not contain natural rubber latex. Use replacement tape 2623.	Box	 A white, Y-shaped adhesive sensor with a central cable labeled 'rainbow'. The cable has a Masimo logo at the top and a small black square at the bottom. The sensor is designed for adult or neonatal use.	
2220	Rainbow® R20-L Infant Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement tape 2624.	Box	 A white, Y-shaped adhesive sensor with a central cable labeled 'rainbow'. The cable has a Masimo logo at the top and a small black square at the bottom. The sensor is designed for infant use.	


Rainbow® ReSposable Sensors (for use with Masimo Rainbow SET Devices with SpHb, SpO₂, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2743	Rainbow® ReSposable™ R2-25 Sensor System Adult R2-25a disposable optical sensors & R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 10-R2-25a /box (DOS) 2-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	Box		
2744	Rainbow® ReSposable™ R2-20 Sensor System Pediatric R2-20a disposable optical sensors & R2-20r reusable sensors (SpO ₂ , SpHb, SpMet) Single patient use 10-R2-20a /box (DOS) 2-R2-20r /box (ROS) Non-sterile, Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2753	Rainbow® ReSposable™ R2-25a Sensors Adult R2-25a disposable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 25-R2-25a /box (DOS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex. For use with R2-25r	Box		
2755	Rainbow® ReSposable™ R2-20a Sensors Pediatric/Slender Digit R2-20a disposable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 25-R2-20a /box (DOS) Non-sterile, Weight 10 – 50 kg This product does not contain natural rubber latex. For use with R2-20r	Box		
2754	Rainbow® ReSposable™ R2-25r Sensors Adult R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex. For use with R2-25a	EA		



Rainbow® ReSposable Sensors (for use with Masimo Rainbow SET Devices with SpHb, SpO₂, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2756	Rainbow® ReSposable™ R2-20r Sensors Pediatric/Slender Digit R2-20r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-20r /box (ROS) Non-sterile, Weight 10 – 50 kg This product does not contain natural rubber latex. For use with R2-20a	EA		
2961	Rainbow® ReSposable™ R2-25 Sample Pack Adult R2-25a disposable optical sensors & R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 3-R2-25a /box (DOS) 1-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	EA		
2962	Rainbow® ReSposable™ R2-20 Sample Pack Pediatric R2-20a disposable optical sensors & R2-20r reusable sensors (SpO ₂ , SpHb, SpMet) Single patient use 3-R2-20a /box (DOS) 1-R2-20r /box (ROS) Non-sterile, Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
2981	Rainbow® ReSposable™ R2-25 Starter Pack Adult R2-25a disposable optical sensors & R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-25a /box (DOS) 1-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	EA		
2982	Rainbow® ReSposable™ R2-20 Starter Pack Pediatric R2-20a disposable optical sensors & R2-20r reusable sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-20a /box (DOS) 1-R2-20r /box (ROS) Non-sterile, Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		


Rainbow® Reusable Spot Check Sensors (for use with Masimo Rainbow SET Devices with SpHb, SpO2, SpMet parameters; SpHb parameter required)

Masimo PN	Description	UM	Product Image	Notes
2333	Rainbow® DCI SC-360 Adult Reusable Sensor, 3ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device	EA		
2334	Rainbow® DCIP SC-360 Pediatric Reusable Sensor, 3 ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device	EA		

Rainbow® Direct Connect Spot Check Sensors (for use with Masimo Rainbow SET Devices with SpHb parameter)

Masimo PN	Description	UM	Product Image	Notes
2241	Rainbow® DC-3 SC-360 Adult Reusable Direct Connect Sensor, 3 ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device and 20 pin instrument connector.	EA		
2242	Rainbow® DC-12 SC-360 Adult Reusable Direct Connect Sensor, 12 ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device and 20 pin instrument connector.			
2402	Rainbow® DCP-3 SC-360 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device and 20 pin instrument connector.	EA		
2403	Rainbow® DCP-12 SC-360 Pediatric Reusable Direct Connect Sensor, 12 ft. (SpHb, SpMet, SpO ₂) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex. 360 Spot checks Requires SpHb parameter on device and 20 pin instrument connector.			
2953 U.S. 2953-OU International	Rainbow® 4D DC SC-500 Adult Reusable Sensor, 3ft. (SpHb, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. 500 Spot checks To be used with Pronto-7 only.	EA		Comes with 500 spot check tests preloaded on the device. Additional Spot Check tests can be loaded in the sensor via the eCommerce Site.
2954 U.S. 2954-OU International	Rainbow® 4D DC SC-1000 Adult Reusable Sensor, 3ft. (SpHb, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. 1000 Spot checks To be used with Pronto-7 only.	EA		Comes with 1000 spot check tests preloaded on the device. Additional Spot Check tests can be loaded in the sensor via the eCommerce Site.

Rainbow® Reusable Sensors (for use with Masimo Rainbow SET Devices with SpCO parameter)

Masimo PN	Description	UM	Product Image	Notes
2696	Rainbow® DCI Adult Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2697	Rainbow® DCIP Pediatric Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.	EA		




Rainbow® Direct Connect Sensors (for use with Masimo Rainbow SET Devices with SpCO, SpO₂, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2201	Rainbow® DCI-dc3 Adult Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2407	Rainbow® DCI-dc8 Adult Reusable Direct Connect Sensor, 8 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			
2202	Rainbow® DCI-dc12 Adult Reusable Direct Connect Sensor, 12 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			
2069	Rainbow® DCIP-dc3 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.	EA		
2640	Rainbow® DCIP-dc8 Pediatric Reusable Direct Connect Sensor, 8 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.			
2070	Rainbow® DCIP-dc12 Pediatric Reusable Direct Connect Sensor, 12 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.			




Rainbow® Patient Cables (for use with Masimo Rainbow SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2405	Rainbow® RC-1 <i>Rainbow 20-pin Patient Cable, 1 ft</i> <i>For use with Rainbow and M-LNCS sensors</i>	EA		
2406	Rainbow® RC-4 <i>Rainbow 20-pin Patient Cable, 4 ft</i> <i>For use with Rainbow and M-LNCS sensors</i>	EA		
2404	Rainbow® RC-12 <i>Rainbow 20-pin Patient Cable, 12 ft</i> <i>For use with Rainbow and M-LNCS sensors</i>			

ReSposable Sensors (for use with Masimo SET Devices with SpO2 and SpMet parameters)				
Masimo PN	Description	UM	Product Image	Notes
2863	ReSposable S2-25a (S-DOS™) Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. For use with M-LNCS or LNCS S-ROS 25/S2-25r	Box		
2865	ReSposable S2-20a (S-DOS™) Pediatric SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex. For use with M-LNCS or LNCS S-ROS 20/S2-20r	Box		
2864	ReSposable M-LNCS™ S2-25r (S-ROS™) Adult SpO ₂ reusable optical sensor (ROS) Reusable 1 /box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. For use with S-DOS 25/S2-25a	EA		
2866	ReSposable M-LNCS™ S2-20r (S-ROS™) Pediatric SpO ₂ reusable optical sensor (ROS) Reusable 1 /box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex. For use with S-DOS 20/S2-20a	EA		
2949	ReSposable LNCS® S2-25r (S-ROS™) Adult SpO ₂ reusable optical sensor (ROS) Reusable 1 /box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex. For use with S-DOS 25/S2-25a	EA		
2950	ReSposable LNCS® S2-20r (S-ROS™) Pediatric SpO ₂ reusable optical sensor (ROS) Reusable 1 /box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex. For use with S-DOS 20/S2-20a	EA		

LNCS® Adhesive Sensors




Masimo PN	Description	UM	Product Image	Notes
1859	LNCS® Adtx Adult SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.	Box		
2317	LNCS® Adtx-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.			
1860	LNCS® Pdtx Pediatric SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2318	LNCS® Pdtx-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.			
1861	LNCS® Inf-L Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex.	Box		
1862	LNCS® Neo-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex.	Box		
1901	LNCS® NeoPt-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex.	Box		


LNCS® Adhesive Sensors with Replaceable Tapes

Masimo PN	Description	UM	Product Image	Notes
2328	LNCS® Inf Infant SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2319	LNCS® Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.			
2329	LNCS® Neo Neonatal SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2320	LNCS® Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.			
2330	LNCS® NeoPt Neonatal SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2321	LNCS® NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.			
2331	LNCS® NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

LNCS® Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting

Masimo PN	Description	UM	Product Image	Notes
2411	LNCS® Trauma Adult Adhesive Sensors Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2412	LNCS® Newborn Neonatal Neonatal Adhesive Sensors Single Patient Use 20/box, Sterilizable Weight < 3 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2413	LNCS® Newborn Infant/Pediatrics Infant / Pediatric Adhesive Sensor Single Patient Use 20/box, Sterilizable Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

LNCS® Sterile Sensors				
Masimo PN	Description	UM	Product Image	Notes
2601	LNCS® Adept-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight > 30 kg This product does not contain natural rubber latex.	Box		Note: The sensors are <u>not</u> intended to be used for evaluations or carried as sales samples. If there is a need for samples of sterile sensors, contact your manager for approval and include the e-mail with FRM-1386.
2602	LNCS® Pdt-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2603	LNCS® Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2604	LNCS® Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2605	LNCS® NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2606	LNCS® NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

LNCS® Reusable Sensors				
Masimo PN	Description	UM	Product Image	Notes
1863	LNCS® DCI Adult Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1864	LNCS® DCIP Pediatric Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
2653	LNCS® DBI Adult Soft Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1895	LNCS® TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2258	LNCS® YI Multisite Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		
1896	LNCS® TF-I Adult Transflectance Forehead, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		

LNCS® Red Patient Cables (for use with Masimo Rainbow SET Devices)

Masimo PN	Description	UM	Product Image	Notes
2365	Red LNC-01 <i>LNCS 20-pin SpO₂ Spot Check Patient Cable.</i>	EA		
2055	Red LNC-04 <i>LNCS 20-pin SpO₂ Patient Cable, 4 ft.</i>	EA		
2056	Red LNC-10 <i>LNCS 20-pin SpO₂ Patient Cable, 10 ft.</i>			
2057	Red LNC-14 <i>LNCS 20-pin SpO₂ Patient Cable, 14 ft.</i>			

LNC™ Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2364	LNC-1 <i>LNCS Series Patient Cable.</i>	EA		
2017	LNC-4 <i>LNCS Series Patient Cable, 4 ft.</i>	EA		
1814	LNC-10 <i>LNCS Series Patient Cable, 10 ft.</i>			
2013	LNC-14 <i>LNCS Series Patient Cable, 14 ft.</i>			
2021	LNC-4-Ext <i>LNC Extension Cable, 4 ft.</i>	EA		
2281	LNC MP10 <i>LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 LNCS Dual Key Patient Cable, 10 ft.</i>	EA		

LNC™ Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
1902	LNC-10 P8 <i>LNCS to IntelliVue Masimo SET Patient Cable, 10 ft.</i>	EA		
2430	LNC-4-GE <i>LNCS to GE Patient Cable, 4ft.</i>	EA		
2016	LNC-10-GE <i>LNCS to GE Patient Cable, 10ft.</i>			
2432	LNC-SL-10 <i>LNCS to Spacelabs Patient Cable, 10 ft.</i>	EA		
2285	LNC-04-DB9 <i>LNCS to DB9 Connector Patient Cable, 4 ft.</i>	EA		
2286	LNC-10-DB9 <i>LNCS to DB9 Connector Patient Cable, 10 ft.</i>			
2287	LNC-14-DB9 <i>LNCS to DB9 Connector Patient Cable, 14 ft.</i>			



LNCS® Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
1816	LNCS® to PC Adapter Cable <i>LNCS Series to LNOP PC Adapter Cable, 18 in.</i>	EA		
2263	LNC™ CMS <i>LNCS Series to Philips 12-pin Instrument Adapter Cable, 3 m.</i>	EA		
2264	LNC™ MAC-GE <i>LNCS Series to GE Adapter Cable, 3 m.</i>	EA		
2266	LNC™ MAC-SL <i>LNCS Series to Spacelabs Adapter Cable, 3 m.</i>	EA		

LNCS® Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
2336	LNC™ MAC-SL2 <i>LNCS Series to Spacelabs Adapter Cable, 3 m.</i>	EA		
2268	LNC™ MAC-180 <i>LNCS Series to N-180 Adapter Cable, 10 ft.</i>	EA		
2269	LNC™ MAC-395 <i>LNCS Series to N-395 Adapter Cable, 10 ft.</i>	EA		
2921	LNC-MAC-4 P8 <i>LNCS Series to Philips FAST-SpO2 Adapter Cable, 4 ft.</i>	EA		
2922	LNC-MAC-10 P8 <i>LNCS Series to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>			

M-LNCS™ Adhesive Sensors

Masimo PN	Description	UM	Product Image	Notes
2508	M-LNCS™ Adtx Adult SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.	Box		
2509	M-LNCS™ Adtx-3 Adult SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.			
2510	M-LNCS™ Pdtx Pediatric SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2511	M-LNCS™ Pdtx-3 Pediatric SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.			

M-LNCS™ Adhesive Sensors with Replaceable Tapes

Masimo PN	Description	UM	Product Image	Notes
2512	M-LNCS™ Inf Infant SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2513	M-LNCS™ Inf-3 Infant SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.			
2514	M-LNCS™ Neo Neonatal/Adult SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2515	M-LNCS™ Neo-3 Neonatal/Adult SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.			
2516	M-LNCS™ NeoPt Neonatal SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2517	M-LNCS™ NeoPt-3 Neonatal SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.			
2518	M-LNCS™ NeoPt-500 Neonatal SpO ₂ Non-Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

M-LNCS™ Adhesive Sensors with Replaceable Tapes

Masimo PN	Description	UM	Product Image	Notes
2522	M-LNCS™ Blue Neonatal/Infant/Pediatric SpO ₂ Adhesive Sensor Single Patient Use 10/box, Non-Sterile Weight 2.5 – 30kg This product does not contain natural rubber latex. Use replacement tape 2822.	Box		

M-LNCS™ Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

Masimo PN	Description	UM	Product Image	Notes
2521	M-LNCS™ Trauma Adult SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2519	M-LNCS™ Newborn Neonatal Neonatal SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight < 3 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2520	M-LNCS™ Newborn Infant/Pediatric Infant/Pediatric SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

M-LNCS™ Reusable Sensors				
Masimo PN	Description	UM	Product Image	Notes
2501	M-LNCS™ DCI Adult Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2502	M-LNCS™ DCIP Pediatric/Slender Digit Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
2503	M-LNCS™ TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2505	M-LNCS™ YI Multisite Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		
2504	M-LNCS™ TF-I Adult Transflectance Forehead, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		

M-LNC™ Patient Cables				
Masimo PN	Description	UM	Product Image	Notes
2523	M-LNC-1 <i>M-LNCS Series Patient Cable, 1 ft.</i>	Box		
2524	M-LNC-4 <i>M-LNCS Series Patient Cable, 4 ft.</i>	Box		
2525	M-LNC-10 <i>M-LNCS Series Patient Cable, 10 ft.</i>			
2526	M-LNC-14 <i>M-LNCS Series Patient Cable, 14 ft.</i>			
2831	M-LNC-MP10 <i>M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 M-LNC Dual Key Patient Cable, 10 ft.</i>	Box		
2532	M-LNC-10 P8 <i>M-LNCS to IntelliVue Masimo SET M-LNC Patient Cable, 10 ft.</i>	Box		
2655	M-LNC-4-GE <i>M-LNCS to GE Patient Cable, 4 ft.</i>	Box		
2533	M-LNC-10-GE <i>M-LNCS to GE Patient Cable, 10 ft.</i>			

M-LNC™ Patient Cables





Masimo PN	Description	UM	Product Image	Notes
2714	M-LNC-SL-10 <i>M-LNCS to Spacelabs Patient Cable, 10 ft.</i>	Box		
2534	M-LNC-04-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 4 ft.</i>	Box		
2535	M-LNC-10-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 10 ft.</i>			
2536	M-LNC-14-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 14 ft.</i>			

M-LNCS™ Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
2542	M-LNC™ MAC-SL2 <i>M-LNC Series to Spacelabs Adapter Cable, 2 ft.</i>	EA		
2541	M-LNC™ MAC-SL <i>M-LNC Series to Spacelabs Adapter Cable, 3 m.</i>	EA		
2537	M-LNC™ MAC-180 <i>M-LNCS Series to N-180 Adapter Cable, 10 ft.</i>	EA		
2538	M-LNC™ MAC-395 <i>M-LNCS Series to N-395 Adapter Cable, 10 ft.</i>	EA		
2748	M-LNC™ MAC P8 <i>M-LNCS to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>	EA		
2539	M-LNC™ MAC CMS <i>M-LNCS Series to Phillips 12-pin Instruments Adapter Cable, 3 m.</i>	EA		

M-LNCS™ Adapter Cables				
Masimo PN	Description	UM	Product Image	Notes
2654	M-LNCS™ to LNC Adapter Cable <i>M-LNCS Series to LNC Patient Cable, 1.5 ft.</i>	EA		
2750	M-LNCS™ to LNC-MAC Adapter Cable <i>M-LNCS Series to LNC-MAC Adapter Cable, 1.5 ft .</i>	EA		

LNOP® Adhesive Sensors

Masimo PN	Description	UM	Product Image	Notes
1001	LNOP® Adt Adult SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	Box		
1025	LNOP® Pdt Pediatric SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
1800	LNOP® Inf-L Infant SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use 1801 replacement tape	Box		
1798	LNOP® Neo-L Neonatal SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use 1799 replacement tape	Box		
1651	LNOP® NeoPt-L Neonatal SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight < 1 kg This product does not contain natural rubber latex. Use 1926 replacement tape	Box		
1970	LNOP® Blue Neonatal /Infant/ Pediatric SpO ₂ Adhesive Sensors Single Patient Use 10/box, Non-sterile Weight 2.5 kg – 30 kg This product does not contain natural rubber latex. Use 1801 replacement tape	Box		

LNOP® Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

Masimo PN	Description	UM	Product Image	Notes
2358	LNOP® Trauma Adult SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex. Use 1926 replacement tape	Box		
2203	LNOP® Newborn Neonatal Neonatal SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight < 3 kg This product does not contain natural rubber latex. Use 1926 replacement tape	Box		
2204	LNOP® Newborn Infant/ Pediatric Infant/ Pediatric SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use 2391 replacement tape	Box		

LNOP® Reusable Sensors

Masimo PN	Description	UM	Product Image	Notes
1269	LNOP® DCI Adult SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1276	LNOP® DCIP Pediatric SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
1794	LNOP® TC-I SpO ₂ Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1793	LNOP® TF-I Adult SpO ₂ Transreflectance Forehead Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1560	LNOP® DC-195 Adult SpO ₂ Reusable Sensor 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1544	LNOP® YI Multisite SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		

LNOP® Red Patient Cables (for use with Masimo Rainbow SET Devices)

Masimo PN	Description	UM	Product Image	Notes
2058	Red PC-04 <i>LNOP Series 20-pin SpO₂ Patient Cable, 4 ft.</i>	EA		
2059	Red PC-08 <i>LNOP Series 20-pin SpO₂ Patient Cable, 8 ft.</i>			
2060	Red PC-12 <i>LNOP Series 20-pin SpO₂ Patient Cable, 12 ft.</i>			

LNOP® Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
1173	PC04 <i>LNOP Series Patient Cable, 4 ft.</i>	EA	A coiled white cable with a white connector on the left and a red connector on the right.	
1005	PC08 <i>LNOP Series Patient Cable, 8 ft.</i>			
1006	PC12 <i>LNOP Series Patient Cable, 12 ft.</i>			
1619	PC04-Ext <i>LNOP Extension Cable, 4 ft.</i>	EA	A coiled white cable with a white connector on the left and a red connector on the right.	
2282	LNOP® MP12 <i>LNOP to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 LNOP Dual Key Patient Cable, 12 ft.</i>	EA	A coiled white cable with a white connector on the left and a blue connector on the right.	
1847	LNOP® PC-12 P8 <i>LNOP to Masimo SET Patient Cable, 12 ft.</i>	EA	A coiled white cable with a white connector on the left and a blue connector on the right.	
1890	PC-12-GE <i>LNOP to GE Patient Cable, 12ft.</i>	EA	A coiled white cable with a white connector on the left and a blue connector on the right.	

LNOP® Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2288	PC-04-DB9 <i>LNOP Series to DB9 Connector Patient Cable, 4ft.</i>	EA		
2289	PC-08-DB9 <i>LNOP Series to DB9 Connector Patient Cable, 10ft.</i>			
2290	PC-12-DB9 <i>LNOP Series to DB9 Connector Patient Cable, 12 ft.</i>			
2433	LNOP® PC-SL-12 <i>LNOP Spacelabs Patient Cable, 12ft.</i>	EA		
1927	MAC-1 <i>LNOP SpO₂ Patient Cable, 12 ft.</i>	EA		

LNOP® Adapter Cables				
Masimo PN	Description	UM	Product Image	Notes
1848	LNOP® MAC P8 <i>LNOP to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>	EA		
1849	LNOP® CMS <i>LNOP Series to Philips 12-pin Instruments Adapter Cable, 12 ft.</i>	EA		
1645	MAC-1 <i>LNOP Series to LNC Adapter Cable, 18 in.</i>	EA		
2858	M-LNCS™ to PC Adapter Cable <i>M-LNCS Series to LNOP PC Patient Cables, 18 in.</i>	EA		
2014	Drager Patient Cable <i>For use with Masimo SET LNOP sensors, 2m</i>	EA		

Direct Connect Sensors (for use with Masimo SET Devices)

Masimo PN	Description	UM	Product Image	Notes
1396	DCSC Adult Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg <i>This product does not contain natural rubber latex.</i>	EA		
2387	DC-8 Adult Reusable Direct Connect SpO ₂ Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg <i>This product does not contain natural rubber latex.</i>			
1969	DC-12 Adult Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile Weight > 30 kg <i>This product does not contain natural rubber latex.</i>			
2651	DBI-dc3 Adult Soft Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg <i>This product does not contain natural rubber latex.</i>	EA		
2652	DBI-dc8 Adult Soft Reusable Direct Connect SpO ₂ Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg <i>This product does not contain natural rubber latex.</i>			



Red™ Direct Connect Sensors (for use with Masimo Rainbow SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2053	Red DCI-dc3 Adult Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.	EA		
2054	Red DCI-dc12 Adult Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.			
2256	Red DCIP-dc3 Pediatric Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight 10kg – 50kg This product does not contain natural rubber latex.	EA		
2257	Red DCIP-dc12 Pediatric Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile Weight 10kg – 50kg This product does not contain natural rubber latex.			
2643	Red DBI-dc3 Adult Soft Reusable Direct Connect Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2644	Red DBI-dc8 Adult Soft Reusable Direct Connect Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			

Mountable Patient Cable Extensions




Masimo PN	Description	UM	Product Image	Notes
1520	PXC16 Mountable Patient Cable Extension, 16 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.	EA		
1555	PXC30 Mountable Patient Cable Extension, 30 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.			
1557	PXC50 Mountable Patient Cable Extension, 50 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.			

Accessories

Sensor Accessories

Masimo PN	Description	UM	Product Image	Notes
2624	Replacement Tapes For Rainbow® R1 20L and R20-L For use with Rainbow R1 20L and R20-L 50/box Non-sterile This product does not contain natural rubber latex.	EA		
2623	Replacement Tapes For Rainbow® R1 25L and R25-L For use with Rainbow R1 25L and R25-L 50/box Non-sterile This product does not contain natural rubber latex.	EA		
1799	Replacement Tapes For LNOP® Neo-L Series Sensors For use with LNOP Neo-L 100/pack Non-sterile This product does not contain natural rubber latex.	EA		
1801	Replacement Tapes For LNOP® Inf-L and LNOP Blue Sensors For use with LNOP Inf-L and LNOP Blue 100/pack Non-sterile This product does not contain natural rubber latex.	EA		
1926	Replacement Wraps for LNOP® NeoPt-L, LNOP® Trauma and Newborn Neonatal Sensors For use with LNOP NeoPt-L, Trauma and Newborn Neonatal 10/pack Non-sterile This product does not contain natural rubber latex.	EA		
2391	Replacement Wraps for LNOP® Newborn Infant / Pediatric Sensors For use with LNOP Newborn Infant / Pediatric 10/pack Non-sterile This product does not contain natural rubber latex.	EA		
2307	Replacement Tapes for LNCS®/M-LNCS™ Inf Series Sensors For use with LNCS/M-LNCS Inf/Inf-3 100/box Non-sterile This product does not contain natural rubber latex.	EA		


Sensor Accessories

Masimo PN	Description	UM	Product Image	Notes
2308	Replacement Tapes For LNCS®/M-LNCS™ Neo Series Sensors For use with LNCS/M-LNCS Neo/Neo-3 100/box Non-sterile This product does not contain natural rubber latex.	EA		
2309	Replacement Wraps for LNCS®/M-LNCS™ NeoPt, NeoPt-3, Trauma and Newborn Neonatal Sensors For use with LNCS/M-LNCS NeoPt/ NeoPt-3/ Trauma and Newborn Neonatal 10/pack Non-sterile This product does not contain natural rubber latex.	EA		
2322	Replacement SoftTouch Wraps for LNCS®/M-LNCS™ NeoPt-500, Newborn Infant / Pediatric For use with LNCS/M-LNCS NeoPt-500/Newborn Infant/Pediatric 10/pack Non-sterile This product does not contain natural rubber latex.	EA		
1602	Foam Wraps for LNCS®/M-LNCS™ YI Sensors For use with YI Sensors 12/pack Washable This product does not contain natural rubber latex.	EA		
1600	Accessory Adhesive Squares For use with YI Sensors 144/pack	EA		
2822	Replacement Tape for M-LNCS™ Blue For use with M-LNCS Blue 50/pack Non-sterile This product does not contain natural rubber latex.	EA		
1596	CleanShield® Multisite Wrap For use with YI Sensors Adult/Pediatric/Neonatal 100/box Non-sterile This product does not contain natural rubber latex.	EA		



Sensor Accessories

Masimo PN	Description	UM	Product Image	Notes
1597	Standard Multisite Wrap For use with YI Sensors Adult/Pediatric/Neonatal 100/box This product does not contain natural rubber latex.	EA		
1598	Standard Petite Wrap For use with YI Sensors Adult slender digit/Ped/Neonatal adhesive 100/box This product does not contain natural rubber latex.	EA		
1601	Accessory Clothing Clips For securing reusable sensor cables to patient clothing 5/pack	EA		
1608	Masimo Headband For use with TF-I Sensors 5/box	EA		
2215	Masimo Headband For use with TF-I Sensors 12/box			
1903	Transflectance Sensor Adhesive Pad For use with TF-I Sensors 20/pack	EA		
21514	LNC™ Replacement Cover For use with all LNC patient cables 25 covers and 50 pins.	EA		



Sensor Demo & Sample Packs

Masimo PN	Description	UM	Product Image	Notes
2625	Rainbow® R1 25 Adult Demo Pack 2 Rainbow R1 25 SpHb, SpO ₂ , SpMet Sensors	EA		
2626	Rainbow® R1 20 Ped Demo Pack 2 Rainbow R1 20 SpHb, SpO ₂ , SpMet Sensors	EA		
2745	Rainbow® R1 25L Adult Demo Pack 2 Rainbow R1 25L SpHb, SpO ₂ , SpMet Sensors	EA		
2746	Rainbow® R1 20L Inf Demo Pack 2 Rainbow R1 20L SpHb, SpO ₂ , SpMet Sensors	EA		
2291	Rainbow® Adult Sample Pack 1 Rainbow R25 and 1 Rainbow R25-L SpO ₂ , SpCO, SpMet Sensors	EA		
2292	Rainbow® Ped/Inf Sample Pack Comes with 1 Rainbow R20 and 1 Rainbow R20-L sensor.	EA		
2656	M-LNCS® Sensor Sample Pack, 3ft. Comes with 1 M-LNCS Adtx-3, 1 M-LNCS Pdtx-3, and 1 M-LNCS Neo-3 adhesive sensors.	EA		


Sensor Demo & Sample Packs

Masimo PN	Description	UM	Product Image	Notes
2657	M-LNCS® Adult Sample Pack <i>Comes with 1 M-LNCS Adtx, 1 M-LNCS Pdtx, and 1 M-LNCS Neo adhesive sensors.</i>	EA		
2658	M-LNCS® Neo Sample Pack <i>Comes with 1 M-LNCS Inf, 1 M-LNCS Neo, and 1 M-LNCS NeoPt adhesive sensors.</i>	EA		
2659	M-LNCS® Sensor Training Kit <i>Includes 1 LNCS Adtx, 1 LNCS Pdtx, 1 LNCS Inf, 1 LNCS Neo, 1 LNCS NeoPt sensor and 1 neonatal foot with application card, direction for use.</i>	EA		
2345	LNCS® Sensor Training Kit <i>Includes 1 LNCS Adtx, 1 LNCS Pdtx, 1 LNCS Inf, 1 LNCS Neo, 1 LNCS NeoPt sensor and 1 neonatal foot with application card, direction for use.</i>	EA		
2346	LNCS® Sensor Sample Pack, 3ft. <i>Comes with 1 LNCS Adtx-3, 1 LNCS Pdtx-3, 1 LNCS Inf-3 and 1 LNCS Neo-3 sensor.</i>	EA		

Sensor Demo & Sample Packs


Masimo PN	Description	UM	Product Image	Notes
2347	LNCS® Adult Sample Pack Comes with 1 LNCS Adtx, 1 LNCS Pdt, 1 LNCS Neo sensor.	EA		
2348	LNCS® Neo Sample Pack Comes with 1 LNCS Inf, 1 LNCS Neo and 1 LNCS NeoPt.	EA		
1898	LNCS® Neo-L/Pdtx/Inf-L Sample Pack Comes with 2 LNCS Neo-L, 1 LNCS Pdt and 1 LNCS Inf-L.	EA		
1228	LNOP® Adt/Pdt Sample Pack Comes with 3 LNOP Adt and 1 LNOP Pdt Adhesive Sensor.	EA		
1283	LNOP® Sample Pack Comes with 2 LNOP Adt, 1 LNOP Pdt, 1 LNOP Inf-L, and 1 LNOP Neo-L Adhesive Sensor.	EA		
2262	LNOP® Sensor Sample Pack Comes with 1 LNOP Adtx, 1 LNOP Pdt, 1 LNOP Inf-L and 1 LNOP Neo-L.	EA		

Sensor Demo & Sample Packs

Masimo PN	Description	UM	Product Image	Notes
1229	LNOP® Neo Sample Pack <i>Comes with 2 LNOP Neo-L, 1 LNOP NeoPt-L and 1 LNOP Inf-L Adhesive Sensor.</i>	EA		
1282	LNOP® Sensors Application Training Packet <i>Comes with 1 LNOP Adt, 1 LNOP Pdt., 1 LNOP Inf-L, 1 LNOP Neo-L and 1 LNOP NeoPt-L, 2 cotton swabs, 2 alcohol swabs and 1 neonatal foot with application card, direction for use.</i>	EA		
1937	MAC-1 Adult Reusable Sensor Kit <i>Comes with 1 MAC-1 Cable #1645 and 1 LNOP DCI.</i>	EA		
1805	MAC-1 Adult Reusable Sensor Kit <i>Comes with 1 MAC-1 Cable #1645 and 1 LNOP DC-195.</i>	EA		

Sensor Demo & Sample Packs

Masimo PN	Description	UM	Product Image	Notes
1806	MAC-1 Adult Adhesive Sensor Kit <i>Comes with a MAC-1 Cable #1645 and 2 each of the LNOP Adt and Pdt sensors.</i>	EA		
1807	MAC-1 Neonatal Adhesive Sensor Kit <i>Comes with a MAC-1 Cable #1645 and 2 each of the LNOP Neo-L and Inf-L sensors</i>	EA		

Bedside Device Accessories				
Masimo PN	Description	UM	Product Image	Notes
1315	Radical® Series Replacement Battery <i>For use with Radical and Radical-7Blue Screen Handhelds.</i>	EA		
2386	Radical-7™ Color Screen Replacement Battery <i>For use with Radical-7 Color Screen only.</i>	EA		
1317	Radical® Series Pole Clamp <i>For use with Radical and Radical-7 Handhelds.</i>	EA		
1395	Radical® Series Handheld Lock <i>To secure Radical and Radical-7 Blue Screen Handhelds in Radical Docking Stations.</i>	EA		
30534	Radical® Series Handheld Lock Key <i>To unlock Radical and Radical-7 Blue Screen handhelds from Docking Stations.</i>	EA		
1584	Radical® Series Power Cord Lock <i>To secure a power cord to Radical Docking Stations 5/pack.</i>	EA		

Bedside Device Accessories				
Masimo PN	Description	UM	Product Image	Notes
2228	Rad-8® Power Cord Lock	EA		
30475	Radical® Series Handheld Case (Shown with Blue Screen Handheld Radical. Blue Radical not included).	EA		
2367	9-pin Masimo SET® Tester Tester for devices containing Masimo SET and connects to LNCS patient cables (Tests only SpO2 and BPM).	EA		
2368	Masimo Rainbow SET® Tester Tester for devices containing Masimo SET or Masimo Rainbow SET, 20 pin instrument connector. (Tests SpHb, SpCO, SpMet, PVI, PI, SpO2 and BPM).	EA		Can be issued with all 20 pin devices
1593	14 Pin Masimo SET® Tester For Oximeter SpO2 and PR Testing Tester for devices containing Masimo SET, 14 pin instrument connector. (Tests only SpO2 and BPM).	EA		

Bedside Device Accessories

Masimo PN	Description	UM	Product Image	Notes
1795	Masimo SET® Tester <i>For Oximeter & Patient Cable SpO₂ and PR Testing. Tester for devices containing Masimo SET and connects to LNOP PC Patient Cables. (Tests only SpO₂ and BPM).</i>	EA		
2227	Rad-8® Mounting Bracket	EA		
2350	Radical-7 Color Handheld Lock <i>Locks Radical-7 Color handheld to the docking station (used with PN 2351).</i>	EA		
2351	Radical-7 Color Handheld Lock Key <i>Used for installation of Radical-7 Color Handheld Lock and unlocking handheld from docking station (Used with PN 2350).</i>	EA		
2389	Mini-Demo Connector - 30 Min <i>Used to demo product (best used to display trending information).</i>	EA		
2390	Mini-Demo Connector - 2 Min <i>Used to demo product (best used for short demonstrations).</i>	EA		

Handheld Device Accessories


Masimo PN	Description	UM	Product Image	Notes
2207	Water Resistant Protective Carrying Case for Battery operated portable handhelds, Black <i>(This can be used with Rad-5/5v).</i>	EA	 Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.	
2208	Water Resistant Protective Carrying Case for Battery operated portable handhelds, Red <i>(This can be used with Rad-5/5v).</i>	EA	 Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.	
13158	Handheld Carrying Case with Shoulder Strap <i>(Shown with 14 pin Rad-5™ Device not included).</i>	EA		Note: Compatible with all Masimo Handheld Devices (Rad-5, Rad-5v, Rad-57, Pronto).

Handheld Device Accessories


Masimo PN	Description	UM	Product Image	Notes
1909	Data Transfer Download Cable <i>Used for connecting Rad-5 (with 14-pin connector) to Personal Computer to enable download of trending data.</i>	EA		
2063	Data Transfer Download Cable, Red – 20 Pin <i>Used for connecting Rad-57's and 20-pin Rad-5's to Personal Computer to enable download of trending data.</i>	EA		
1908	TrendCom™, Trend Download Software <i>Used with appropriate Data Transfer Download cable to initiate download and de-compression of Rad-5 and Rad-57 trending data to a Personal Computer. Used with serial cable to initiate download and decompression of Radical-7, Rad-87 and Rad-8 trending data.</i>	EA		
1842	Handheld Boot, Gray <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
1980	Handheld Boot, Yellow <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
1981	Handheld Boot, Red <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
1982	Handheld Boot, Orange <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		

Handheld Device Accessories

Masimo PN	Description	UM	Product Image	Notes
2097	Handheld Boot, Royal Blue <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2098	Handheld Boot, Light Blue <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2099	Handheld Boot, Pink <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2357	Light Shield <i>For use in strobing environment and high sunlight intensity areas. 5/Pack</i>	EA		
2815	Masimo Ambient Shield Accessory <i>For use with Rainbow Sensors. 10/pack</i>	EA		

SatShare® Cables					
Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1595	AT01	Airshields Athena Monitor with Nellcor SpO ₂	EA		
1324	DO01	Datex-Ohmeda AS/3 Patient Monitor [M-NSAT] System 5 Compact Monitor [M-NSAT]	EA		
1326	DO02	Datex-Ohmeda AS/3 Patient Monitor [M-NESTPR-00-00 / M-NESTR-00-00] CardioCap/5 Monitor Light Monitor System 5 Compact Monitor [M-NESTPR-00-00 / M-NESTR-00-00]	EA		
1528	DO03	Datex-Ohmeda AS/3 Patient Monitor [M-EST]	EA		
1533	DO04	Datex-Ohmeda AS/3 Patient Monitor [M-ESTP] System 5 Compact Monitor [M-EST / M-ESTP]	EA		
2023	DO05	Datex-Ohmeda System 5 Compact Monitor [M-NESTPR-00-01]	EA		
2223	DO06	Datex-Ohmeda System 5 Compact Monitor [M-PRESTN..01]	EA		
2969	DO07	Datex-Ohmeda System 5 Anesthesia Monitor [E-PRESTN]	EA		
1325	DS01	Corometrics Maternal / Fetal Monitor Model 118 Maternal / Fetal Monitor Model 120 Critikon Dinamap MPS Select Monitor Datascope Expert Monitor with Nellcor SpO ₂ [HS-500] Accutor Monitor with Nellcor SpO ₂ Passport Monitor with Nellcor SpO ₂ North American Draeger Evita 2 Dura, Evita 4 Anesthesia Machine Julian Anesthesia Machine Narcomed 4 Anesthesia Machine Vitalert 1000, 3100 Monitor Edentec / Nellcor Puritan Bennett Edentrace II Plus (Model 3711) Medical Data Electronics Escort Prism Monitor [20001] Escort / Escort II Monitor Welch Allyn Protocol Propaq CS Monitor Propaq Encore Monitor Propaq Series 100 Monitor Witt Biomedical Frontend Patient Monitor	EA		
1539	DS02	Datascope Passport Monitor with Datascope SpO ₂ Accutor Monitor with Datascope SpO ₂	EA		

SatShare® Cables

Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1365	DS03	Datascope Expert Fukuda Denshi Dynascope DS-5300W [HS-500]	EA		
2979	DS04	Datascope Fukuda Denshi Dynascope DS-5400 [HS-500]	EA		
1357	HP03	Philips / Agilent / HP Viridia CMS Model [M1020A] M1165, M1166, M1167, M1175, M1176, M1177 Monitors [M1020A] M1092, M1094, M1095, M1097 Display Modules, using the M1046 Computer Module [M1020A] Omnicare / Viridia Model 24/26 Series Monitor System M1204 Model 78352 Monitor [M1020A] Series 50XM Fetal Monitor	EA		
1321	MQ01	GE / Marquette Solar 7000 / 8000 Series Monitor, Tramscope 12 / 12C Monitor, Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor, Smart- Pac Transport Display Monitor [Tram 250SL / Tram 450SL / Tram 650SL / Tram 850A&SL / SpO2 Solar Parameter Module] Eagle 3000 Series Monitor Eagle 4000 Series Monitor Dash Series Monitor North American Draeger Narkomed 6100 Anesthesia Systems with Integrated Patient Monitoring	EA		
1789	MS01	Datex-Ohmeda S5 [E-Masimo-00] Draeger Infinity Delta [SmartPOD with Masimo SET] GE/Marquette Solar 8000 [Tram 451M / Tram 851M / Masimo SET single parameter module] Philips IntelliVue Masimo SET MMS SpO2 module	EA		
1367	NK01	Nihon Kohden Model BMS-8301A/MV-831RA and Model BMS-8502A [AL-800PA]	EA		
1368	NK02	Nihon Kohden BSM-2301A and BSM-5105K	EA		

SatShare® Cables

Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1366	OH01	Airshields Athena Monitor with Ohmeda SpO ₂ Ohmeda RGM Monitor (Respiratory Gas Machine) GE / Marquette Solar 7000 / 8000 Series Monitor, Tramscope 12 / 12C Monitor, Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor, Smart-Pac Transport Display Monitor [Tram 100A / Tram 200A&SL / Tram 400A / Tram 500A / Tram 600A&SL / Tram 800A&SL]	EA		
1646	OH02	Hellige SMU 611 Monitor <i>(Not available for sale in the United States)</i>	EA		
1840	PM01	Philips Medical M3/M4 [MMS Module] Intellivue Monitor [MMS Module]	EA		
1323	SL01	SpaceLabs Medical PC Express Monitor (90308), PC Scout Monitor (90309), PC Bedside / Central Monitor (90303, 90311, 90312), PC2 Bedside / Central Monitor (90305, 90315, 90316) [90465 / 90466 / 90467 / 90489] 600 Series Monitor	EA		
1362	SL02	SpaceLabs Medical Ultraview 1030 / 1050 Portable Bedside Monitor (90367, 90369) [90496] Ultraview 1500 / 1600 / 1700 Bedside Monitor (90363, 90364, 90387) [90496] Ultraview UCW Bedside Monitor (90385) [90496]	EA		
1322	SM01	Siemens Medical Systems Sirecust 1281 Monitor SC 6002XL, SC7000, SC9000XL	EA		
1605	SatShare Extension Cable		EA		

* Module compatibility information is subject to change without notice. For more current information, please contact Masimo Corporation.