

# FandPLimitTool

For single-molecule localization and resolution

# MUMDesignTool

Designing the focal plane spacing for multifocal plane microscopy



## Ward Ober Lab

This document is a user's manual for the Fundamental and Practical Limit Tool (FandPLimitTool) and the MUMDesignTool. It explains how to calculate 2D and 3D localization accuracy limits as well as 2D and 3D resolution limits for single-molecule microscopy data. In addition, it illustrates how to design the focal plane spacing for a multifocal plane microscopy setup with up to 10 focal planes.

**UT SOUTHWESTERN**

**Room ND6.310  
6001 Forest Park Rd  
Dallas, TX 75390-8576**

**+1 (214) 648-1922**

**+1 (214) 648-1259**

# **FandPLimitTool**

For single molecule localization  
and resolution

# **MUMDesignTool**

Designing the focal plane spacing  
for multifocal plane microscopy

## **User's Manual**

Copyright © 2004 – 2014 Ward Ober Lab

<http://www.wardoberlab.com/>

[wardlab@utsouthwestern.edu](mailto:wardlab@utsouthwestern.edu)

Updated on June 30, 2014 by Amir Tahmasbi

Development Team:

Raimund J. Ober, Amir Tahmasbi, Anish V. Abraham, Jerry Chao, Sripad Ram, Felix W. Tang, and E. Sally Ward

## Contents

List of figures.....	4
List of tables.....	6
Introduction.....	7
What is new? .....	7
Software overview.....	7
Before you begin.....	9
(WARNING) User defined functions .....	9
1. 2D-localization accuracy measure .....	10
1.1. Summary of imaging conditions .....	10
1.2. Providing required parameters.....	11
1.3. Providing advanced parameters .....	13
1.4. Viewing and saving the parameter summary .....	16
1.5. Executing the task and viewing results .....	16
2. Varying a parameter (multivalue mode).....	18
2.1. Providing inputs for the multivalue mode .....	18
2.2. Plotting the limit of accuracy versus a parameter value .....	19
3. Simulating and visualizing model and data images .....	20
3.1. Calculate a model image.....	20
3.2. Simulate data images.....	22
4. Saving and loading the settings .....	25
4.1. Saving the settings .....	25
4.2. Loading the settings .....	25
5. Exporting results, parameters, images and report .....	26
5.1. Exporting the results and parameters .....	26
5.2. Exporting the model and data images .....	28
5.3. Exporting a LaTeX report.....	28
5.4. Exporting everything simultaneously.....	29
6. Console window, log file and bug reporting .....	30
6.1. The console window .....	30

6.2.	The log file and bug reporting.....	30
7.	The MUMDesignTool .....	32
7.1.	Summary of simulation conditions .....	32
7.2.	Providing the inputs and configuring the MUMDesignTool .....	33
7.3.	Designing the plane spacing using the Rapid Mode .....	34
7.4.	Precise Mode .....	38
8.	3D-Localization accuracy measure (Born-Wolf Model) .....	40
8.1.	Summary of imaging conditions .....	40
8.2.	Providing required parameters.....	41
8.3.	Providing advanced parameters .....	42
8.4.	Executing the task and viewing results .....	45
9.	2D-Resolution measure.....	47
9.1.	Summary of imaging conditions .....	47
9.2.	Providing required parameters.....	48
9.3.	Providing advanced parameters .....	49
9.4.	Executing the task and viewing results .....	51
10.	3D-Resolution measure (Born-Wolf Model) .....	53
10.1.	Summary of imaging conditions .....	53
10.2.	Providing required parameters.....	54
10.3.	Providing advanced parameters .....	55
10.4.	Executing the task and viewing results .....	58
11.	3D-Localization accuracy measure - Multifocal Plane .....	59
11.1.	Summary of imaging conditions .....	59
11.2.	Providing required parameters.....	60
11.3.	Providing advanced parameters .....	61
11.4.	Executing the task and viewing results .....	64
	References .....	65

## List of figures

Figure 1.2.1: Selecting the 2D Localization Accuracy calculate type .....	11
Figure 1.2.2: Required parameters window (2D Localization) .....	12
Figure 1.2.3: Electron multiplication parameters .....	13
Figure 1.3.1: Advanced parameters window (2D Localization).....	14
Figure 1.3.2: Pixel Integration Trapezoidal options .....	14
Figure 1.3.3: FIM Single Integration Trapezoidal options .....	15
Figure 1.3.4: FIM Double Integration Trapezoidal options.....	15
Figure 1.3.5: Noise factor calculation parameters .....	15
Figure 1.4.1: Viewing and saving the parameter summary (2D Localization) .....	16
Figure 1.5.1: Calculation and Final confirmation .....	17
Figure 1.5.2: 2D localization accuracy Results window .....	17
Figure 2.1.1: Multivalue Mode setting window .....	18
Figure 2.2.1: Multivalue results and plot windows.....	19
Figure 3.1.1: Unavailable model and data images .....	20
Figure 3.1.2: Calculate a model image .....	21
Figure 3.1.3: Model image viewer .....	21
Figure 3.1.4: Model image viewer in Multivalue Mode.....	22
Figure 3.2.1: Enter the number of data images .....	23
Figure 3.2.2: Simulating data images.....	23
Figure 3.2.3: Simulated image viewer in Multivalue Mode.....	24
Figure 4.1.1: Saving and loading the settings.....	25
Figure 5.1.1: Exporting is not available .....	26
Figure 5.1.2: Exporting the results and parameters .....	27
Figure 5.1.3: An example for the exported results and parameters .....	27
Figure 5.2.1: Exporting the model and data images .....	28
Figure 5.3.1: Exporting a LaTeX report .....	29
Figure 6.1.1: Console window .....	30
Figure 6.2.1: Log file .....	31
Figure 6.2.2: About the FandPLimitTool .....	31

Figure 7.2.1: The MUMDesignTool configuration window.....	33
Figure 7.2.2: MUMDesignTool calculation in progress.....	34
Figure 7.3.1: MUMDesignTool Rapid Mode (main window).....	35
Figure 7.3.2: Changing the number of focal planes in real-time .....	36
Figure 7.3.3: The lateral elements of FIM and PLAM (MUMDesignTool Rapid Mode) .....	37
Figure 7.3.4: Exporting the MUMDesignTool results .....	37
Figure 7.4.1: Entering the noise parameters.....	38
Figure 7.4.2: Calculation procedure in the Precise Mode.....	39
Figure 7.4.3: MUMDesignTool Precise Mode .....	39
Figure 8.2.1: 3D Localization Accuracy problem scenario and its required parameters window .....	41
Figure 8.3.1: 3D Localization Accuracy advanced parameters window .....	42
Figure 8.3.2: Pixel Integration Trapezoidal options .....	43
Figure 8.3.3: Model Integration Trapezoidal options.....	43
Figure 8.3.4: FIM Single Integration Trapezoidal options .....	44
Figure 8.3.5: FIM Double Integration Trapezoidal options.....	44
Figure 8.4.1: Calculation and Final confirmation .....	45
Figure 8.4.2: 3D localization Accuracy results window.....	46
Figure 9.2.1: 2D Resolution Accuracy calculate type and its required parameters window .....	48
Figure 9.3.1: 2D Resolution Accuracy advanced parameters window.....	49
Figure 9.3.2: Pixel Integration Trapezoidal options .....	50
Figure 9.3.3: FIM Single Integration Trapezoidal options .....	50
Figure 9.3.4: FIM Double Integration Trapezoidal options.....	51
Figure 9.4.1: Calculation and Final confirmation .....	51
Figure 9.4.2: Results window (2D Resolution).....	52
Figure 10.2.1: 3D Resolution Accuracy calculate type and its required parameters window.....	54
Figure 10.3.1: 3D Resolution Accuracy advanced parameters window .....	55
Figure 10.3.2: Pixel Integration Trapezoidal options.....	56
Figure 10.3.3: Model Integration Trapezoidal options .....	56
Figure 10.3.4: FIM Single Integration Trapezoidal options.....	57
Figure 10.3.5: FIM Double Integration Trapezoidal options .....	57

Figure 10.4.1: Calculation and Final confirmation (3D Resolution) .....	58
Figure 10.4.2: Results window (3D Resolution) .....	58
Figure 11.2.1: 3D Localization Accuracy (MUM) calculate type and its required parameters window	60
Figure 11.3.1: 3D Localization Accuracy (MUM) advanced parameters window .....	61
Figure 11.3.2: Pixel Integration Trapezoidal options.....	62
Figure 11.3.3: Model Integration Trapezoidal options .....	62
Figure 11.3.4: FIM Single Integration Trapezoidal options.....	63
Figure 11.3.5: FIM Double Integration Trapezoidal options .....	63
Figure 11.4.1: Calculation and Final confirmation (3D localization MUM) .....	64
Figure 11.4.2: Results window (3D localization MUM).....	64

## List of tables

Table 1.1.1: 2D Localization Accuracy - Summary of imaging conditions.....	10
Table 7.1.1: Summary of simulation parameters for the MUMDesignTool.....	32
Table 8.1.1: 3D Localization Accuracy - Summary of imaging conditions.....	40
Table 9.1.1: 2D Resolution Accuracy - Summary of imaging conditions .....	47
Table 10.1.1: 3D Resolution Accuracy - Summary of imaging conditions.....	53
Table 11.1.1: 3D Localization Accuracy (MUM) - Summary of imaging conditions.....	59

## Introduction

This document is a user's manual for the FandPLimitTool V1.2 and MUMDesignTool V1.0. It provides detailed instructions for performing five specific calculations which are representative of the types of calculations that can be performed with this application. Moreover, it provides detailed information about how to use the MUMDesignTool to design the focal plane spacing for a multifocal plane microscopy (MUM) setup [1].

Consider the following terminology which will be used multiple times throughout this document.

**FLAM:** fundamental localization accuracy measure

**PLAM:** practical localization accuracy measure

**FREM:** fundamental resolution measure

**PREM:** practical resolution measure

## What is new?

We have added several new features to the FandPLimitTool V1.2 which make it more powerful and faster than its previous version (V1.1). The new features are listed as follows and each option will be discussed in detail later.

- A new software module 'the MUMDesignTool' is incorporated with the package which provides tools to design the focal plane spacing for multifocal plane microscopy.
- New functionality that allows the calculation of the 3D resolution measure.
- A new function called "vary a parameter" is added to the software which allows users to calculate and plot the limits of accuracy for a range of different experimental parameters.
- New functionality that allows simulating single molecule model and data images.
- New functionality that allows determining limits of accuracy in the presence of stochastic signal amplification, i.e. when using an electron multiplying CCD (EMCCD) detector.
- The code is significantly optimized. The software therefore runs much faster than the previous version.
- Exporting tools are enhanced and user can now export the results, parameters, simulated model and data images and graphs as text, image and LaTeX report files, respectively.
- Summary of all input parameters can be explored within the software.
- Enhanced user interface design.
- A console window and a log file to facilitate bug reporting.


## Software overview

The FandPLimitTool allows users to determine the fundamental and practical limits of 2D and 3D localization accuracy and resolution for a variety of imaging conditions. The software has been

built using object oriented programming methodologies and has been designed for flexibility to accommodate new models and numerical calculation routines in the future.

Each calculation supported by this application requires a different set of parameters. The parameters are divided into two categories: *Required parameters*, and *Advanced parameters*.

Required parameters are those parameters that specify the imaging conditions under which the calculations are to be performed. Users must provide values for these parameters before a calculation can proceed. These parameters by default are visible to the user when software starts. They can also be accessed by clicking the **Required Parameters** option under the **Parameters Menu** or through its shortcut **Ctrl+U**. Although default values are always provided for all required parameters, users must modify the fields to suit their desired imaging conditions.



Advanced parameters, on the other hand, are optional inputs that govern the numerical routines involved in performing the required calculations. These parameters are optional; the default values supplied for these inputs will work for typical imaging conditions. Advanced parameters can be accessed by clicking the **Advanced Parameters** option under the **Parameters Menu** or through its shortcut **Ctrl+N**. For details about specifying values for advanced parameters, please refer to Section 1 to Section 5.

Clicking the **Parameter Summary** option under the **View Menu** or hitting **Ctrl+P** provides a summary of the required and advanced parameters in a plain text format. This helps to conveniently review and save the parameters as a text file. In addition, clicking the **Console** option under the **View Menu** shows a console window in which the details of calculations will be listed. This console window is also helpful for bug reporting purposes.

To calculate the limits of accuracy after specifying all desired parameters, click the **Calculate Limits of Accuracy** button. When the calculations are completed, a **Results window** will be brought into view displaying the results. The **Results window** can be accessed at any time by clicking the **Results** option under the **View Menu** or through its shortcut **Ctrl+R**.

Similarly, to calculate a model image after specifying all desired parameters press the **Calculate Model Image**. A model image is an image profile that is devoid of stochasticity and noise. To simulate realizations of this model image, press the **Simulate Data Images** button. Note that the simulated data images are corrupted by stochasticity and noise. For instance, depending on the

desired set of parameters, the image can be corrupted by extraneous noise sources such as background, readout or stochastic signal amplification.

Finally, clicking the **MUMDesignTool** option under the **Tools Menu** or hitting **Ctrl+M** calls the *MUMDesignTool* application which will be discussed in detail later in this tutorial.

## Before you begin

Before following this tutorial, please make sure that the software has been correctly installed and started. Please refer to the *FandPLimitTool – Installation Guide* for instructions on installing and starting the application.

## (WARNING) User defined functions

The **Advanced Inputs Window** has fields where the user is required to select which function to use to perform a certain part of the calculation. In the options available for such fields, the “User Defined” option is also often available. This option allows the user to specify the names of their own custom function for performing the associated part of the calculation.

However, in the executable version of the program, this option is not fully supported. Users cannot specify the names of their own functions for the user defined options as such functions will not have been compiled along with the original program to generate the executable and therefore will not be recognized by the executable version of the program.

## 1. 2D-localization accuracy measure

This section provides an example of how to use the FandPLimitTool to calculate the accuracy with which the 2D location of an object can be estimated from its image (PLAM). The object is assumed to be in focus with respect to the objective lens and the image of the object is assumed to be an Airy profile.

The calculation takes into account the stochastic nature of the photon detection processes and assumes the photons detected from the object to be Poisson distributed. In addition, the calculation also takes into account extraneous noise sources - in particular the background component which is modeled as additive Poisson noise, the readout noise of the detector which is modeled as independent, additive Gaussian noise and the stochastic signal amplification noise which is present in EMCCD detectors [3].

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the region of interest (ROI) containing the image of the object is assumed to be a rectangular array of pixels. The object is assumed to be in the center of this rectangular pixel array and the location of the object is specified with respect to the object space. Refer to [2] for details on determining the 2D localization accuracy in single molecule microscopy.

### 1.1. Summary of imaging conditions

The following table summarizes the imaging conditions under which the calculations are performed.

Property	Value
Numerical aperture	1.45
Wavelength	0.69 $\mu\text{m}$
Photon detection rate	10,000 <i>photons/second</i>
Exposure time	0.13 <i>seconds</i>
Magnification	100 X
Location of object in ROI	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$
Pixel size	13 $\mu\text{m}$ $\times$ 13 $\mu\text{m}$
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>
Background level (Poisson rate)	1000 <i>photons/pixel/second</i>
Gaussian noise mean	0 $e^-$
Gaussian noise standard deviation	8 $e^-$ <i>rms</i>
Use EM	<i>unchecked</i>
EM gain	100
EM model	<i>Geometric</i>

**Table 1.1.1: 2D Localization Accuracy - Summary of imaging conditions**

## 1.2. Providing required parameters

**Step 1:** Select the problem scenario of “2D Localization Accuracy - Airy: Pixelated with Poisson + Gaussian Noise” from the Calculate pull-down menu in the main application window.





Figure 1.2.1: Selecting the 2D Localization Accuracy calculate type

**Step 2:** The required parameters window, which is visible by default, will be updated accordingly. If this window is not visible, access it by clicking the **Required Parameters** option under the **Parameters Menu** in the main application window.



**Figure 1.2.2: Required parameters window (2D Localization)**


**Step 3:** In the - Fundamental inputs - section of the window, enter the values for the various fields as follows: *Numerical aperture = 1.45*, *Wavelength = 0.69*, *Photon detection rate = 10000*, and *Exposure time = 0.13*.

**Step 4:** In the - Additional inputs - section of the window, enter the values for the various fields as follows: *Magnification = 100*,  *$x_0 = 13 \times 4.5 / 100$* ,  *$y_0 = 13 \times 4.5 / 100$* , *Pixel size = [13, 13]*, and *ROI size = [9, 9]*.

**Remark 1:** The center of a *9 pixels × 9 pixels* pixel array where the dimensions of each pixel is *13 μm × 13 μm*, is *(4.5 × 13) μm × (4.5 × 13) μm*. To convert this location coordinates to the object space, the coordinate values are divided by the magnification. Hence the location coordinates are specified as *13 × 4.5 / 100 μm × 13 × 4.5 / 100 μm*.

**Step 5:** In the - Extraneous noise sources - section of the window, enter the values for the various fields as follows: *Background level = 1000* and *Standard deviation = 8*.

**Remark 2:** In case that the localization accuracy is required for an EMCCD detector, the *Use EM* checkbox should be checked and the desired value for the electron multiplication gain should be entered in the *EM gain* field. The default electron multiplication model is Geometric. However, different electron multiplication models, e.g. High gain and Gaussian, can be chosen from the *EM model* combo box based on the requirements. By holding the mouse pointer on the *EM model* combo box, a brief description of each option is shown on a tool tip text. For more information about different electron multiplication models consult [3].



**Figure 1.2.3: Electron multiplication parameters**

**Step 6:** In the - Parameters to be estimated - section of the window, ensure that all check-boxes for all parameters are checked.


**Step 7:** To improve screen clarity, close the **Required parameters window**.

**Remark 3:** Closing the window is not in any way required by the application for the calculations to be performed and is only a suggestion in order to reduce the number of open windows and hence improve screen clarity. You can access the window again at any time by clicking the **Required Parameters** option under the **Parameters Menu** in the main application window.

### 1.3. Providing advanced parameters

**NOTE:** All inputs in the Advanced Parameters window are optional. The default values provided have been found to work adequately for typical imaging conditions.


**Step 1:** Access the advanced parameters window by clicking the **Advanced Parameters** option under the **Parameters Menu** in the main window.



**Figure 1.3.1: Advanced parameters window (2D Localization)**

**Step 2:** In the - Model profile calculations - section of window, access the Pixel integration method pull-down menu and select “Trapezoidal”.

**Step 3:** Click the **Additional parameters** button next to the Pixel integration method pull-down menu. In the dialog box that appears, enter values for fields as follows:  $XGridding = 13$ ,  $YGridding = 13$ . Click OK when done.



**Figure 1.3.2: Pixel Integration Trapezoidal options**

**Step 4:** Access the Background type pull-down menu and select “Constant”. There are no additional parameters to be accessed via the **Additional parameters** button for the “Constant” option.

**Step 5:** In the - Fisher information matrix calculations - section of the window, access the Single integration method and select “Trapezoidal”.

**Step 6:** Click on the **Additional parameters** button next to the Single integration method pull-down menu. In the dialog box that appears, enter  $Step\ size = 0.05$ .




Figure 1.3.3: FIM Single Integration Trapezoidal options

**Step 7:** Access the Double integration method pull-down menu in the - Fisher information matrix calculations - section of the window and select "Trapezoidal".




Figure 1.3.4: FIM Double Integration Trapezoidal options

**Step 8:** Click on the **Additional parameters** button next to the Double integration method pull-down menu. In the dialog box that appears, enter the values  $XGridding = 13$  and  $YGridding = 13$ .

**Step 9:** The - Fisher information matrix calculations - section has a subsection titled - Noise factor calculation parameters - . Enter values  $Limits\ of\ summation = [1, 3000]$  and  $Limits\ of\ integration = [-5000, 5000]$  in this subsection.

**Remark 4:** The limits of integration and summation are initially set to default values which are automatically calculated based on the imaging conditions. It is possible specify these limits by changing the value of the combo box to the right of the limits of integration and summation fields.




Figure 1.3.5: Noise factor calculation parameters


**Step 10:** To improve screen clarity, close the **Advanced parameters window**.

**Remark 5:** Closing the window is not in any way required by the application for the calculations to be performed and is only a suggestion in order to reduce the number of open windows and hence improve screen clarity. You can access the window again at any time by clicking the **Advanced Parameters** option under the **Parameters Menu** in the main application window.

## 1.4. Viewing and saving the parameter summary

**NOTE:** This is a new feature of the FandPlimitTool V1.2. Typically it is necessary to save the required and advanced parameters for later developments and to keep track of calculations. This new feature provides a summary of the required and advanced parameters in a plain text format.

**Step 1:** Access the required inputs summary window by hitting the **Parameter Summary** option under **View Menu**.



**Figure 1.4.1: Viewing and saving the parameter summary (2D Localization)**


**Step 2:** You will see the summary of the required and advanced parameters in a plain text format. This summary can be saved on the disc as a .txt file by hitting the **Save as text** button.

## 1.5. Executing the task and viewing results

**Step 1:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed. During the calculation, a console


window will be shown to inform the user about the different steps of the calculation. This console window will also display errors and warnings, if there is any, which is very useful for bug reporting purposes.

**NOTE:** The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.



**Figure 1.5.1: Calculation and Final confirmation**

**Step 2:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy** and the **Results window** will be brought in focus. In addition, the console window will be hidden when the calculation is over. The console window can be accessed any time by selecting the **Console** option under the **View Menu**.



**Figure 1.5.2: 2D localization accuracy Results window**

## 2. Varying a parameter (multivalue mode)

Multivalue Mode is a new functionality designed for the calculation and visualization of the fundamental and practical limits on the 2D and 3D localization accuracy and resolution for a range of parameter values. This section provides an example of how to use the Multivalue Mode to calculate the 2D PLAM for a range of an input parameter. We assume that user has applied all the steps in Section 1 (2D-Localization Accuracy Limit).

### 2.1. Providing inputs for the multivalue mode

**Step 1:** Access the multivalue mode window by toggling the **Vary a Parameter** button at the bottom of the Required Parameters window.




Figure 2.1.1: Multivalue Mode setting window

**Step 2:** In the – Multivalue parameter - section of window, access the Parameter pull-down menu and select “Photon detection rate”.

**Step 3:** In the Values field, enter 100:100:1000. This will automatically define a vector containing [100 200 300 400 500 600 700 800 900 1000].

**Step 4:** check the Auto-center image for each calculation checkbox. This will automatically center the particle in the ROI.

**Step 5:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate Limits of Accuracy** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed. In addition the console window will become visible.


**NOTE:** The calculation could indeed take a couple of hours to complete depending on the number of data points and hardware capabilities of the platform on which the application is being executed.

**Step 6:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy**.

**NOTE:** Steps 5 and 6 are explained in more details in Section 1.5 (Executing the task and viewing results).

## 2.2. Plotting the limit of accuracy versus a parameter value

**Step 1:** In the **Multivalue Mode window**, click the **Plot** button. This will open the **Multivalue Plot Window** shown in the following figure.



**Figure 2.2.1: Multivalue results and plot windows**

**Step 2:** In the **Multivalue Plot window**, access the **Select Plot** pull-down menu and select “Limit of accuracy of  $x_0$  vs Photon detection rate”. A plot of the limit of accuracy of  $x_0$  versus the photon detection rate will be shown as illustrated in the above figure.


**Step 3:** In the **Multivalue Mode window**, click the **Results** button. This will open the **Multivalue Results Window** shown in the above figure. The selected values for the parameter will then be shown in a list box. Clicking different lines of this list box will show the corresponding results in the **Results Window**.

### 3. Simulating and visualizing model and data images

Another new feature of the FandPLimitTool V1.2 is that user can now simulate and visualize model and data images of a single molecule associated with a desired problem scenario and a specific set of parameters. A model image is in fact the image profile of a single molecule (or two single molecules in case of a resolution problem) in the absence of stochasticity and extraneous noise. A simulated data image, on the other hand, is a realization of a model image and is corrupted by stochasticity and extraneous noise sources such as background, readout and stochastic signal amplification. The following sections provide examples for calculating and visualizing model and data images, respectively. We assume that user has applied all the steps in Section 1 (2D-Localization Accuracy Limit).


#### 3.1. Calculate a model image

**NOTE:** For some of the problem scenarios (mostly fundamental limits), a model image may not be available. In such cases, the **Calculate Model Image** button will be grayed out and as such a model image cannot be calculated.




**Figure 3.1.1: Unavailable model and data images**

**Step 1:** If a model image is available for current problem scenario, in the main application window, click the **Calculate Model Image** button. The **Calculate Model Image** button text will change to **Calculating Model Image ...** and the application will temporarily become unresponsive while the calculations are being performed.



**Figure 3.1.2: Calculate a model image**

**Step 2:** When the calculation is completed, the **Calculating Model Image ...** button text will change back to **Calculate Model Image** and the display tool will be brought in focus.




**Figure 3.1.3: Model image viewer**

**Step 3:** In the viewer window, you can visualize the model image and change the visual effects.

**NOTE:** The viewer belongs to our previously developed Microscopy Image Analysis Tool (MIATool). For more information please see [4].

**Remark 6:** In case that the **Show Model** functionality is used together with the **Multivalued Mode**, the viewer will show a stack of images associated with the different values of the Multivalued parameter. In this case, user can visualize the different model images by scrolling the slider to the left in the **Viewer Slider window**, as shown below. The current value of the parameter (e.g. Magnification = 150) will also be shown in the viewer window (see below).




**Figure 3.1.4: Model image viewer in Multivalued Mode**

**Step 4:** To improve screen clarity, close the **Model Image Viewer** using the standard exit buttons on any of the viewer windows.

## 3.2. Simulate data images


**NOTE:** For some of the problem scenarios (mostly fundamental limits), a data image may not be available. In such cases, the **Simulate Data Images** button will be grayed out and as such a data image cannot be simulated (this is shown in Figure 3.1.1).

**Step 1:** If a data image is available for current problem scenario, in the main application window, click the **Simulate Data Images** button. This will present a dialog box to enter the desired number of data images (i.e. realizations) of current model.



**Figure 3.2.1: Enter the number of data images**

**Step 2:** Enter the desired number of data images. Click **OK** to proceed. The **Simulate Data Images** button text will change to **Simulating X Data Images ...**, where X is the number of desired images specified by user. The application will temporarily become unresponsive while the calculations are being performed.



**Figure 3.2.2: Simulating data images**


**Step 3:** When the calculation is completed, the **Simulating X Data Images ...** button text will change back to **Simulate Data Images** and the display tool will be brought in focus.

**Step 4:** In the viewer window, you can visualize the data images and change the visual effects. You can also visualize different realization by scrolling the slider to the right in the **Viewer Slider window**.

**NOTE:** The viewer belongs to our previously developed Microscopy Image Analysis Tool (MIATool). For more information please see [4].

**Remark 7:** In case that the **Simulate Data Images** functionality is used together with the **Multivalued Mode**, the viewer will show a stack of data images associated with the different values of the Multivalued parameter. In this case, user can visualize images that correspond to different model by scrolling the slider to

the left in the **Viewer Slider window**, as shown below. On the other hand, different realizations (data images) of a specific model can be accessed using the slider to the right (see below).



**Figure 3.2.3: Simulated image viewer in Multivalue Mode**

## 4. Saving and loading the settings

User can conveniently save and load the software settings which include the current problem scenario, required and advanced parameters, etc., into and from the disk, respectively. This functionality is useful for retrieving the adjustments and results of a specific calculation after closing and reopening the software. We next briefly describe how to make use of this important functionality.

### 4.1. Saving the settings

**Step 1:** To save the settings, click the **Save Settings** option under the **File Menu** or use the shortcut **Ctrl+S** (see below). This will present a dialog box to enter a file name. The file extension is **.fim**. Enter your desired name and click **Save**.




Figure 4.1.1: Saving and loading the settings

### 4.2. Loading the settings


**Step 1:** To load the settings, click the **Load Settings** option under the **File Menu** or use the shortcut **Ctrl+L** (see Figure 4.1.1). This will present a dialog box to choose the desired file with extension **.fim**. Select a previously saved setting and press **Open**.

## 5. Exporting results, parameters, images and report

The FandPLimitTool V1.2 provides very flexible tools for exporting the results of calculations, simulated images and even a report file. More specifically, user can export the results and parameters of a desired problem scenario, the simulated model and data images and a LaTeX report that can be generated when performing a multivale calculation (i.e. when varying a parameter). The simulated model and data images can also be exported when varying a parameter. In such cases, a sequence of images will be exported. In this section, we discuss different steps of exporting various results.

### 5.1. Exporting the results and parameters

**NOTE:** To export the desired outputs such as results, images and a report, appropriate calculations has to be performed prior to exporting. If the results, model and data images and a LaTeX report are not available, the corresponding menu options under the **Export Menu** will be grayed out (see below). After the appropriate calculation is performed the corresponding menu option becomes accessible.



**Figure 5.1.1: Exporting is not available**

**Step 1:** Perform either a standard (see e.g. Section 1) or a multivalue (see Section 2) calculation for the limits of accuracy. Once the calculation is performed and the results are available within the software, the **Results and Parameters** option under the **Export Menu** becomes accessible.

**Step 2:** Select the **Results and Parameters** option under the **Export Menu** (see figure below). This will present a folder browser.

**Step 3:** Select a folder name in which you would like to save the results and parameters and click OK. The results and parameters will be saved as two independent text files in the specified folder.




Figure 5.1.2: Exporting the results and parameters

**NOTE:** Depending on the type of calculation (i.e. singlevalue and multivalued), the exported results text file may include one value or multiple values for the limit of accuracy of each parameter. An example for the exported text file for a multivalued calculation is shown below.




Figure 5.1.3: An example for the exported results and parameters

## 5.2. Exporting the model and data images

**Step 1:** Exporting the model images is available upon performing at least one of the three different actions which are the calculation of the limits of accuracy (see e.g. Section 1), the calculation of a model image (see Section 3.1) and the simulation of data images (see Section 3.2). On the other hand, both model and data images can be exported only after simulating data images. Once one of the mentioned calculations is performed, the **Model and Data Images** option under the **Export Menu** becomes accessible.

**Step 2:** Select the **Model and Data Images** option under the **Export Menu** (see figure below). This will present a folder browser.

**Step 3:** Select a folder name in which you would like to save the simulated model and data images and click OK. The model and data images will be saved within two separate folders under the specified folder.

**NOTE:** Depending on the type of calculation (i.e. singlevalue or multivalue), the exported model and data image folders may contain a sequence of images.




Figure 5.2.1: Exporting the model and data images

## 5.3. Exporting a LeTeX report

**NOTE:** A LeTeX report can be generated only after performing a multivale calculation (i.e. calculating the limits of accuracy when varying a parameter).

**Step 1:** Perform a multivalue calculation for the limits of accuracy (see Section 2). Once the calculation is performed and the results are available within the software, the **LaTeX Report** option under the **Export Menu** becomes accessible.

**Step 2:** Select the **LaTeX Report** option under the **Export Menu** (see figure below). This will present a folder browser.




Figure 5.3.1: Exporting a LaTeX report

**Step 3:** Select a folder in which you would like to save the LaTeX report and click OK. A .tex file together with several .eps images will be saved within a folder named LaTeX Report under the specified folder.

## 5.4. Exporting everything simultaneously

**NOTE:** If all the required calculations are performed, the export **All...** option becomes available and user can export everything simultaneously.

**Step 1:** Select the **All...** option under the **Export Menu**. This will present a folder browser.

**Step 2:** Select a folder name in which you would like to export all the results including the limits of accuracy, the simulated model and data images and the LaTeX report, and then click OK.

## 6. Console window, log file and bug reporting

The FandPLimitTool V1.2 is equipped with a new console window and a log file. The console window, which can be accessed directly within the software, not only can provide users with the details of calculations in real-time, but also is very useful for identifying warning, errors and reporting bugs. The information that are shown in the console window are also stored in a log file.

### 6.1. The console window

As mentioned earlier, during the long calculations a console window is automatically displayed to inform the user about different steps of the calculation. This console window can also display errors and warnings, if there is any, which is very useful for bug reporting purposes.

**Step 1:** To access the console window at any time, select the **Console** option under the **View** **Menu**.

**Step 2:** To hide the console window, simply close it using the standard close button. Note that closing the console window will not stop storing the information in the log file.




Figure 6.1.1: Console window

### 6.2. The log file and bug reporting

The information that are displayed in the console window are also stored in a log file "fandp\_log.mia" in the following path:

C:\Users\{username}\AppData\Roaming\FandPLimitTool




Figure 6.2.1: Log file

To report a bug (that leads to an error or warning), please send this log file to us. More information about the current versions of the FandPLimitTool and the MUMDesignTool, their homepages and licensing can be found in the corresponding **About windows** (e.g. see below).




Figure 6.2.2: About the FandPLimitTool

## 7. The MUMDesignTool

This section describes how to use the MUMDesignTool, a new software module in the new version of the FandPLimitTool, to calculate and plot the Fisher information matrix (FIM) and the PLAM [1] along the z-axis. The MUMDesignTool is developed to help in designing the plane spacing for a multifocal plane microscopy (MUM) setup including up to 10 focal planes.

The MUMDesignTool provides two working modes. The *Rapid Mode* takes into account the stochastic nature of the photon detection processes and assumes that the observed data is Poisson distributed [1]. However, this mode does not take into consideration other extraneous noise sources such as the background, readout and stochastic signal amplification. This mode yields a fast approach for the calculation of the FIM and the PLAM, and allows changing the number of focal planes and their spacings in real-time. In addition, the *Precise Mode* takes into account the extraneous noise sources, such as the background noise which is independently Poisson distributed, the readout noise which is modeled as independent additive Gaussian noise [2], and in case of using an EMCCD detector, the stochastic electron multiplication noise that is modeled as a branching process [3].

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the ROIs containing the images of the object are assumed to be rectangular arrays of pixels. The object is assumed to be in the center of this rectangular pixel array and the location of the object is specified with respect to the object space. Refer to [1] for more details regarding the design of focal plane spacing for a MUM setup.

### 7.1. Summary of simulation conditions


The following table summarizes the imaging conditions under which the calculations are performed.

Property	Value
Number of focal planes	2
Photon detection rate	10,000 <i>photons/second</i>
Photon count percentage	[50, 50] %
Object medium refractive index	1.515
Numerical aperture	1.45
Wavelength	0.69 $\mu\text{m}$
Exposure time	0.1 <i>seconds</i>
Magnification	100 X
Location of object in ROI	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$
Pixel size	13 $\mu\text{m}$ $\times$ 13 $\mu\text{m}$
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>
Defocus start	-2 $\mu\text{m}$
Defocus increment	0.05 $\mu\text{m}$
Defocus end	2 $\mu\text{m}$

**Table 7.1.1: Summary of simulation parameters for the MUMDesignTool**

## 7.2. Providing the inputs and configuring the MUMDesignTool

**Step 1:** Access the **Configure Setting - MUM Design Tool** window by clicking the **MUMDesignTool** option under **Tools Menu** in the main application window or through its shortcut **Ctrl+M**.



**Figure 7.2.1: The MUMDesignTool configuration window**

**Step 2:** In the – Simulation parameters - section of the window, enter the values for the various fields as follows: *Number of focal planes = 2, Total photon count = 10000, Photon count percentage = [50, 50], Immersion medium refractive index = 1.515, Numerical aperture = 1.45, Wavelength = 0.69, Exposure time = 0.1, Magnification = 100,  $x_0 = 13 \times 4.5 / 100$ ,  $y_0 = 13 \times 4.5 / 100$ , Pixel size = [13, 13], and ROI size = [9, 9].*

**Note:** The center of a  $9 \text{ pixels} \times 9 \text{ pixels}$  pixel array where the dimensions of each pixel is  $13 \mu\text{m} \times 13 \mu\text{m}$ , is  $(4.5 \times 13) \mu\text{m} \times (4.5 \times 13) \mu\text{m}$ . To convert this location coordinates to the object space, the coordinate values are divided by the magnification. Hence the location coordinates are specified as  $13 \times 4.5 / 100 \mu\text{m} \times 13 \times 4.5 / 100 \mu\text{m}$ .

**Step 3:** In the – Defocus range - section of the window, enter the values for the various fields as follows: Defocus start = -2, Defocus increment = 0.05, Defocus end = 2.

**Step 4:** In the configure setting – MUM Design Tool window, click the **Done** button. The **Done** button text will change to **Calculation in progress ...** and the application will temporarily become

unresponsive while the calculations are being performed. During the calculation, a console window will be shown to inform the user about the different steps of the calculation. This console window will also display errors and warnings, if there is any, which is very useful for bug reporting purposes.

**NOTE:** The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.




Figure 7.2.2: MUMDesignTool calculation in progress

**Step 5:** When the calculation is completed, the **MUMDesignTool (Rapid Mode) window** will be brought in focus (more details about this window will be discussed in the next section).

### 7.3. Designing the plane spacing using the Rapid Mode

Figure 7.3.1 shows the main window of the MUMDesignTool. The - Simulation summary - section of the window represents a summary of the simulation parameters. The - Focal plane offset - section of the window includes a pull-down menu for selecting the number of focal planes as well as 10 slides to change the plane spacings. Note that these settings can be performed in real-time. The right-hand side panels, which show the FIM and PLAM along the z-axis for the MUM setup, will be updated automatically by changing the plane spacing and other parameters.

The - PDR - section of the window assists us with choosing the total photon detection rate of the setup as well as the photon percentage for each focal plane. Changing the number of planes will automatically update the photon percentages to be consistent with the total photon detection rate.

The – Display options – section of the window provides controls for better calculation and visualization of the FIM and PLAM, and also provides tools for exporting the results. For instance, using radio buttons  $z_0$ ,  $x_0$  and  $y_0$ , user can switch the FIM and PLAM plots to visualize the behavior of the corresponding coordinate along the  $z$ -axis. The embedded checkboxes help to hide or show different curves within the figures and enhance the visibility. Moreover, a desired region in the plot can be highlighted using other controls in this panel.

The **Reset** button in this panel resets the MUMDesignTool by recalling the **Configure Setting – MUMDesignTool Window**. The **Split windows** button shows the FIM and PLAM plots in new decoupled windows. Also, the **Export** button exports both the results as .txt files and the plots as .pdf files. Finally, the **Precise Mode** button runs the precise mode calculation of the MUMDesignTool to provide accurate results for the behavior of the FIM and PLAM curves (This function will be discussed in detail later).

We next give an example for designing the focal plane spacing for a 3-plane MUM setup using the strong coupling spacing scenario (for more information see [1]). We assume that user has applied the steps in Section 7.2 (Providing the inputs and configuring the MUMDesignTool). Here are the steps:





Figure 7.3.1: MUMDesignTool Rapid Mode (main window)

**Step 1:** In the – Focal plane offset - section of the window select “3” from the Number of focal planes pull-down menu.



**Figure 7.3.2: Changing the number of focal planes in real-time**

**Step 2:** You will see an alert stating that “Photon count percentages have been reset”. Click **OK** to proceed.

**Step 3:** In the – Focal plane offset - section of the window set the value of the plane spacing slider 2 to be equal to “0.6 microns”. Also, set the value of the plane spacing slider 3 to be equal to “1.2 microns”. This will overlap the zero of the  $z_0$ -FIM of each focal plane with the local maximum of that of the adjacent focal plane (for more information see [1]). As a result, the FIM and PLAM curves will be updated.

**Step 4:** In the – Display options - section of the window, select “Sum line” from the View mode pull-down menu. This will update the FIM such that only the sum line data points will be visible.

**Step 5:** In the – Display options - section of the window, change the right-hand side field of Highlight region to “0.6” and check the “Show” checkbox. This will highlight the region between plane 1 and plane 2 for better visualization.

**Step 6:** In the – Limit ruler - section of the window, set the  $z_0$  slider to 19 nm. At this step, you will see the  $z_0$  elements of the FIM and PLAM curves.

**Step 7:** To see the  $x_0$  elements of the FIM and PLAM curves, in the – Display options - section, select the  $x_0$  radio button. The FIM and PLAM curves will be updated accordingly as can be seen in the following figure.




Figure 7.3.3: The lateral elements of FIM and PLAM (MUMDesignTool Rapid Mode)

**Step 8:** To save the calculated data points and the plots click the **Export** button in the – Display options – section. You will be asked to select a folder.





Figure 7.3.4: Exporting the MUMDesignTool results

**Step 9:** In the **Browse For Folder** window, select a folder and click the **OK** button. This will export the results as text files and the plots as pdf files.

## 7.4. Precise Mode

Once user has designed the plane spacing using the Rapid Mode, the next step is to verify whether the designed spacing is appropriate in the presence of extraneous noise sources. This verification is performed using the Precise Mode of the MUMDesignTool. We next provide an example for this verification step.

**Step 1:** Click the **Precise Mode** button in the – Display options – section to calculate the accurate FIM and PLAMs along the z-axis considering the extraneous noise parameters. This will load the **Noise Parameter Setting (Precise Mode)** window as shown below.



**Figure 7.4.1: Entering the noise parameters**

**Step 2:** In **Noise Parameter Setting (Precise Mode)** window, select the “Poisson + Gaussian + Branching” option in the **Noise Type** drop down menu. This will consider background, readout and branching noise for the calculations. This option is suitable when using an EMCCD detector.

**NOTE:** Holding the mouse pointer on this drop down menu will give a brief explanation (as a tool tip text) about the suitable option for different types of detectors (see Figure 7.4.1).

**Step 3:** Enter the values *Background noise* = [35 60 85], the *Readout noise std* = [8 5 7], the *EM gain* = [250 100 110] and the *EM model* = [1 1 1]. Click **OK** when done. The **Precise Mode** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed (see figure below). The console window will also be displayed.

**NOTE:** The calculation could indeed take an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.




Figure 7.4.2: Calculation procedure in the Precise Mode

**Remark 8:** Each element in the *Background noise*, *Readout noise std*, *EM gain* and *EM model* vectors corresponds to one of the focal planes. For example, in the above case, the vectors have 3 elements associated with the 3 focal planes. For more information about the different EM models see [3].

**Step 4:** When the calculation is completed, the **Calculating ...** button text will change back to **Precise Mode** and the **MUMDesignTool (Precise Mode) window** will be brought in focus (see figure below).




Figure 7.4.3: MUMDesignTool Precise Mode

## 8. 3D-Localization accuracy measure (Born-Wolf Model)

This section provides an example of how to use the FandPLimitTool to calculate the accuracy with which the 3D location of an object can be estimated from its image (PLAM). The image of the object is assumed to be a 3D Point Spread Function profile.

The calculation takes into account the stochastic nature of the photon detection processes and assumes the photons detected from the object to be Poisson distributed. In addition, the calculation also takes into account extraneous noise sources - in particular the background component which is modeled as Poisson noise [1], the readout noise of the detector which is modeled as independent, additive Gaussian noise [2] and the electron multiplication noise of an EMCCD detector which is modeled as a branching process [3].

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the ROI containing the image of the object is assumed to be a rectangular array of pixels. The object is assumed to be in the center of this rectangular pixel array and the location of the object is specified with respect to the object space.

Refer to [5] for details on determining the 3D localization accuracy in single molecule microscopy.

### 8.1. Summary of imaging conditions

The following table summarizes the imaging conditions under which the calculations are performed.


Property	Value
$Z_0$	0.5 $\mu\text{m}$
Object medium refractive index	1.515
Numerical aperture	1.3
Wavelength	0.52 $\mu\text{m}$
Photon detection rate	10,000 <i>photons/second</i>
Exposure time	0.13 <i>seconds</i>
Magnification	100 X
Location of object in ROI	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$
Pixel size	13 $\mu\text{m}$ $\times$ 13 $\mu\text{m}$
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>
Background level (Poisson rate)	1000 <i>photons/pixel/second</i>
Gaussian noise mean	0 $e^-$
Gaussian noise standard deviation	8 $e^-$ <i>rms</i>
Use EM	<i>unchecked</i>
EM gain	100
EM model	<i>Geometric</i>

**Table 8.1.1: 3D Localization Accuracy - Summary of imaging conditions**

## 8.2. Providing required parameters

**Step 1:** Select the calculation option “3D Localization Accuracy - Airy: Pixelated with Poisson + Gaussian Noise” from the Calculate pull-down menu in the main application window.

**Step 2:** The required parameters window should be visible by default. If it is not visible, access it by clicking the **Required Parameters** option under the **Parameters Menu**.



**Figure 8.2.1:** 3D Localization Accuracy problem scenario and its required parameters window

**Step 3:** In the - Fundamental inputs - section of the window, enter the values for the various fields as follows:  $z_0 = 0.5$ , *Object medium refractive index* = 1.515, *Numerical aperture* = 1.3, *Wavelength* = 0.52, *Photon detection rate* = 10000, and *Exposure time* = 0.13.

**Step 4:** In the - Additional inputs - section of the window enter the values for the various fields as follows: *Magnification* = 100,  $x_0 = 13 \times 4.5 / 100$ ,  $y_0 = 13 \times 4.5 / 100$ , *Pixel size* = [13, 13], and *ROI size* = [9, 9].

**Remark 9:** The center of a  $9 \text{ pixels} \times 9 \text{ pixels}$  pixel array where the dimensions of each pixel is  $13 \mu\text{m} \times 13 \mu\text{m}$ , is  $(4.5 \times 13)\mu\text{m} \times (4.5 \times 13)\mu\text{m}$ . To convert this location coordinates to the object space, the coordinate

values are divided by the magnification. Hence the location coordinates are specified as  $13 \times 4.5 / 100 \mu\text{m} \times 13 \times 4.5 / 100 \mu\text{m}$ .

**Step 5:** In the - Extraneous noise sources - section of the window, enter the values for the various fields as follows: *Background level = 1000* and *Standard deviation = 8*.


**Remark 10:** In case that the localization accuracy is required for an EMCCD detector, the *Use EM* checkbox should be checked and the desired value for the electron multiplication gain should be entered in the *EM gain* field. For more information see Remark 2.

**Step 6:** In the - Parameters to be estimated - section of the window, ensure that all check-boxes for all parameters are checked.

### 8.3. Providing advanced parameters

**NOTE:** All inputs in the **Advanced Parameters window** are optional. The default values provided have been found to work adequately for typical imaging conditions.


**Step 1:** Access the advanced parameters window by clicking the **Advanced Parameters** option under the **Parameters Menu**.



**Figure 8.3.1: 3D Localization Accuracy advanced parameters window**

**Step 2:** The Advanced parameters window contains an - Advanced fundamental inputs - section. The default value of the field in this section can be left unchanged.


**Step 3:** In the - Model profile calculations - section of the window, access the Pixel integration method pull-down menu and select “Trapezoidal”.



**Figure 8.3.2: Pixel Integration Trapezoidal options**

**Step 4:** Click the **Additional parameters** button next to the Pixel integration method pull-down menu. In the dialog box that appears, enter the values  $XGridding = 13$  and  $YGridding = 13$ . Click OK when done.

**Step 5:** Access the Model integration method pull-down menu in the same section of the window and select “Trapezoidal”.




**Figure 8.3.3: Model Integration Trapezoidal options**

**Step 6:** Click the **Additional parameters** button next to the Model integration method pull-down menu. In the dialog box that appears, enter the values  $Step\ size = 0.0005$ ,  $Lower\ limit = 0$ , and  $Upper\ limit = 1$ . Click OK when done.

**Step 7:** Access the Background type pull-down menu and select “Constant”. There are no additional parameters to be accessed via the Additional parameters button for the “Constant” option.


**Step 8:** In the - Fisher information matrix calculations - section of the window, access the Single integration method pull-down menu and select “Trapezoidal”.



**Figure 8.3.4: FIM Single Integration Trapezoidal options**

**Step 9:** Click the **Additional parameters** button next to the Single integration method pull-down menu. In the dialog box that appears, enter the values *Step size = 0.0005*, *Lower limit = 0*, and *Upper limit = 1*. Click OK when done.

**Step 10:** Access the Double integration method pull-down menu in the same section of the window and select “Trapezoidal”.



**Figure 8.3.5: FIM Double Integration Trapezoidal options**

**Step 11:** Click on the **Additional parameters** button next to the Double integration method pull-down menu. In the dialog box that appears, enter the values *XGridding = 13* and *YGridding = 13*.


**Step 12:** The - Fisher information matrix calculations - section has a subsection titled - Noise factor calculation parameters - . Enter the following values in this subsection: *Limits of summation = [1, 3000]*, *Limits of integration = [-5000, 5000]* and *Trapezoidal integration step size = 0.05*.

**Remark 11:** The limits of integration and summation are initially set to default values which are automatically calculated based on the imaging conditions. It is possible specify these limits by changing the value of the combo box to the right of the limits of integration and summation fields.

## 8.4. Executing the task and viewing results

**Step 1:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed. During the calculation, a console window will be shown to inform the user about the different steps of the calculation.

**NOTE:** The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.



**Figure 8.4.1: Calculation and Final confirmation**

**Step 2:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy** and the **Results window** will be brought in focus. In addition, the console window will be hidden when the calculation is over. The console window can be accessed any time by selecting the **Console** option under the **View Menu**.




Figure 8.4.2: 3D localization Accuracy results window

## 9. 2D-Resolution measure

This section provides an example of how to use the FandPLimitTool to calculate the accuracy with which the 2D location of two objects can be resolved from their image (PREM). The objects are assumed to be in focus with respect to the objective lens and the image of each object is assumed to be an Airy profile.

The calculation takes into account the stochastic nature of the photon detection processes and assumes the photons detected from the object to be Poisson distributed. In addition, the calculation also takes into account extraneous noise sources - in particular the background component which is modeled as additive Poisson noise, the readout noise of the detector which is modeled as independent, additive Gaussian noise [2] and the electron multiplication noise of an EMCCD detector which is modeled as a branching process [3].

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the ROI containing the image of the objects is assumed to be a rectangular array of pixels. The midpoint between the objects is assumed to be at the center of this rectangular pixel array and the location of the midpoint between the objects is specified with respect to the object space.

Refer to [6] for details or determining the 2D resolution accuracy in single molecule microscopy.

### 9.1. Summary of imaging conditions


The following table summarizes the imaging conditions under which the calculations are performed.

Property	Object 1	Object 2
Numerical aperture	1.3	1.3
Wavelength	0.52 $\mu\text{m}$	0.52 $\mu\text{m}$
Photon detection rate	10,000 <i>photons/second</i>	10,000 <i>photons/second</i>
d	0.1 $\mu\text{m}$	
phi	$\pi/4$ radians	
Exposure time	0.13 <i>seconds</i>	
Magnification	100 X	
Location of object in ROI	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$	
Pixel size	13 $\mu\text{m} \times 13 \mu\text{m}$	
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>	
Background level (Poisson rate)	1000 <i>photons/pixel/second</i>	1000 <i>photons/pixel/second</i>
Gaussian noise mean	0 $e^-$	
Gaussian noise standard deviation	8 $e^-$ rms	
Use EM	<i>unchecked</i>	
EM gain	100	
EM model	<i>Geometric</i>	

**Table 9.1.1: 2D Resolution Accuracy - Summary of imaging conditions**

## 9.2. Providing required parameters

**Step 1:** Select the calculation option “2D Resolution Accuracy - Airy: Pixelated with Poisson + Gaussian Noise (Symmetric Parameterization)” from the Calculate pull-down menu in the main application window.



**Figure 9.2.1: 2D Resolution Accuracy calculate type and its required parameters window**

**Step 2:** The required parameters window, which is visible by default, will be updated accordingly. If this window is not visible, access it by clicking the **Required Parameters** option under the **Parameters Menu** in the main application window.

**Step 3:** In the - Object selection - section of the **Required parameters window**, select “Object 1”.

**Step 4:** In the - Fundamental inputs – Object 1 - section of the window, enter the values for the various fields as follows: *Numerical aperture = 1.3, Wavelength = 0.52, Photon detection rate = 10000, d = 0.1, phi = pi/4 and Exposure time = 0.13.*

**Step 5:** In the - Additional inputs - section of the window, enter the values for the various fields as follows: *Magnification = 100*, *sx = 13×4.5/100*, *sy = 13×4.5/100*, *Pixel size = [13, 13]*, and *ROI size = [9, 9]*.

**Remark 12:** The center of a *9 pixels × 9 pixels* pixel array where the dimensions of each pixel is *13 μm × 13 μm*, is *(4.5 × 13)μm × (4.5 × 13)μm*. To convert this location coordinates to the object space, the coordinate values are divided by the magnification. Hence the location coordinates are specified as *13×4.5/100 μm × 13×4.5/100 μm*.

**Step 6:** In the - Extraneous noise sources – Object 1 - section of the window, enter the values for the various fields as follows: *Background level = 1000* and *Standard deviation = 8*.

**Remark 13:** In case that the localization accuracy is required for an EMCCD detector, the *Use EM* checkbox should be checked and the desired value for the electron multiplication gain should be entered in the *EM gain* field. For more information see Remark 2.


**Step 7:** In the - Object selection - section of the **Required parameters window**, select “Object 2”. Then repeat steps 4 and 6 for the Object 2.

**Step 8:** In the - Parameters to be estimated - section of the window, ensure that all check-boxes for all parameters are checked.

### 9.3. Providing advanced parameters


**NOTE:** All inputs in the **Advanced Parameters window** are optional. The default values provided have been found to work adequately for typical imaging conditions.

**Step 1:** Access the advanced parameters window by clicking the **Advanced Parameters** option under the **Parameters Menu**.



**Figure 9.3.1: 2D Resolution Accuracy advanced parameters window**

**Step 2:** In the - Model profile calculations - section of window, access the Pixel integration method pull-down menu and select “Trapezoidal”.




**Figure 9.3.2: Pixel Integration Trapezoidal options**

**Step 3:** Click the **Additional parameters** button next to the Pixel integration method pull-down menu. In the dialog box that appears, enter values for fields as follows:  $XGridding = 13$ ,  $YGridding = 13$ . Click **OK** when done.

**Step 4:** Access the Background type pull-down menu and select “Constant”. There are no additional parameters to be accessed via the **Additional parameters** button for the “Constant” option.

**Step 5:** In the - Fisher information matrix calculations - section of the window, access the Single integration method and select “Trapezoidal”.



**Figure 9.3.3: FIM Single Integration Trapezoidal options**

**Step 6:** Click on the **Additional parameters** button next to the Single integration method pull-down menu. In the dialog box that appears, enter  $Step\ size = 0.05$ .

**Step 7:** Access the Double integration method pull-down menu in the - Fisher information matrix calculations - section of the window and select “Trapezoidal”.

**Step 8:** Click on the **Additional parameters** button next to the Double integration method pull-down menu. In the dialog box that appears, enter the values  $XGridding = 13$  and  $YGridding = 13$ .




Figure 9.3.4: FIM Double Integration Trapezoidal options

**Step 9:** The - Fisher information matrix calculations - section has a subsection titled - Noise factor calculation parameters - . Enter values *Limits of summation* = [1, 3000] and *Limits of integration* = [-5000, 5000] in this subsection.

**Remark 14:** The limits of integration and summation are initially set to default values which are automatically calculated based on the imaging conditions. It is possible specify these limits by changing the value of the combo box to the right of the limits of integration and summation fields.

## 9.4. Executing the task and viewing results

**Step 1:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed.





Figure 9.4.1: Calculation and Final confirmation

**NOTE:** During the calculation, a console window will be shown to inform the user about the different steps of the calculation. The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.

**Step 2:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy** and the **Results window** will be brought in focus. In addition, the console window will be hidden when the calculation is over. The console window can be accessed any time by selecting the **Console** option under the **View Menu**.



**Figure 9.4.2: Results window (2D Resolution)**

## 10. 3D-Resolution measure (Born-Wolf Model)

This section provides an example of how to use the FandPLimitTool to calculate the accuracy with which the 3D location of two objects can be resolved from their image (PREM). The image of each object is assumed to be 3D Point Spread Function.

The calculation takes into account the stochastic nature of the photon detection processes and assumes the photons detected from the object to be Poisson distributed. In addition, the calculation also takes into account extraneous noise sources - in particular the background component which is modeled as additive Poisson noise, the readout noise of the detector which is modeled as independent, additive Gaussian noise [2] and the electron multiplication noise of an EMCCD detector which is modeled as a branching process [3].

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the ROIs containing the images of the object are assumed to be rectangular arrays of pixels. The object is assumed to be in the center of this rectangular pixel array and the location of the object is specified with respect to the object space.

Refer to [6] for details on 3D resolution accuracy in single molecule microscopy.

### 10.1. Summary of imaging conditions

The following table summarizes the imaging conditions under which the calculations are performed.

Property	Object 1	Object 2
Numerical aperture	1.3	1.3
Wavelength	0.52 $\mu\text{m}$	0.52 $\mu\text{m}$
Photon detection rate	10,000 <i>photons/second</i>	10,000 <i>photons/second</i>
Object medium refractive index	1.515	
d	0.1 $\mu\text{m}$	
phi	$\pi/4$ radians	
omega	$\pi/3$ radians	
Exposure time	0.13 <i>seconds</i>	
Magnification	100 X	
Location of object in ROI (sx, sy)	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$	
Axial location of object (sz)	0 $\mu\text{m}$	
Pixel size	13 $\mu\text{m}$ $\times$ 13 $\mu\text{m}$	
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>	
Background level (Poisson rate)	1000 <i>photons/pixel/second</i>	1000 <i>photons/pixel/second</i>
Gaussian noise mean	0 $e^-$	
Gaussian noise standard deviation	8 $e^-$ <i>rms</i>	
Use EM	<i>unchecked</i>	
EM gain	100	
EM model	<i>Geometric</i>	


**Table 10.1.1: 3D Resolution Accuracy - Summary of imaging conditions**

## 10.2. Providing required parameters

**Step 1:** Select the calculation option “3D Resolution Accuracy – Point Spread Function (Born-Wolf Model): Pixelated with Poisson + Gaussian Noise” from the Calculate pull-down menu in the main application window.

**Step 2:** The required parameters window, which is visible by default, will be updated accordingly. If this window is not visible, access it by clicking the **Required Parameters** option under the **Parameters Menu** in the main application window.

**Step 3:** In the - Object selection - section of the **Required parameters window**, select “Object 1”.



**Figure 10.2.1: 3D Resolution Accuracy calculate type and its required parameters window**

**Step 4:** In the - Fundamental inputs – Object 1 - section of the window, enter the values for the various fields as follows: *Numerical aperture = 1.3*, *Wavelength = 0.52*, *Photon detection rate = 10000*, *Object medium refractive index = 1.515*, *d = 0.1*, *phi = pi/4*, *omega = pi/3* and *Exposure time = 0.13*.

**Step 5:** In the - Additional inputs - section of the window, enter the values for the various fields as follows: *Magnification* = 100,  $s_x = 13 \times 4.5 / 100$ ,  $s_y = 13 \times 4.5 / 100$ ,  $s_z = 0$ , *Pixel size* = [13, 13], and *ROI size* = [9, 9].


**Remark 15:** The center of a  $9 \text{ pixels} \times 9 \text{ pixels}$  pixel array where the dimensions of each pixel is  $13 \mu\text{m} \times 13 \mu\text{m}$ , is  $(4.5 \times 13) \mu\text{m} \times (4.5 \times 13) \mu\text{m}$ . To convert this location coordinates to the object space, the coordinate values are divided by the magnification. Hence the location coordinates are specified as  $13 \times 4.5 / 100 \mu\text{m} \times 13 \times 4.5 / 100 \mu\text{m}$ .

**Step 6:** In the - Extraneous noise sources – Object 1 - section of the window, enter the values for the various fields as follows: *Background level* = 1000 and *Standard deviation* = 8.

**Step 7:** In the - Plane selection - section of the **Required parameters window**, select “Object 2”. Then repeat steps 4 and 6 for the Object 2.

### 10.3. Providing advanced parameters

**NOTE:** All inputs in the **Advanced Parameters window** are optional. The default values provided have been found to work adequately for typical imaging conditions.




**Figure 10.3.1: 3D Resolution Accuracy advanced parameters window**

**Step 1:** Access the advanced parameters window by clicking the **Advanced Parameters** option under the **Parameters Menu**.

**Step 2:** Skip the - Object selection - and - Advanced fundamental inputs (Object 1) - sections of the window. Their fields do not need to be changed as they are automatically calculated from the values entered in the Required Inputs window.


**Step 3:** In the - Model profile calculations - section of window, access the Pixel integration method pull-down menu and select “Trapezoidal”.



**Figure 10.3.2: Pixel Integration Trapezoidal options**

**Step 4:** Click the **Additional parameters** button next to the Pixel integration method pull-down menu. In the dialog box that appears, enter values for fields as follows: *XGridding* = 13, *YGridding* = 13. Click **OK** when done.

**Step 5:** Access the Model integration method pull-down menu in the same section of the window and select “Trapezoidal”.




**Figure 10.3.3: Model Integration Trapezoidal options**

**Step 6:** Click the **Additional parameters** button next to the Model integration method pull-down menu. In the dialog box that appears, enter the values *Step size* = 0.0005, *Lower limit* = 0, and *Upper limit* = 1. Click **OK** when done.

**Step 7:** Access the Background type pull-down menu and select “Constant”. There are no additional parameters to be accessed via the **Additional parameters** button for the “Constant” option.


**Step 8:** In the - Fisher information matrix calculations - section of the window, access the Single integration method pull-down menu and select “Trapezoidal”.



**Figure 10.3.4: FIM Single Integration Trapezoidal options**

**Step 9:** Click the **Additional parameters** button next to the Single integration method pull-down menu. In the dialog box that appears, enter the values *Step size* = 0.0005, *Lower limit* = 0, and *Upper limit* = 1. Click **OK** when done.

**Step 10:** Access the Double integration method pull-down menu in the same section of the window and select “Trapezoidal”.



**Figure 10.3.5: FIM Double Integration Trapezoidal options**

**Step 11:** Click on the **Additional parameters** button next to the Double integration method pull-down menu. In the dialog box that appears, enter the values *XGridding* = 13 and *YGridding* = 13.


**Step 12:** The - Fisher information matrix calculations - section has a subsection titled - Noise factor calculation parameters -. Enter the following values in this subsection: *Limits of summation* = [1, 3000], *Limits of integration* = [-5000, 5000] and *Trapezoidal integration step size* = 0.05.

**Remark 16:** The limits of integration and summation are initially set to default values which are automatically calculated based on the imaging conditions. It is possible to specify these limits by changing the value of the combo box to the right of the limits of integration and summation fields.

## 10.4. Executing the task and viewing results


**Step 1:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed.

**NOTE:** The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.



**Figure 10.4.1: Calculation and Final confirmation (3D Resolution)**

**Step 2:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy** and the **Results window** will be brought in focus.



**Figure 10.4.2: Results window (3D Resolution)**

## 11. 3D-Localization accuracy measure - Multifocal Plane

This section provides an example of how to use the FandPLimitTool to calculate the accuracy with which the 3D location of an object can be estimated from the images acquired by a 2-plane MUM setup (PLAM). The images of the object are assumed to be 3D Point Spread Function profiles.

The calculation takes into account the stochastic nature of the photon detection processes and assumes the photons detected from the object to be Poisson distributed. In addition, the calculation also takes into account extraneous noise sources - in particular the background component which is modeled as additive Poisson noise, the readout noise of the detector which is modeled as independent, additive Gaussian noise and the electron multiplication noise which is modeled as a branching process.

It is assumed that the pixels are rectangular in shape with no dead space between any two pixels. Further, the ROIs containing the images of the object are assumed to be rectangular arrays of pixels. The object is assumed to be in the center of this rectangular pixel array and the location of the object is specified with respect to the object space. Moreover, it is assumed that the photon detection process at each detector is independent of that of other detectors. Refer to [1] and [7] for details on determining the 3D localization accuracy in multifocal plane microscopy (MUM).

### 11.1. Summary of imaging conditions

The following table summarizes the imaging conditions under which the calculations are performed.

Property	Plane 1	Plane 2
$Z_0$	0.5 $\mu\text{m}$	
Alpha	15.708	
Object medium refractive index	1.515	
Numerical aperture	1.3	
Wavelength	0.52 $\mu\text{m}$	
Photon detection rate	10,000 <i>photons/second</i>	
Exposure time	0.13 <i>seconds</i>	
Magnification	100 X	
Photon detection rate ratio	[0.5, 0.5]	
Focal plane distance	[0, 0.5] $\mu\text{m}$	
Location of object in ROI	0.585 $\mu\text{m}$ , 0.585 $\mu\text{m}$	
Pixel size	13 $\mu\text{m}$ $\times$ 13 $\mu\text{m}$	
ROI Size	9 <i>pixels</i> $\times$ 9 <i>pixels</i>	
Background level (Poisson rate)	25 <i>photons/pixel/second</i>	25 <i>photons/pixel/second</i>
Gaussian noise standard deviation	8 $e^-$ <i>rms</i>	
Use EM	<i>unchecked</i>	
EM gain	100	
EM model	<i>Geometric</i>	

**Table 11.1.1: 3D Localization Accuracy (MUM) - Summary of imaging conditions**


## 11.2. Providing required parameters

**Step 1:** Select the calculation option “3D Localization Accuracy - Multifocal: Pixelated with Poisson + Gaussian Noise” from the Calculate pull-down menu in the main application window.

**Step 2:** The required parameters window, which is visible by default, will be updated accordingly. If this window is not visible, access it by clicking the **Required Parameters** option under the **Parameters Menu** in the main application window.

**Step 3:** In the - Plane selection - section of the **Required parameters window**, select “Plane 1”.

**Step 4:** In the - Fundamental inputs - Plane 1 - section of the window, enter the values for the various fields as follows:  $z_0 = 0.5$ ,  $\text{Alpha} = 15.708$ ,  $\text{Object medium refractive index} = 1.515$ ,  $\text{Numerical aperture} = 1.3$ ,  $\text{Wavelength} = 0.52$ ,  $\text{Photon detection rate} = 10000$ , and  $\text{Exposure time} = 0.13$ .



**Figure 11.2.1: 3D Localization Accuracy (MUM) calculate type and its required parameters window**

**Step 5:** In the - Additional inputs - section of the window, enter the values for the various fields as follows:  $\text{Magnification} = 100$ ,  $\text{Photon detection rate ratio} = [0.5, 0.5]$ ,  $\text{Focal plane distance} = [0, 0.5]$ ,  $x_0 = 13 \times 4.5 / 100$ ,  $y_0 = 13 \times 4.5 / 100$ ,  $\text{Pixel size} = [13, 13]$ , and  $\text{ROI size} = [9, 9]$ .

**Remark 17:** The center of a  $9 \text{ pixels} \times 9 \text{ pixels}$  pixel array where the dimensions of each pixel is  $13 \mu\text{m} \times 13 \mu\text{m}$ , is  $(4.5 \times 13)\mu\text{m} \times (4.5 \times 13)\mu\text{m}$ . To convert this location coordinates to the object space, the coordinate values are divided by the magnification. Hence the location coordinates are specified as  $13 \times 4.5 / 100 \mu\text{m} \times 13 \times 4.5 / 100 \mu\text{m}$ .

**Step 6:** In the - Extraneous noise sources – Plane 1 - section of the window, enter the values for the various fields as follows: *Background level* = 25 and *Standard deviation* = 8.

**Step 7:** In the - Plane selection - section of the **Required parameters window**, select “Plane 2”. Then repeat step 6 for Plane 2.


**Remark 18:** In case that the localization accuracy is required for an EMCCD detector, the *Use EM* checkbox should be checked and the desired value for the electron multiplication gain should be entered in the *EM gain* field. For more information see Remark 2.

**Step 8:** In the - Parameters to be estimated - section of the window, ensure that all check-boxes for all parameters are checked.

### 11.3. Providing advanced parameters

**NOTE:** All inputs in the **Advanced Parameters window** are optional. The default values provided have been found to work adequately for typical imaging conditions.


**Step 1:** Access the advanced parameters window by clicking the **Advanced Parameters** option under the **Parameters Menu**.



**Figure 11.3.1: 3D Localization Accuracy (MUM) advanced parameters window**

**Step 2:** The Advanced parameters window contains an - Advanced fundamental inputs - section. The default value of the field in this section can be left unchanged.


**Step 3:** In the - Model profile calculations - section of window, access the Pixel integration method pull-down menu and select “Trapezoidal”.



**Figure 11.3.2: Pixel Integration Trapezoidal options**

**Step 4:** Click the **Additional parameters** button next to the Pixel integration method pull-down menu. In the dialog box that appears, enter values for fields as follows:  $XGridding = 13$ ,  $YGridding = 13$ . Click **OK** when done.

**Step 5:** Access the Model integration method pull-down menu in the same section of the window and select “Trapezoidal”.




**Figure 11.3.3: Model Integration Trapezoidal options**

**Step 6:** Click the **Additional parameters** button next to the Model integration method pull-down menu. In the dialog box that appears, enter the values  $Step\ size = 0.0005$ ,  $Lower\ limit = 0$ , and  $Upper\ limit = 1$ . Click **OK** when done.

**Step 7:** Access the Background type pull-down menu and select “Constant”. There are no additional parameters to be accessed via the **Additional parameters** button for the “Constant” option.


**Step 8:** In the - Fisher information matrix calculations - section of the window, access the Single integration method pull-down menu and select “Trapezoidal”.



**Figure 11.3.4: FIM Single Integration Trapezoidal options**

**Step 9:** Click the **Additional parameters** button next to the Single integration method pull-down menu. In the dialog box that appears, enter the values *Step size = 0.0005*, *Lower limit = 0*, and *Upper limit = 1*. Click OK when done.

**Step 10:** Access the Double integration method pull-down menu in the same section of the window and select “Trapezoidal”.



**Figure 11.3.5: FIM Double Integration Trapezoidal options**

**Step 11:** Click on the **Additional parameters** button next to the Double integration method pull-down menu. In the dialog box that appears, enter the values *XGridding = 13* and *YGridding = 13*.

**Step 12:** The - Fisher information matrix calculations - section has a subsection titled - Noise factor calculation parameters - . Enter the following values in this subsection: *Limits of summation = [1, 3000]*, *Limits of integration = [-5000, 5000]* and *Trapezoidal integration step size = 0.05*.

**Remark 20:** The limits of integration and summation are initially set to default values which are automatically calculated based on the imaging conditions. It is possible specify these limits by changing the value of the combo box to the right of the limits of integration and summation fields.

## 11.4. Executing the task and viewing results

**Step 1:** In the main application window, click the **Calculate Limits of Accuracy** button. This will present a dialog box warning that the calculation may take a long time to complete. Click **Yes** to proceed. The **Calculate** button text will change to **Calculating ...** and the application will temporarily become unresponsive while the calculations are being performed. During the calculation, a console window will be shown to inform the user about the different steps of the calculation.

**NOTE:** The calculation could indeed take half an hour to complete depending on the hardware capabilities of the platform on which the application is being executed.




Figure 11.4.1: Calculation and Final confirmation (3D localization MUM)

**Step 2:** When the calculation is completed, the **Calculating ...** button text will change back to **Calculate Limits of Accuracy** and the **Results window** will be brought in focus. In addition, the console window will be hidden when the calculation is over.




Figure 11.4.2: Results window (3D localization MUM)

## References

- [1] A. Tahmasbi, S. Ram, J. Chao, A. V. Abraham, F. W. Tang, E. S. Ward, and R. J. Ober, "Designing the focal plane spacing for multifocal plane microscopy," *Opt. Express*, vol. 22, no. 14, pp. 16706-16721, 2014.
- [2] R. J. Ober, S. Ram, and E. S. Ward, "Localization accuracy in single molecule microscopy," *Biophys. J.*, vol. 86, pp. 1185-1200, 2004.
- [3] J. Chao, E. S. Ward, and R. J. Ober, "Fisher information matrix for branching processes with application to electron-multiplying charge-coupled devices," *Multidim. Sys. Sig. Proc.*, vol. 23, no. 3, pp. 349-379, 2012.
- [4] J. Chao, E. S. Ward, E., and R. J. Ober, "A software framework for the analysis of complex microscopy image data," *IEEE Trans. Inform. Technol. Biomed.*, vol. 14, no. 4, pp. 1075-1087, 2010.
- [5] S. Ram, R. J. Ober, and E. S. Ward, "How accurately can a single molecule be localized in three dimensions when imaged through a fluorescence microscope?" *Proc. of SPIE*, vol. 5699, pp. 426-435, 2005.
- [6] S. Ram, E. S. Ward, and R. J. Ober, "Beyond Rayleigh's criterion: A resolution measure with application to single-molecule microscopy," *Proc. Natl. Acad. Sci. USA*, vol. 103, no. 12, pp. 4457-4462, 2006.
- [7] S. Ram, A. V. Abraham, E. S. Ward, and R. J. Ober, "A novel 3d resolution measure for optical microscopes with applications to single molecule imaging," *Proc. of SPIE*, vol. 6444, 2006.