



Products 2012-2013

AURICULAR THERAPY

AURICULO MEDECINE

ACUPUNCTURE AND TCM

REFLEXO, PHYSIOTHERAPY

PAIN TREATMENT

LIGHT AND LASER THERAPY

TRAINING MATERIAL ...





An **efficient and dynamic** team in **Lyon, France**, to answer to all your questions in french, english, german, spanish, portuguese or italian language !

Address :

Chemin des Mûriers
F - 69540 Irigny – France

Tél. : +33 (0) 472 663 322
Fax. : +33 (0) 478 508 903

Email : contact@sedatelec.com

Our website : <http://www.sedatelec.com>

Opening hours : from Monday to friday 9h00-12h30 and 14h00-17h30

Sedatelec customers service :

Christophe Suffissais, *France, Italy Commercial Manager*
Marketing and Business developper

Tél.: **+33 (0) 472 663 328**
Email : christophesuffissais@sedatelec.com

Martine Ribes, *Export Commercial Manager*

Tél.: **+33 (0) 472 663 324**
Email : martineribes@sedatelec.com

Mickaël Renard, *France Commercial*

Tél.: **+33 (0) 472 663 321**
Email : mickaelrenard@sedatelec.com

Bettina Samitsch, *France Assistante Commercial*

Tél.: **+33 (0) 472 663 322**
Email : bettinasamitsch@sedatelec.com

Our vocation:

To serve and accompany the expansion of Complementary Medicines in the field of **Auriculotherapy, Auriculomedicine, Acupuncture.**

To be **close to our customers**, Sedatelec has a significant network of specialised dealers **worldwide.**

All the team is at your disposal to inform, to advise and help you to choose the product you need in your practice.

Our motto: "only what you need for your patients".

Summary

| | |
|--|-------------|
| Sedatelec, our fields of expertise..... | p.6 |
| AURICULOTHERAPY | p.7 - p.16 |
| Instruments and didactic support | p.7 |
| Rubber ears, ear stamps | |
| Pressure feelers, Glass rods, Angle protractor, anti tobacco protractor | |
| Patient grounding clip | |
| CD Audio Dr Paul NOGIER, software and CD-Rom | p.9 |
| ASP [®] needles | p.10 |
| Ear Semi-Permanent needles - ASP [®] Original | p.10 |
| ASP [®] Original, Classic, Gold, Titanium, Adhesive plasts and Test box | |
| Ear Permanent needles – ASP [®] Perma | p.13 |
| ASP [®] Perma Regular (Ø 0,8 mm), ASP [®] Perma Small (Ø 0,6 mm) | |
| Electrical point-finders and stimulators | p.14 |
| AGISCOP D [®] , AGISCOP DT [®] , AGISCOP D [®] Gold Silber, AGISCOP DT [®] Gold Silber | |
| Electrical point-finders and stimulators Modulo [®] 100 | p.15 |
| Modulo [®] 100 Auriculo GS, Modulo [®] 100 Auriculo GS Automatic, Stimulators in Auriculotherapy micro-currents, luminous and Laser | |
| Electronic device for Moxibustion (Auriculotherapy and Acupuncture) | p.26 |
| Premio [®] 10 moxa | |
| AURICULOMEDICINE | p.17 - p.21 |
| Instruments | p.17 |
| Black-White, North-South, Gold-Silver, Gold-Aluminium detectors, Dim bright lamp (Heine lamp or DBL), Polartrons North, South, Unidirectional North, Unidirectional South, large magnet, VAS clip with or without cable | |
| Luminous detectors and stimulators | p.19 |
| DB 165 (detector with beating), GIR 30 + | |
| Detectors and stimulators Modulo [®] 200 | p.21 |
| Modulo [®] 200 Auriculomedicine, Detectors and stimulators in Auriculomedicine micro laser, luminous and micro-currents | |
| FILTERS AURICULOMEDICINE, CHROMOTHERAPY | p.22 - p.25 |
| NOGIER Colour filters | p.22 |
| Colour filters, Frequential Colour Programme Paul Nogier, filters L, Specific Filters, Polaroid plates, Test-rings for filters, storage box for filters | |

Summary (following)

| | |
|--|-------------|
| ACUPUNCTURE | p.26 - p.33 |
| Electronic device for Moxibustion (Auriculotherapy and Acupuncture) | p.26 |
| Premio [®] 10 moxa | |
| Detectors and stimulators micro-currents, Laser Modulo [®] 100 | p.28 |
| Modulo [®] 100 Acupuncture | |
| Detectors and stimulators micro-currents, Laser Modulo [®] 100 for Reflexo – Acupuncture | p.29 |
| Modulo [®] 100 Micro-currents Reflexo Acutrodes, Modulo [®] 100 Micro-currents Reflexo Needles | |
| Modulo [®] 100 Laser stimulators | p.29 |
| Modulo [®] 100 Acu 15 W Pulsed IR LASER, Modulo [®] 100 Acu 40 W Pulsed IR LASER, Modulo [®] 100 Acu 75 W Pulsed IR LASER, Modulo [®] 100 Acu 40 W Superpulsed IR LASER | |
| Dual Electric Stimulator, Electro-Acupuncture | p.31 |
| AGISTIM duo, Styl Modul, kit 4 needles module, reusable electrodes for stimulation, Electrolipolisis sterile needles ASE, needles or electrodes connectors | |
| LASERTHERAPY, LOCAL TREATMENT | p.34 - p.37 |
| Modulo [®] 200 for LASER Therapy | p.34 |
| Modulo [®] 200 Pulsed IR LASER 15 W, Modulo [®] 200 Pulsed IR LASER 40 W, Modulo [®] 200 Pulsed IR LASER 75 W, Modulo [®] 200 Super Pulsed IR LASER 40 W, | |
| LASER for Physiotherapy | p.36 |
| Modulo [®] 200 Physio Superpulsed IR LASER | |
| LASER for Dental | p.36 |
| Modulo [®] 200 Dental Fiber Superpulsed IR LASER | |
| FREQUENTIAL THERAPY MICRO-CURRENTS | p. 37 |
| Modulo [®] 200 Micro-currents Acutrodes, Modulo [®] 200 Micro-currents Acuneedles | |
| MODULO [®] (retail and examples of Modulo [®] 's personalisation) | p.38 - p.41 |
| Spare parts for Modulo [®] | p.38 |
| Examples of combinations with Modulo [®] range | p.42 |
| IRIDOLOGY, SYMPATHICOTHERAPY | p.44 |
| Iriscope, endo-nasal stylets and tipped applicators | |
| CHARTS AND BOOKS | p.45 - p.49 |
| INDEX | p.50- p.51 |
| SALES CONDITIONS | p.52 |
| ORDER | p.53 |



About us :

Your **French medical equipment manufacturer** is located in the immediate vicinity of **Lyon, France**, cradle of **NOGIER Auriculo**.

We are a company specializing in the **design** and **distribution of medical equipment**. Our staff about 20 people, exporting about 80% of our production in more than forty countries, in Europe, the American continent and elsewhere in the world. Since lot of years, Sedatelec **designs, develops, manufactures** and distributes mainly three product lines:

- **Sterile semi-permanent needles for single use** for ear Acupuncture,
- **Electronic devices for diagnosis and treatment**: measure of body's impedance, electrical stimulators, luminous stimulators of different wavelengths, generators of electric and magnetic fields, laser biostimulation, electronic "moxa-like", device for frequencial phototherapy,
- **Instruments and small devices for Auriculo and Acupuncture** (feelers, filters...),

In complement with educational materials, works, charts and books relating to those techniques.

Our manufacturing Department makes every effort to provide you with quality products in the shortest time.

Quality:

In order to guarantee constant product's quality and to meet the requirements of materiovigilance and traceability, Sedatelec has set up a quality organization in accordance with the reference standard ISO 13485. An independent notified body inspects this organization annually: the TÜV SÜD Product Service.



Sedatelec products carry the EC or CE₀₁₂₃ mark. Most of our devices are approved for marketing in Canada. ASP needles and Premio 10 Moxa are also authorised for sale in the United States.

Guarantee:

Our appliances have 2 years guarrantee. Small instruments 1 year.

After Market Service:

We repair all our appliances up to 10 years after ceasing manufacture. We can provide you with a large number of spare parts.

Special quotation:

Orders are considered by our Customers Relationships Department; Feel free to consult us...

Fast Service:

Our products are generally available from stock and we have a network of distributors present in many countries.

Our priority remains the same: the confidence and satisfaction of our customers and partners.

Instruments and didactic support for Auriculotherapy

Original Rubber ears (left and right):

Very realistic and incomparable mould of the ear, perfect to show the points to patients or students. Usefull for patient education purposes, these ears can also be used by inexperienced therapists to master their knowledge of the external ear and anatomical identification of the points. The excellent quality rubber allow multiple needle punctures without deteriorating the ear. They are available left or right versions.

Size: 85 x 55 x 40 mm

Weight: 78 g




| Ref. | Denomination |
|-------------|------------------|
| OREG | Left rubber ear |
| ORED | Right rubber ear |

Set of 1, supplied under plastic bag

Ear stamps (wood or self-inking):

Useful during consultations to establish and follow-up each patient's medical file, its record the points treated. The quality of impression highly is appreciate by therapists, 2 models of imprints allowing, wood model or self-inking model:

N° 1



| Ref. | Denomination |
|--------------|--------------------|
| TAMB1 | Wood ear stamp N°1 |
| TAMB2 | Wood ear stamp N°2 |

Size: 116 x 61 x 52 mm, weight: 153 g

Set of 1, supplied under plastic bag

N° 2



| Ref. | Denomination |
|--------------|---------------------------|
| TAME1 | Self-inking ear stamp N°1 |
| TAME2 | Self-inking ear stamp N°2 |

Size: 115 x 55 x 66 mm, weight: 100 g



N° 1



N° 2

Set of 1, supplied under plastic bag

Instruments and didactic support for Auriculotherapy

Pressure feelers:

Used to **localize painful points on the ear** in Auriculotherapy (250g) or for the investigation of middle tissue in Auriculomedicine (130g). These are essential instruments to begin in auriculo therapy. 400g (beige) and 600g (orange) for trained therapists.

4 pens are available, from low (130g) to high pressure (650g)

Pressure feelers black, blue, orange and beige:

- ❖ **Black** : 130 grams, to investigate on middle tissues
- ❖ **Blue** : 250 grams, to detect acute zones
- ❖ **Orange** : 400 grams, to detect insensitive zones
- ❖ **Beige** : 650 grams, to detect insensitive zones



| Ref. | Denomination |
|---------------|---|
| PAL250 | Pressure feeler 250 g (blue) |
| PAL130 | Pressure feeler 130 g (black) |
| PAL400 | Pressure feeler 400 g (orange) |
| PAL650 | Pressure feeler 650 g (beige) |

Set of 1, supplied under plastic bag

Size: 145 x 9 mm Weight: 25 g

Glass rods:

Thanks to their fine, smooth tip, they are used to massage very precise auricular points. Made of polished glass they offer excellent quality stimulation without risk of irritating the skin.



| Ref. | Denomination |
|--------------|-----------------------|
| BAT02 | Glass rods (set of 2) |

Set of 2, supplied under plastic bag

Size: 134 x 6 mm Weight: 8 g

Protractors (30°, anti-tobacco):

Diagnostic support material :

One transparent protractor with angles traced at 30° to apply ear geometry

The transparent protractor with 7 grooves designed for the anti- tobacco programme



| Ref. | Denomination |
|--------------|-------------------------|
| RAP30 | Angle protractor 30° |
| RAPAT | Anti-tobacco protractor |

Set of 1, supplied under plastic bag
Size: 140 x 110 mm, weight: 48 g

Instruments and didactic support for Auriculotherapy

Patient grounding clip:

Accessory used **to flow electrostatic charges from the patient to earth**. Contains a fused safety device.



| Ref. | Denomination |
|------------|------------------------|
| MAT | Patient grounding clip |

| Ref. | Spare part Denomination |
|---------------|-------------------------|
| MAT-01 | Grounding clip fuse |

Wire length: 2,8 metres, weigh: 143,5 g
Set of 1, supplied with one additional fuse under plastic bag

Software, CD-ROM NATOM:

NATOM Auriculo: Anatomy and Auriculotherapy software (Dr Rouxeville, Dr Vulliez)
1300 anatomic charts, 70 Auriculotherapy charts with possibility to modify, save, send pictures, insert in patients data base.

French, English, German / PC version




| Ref. | Auteurs - Titres - Editeurs |
|------------------------------|--|
| * AVNATOM01-PC-FR | Logiciel d'Anatomie et d'Auriculothérapie Ch. Vulliez - Y. Rouxeville Edition NATOM 1.2 /PC – en français |
| * AVNATOM01-PC-AN | Anatomy and Auriculotherapy software Ch. Vulliez - Y. Rouxeville Edition NATOM 1.2 /PC – in english |
| * AVNATOM01-PC-AL | Anatomy and Auriculotherapy software Ch. Vulliez - Y. Rouxeville Edition NATOM 1.2 /PC – in german |
| * AVNATOM01-MAJ-PC-FR | Up dated software of Anatomy and Auriculotherapy Edition NATOM |

(*) Limited to our stock
Possibility of multiple settings at the request

CD Audio Dr Paul NOGIER:

Extracts of seminary between 1976 and 1981 in MP3 format, published by Dr André Lentz, **10 hours** registration **Dr Paul NOGIER** historical seminars (in French only).



| Ref. | Denomination |
|--------------------|--|
| AVPNOG01-FR | Dr Paul NOGIER - Extraits séminaires 1976-1981 Edité par Dr André LENTZ |

Ear Semi-Permanent Needles - ASP[®] Original

ASP[®] Original Range:

The reference in ear acupuncture by Dr. Paul NOGIER

Daily used by many practitioners, hospitals and pain-relief clinics for more than 40 years in the world. Created in France in collaboration with Dr. Paul NOGIER, Ear Semi-Permanent needle by Sedatelec became the reference for Auriculotherapy. With many years of research and improvement, **ASP[®] needle has reached a level of incomparable quality.** This small needle ASP[®] Original is housed in a plastic applicator to facilitate insertion, exactly on the point and in the best conditions of sterility. All laboratory tests are performed by an independent laboratory. They are CE 0123 marked.

The **legendary "click"** allows a precise, quick and perfectly calibrated implantation. The practitioner can have full confidence in the ASP[®] allowing full attention to the patient.

- ❖ A choice of 3 packages : box of 8 (Packsolo), easier to carry in displacement, 80, 200 needles, with adhesive plasts, respectively 9, 120, 240 plasts.
- ❖ Single use disposable needles, they are sterilised by gamma rays (25 kGy).
- ❖ **Excellent reproductibility** and effects
- ❖ Blisters pack of 8 needles with batch number and expiry date.
- ❖ A choice of 3 medical grade metals: **Steel Classic, Gold or Titanium.**
- ❖ **Stimulation magnet** designed to activate periodically the action of the ASP[®] needle only Classic or Gold version.
- ❖ The "click" of Sedatelec guarantee of injector high quality
- ❖ **A calibrated pressure** generates necessary impulse to implant.

Implantation is fast, very precise and without risk of infection, if the practitioner has taken care to correctly disinfect the patient's ear.
70° alcohol is a very easy and efficient mean subject to wait for 30 seconds to allow evaporation.

Some use indications:

Permanent algies
Psychological disorders (depression, anxiety, etc.)

Chronic functional disorders
Addiction disorders (alcohol, tobacco cessation, etc.)

Obesity, weight loss
Acute pain (migraine, neuralgia, sciatic, etc.)
Skin and respiratory allergy, chronic pain (zona, neuropathy, bone metastases)
Side effects of chemotherapy

Some heart disease (angina, heart failure, etc.)
Certain dermatological diseases, ENT, ophthalmic
Stress-related disorders (insomnia, tiredness, stomach pains, hypertension, ...)



Ear Semi-Permanent needles - ASP® Original

ASP® Original Classic:

Ear Semi Permanent needle in stainless steel without nickel that has demonstrated its effectiveness for more than 40 years. Leaflet instruction for use and respectively 9, 120, 240 adhesive plasts are provided in each box.



Box of 80

| Ref. | Denomination |
|---------------|------------------------------|
| ASPA80 | ASP® Original Classic (1x80) |

Size: 170x57x45mm Weight: 95 g

Box of 200



| Ref. | Denomination |
|----------------|-------------------------------|
| ASPA200 | ASP® Original Classic (1x200) |

Size: 172x117x56mm Weight: 225 g



Box of 8

| Ref. | Denomination |
|-----------------|-----------------------------|
| ASPASOLO | ASP® Original Classic (1x8) |

Size: 161x57x9mm Weight: 25 g

| Ref. | Denomination |
|-------------------|--------------------------|
| ASPASOLO25 | ASP® SOLO Classic (25X8) |

Size: 255x170x65mm Weight: 700g

ASP® Original Gold:

Ear Semi Permanent needle in stainless steel covered with pure gold, used to stimulate « the point in gold » and some allergic patient. Leaflet instruction for use and respectively 9, 120, 240 adhesive plasts are provided in each box.



Box of 80

| Ref. | Denomination |
|---------------|--------------------------|
| ASPG80 | ASP Original Gold (1x80) |

Size: 170x57x45mm Weight: 95 g

Box of 200



| Ref. | Denomination |
|----------------|---------------------------|
| ASPG200 | ASP Original Gold (1x200) |

Size: 172x117x56mm Weight: 225 g



Box of 8

| Ref. | Denomination |
|-----------------|-------------------------|
| ASPGSOLO | ASP Original Gold (1x8) |

Size: 161x57x9 mm Weight: 25 g

| Ref. | Denomination |
|-------------------|----------------------|
| ASPGSOLO25 | ASP SOLO Gold (25X8) |

Size: 255x170x65mm Weight: 700g

Ear Semi-Permanent needles - ASP[®] Original

ASP[®] Original Titanium:

Ear Semi Permanent needle of medical character, its use is mainly oriented to clinics and hospitals. Produced in Titanium medical grade for a maximum level of biocompatibility, for the treatment of patient with heavy health conditions (diabete, immunocompromised or with multiple allergies). **It can remain in place under MRI examination.**

Leaflet instruction for use and respectively 9, 120, 240 adhesive plasts are provided in each box.



Box of 80

| Ref. | Denomination |
|---------------|------------------------------|
| ASPT80 | ASP Original Titanium (1x80) |

Size: 170x57x45mm Weight: 85 g

Box of 200



| Ref. | Denomination |
|----------------|-------------------------------|
| ASPT200 | ASP Original Titanium (1x200) |

Size: 172x117x56mm Weight: 220 g



Box of 8

| Ref. | Denomination |
|-----------------|-----------------------------|
| ASPTSOLO | ASP Original Titanium (1x8) |

Size: 161x57x9mm Weight: 25 g


25 boxes of 8 :

| Ref. | Denomination |
|-------------------|--------------------------|
| ASPTSOLO25 | ASP SOLO Titanium (25X8) |

Size: 255x170x65mm Weight: 680g

Adhesive plasters ASP[®] Original:

ASP plasts help to maintain the needle in place after the implantation. They come into contact with the patient's skin and have undergone biocompatibility tests, which have demonstrated their hypoallergenic properties. Adhesive plasts are also provided separately, packaged in units of 600 (5*120).



| Ref. | Denomination |
|-------------------|--|
| ASPLAST600 | Adhesive plasters ASP Original (5x120) |

Size: 165x90 mm

Weight: 15 g

ASP[®] Original Test box:

Samples to introduce, discover, or try all of **ASP[®] Original** line (ASP[®] Original Classic, ASP[®] Original Gold, ASP[®] Original Titanium). Available at Sedatelec or Sedatelec distributors.



| Réf. | Denomination |
|--------------------|--------------------------------|
| ASPTTEST-LL | Test box discover ASP Original |

Size: 161x57x9mm

Weight: 25 g

The ASP[®] Original box test unit, ref ASPTTEST-LL, contains 4 needles ASP[®] Original Classic + 2 ASP[®] Original Gold + 2 ASP[®] Original Titanium

Ear Permanent Needle – ASP Perma

The ASP® Perma overrides a renewed pose of ASP® Original (Classic, Gold or Titanium).

Indications for the ASP Perma are **chronic pathologies** such as permanent pain, mental disorders, persistent neurological problems, addiction and chronic functional disorders.

ASP Perma needles must be used by medical practitioners familiar to the techniques of auricular acupuncture using implants.

This small needle is housed in a plastic applicator to facilitate insertion, exactly on the point and in the best conditions of sterility. the ASP Perma are sterilized by radiation Gamma. All biological tests are performed by an independent laboratory. They are CE marked 0123.

- ❖ **Sustainable auricular stimulation for diseases or chronic pains.**
- ❖ 0,8 or 0,6 mm needles diameters produced in titanium medical grade, 100% biocompatible.
- ❖ Subcutaneous implantation with the unit **piston metal** applicator, **precise and sterile.**
- ❖ Possible withdrawal by small surgery gesture. Respect the precaution in use.
- ❖ 2 packages : **8 or 80 needles.**

Implantation is fast, very precise and without risk of infection, if the practitioner has taken care to correctly disinfects the patient's ear. 70° alcohol is a very easy and efficient mean subject to wait for 30 seconds to allow evaporation. In addition, cleaning and disinfecting the pinna of the ear with alcohol, greatly facilitates the subsequent use of any electrical detector to precisely locate the points. The healing is very rapid.

Two diameters are available: Regular 0,8 mm and Small 0,6 mm

ASP® Perma Regular Titanium:

Diameter: Ø 0,8 mm (dark blue colour) with a choice of two packages in box 8 or 80 needles




Box of 8

| Ref. | Denomination |
|---------------|--------------------------|
| ASPP8 | ASP Perma Regular (1x8) |
| ASPP80 | ASP Perma Regular (1x80) |

Size: 161x57x9mm Weight: 23 g

ASP® Perma Small Titanium:

Diameter: Ø 0,6 mm (light blue colour) with a choice of two packages in box 8 or 80 needles



| Ref. | Denomination |
|----------------|------------------------|
| ASPPS8 | ASP Perma Small (1x8) |
| ASPPS80 | ASP Perma Small (1x80) |

Size: 161x57x9mm Weight: 22 g

Electrical point-finders and stimulators

AGISCOP D[®], AGISCOP DT[®]:

AGISCOP D[®] or AGISCOP D Gold - Silber[®] are able to measure a wide impedance range.

These are devices equipped with dual impedance measurement systems allowing **an ear point** or **acupuncture point** to be localised precisely.

Differential detection remains imperative for the measurement and the localisation of the point in the pinna of the ear.

The ERT nozzle allows simplified detection, which is particularly suitable for searching acupuncture points.

Two positions are used in electrical detection:

- ✓ **Position +** is the most common detection position. It allows points to be localised by measuring falls in impedance compared to the environment.
- ✓ **Position -** little used, allows increases in impedance to be measured. It is not recommended for beginners.

A single sensor combined with very fine sensitivity adjustment allows any point to be detected: acupuncture point with high sensitivity setting (5-6 or ERT nozzle) or pathological ear points (1-3).

When a point is detected, the instrument emits an audible beep and the pilot light indicator flashes.

AGISCOP DT[®] or AGISCOP DT Gold - Silber[®] version : the device has **4 treatment positions** with the probe: **2 frequencies for treatment and 2 level of energy.**

- ✓ The high frequency / high-energy position is little used on the ear. The electrical energy is occasionally painful on some sensitive points of the pinna of the ear. In this case, move to a lower energy level.
- ✓ Positioned to the right: **frequency F** (73 Hz)
- ✓ Positioned to the left: **frequency U** (1.14 Hz)

Treatment 30 seconds to 1 minute per point



| Ref. | Denomination |
|--------------|---------------------------|
| AGID | AGISCOP D (complete kit) |
| AGIDT | AGISCOP DT (complete kit) |

| Technical characteristics | |
|---------------------------|---|
| Detection | Current intensity: 5 µA |
| Treatment | Frequencies: 1,14 et 73 Hz |
| Energy per pulse | Low level: 0.15µJ High level: 1.5 µJ |
| Power supply | 2 alkaline batteries or accu 9V |
| Consumption | 6-8 mA |
| Size | 120 x 68 x 45 mm |
| Weigh | 347 g |

Delivered in a carrying case with a special acupuncture nozzle (ERT), 2 batteries and manual.

Electrical point-finders and stimulators


AGISCOP D[®] and AGISCOP DT[®] version Gold Silber:

An AGISCOP D Gold - Silber[®] and AGISCOP DT Gold - Silber[®] version exists specially created for the German school to localise points with High or low impedance compared to the environment.

Differential detection remains essential for the measurement and the localisation of points in the pinna of the ear.

| Ref. | Denomination |
|----------------|--------------------------|
| AGIDGS | AGISCOP D Gold - Silber |
| AGIDTGS | AGISCOP DT Gold - Silber |

Spare for AGISCOP D[®], DT[®], D[®] GS, DT[®] GS :



| Ref. | Spare parts denomination |
|---------------|----------------------------|
| AGI-01 | White sensor AGISCOP |
| AGI-02 | Sensor wire AGISCOP |
| AGI-03 | Grounding electrode + wire |
| AGI-04 | ERT nozzle AGISCOP |

Modulo[®] 100 Auriculo:

NEW

Essential for Auriculotherapy.

Excellent quality of electrical differential detection, identification, measure and qualification of Auriculotherapy points, by selecting either fall or increase of impedance differential (differential electric module).

Precise indication of "pathological" points to be treated.



Visualisation of impedance variation through a useful LED front display.

A variable threshold (red LED), pre-calibrated that can be adapted to each patient, adds a sound signal to the display. It is thus possible to select the points that truly meet the search criteria.

The **electrical stimulation** of points occurs preferentially with "**Scanning**" mode rising and falling through the 7 NOGIER frequencies (↔).

The **stimulation** is also possible in low Frequency 1.14 Hz (LF) or High Frequency 73 Hz (HF).

2 levels of energy are possible in micro-currents stimulation, normal or weak: Energy+ or -
A sound signal is emitted every 30 seconds to inform the time of stimulation.

Electrical point-finders and stimulators

Modulo® 100 Auriculo (following)

Essential for Auriculotherapy.

Stimulation by successive scans↔ of the 7 fundamental Frequencies NOGIER.



- ❖ **Micro – impulsions treatment**
- ❖ **Stimulation with low frequency NOGIER (1,14Hz)**
- ❖ **Stimulation with high frequency NOGIER (73Hz)**
- ❖ **Quick treatment (around 30 seconds)**
- ❖ Non invasive
- ❖ Excellent alternative to needles

| Ref. | Denomination |
|--------------------|--------------------------------|
| Modulo100AT | Modulo® 100 Auriculotherapy GS |

Supplied with Base Modulo 100 (ref. MODULO 100), battery charger (ref. MOD-01), differential electric Auriculotherapy module (ref. MOD-12) and operating instructions.

Modulo® 100 Auriculo Automatic version

Automatic Gold Silber electric **differential detection** allows to identify, to measure and to **qualify automatically the points** in Auriculotherapy that are either rising or dropping with respect to the area around.



| Ref. | Denomination |
|------------------------|-------------------------------------|
| Modulo100ATAuto | Modulo® 100 Auriculotherapy GS Auto |

Supplied with Base Modulo 100 (ref. MODULO 100), battery charger (ref. MOD-01), auto GS differential electric Auriculotherapy module (ref. MOD-13) and operating instructions.

Our advise to complete these Modulo® 100 Auriculo or Modulo® 100 Auriculo Automatic version is **to choose optional:**

- ⇒ **External grounding handle** which allows more freedom for the detection and the electrical stimulation and offers a greater comfort depending on your environment
- ⇒ **15 W Pulsed IR LASER Module:** for infrared LASER stimulation, painless alternative to needles or to electrical stimulation.
- ⇒ **LASER safety goggles,** indispensable to use laser module.
- ⇒ **Acupuncture direct electric module:** accuracy of points and stimulation in tonification, dispersion or harmonisation for acupuncture.

To equip and enrich your Modulo® 100 Auriculo or your Modulo® 100 Auriculo GS Auto, other options and interchangeable accessories are available in the **Modulo® range. See the spare parts of Modulo® range page 39.**

Instruments for Auriculomedicine

Detectors:

Essential tools for Auriculomedicine.

Used to localize pathological points with the VAS.

4 models available:

Black – White, North – South, Gold – Silver, Gold – Aluminium:



| Ref. | Denomination |
|---------------|---------------------------|
| DETNB | Black – White detector |
| DETNS | North – South detector |
| DETOA | Gold – Silver detector |
| DETOAL | Gold – Aluminium detector |

Set of 1, supplied under plastic bag
Size: 128 x 70 mm Weight: 16 g

Dim Bright Lamp (Heine Lamp or DBL):

Also called Heine lamps, these lamps emit an adjustable beam of white light designed to provoke a Nogier pulse reaction.

They are especially used to test transfers in Auriculomedicine.

The DBL can be used with a short or long nozzle. The version short nozzle is recommended auriculo therapy



| Ref. | Denomination |
|--------------|------------------|
| LEPCC | Short nozzle DBL |
| LEPCL | Long nozzle DBL |

Size: 195 x 35 x 35 mm Weight: 274 g

| Ref. | Spare parts Denomination |
|----------------|-------------------------------------|
| LEP-01 | DBL handle |
| LEP-02 | Short nozzle + bulb (077) |
| LEP-03 | Long nozzle + bulb (037) |
| LEP-021 | Bulb (077) (short nozzle) |
| LEP-022 | Bulb (035) (former short nozzle) |
| LEP-031 | Bulb DBL (037) (former long nozzle) |
| LEP-04 | Leather case |

The Heine lamps are supplied in a leather case with a spare bulb.

Instruments for Auriculomedicine

Polartrons North, South, Unidirectional North, Unidirectional South:

Used for treatment, they permit to do a magnetic stimulation on the body or ear points to be treated.

Very simple to use, they are particularly efficient in the treatment of different pathologies (haematomas ...).

4 models available:



Set of 1, supplied under plastic bag

| Ref. | Denomination |
|--------------|----------------------------------|
| POLN | Polartron North Multidirectional |
| POLS | Polartron South Multidirectional |
| POLUN | Polartron North Unidirectional |
| POLUS | Polartron South Unidirectional |

Size: 245x85x25mm Weight: 340 g

Large magnet:

Set of 1, supplied under plastic bag

| Ref. | Denomination |
|--------------|--------------|
| POLGA | Large magnet |

Size: 103x20x20mm Weight: 140 g

VAS clip with ou without cable:



Set of 1, supplied under plastic bag

| Ref. | Denomination |
|---------------|-----------------------|
| MAT-02 | VAS clip with wire |
| MAT-03 | VAS clip without wire |

Size: 190x50x25mm Weight: 102 g

Notes

| |
|--|
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |

Luminous detectors, stimulators for Auriculomedicine

DB 165, detector with beating:

The DB 165 is a luminous detector for Auriculomedicine developed from works by Dr. Paul NOGIER.

It has 1 high frequency, 1 low frequency and 1 so-called beat frequency emitted when the two keys are pressed simultaneously.

The DB 165 is equipped with 2 removable nozzles: pulsed red light emission and pulsed white light emission.

The red light allows the laterality of patients to be tested and pathological areas to be investigated. This use is discussed in the "**Introduction pratique à l'Auriculomédecine-La Photoperception cutanée**" ("Practical introduction to Auriculomedicine - Skin photoperception") book by Dr. Raphaël NOGIER.

Use of the white light is less well known. It appears that it allows finding more points, especially psychological.

If all three tests of the frequencies on a point cause a reaction to the pulse (VAS), it will be interesting to treat.



| Ref. | Denomination |
|--------------|--|
| DB165 | DB 165 (complete kit, red + white light) |

| Ref. | Spare parts Denomination |
|-----------------|---------------------------|
| DB165-01 | Red light nozzle DB 165 |
| DB165-02 | White light nozzle DB 165 |

| Technical characteristics | |
|--|---|
| Red light nozzle: Wavelength: 630 nm Maximum power: 0,94 mW | White light nozzle: Wavelength: 450 – 750 nm Maximum power: 2,6 mW |
| Pulse frequencies: Black key: high frequency FH = 8,74 Hz Red key: low frequency FB = 3,79 Hz Simultaneous pressing on the 2 keys: beating frequency | Power supply: 1 battery or accu 9V Consumption: 33 mA Size: 125 x 60 x 25 Weight: 141,80 g |
| Supplied in a leather case with: - 2 nozzles: red and white light - 1 battery and 1 user guide | |

Luminous detectors, stimulators for Auriculomedicine

GIR 30 +:

Essential, reliable, easy to handle and perfectly suited to the practice of Auriculomedicine with VAS (from now on it replace Girlase).

GIR 30 + is an emission device of **red and infrared pulsed light** (**white** light in option) designed according **work of Dr Paul NOGIER**.

With its detection key, it allows to search the pathological points on the ear using the VAS. Moving to treatment is achieved simply by pressing the T key. The frequency selected will flash, allowing confirmation that the instrument is in treatment mode.

Local treatment on small skin surface areas is also possible.

It offers the ability of being able to **vary the frequency – 5%, +/- 10%, 20% ou 30% compared to the baseline frequency**. These changes allow the diagnosis and treatment of diseases, which are well known, to auricular doctors. Laterality frequency is given at – 5%.

Thanks to its interchangeable nozzles, the GIR 30+ offers different types of irradiation, which are used depending on the localisations of areas or points to be stimulated. The red is effective in stimulating superficial points on the skin covering and the infra-red which is better absorbed by the tissues is designed for points located deeper.

Pocket and compact GIR 30 + is equipped with an independent DIR test unit to control infrared emission.

Option: the use of the white light allows finding more points, especially psychological.



| Ref. | Denomination |
|-----------------|------------------------------------|
| GIR30P | GIR 30 + (complete kit, red + IR) |
| GIR30-06 | Option white light nozzle GIR 30 + |

| Ref. | Spare parts Denomination |
|-----------------|--------------------------|
| GIR30-01 | Infrared light nozzle |
| GIR30-02 | Red light nozzle |
| GIR30-03 | DIR, Infrared detector |
| GIR30-04 | Leather case |

| Technical characteristics | |
|---|---|
| Red light nozzle: Wavelength: 630 nm Maximum power: 0,94 mW | Infra-red light nozzle: Wavelength: 880 nm Maximum power: 0,88 mW |
| White light nozzle: Wavelength: 450 – 750 nm Maximum power: 2,6 mW | Power supply: 1 battery or accu 9V Consumption: 25 à 35 mA Size: 170 x 78 x 28 mm Weight: 267,70 g |
| Supplied in a carrying case with: <ul style="list-style-type: none"> - 1 infra-red and 1 red nozzle - 1 test unit "DIR" (Infra-Red Detection) - 1 battery and 1 user manual | |

Luminous detectors, stimulators for Auriculomedicine

Modulo[®] 200 Auriculomedicine:

Essential for the regular practise of Auriculomedicine.

Detection by taking the **Dr Paul NOGIER Pulse** (RAC or VAS) and luminous **stimulation** or Laser **stimulation** of the body skin zones **according the 7 frequencies of Dr Paul NOGIER (A, B, C, D, E, F et G).**

Possible modulation of frequencies by variations of ± 5 , ± 10 , ± 20 , $\pm 30\%$ et -35% .

Advantages: every **light and Laser Module[®]** has the property of low emission for the diagnosis with RAC and normal emission for the treatment.

Red light module: classically used in Auriculomedicine for the superficial tissues.

White light module: larger spectrum mobilizes more receptors.

Infra Red Micro Laser module: to stimulate reactive points and zones by the Laser IR coherent light at 905 nm at very low power, invisible laser emission, equivalent to the GirLase 's one.

Classification: class 3R according IEC 60825-1.

Thanks to the ergonomic handpiece, which allows scrolling the Nogier frequencies and their variations, you will evaluate the various physiological disorders.

A sound signal is emitted every 30 seconds to inform the time of stimulation. Stimulation with light, or Laser **according frequencies and variations of frequencies of Dr Paul NOGIER frequencies.**



- ❖ Diagnosis, treatment (individual NOGIER Frequencies)
- ❖ Universal Frequency P.NOGIER (1,14Hz)
- ❖ Modulation of every frequency Dr Paul NOGIER from 0 to + 30% and - 30%
- ❖ Choice between red, white light or Infrared Micro Laser
- ❖ Programming free associations NOGIER frequencies

| Ref. | Denomination |
|--------------------|--|
| Modulo200AM | Modulo [®] 200 Auriculomedicine |

Supplied with Base Modulo[®] 200 (ref. MODULO 200), battery charger (ref. MOD-01), Micro Laser Infrared module (ref. MOD-21), red light Module (ref. MOD-23), white light Module (ref. MOD-24) and manual.

Our advise to complete Modulo[®] 200 Auriculomedicine is to **add optional:**

- ⇒ **Infrared light module**, which penetrates deeper than the other light modules the skin tissues.
- ⇒ **Laser: Pulsed Laser IR 15 W module** which penetrates till 1cm deep the skin and the **safety laser goggles**, indispensable to use laser module.
- ⇒ « **Micro currents** »: **direct electric module or Acutrode electric module** to stimulate reactive points and zones by micro-currents thanks to its electrode and **External grounding mass**, indispensable to use electric module.

To equip and enrich your Modulo[®] 200 Auriculomedicine, other options and interchangeable accessories are available in the **Modulo[®] range**. **See the spare parts of Modulo[®] range page 39.**

Filters Auriculomedicine, Chromotherapy

Nogier colour filters for colour therapy and Auriculomedicine:

Close, away or placed on the skin of the body or the ear, filters are intended to be used to diagnose a pathological state by measurement of the VAS or NOGIER pulse.

Materials used in the manufacture of the rings were specially designed to be neutral to the pulse reaction. It exists many colour filters, the frequential colour programme, metal filters, and magnetic or polarised filters.



| Ref. | Denomination |
|------------------|-------------------------------|
| FIL1C | Filter 1 colour |
| FIL2S | Filter 2 colours (superposed) |
| FIL3S | Filter 3 colours (superposed) |
| FIL3E | Filter 3 colours (star) |
| FIL7E-PCF | Filter 7 colours (star) |

Size: 70 x 46 x 3 mm Weight: 12 g

2 colours filters allow testing two effects simultaneously (2 superposed colour gelatines in the same filter)

3 colours filters offer two possibilities: superposed colours or star-shaped. It is up to you to choose the type of filters, which are suitable for you.

Filters 1 colour (choose your colour, ref : FIL1C - n°):

| Colour | N° | Colour | N° | Colour | N° |
|--------|----------------|---------|----------------------------|-----------|---|
| Yellow | 3 frequency C | Red | 23A | Blue | 44 frequency E and phase 3 |
| Yellow | 4 | Red | 24 frequency D | Blue | 44A |
| Yellow | 6 | Red | 25 frequency B and phase 1 | Blue | 45A  |
| Yellow | 8 | Red | | Blue | |
| Yellow | 8N5 (Relais) | Red | | Blue | |
| Yellow | 9 | Red | 26* | Blue | 46* |
| Yellow | 11 | Red | 29 | Blue | 47 |
| Yellow | 12 | Magenta | 30 frequency G | Blue | 47B |
| Yellow | 13 | Magenta | 31 | Blue | 48 |
| Yellow | 15 | Magenta | 32 | Blue | 48A |
| Yellow | 16 | Magenta | 33 | Blue | 50 |
| Yellow | 72B | Red | 70 | Blue | 75 |
| Yellow | 73 | Red | 92 | Blue | 78 |
| Yellow | 90 | Green | 52 | Blue | 80B |
| | | Green | 53 | Blue | 80C |
| | | Green | 54* | Blue | 82 |
| Orange | 21 frequency A | Green | 55 | Blue | 82C |
| Orange | 22 | Green | 57 | Blue | 94 |
| Orange | 85 | Green | 58 phase 2 | Blue | 98 frequency F |
| Orange | 106 | Green | 59 | Turquoise | 38 |
| | | Green | 60* | Cyan | 38A |
| | | Green | 61 | | |
| Violet | 34 | Green | 64 | Brown | 87 |
| | | Green | 65 | | |
| Violet | 34A | Green | 65A | Brown | 87C |
| Violet | 35 | Green | 74 | Brown | 93 |
| Violet | 97 | Green | 99 | Black | 89B |

(*) limited to our stocks

Filters Auriculomedicine, Chromotherapy

Filters 2 colours, 3 colours, 7 colours (choose your colours, ref : FIL2S - n°, ref : FIL3S - n°, ref : FIL3E - n° or ref : FIL7E - n°):

| Filters (2 colours superposed) |
|--------------------------------|
| 12 / 44A |
| 15 / 12 |
| 21 / 24 |
| 23A / 44 |
| 30 / 16 |
| 30 / 32 |
| 33 / 90 |
| 64 / 44 |
| 64 / 44A |
| 98 / 47B |

| Filters (3 colours superposed) |
|---------------------------------|
| 25 / 58 / 24 |
| 21 / 25 / 44 |
| 21/ 65A / 58 |
| 3 / 12 / 15 |
| 25 / 44 / 58 (3 phases) |
| 30 / 58 / 99 (sleep) |
| Filters star-shaped (3 colours) |
| 21/65A/24 |
| 24 / 44 / 58 |
| 24 / 58 / 65A |
| 25 / 44 / 58 (3 phases) |
| Filters star-shaped (7 colours) |
| 21 / 25 / 3 / 24 / 44 / 98 / 30 |

We can make special filters just on request: superposed colours, half moon, star shaped.

Frequential Colour Programme Paul Nogier (PCFPN):

The Frequential Colour Programme Paul Nogier, **PCFPN**, represents fundamental kit for use of filters in Auriculomedicine.

It includes 7 colours filters representing the 7 frequencies NOGIER 21 (A) - 25 (B) - 3 (C) - 24 (D) - 44 (E) - 98 (F) - 30 (G) plus one filter 7 star-shaped colours including all these colours in the same filter.

| Ref. | Denomination |
|--------------|---------------------------------|
| PCFPN | Frequential Colour Programme PN |

L Filters (choose your colour, ref : FILL1C - n°):

The L colour filter type offers a useful alternative to test the equivalencies of products or substances. If you are interested in these filters, we can also have other products.

| Ref. | Denomination |
|---------------|-------------------|
| FILL1C | Filter L 1 colour |

Filters L1 colour

| N° | N° | N° | N° |
|-------------------------|------------------------------|----------------------------|-----------------------------|
| L 116 (ectoderm) | L 328 (acetylcholine) | L 196 (=45) (blue) | L 124 (= 56) (green) |
| L 410 (endoderm) | L 148 (noradrenalin) | L 071 (=49) (blue) | |
| L 248 (mesoderm) | | L 198 (=49B) (blue) | |

Auriculomedicine filters, Chromotherapy

Specific filters:

Used in Auriculomedicine **to test patient reactivity, energy levels**, blockade of junctions or **to calibrate the pulse**.

Different filters are available to you to help your diagnosis:

- ❖ Gold filter
- ❖ Silver filter
- ❖ Gold/silver filter
- ❖ Black/white filter
- ❖ Polarised Black/white filter
- ❖ Magnetic filter
- ❖ Polarised Magnetic filter
- ❖ Simple Polaroid filter
- ❖ Crossed Polaroid filter
- ❖ Clay filter
- ❖ Space-Time filter ...



| Ref. | Denomination |
|---------------------|---------------------------------------|
| FILOR | Gold filter |
| FILAR | Silver filter |
| FILOA | Gold/silver filter |
| FILNB | Black/white filter |
| FILNBPOL | Black/white Polarised filter |
| FILMAG | Magnetic South filter |
| FILMAGPOL | Magnetic Polarised filter |
| FILMAGN/S | North-South Magnetic Filter |
| FILMAGN/SPOL | North-South Magnetic Polarised Filter |
| FILPOLS | Simple Polaroid filter |
| FILPOLC | Crossed Polaroid filter |
| FILARG | Clay filter |
| FILET04 | Space-Time filter (set of 4) |

Polaroid plates:

Polaroid plates of different size are available:

6 x 12 cm / 10 x 15 cm / 20 x 30 cm

| Ref. | Denomination |
|----------------|----------------------------|
| POL612 | Plaque Polaroid 6 x 12 cm |
| POL1015 | Plaque Polaroid 10 x 15 cm |
| POL2030 | Plaque Polaroid 20 x 30 cm |

Size: 10 x 15 mm

Weight: 18 g

Auriculomedicine filters, Chromotherapy

Test-Rings for filters (anotest):

Male and female transparent discs, absorbent neutral rubber discs with respect to the NOGIER pulse.

Material tested by Auriculomedicine practitioners who have been using these techniques for many years.

All you need to do your own test-rings...(foods or drugs allergies, etc ...)



| Ref. | Denomination |
|----------------|---------------------------|
| ATEST10 | 10 test-rings for filters |

Supplied in set of 10 test-rings.

Possible to purchase separately spare rings and discs (ATEST-01, ATEST-02)

| Ref. | Spare parts Denomination |
|-----------------|---|
| ATEST-01 | Set of 10 rubber rings |
| ATEST-02 | Set of 20 discs (10 male discs plus 10 female discs) |

Storage box to store 20 filters or anotest test-rings (personal or colours):

It allows to organize your various filters, up to 20 filters.



| Ref. | Denomination |
|-----------------|-------------------------------|
| ATEST-03 | Empty storage box for filters |

Notes ...

| |
|--|
| |
| |
| |
| |
| |
| |

Electronic moxibustion device, Acupuncture & Ear Acupuncture



Premio® 10 moxa

Premio® 10 moxa is a device very innovating, patented, which was especially studied to **ensure the same effectiveness as the traditional moxibustion**, whose implementation is **very easy, very fast**, without danger and corresponds to the current standards of the practice in medical cabinet.

The emitter designed for the Premio 10 moxa has **spectral emission characteristics itself replicating those of burning Artemisia**:

The **same spectral band**, and the **same distribution between short, medium and long infrared**, ensures an **mobilization** of the physiological receptors, whatever the **type** and **depth**, equivalent to moxa rolls.

It really is a **"Moxa-Like"** that finds the **practice and effectiveness of the moxibustion, with the same spectrum, the same applications, same results**, without drawbacks of toxic fumes, strong and lingering, ignition and extinction difficult, burning ash or cool the moxa, difficulty to use in a medical environment.



The distribution of wavelengths across such a wide range is a basic characteristic of Artemisia, and largely explains its effectiveness, due to its being appropriate for the physiological, thermo- and photo-receptors of the organism.

Electronic moxibustion device, Acupuncture & Ear Acupuncture

Premio® 10 moxa

Essential for regular practice of moxibustion with a quality heat same as Artemisia.
Lasting effect and sustainable investment for a sustainable moxibustion.



| Ref. | Denomination |
|-----------------|----------------|
| PREMIO10 | Premio 10 moxa |

Size : 180 mm
Weight : 230 g

| Ref. | Spare part denomination |
|--------------------|-----------------------------|
| PREMIO10-01 | Premio 10 moxa concentrator |

Premio® 10 moxa applications :



In Acupuncture
Familiar gestures
Known results



The "heat needle" effect :
A tingling penetrating heat
targeting on a specific point



In Auriculartherapy
A new precise stimulation technique in seconds

In local heat treatment
A precise, safe, effective application.



Electric detector and stimulator, Laser Acupuncture



Modulo® 100 Acupuncture :

Perfectly adaptable for micro-current stimulation and Laser Acupuncture, localization, detection, precision and electric qualification of acupuncture points.

Representation through the diodes, impedance falls, ensure timeless and accuracy of detection.

Timeless and accuracy of treatment options.

Stimulation points by electrical or LASER is performed in three different modes :

Harmonisation (↔ balance by sending small electrical impulses frequencies), **Tonifie** (LF) and **Disperse** (HF).

Two possible micro-current energy levels in power, normal (ENERGY +) or low (ENERGY -).

Treatment with the 15 W Pulsed IR **LASER Module**® ensures optimum penetration (infrared 905 nm) and the contribution of energy required by a point in time is from 30 to 60s. Painless LASER stimulation is particularly suitable for children and apprehensive patients.

Invisible LASER emission. Classification : class 3B according to IEC 60825-1.

A beep sounds every 30 seconds to indicate the stimulation time.

- ❖ **Dispersion**
- ❖ **Tonification**
- ❖ **Harmonisation**
- ❖ 7 fundamental NOGIER frequencies
- ❖ Aseptic
- ❖ Painless in LASER treatment
- ❖ **Non invasive**
- ❖ Great alternative to needles
- ❖ **Short time stimulation (about 30 to 60 seconds)**
- ❖ Stimulate photoreceptors respecting tissues



| Ref. | Denomination |
|---------------------|-------------------------|
| Modulo100Acu | Modulo® 100 Acupuncture |

Delivered with Basis Modulo® 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), Direct Electric Acupuncture Module (ref. MOD-11), 15 W Pulsed IR LASER Module (ref. MOD-32), LASER Goggles (ref. MOD-52) and manual.

Our advice to complete the Modulo® 100 Acupuncture is an option **to add the 40 W Pulsed IR LASER Module** with a mean penetration of the skin surface up to 2 to 3 cm deep and add the **75 W Pulsed IR LASER Module** with a mean penetration of the skin surface up to 4 cm deep.

In « micro-currents » : Acuneedles Module, electrical stimulation by needles, help to stimulate the points and hot spots by micro-currents, on 1 to 4 needles simultaneously.

Complete and extend your Modulo® 100 Acupuncture, others options are available and interchangeable accessories in the **range Modulo®**. **View the options parts of the range Modulo® page 39.**

Electric detector & stimulator, Reflexo and LaserAcupuncture

Modulo® 100 Micro-currents Reflexo Acutrodes :

Localization, detection, precision and electric qualification of acupuncture points and micro-currents stimulation points (body, face ...), without needles, with electrodes, **from 1 to 4 points.**




| Ref. | Denomination |
|----------------------------------|-------------------------------|
| Modulo100ReflexoAcutrodes | Modulo® 100 Reflexo Acutrodes |

Delivered with Basis Modulo® 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), Direct Electric Acupuncture Module (ref. MOD-11), Acutrodes Module (ref. MOD-41), External grounding electrode (ref. MOD-51) and manual.

Modulo® 100 Micro-currents Reflexo Needles :

NEW

Localization, detection, precision and electric qualification of acupuncture points and micro-currents stimulation points (body, face ...), with connected needles **from 1 to 4.** The connectors provide a perfect and constant contact with the needles while ensuring total isolation of the skin.



| Ref. | Denomination |
|-----------------------------------|-----------------------------|
| Modulo100ReflexoAcuneedles | Modulo® 100 Reflexo Needles |

Delivered with Basis Modulo® 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), Direct Electric Acupuncture Module (ref. MOD-11), Acuneedles Module (ref. MOD-42), External grounding electrode (ref. MOD-51) and manual.

Modulo® 100 Acu 15 W Pulsed IR LASER:



Average penetration of the skin surface to about 1 cm deep with the 15 W Pulsed IR LASER Module.

Invisible LASER emission. Classification : class 3B according to IEC 60825-1.



| Ref. | Denomination |
|----------------------------|--------------------------------------|
| Modulo100AcuLaser15 | Modulo® 100 Acu 15 W Pulsed IR LASER |

Delivered with Basis Modulo® 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), 15 W Pulsed IR LASER Module (ref. MOD-32), LASER goggles (ref. Mod-52) and manual.

LASER Acupuncture stimulation

Modulo[®] 100 Acu 40 W Pulsed IR LASER:



Average penetration of the skin surface to about 2 cm deep with the 40 W Pulsed IR LASER Module.

Invisible LASER emission. Classification : class 3B according to IEC 60825-1.



| Ref. | Denomination |
|----------------------------|---|
| Modulo100AcuLaser40 | Modulo 100 [®] Acu Laser IR Pulsed 40W |

Delivered with Basis Modulo[®] 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), 40 W Pulsed IR LASER Module (ref. MOD-33), LASER goggles (ref. Mod-52) and manual.

Modulo[®] 100 Acu 75 W Pulsed IR LASER:



Average penetration of the skin surface to about 4 cm deep with the 75 W Pulsed IR LASER Module.

Invisible LASER emission. Classification : class 3B according to IEC 60825-1.



| Ref. | Denomination |
|----------------------------|--|
| Modulo100AcuLaser75 | Modulo 100 [®] Acu 75 W Pulsed IR LASER |

Delivered with Basis Modulo[®] 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), 75 W Pulsed IR LASER Module (ref. MOD-35), LASER goggles (ref. Mod-52) and manual.

Modulo[®] 100 Acu 40 W Superpulsed IR LASER:



Performance in Joule/seconds ! The Super Pulsed LASER provides energy up to 1 Joule/30 seconds.

Average penetration of the skin surface to about 2 cm deep with the 40 W Super Pulsed IR LASER Module. Invisible LASER emission. Classification : class 3B according to IEC 60825-1.



| Ref. | Denomination |
|-------------------------------|---|
| Modulo100AcuLaser40 SP | Modulo [®] 100 Acu 40 W Superpulsed IR LASER |

Delivered with Basis Modulo[®] 100 (ref. MODULO 100), Stand and battery charger (ref. MOD-01), 40 W SuperPulsed IR LASER Module (ref. MOD-37), LASER goggles (ref. Mod-52) and manual.

Consider your options of choice, many interchangeable modules : electric, LASER or light can be added and also accessories

Double electrical stimulator, electro acupuncture

NEW

AGISTIM ^{DUO}

Essential electrical stimulation : excellent reliability and performance.

2 x 2 channel stimulator for electro-acupuncture can stimulate up to 8 needles or 8 electrodes, 16 needles optional.

"Duo" means that the l'AGISTIM ^{DUO} is equivalent to 2 stimulators in a single device. The four outputs run 2 by 2, providing opportunities for different settings.

It also allows treatment of electrolipolysis.

Opportunity to choose **frequencies and emission modes** different on the two groups of the two channels with adjustable intensity on each channel.

The frequency ranges are to be selected **between 1 and 99 Hz**. They are accurate to about 0,1 or to about 1 hertz closely the chosen scale.

It is also possible to directly select the **7 NOGIER frequencies**.

8 transmission modes or different processing : continuous stimulation with rest periods of varying length, one frequencies scan function, the AGISTIM ^{DUO} can adapt the type of stimulation for each patient and the condition being treated.

The electrical stimulation enhance the action of needles. In some cases it can replace many manipulations and even replace needles for children and people fearful.

The AGISTIM ^{DUO} is used in :

- ❖ **Electroacupuncture** : treatment faster than the needles used alone
- ❖ **Neurostimulation** : stimulation of nerve fibers with the electrodes
- ❖ **Analgesia** : release of endorphins based on selected frequencies
- ❖ **Electrolipolysis** : effective reduction of localized lipodystrophy with appropriate needles (ASE)



Reusable stimulation electrodes

Of round shape Ø 30 mm
Or of square-shaped (50 mm)
Are also available (page 33)

L'AGISTIM ^{DUO} offers maximum security. It runs on sector power by security adaptable to all sectors.

An **integrated security system** independent stops the unit if the frequency or intensity of currents exceed the maximum values. In case of error or omission, the device reports explicitly the problem still ensuring proper operation.

Processing parameters are always visible on the display screen. The processing time is adjustable from 1 to 60 min.

Double electrical stimulator, electro acupuncture

AGISTIM ^{DUO}

NEW


| Ref. | Denomination |
|---------------|---------------------------------------|
| AGIDUO | AGISTIM ^{DUO} (complete kit) |

Delivered in a carrying case with :

- 1 test box
- 4 needles connector
- 2 electrodes connector (white - grey)
- 4 stimulation electrodes Ø 30 mm
- 1 manual

| Specifications | |
|--|---|
| Current | Alternating, pulsed with a mean value of zero |
| Impulse time | Asymmetrical rectangular |
| Impulse time | 0,4 ms |
| Amplitude of pulse | Adjusted independently on each channel between 0 and 12 mA |
| Maximum RMS intensity per channel | 3,4 mA |
| Impulse frequencies | 1 to 99 Hz or NOGIER frequencies : A, B, C, D, E, F, G, U |
| Stimulation mode | continuous - short rest time, long rest time - frequency scanning |
| Power supply | Sector supply |
| Consumption | 6,0 VA |
| Measurements | 60 x 230 x 300 mm |
| Weight | 1416,30 g |

To equip and enrich your **AGISTIM ^{DUO}** options and accessories :

Styl Modul for Acupuncture, Auriculartherapy :

Electrical frequencies stimulation very accurate on a point raised without needles or electrodes sans pose d'aiguille ou d'électrodes. Very useful especially on ear points.

Possibility of dealing with NOGIER frequencies and any frequency between 1 and 99 Hz. This device connects directly to the outputs of AGISTIM ^{DUO}.



| Ref. | Denomination |
|------------------|-------------------------------------|
| AGIDUO-09 | Styl Module (point stimulation pen) |

Spare: Test box, power supply :

| Ref. | Spare parts denomination |
|-------------------|--|
| AGIDUO-05 | AGISTIM DUO test box |
| AGIDUO-063 | Main power supply 110-240V AGISTIM Duo |

Options and accessories for electro-acupuncture

To equip and enrich your **AGISTIM DUO** options and accessories :

4 Reusable stimulation electrodes :

Ø 30 mm round or square (50 mm)

| Ref. | Spare parts denomination |
|-------------------|---|
| AGIDUO-041 | 4 reusable stimulation electrodes (Ø 30 mm) |
| AGIDUO-042 | 4 stimulation electrodes (50 x 50 mm) |

Kit 4 needles module :

Double the number of needles connected to an AGISTIM DUO easy attachment system : attachment of the needle with one hand by pulling on the wire.

Delivered with 4 colors of cables to identify the needles during stimulation:

| Ref. | denomination |
|---------------------|---|
| AGIDUO-03 | 4 needles kit (the four 4 needle modules) |
| AGIDUO-03-S1 | Kit 4 needles double length (the four 4 needle modules) |

With spare parts, choose the color of the cable : **grey, pink, purple or green**:

| Ref. | denomination | Ref. | denomination |
|-------------------|-------------------------------|----------------------|---|
| AGIDUO-031 | 4 Needle Module grey | AGIDUO-031-S1 | 4 Needle Module grey , double length |
| AGIDUO-032 | 4 Needle Module pink | AGIDUO-032-S1 | 4 Needle Module pink , double length |
| AGIDUO-033 | 4 Needle Module purple | AGIDUO-033-S1 | 4 Needle Module purple , double length |
| AGIDUO-034 | 4 Needle Module green | AGIDUO-034-S1 | 4 Needle Module green , double length |

Electrolipolysis sterile needles for single use: ASE

Disposables used to practice electrolipolysis (a full stainless steel handle, smooth).

Sterilization with ethylene oxide, equipped with a **diamond shape tip**, they become established easily in fat mass, no bruising and no pain for the patient. Related to AGISTIM DUO, they can treat cellulite, reduce the orange peel effect and reduce localized lipodystrophy.

Delivered packaged in packs of 8 and needles available in boxes of 120 needles
2 lengths: 100 and 125 mm / Diameter: 0.30 mm



| Réf. | Libellé |
|------------------|---------------------|
| ASE100120 | ASE PD 100 (15 x 8) |
| ASE125120 | ASE PD 125 (15 x 8) |

Needles connectors:

Ensure **proper and constant contact** with needles while ensuring **total isolation of the skin**.



| Réf. | Libellé pièces détachées |
|-------------------|--|
| AGIDUO-011 | Needles connector (1 white - 1 grey) |
| AGIDUO-012 | Needles connector (1 white - 1 pink) |
| AGIDUO-013 | Needles connector (1 white - 1 purple) |
| AGIDUO-014 | Needles connector (1 white - 1 green) |
| AGIDUO-02 | Connecteur électrodes (1 blanc - 1 gris) |

LASER Therapy, local treatment

Modulo® 200:

NEW



The major, essential LASER for local treatment

Each choice of LASER Module has the property :

- ❖ **15 W Pulsed IR LASER Module**, with average penetration from the skin surface to a depth of about 1 cm
- ❖ **40 W Pulsed IR LASER Module**, with average penetration from the skin surface to a depth of about 2 to 3 cm
- ❖ **75 W Pulsed IR LASER Module** with average penetration from the skin surface to a depth of about 4 to 5 cm
- ❖ **40 W Super Pulsed IR LASER Module**, with average penetration from the skin surface to a depth of about 2 - 3 cm with an energy of 1 Joule/30 seconds

Associations "Frequency / Percentage" can be programmed for special effects.

Analgesic, Regenerating and Relaxing programs are pre recorded.

LASER stimulation frequencies according to Dr. Paul NOGIER.

Treatment with the 7 individuals NOGIER frequencies

NOGIER Universel frequency (1,14Hz)

Modulating each frequency of Dr Paul NOGIER from 0 to + 30% and from 0 to - 30%

Programming free of associations Nogier frequencies

Modulo® 200 Pulsed IR LASER 15 W:



Average penetration from the skin surface up to about 1 cm deep with 15 W Pulsed IR LASER Module. Invisible LASER emission. Classification : class 3B according IEC 60825-1.



Ref.

Denomination

Modulo200Laser15

Modulo® 200 Pulsed IR LASER 15 W

Delivered with :

- 1 Basis Modulo® 200 (ref. MODULO 200)
- 1 Stand and battery charger (ref . MOD-01)
- 1 15W Pulsed IR LASER Module (ref. MOD-32)
- 1 LASER safety goggles (ref. MOD-52)
- 1 manual

LASER Therapy, local treatment

Modulo® 200 Pulsed IR LASER 40 W:



Average penetration from the skin surface up to about 2 to 3 cm deep with 40 W Pulsed IR LASER Module. Invisible LASER emission.
Classification : class 3B according IEC 60825-1.



| Ref. | Denomination |
|-------------------------|----------------------------------|
| Modulo200Laser40 | Modulo® 200 Pulsed IR LASER 40 W |

Delivered with :

- 1 Basis Modulo® 200 (ref. MODULO 200)
- 1 Stand and battery charger (ref . MOD-01)
- 1 40W Pulsed IR LASER Module (ref. MOD-33)
- 1 LASER safety goggles (ref. MOD-52)
- 1 manual

Modulo® 200 Pulsed IR LASER 75 W:



Average penetration from the skin surface up to about 4 to 5 cm deep with 75 W Pulsed IR LASER Module. Invisible LASER emission.
Classification : class 3B according IEC 60825-1.



| Ref. | Denomination |
|-------------------------|----------------------------------|
| Modulo200Laser75 | Modulo 200® Pulsed IR LASER 75 W |

Delivered with :

- 1 Basis Modulo® 200 (ref. MODULO 200)
- 1 Stand and battery charger (ref . MOD-01)
- 1 75W Pulsed IR LASER Module (ref. MOD-35),
- 1 LASER safety goggles (ref. MOD-52)
- 1 manual

Modulo® 200 IR LASER 40 W Superpulsed:



Performance in Joule/seconds ! The Super Pulsed LASER provides energy up to 1 Joule/30 seconds.
Average penetration from the skin surface up to about 2 to 3 cm deep with 40 W Pulsed IR LASER Module. Invisible LASER emission.
Classification : class 3B according IEC 60825-1.



| Ref. | Denomination |
|---------------------------|---------------------------------------|
| Modulo200Laser40SP | Modulo® 200 Superpulsed IR LASER 40 W |

Delivered with :

- 1 Basis Modulo® 200 (ref. MODULO 200)
- 1 Stand and battery charger (ref . MOD-01)
- 1 40W Super Pulsed IR LASER Module (ref. MOD-37)
- 1 LASER safety goggles (ref. MOD-52)
- 1 manual

Physiotherapy LASER, Dental LASER

Modulo® 200 Physio IR LASER 40 W Superpulsed:



LASER Stimulation according to Dr Paul NOGIER frequencies, is particularly appropriated for several **physiotherapy applications with 40 W Superpulsed IR LASER module**.
This superpulsed LASER brings an energy of 1 Joule every 30 seconds.



Three basic programmes included: "**REGENERATING**", "**ANALGESIC**" and "**MUSCLE RELAXING**" for common indications (tissue alterations, pain, contractures) easily selected

Associations "Frequency / Percent" can be programmed for special effects.

Invisible LASER emission.

Classification: class 3B according to IEC 60825-1.

| Ref. | Denomination |
|--------------------------|--------------------------------|
| Modulo200PhysioSP | Modulo® 200 Physio SP IR LASER |

Supplied with Base Modulo® 200 (ref. MODULO 200), battery charger (ref. MOD-01), 40W Superpulsed IR Laser Module (ref. MOD-37), safety Laser goggles (ref. MOD-52) and operating instructions.

Modulo® 200 Dental LASER 30 W Superpulsed:

NEW



LASER Stimulation according to Dr Paul NOGIER frequencies, specially designed for work **inside the mouth (dental, stomato...)** with **30W Superpulsed IR LASER module with optic fiber (Ø 6mm)** and **any application where the sterilisation of the tip is necessary**.

It brings an energy of 0,7 Joule every 30 seconds. Supplied with 2 removable Optic Fibers for sterilization.

Treatment with the 7 individual Frequencies NOGIER and Universal NOGIER frequency U (1,14Hz), modulation of each frequency Dr Paul NOGIER from 0 to + 30% and from 0 to - 35%, Store up to 5 programs for automatic application of personalised series of specific Dr Paul NOGIER's frequencies.

Invisible LASER emission. Classification: class 3B according to IEC 60825-1.



| Ref. | Denomination |
|-----------------------------|--------------------------------------|
| Modulo200DentalLaser | Modulo 200® Dental Fiber SP IR LASER |

Supplied with Base Modulo® 200 (ref. MODULO 200), battery charger (ref. MOD-01), 30W SuperPulsed IR Laser module with optic fiber (ref. MOD-37F), safety Laser goggles (ref. MOD-52) and operating instructions.

Optical fiber:

Alone Optical Fiber, removable for sterilization.

| Ref. | Denomination |
|-------------------|---------------------------------------|
| MOD-37F-01 | Optic fiber for MOD-37F (fiber Ø 6mm) |

Frequential Therapy and micro-currents

The unavoidable and indispensable for frequential therapy

Stimulation with micro-currents of points and cutaneous zones of the body according to **Universal frequency U** and the **7 individual frequencies of Dr Paul NOGIER** (A, B, C, D, E, F et G), possible **modulation** of each frequency by **variation of ± 5 , ± 10 , ± 20 , $\pm 30\%$ and -35%** . Associations of "Frequency / Percentage" can be programmed for specific effects. Programs for "**REGENERATING**", "**ANALGESIC**" and "**MUSCLE RELAXING**" are already stored.

The application of short electrical pulses of low intensity to stimulate biological receptors of detected points, in a similar manner to the application of needle, but without puncture or pain.

The use of electrical modules requires to connect the external grounding handle and to let the patient hold the half metallic sphere in the hollow of his hand, for an optimal electric contact. A beep is emitted every 30 seconds to tell you the time of stimulation.

Modulo 200[®] Micro-currents Acutrodes:

Each module has its property:

- ❖ **Direct** Electrical module allows to **stimulate a point and a reactive zone**
- ❖ **Acutrodes** module allows stimulating points and reactive zones with 1 to 4 electrodes simultaneously.



| Ref. | Denomination |
|---------------------------|--|
| Modulo200Acutrodes | Modulo [®] Micro-currents Acutrodes |

Supplied with Base Modulo[®] 200 (ref. MODULO 200), battery charger (ref. MOD-01), Direct electrical Acupuncture module (ref. MOD-11), micro-current Acutrodes module (ref. MOD-41), external grounding handle (ref. MOD-51) and operating instruction.

Modulo[®] 200 Micro-currents Acuneedles:

Micro-currents Acuneedles module, allows to stimulate points and reactive zones with 1 to 4 acupuncture needles simultaneously.



| Ref. | Denomination |
|----------------------------|---|
| Modulo200Acuneedles | Modulo [®] Micro-currents Acuneedles |

Supplied with Base Modulo[®] 200 (ref. MODULO 200), battery charger (ref. MOD-01), Direct electrical Acupuncture module (ref. MOD-11), micro-current Acuneedles module (ref. MOD-42), external grounding handle (ref. MOD-51) and operating instruction.

Modularity, evolutivity of Modulo[®] 100 and Modulo[®] 200

You can create and customize your Modulo[®] by selecting modules, accessories, a base and battery charger:

Base Modulo[®] 100

The Base Modulo[®] 100 (main body and hand piece) must be associated with a battery charger. The hand piece receives **differential or direct electric modules for detection of excellent quality** and receives the modules and accessories of your choice for **multiple treatments**.



- ❖ Simplicity and ease of use, ergonomics and efficiency
- ❖ Any type of emission (electric, luminous, LASER...) according to selected module
- ❖ Viewing and sound when detection of the reactive point with electrical modules
- ❖ Checking of grounding electrode handling for an optimal detection
- ❖ Stimulation according to three modes of frequencies: LF (Low Frequency), HF (High Frequency), ↔ (scanning)

| Ref. | Denomination |
|------------------|--------------------------------------|
| MODULO100 | Base Modulo [®] 100 (alone) |

Size: 120 x 120 x 185 mm
Weight: 415 g
Supplied by one in carton box

Base Modulo[®] 200

The Base Modulo[®] 200 (main body and hand piece) must be associated with a battery charger. The hand piece receives **the various modules and accessories** of your choice for **multiple treatments** and allows a possible detection with the RAC or NOGIER Pulse.



- ❖ Simplicity and ease of use, ergonomics and efficiency
- ❖ Any type of emission (electric, luminous, LASER...) according to selected module
- ❖ Detection with RAC or NOGIER Pulse: low level stimuli and scrolling through 7 fundamental NOGIER Frequencies (from A to G) and their variations ± 5 , ± 10 , ± 20 , $\pm 30\%$ including dental -35% and Universal frequency U
- ❖ Treatment: normal level stimuli at fundamental NOGIER frequencies and their variations defined
- ❖ Three basic programs included: "Analgesic" (E and G frequencies), "Regenerating" (A, B and F frequencies), and "Muscle Relaxing" (C, D and G frequencies) for common indications (tissue alterations, pain, contractures) easily selected
- ❖ Store up to 5 programs for automatic application of personalised series of specific frequencies.

| Ref. | Denomination |
|------------------|--------------------------------------|
| MODULO200 | Base Modulo [®] 200 (alone) |

Size: 120 x 120 x 185 mm
Weight: 415 g
Supplied by one in carton box

Modularity, evolutivity of Modulo[®] 100 and Modulo[®] 200

You can create and customize your Modulo[®] by selecting modules, accessories, a base and battery charger:

Electrical Acupuncture, Auriculotherapy Modules:

Adjusted to each detection and stimulation type:

MOD-11 for Acupuncture,

MOD-12 or MOD-13 for Auriculotherapy



| Ref. | Denomination |
|---------------|--|
| MOD-11 | Direct Electric Module Acupuncture |
| MOD-12 | Differential Electric Module Auriculotherapy Gold Silber |
| MOD-13 | Automatic Differential Electric Module Auriculotherapy Gold Silber |

Size: 110 x 20 mm

Weight: 28 g, Supplied by one in carton box

Micro-currents Acutrodes Module:

Stimulates in micro-currents, without needles, **1 to 4 points** in Acupuncture.

External grounding handle (MOD-51) is strongly recommended to stimulate with Acutrode Module and **indispensable for use with Modulo[®] 200**.



| Ref. | Denomination |
|-----------------|---|
| MOD-41 | Acutrodes Module (+ 100 adhesives plasters) |
| MOD-41-1 | 500 adhesives plasters pour Acutrodes |

Wire length: 700 mm

Module size: 90 x 20 mm

Weight: 55 g, Supplied by one in carton box


Micro-currents Acuneedles Module:

NEW

Stimulates in micro-currents, with needles, **1 to 4 points** in Auriculotherapy and/or in Acupuncture.

Connectors provide a perfect and constant contact with the needles and guarantee a total isolation with the skin.

External grounding handle (MOD-51) is strongly recommended to stimulate with Acuneedle Module and **indispensable for use with Modulo[®] 200**.



| Ref. | Denomination |
|---------------|-------------------|
| MOD-42 | Acuneedles Module |

Wire length: 700 mm

Module size: 90 x 20 mm

Weight: 55 g, Supplied by one in carton box

Modularity, evolutivity of Modulo[®] 100 and Modulo[®] 200

You can create and customize your Modulo[®] by selecting modules, accessories, a base and battery charger:

External Grounding Handle:

External grounding handle gives more freedom with Modulo[®] 100. It is **indispensable for electrical stimulation with Modulo[®] 200**. It is highly recommended for electrical stimulation with Module Acutrodes and essential for electrical stimulation with this module Acuneedles.



| Réf. | Denomination |
|---------------|---------------------------|
| MOD-51 | External grounding handle |

Size: 45 x 30 mm

Wire size: 700 mm

Weight: 87 g, Supplied by one in carton box

Pulsed IR LASER Modules:

According to depth penetration required (1, 3 or 5 cm), three powers are adaptable 15, 40 or 75W. Invisible LASER emission. Classification: class 3B according to IEC 60825-1.



| Ref. | Denomination |
|---------------|-----------------------------|
| MOD-32 | 15 W Pulsed IR LASER Module |
| MOD-33 | 40 W Pulsed IR LASER Module |
| MOD-35 | 75 W Pulsed IR LASER Module |

Size: 90 x 20 mm

Weight: 28 g, Supplied by one in carton box

Superpulsed IR LASER Module:

When high power is required: **1 Joule/30 seconds**.

Invisible LASER emission. Classification: class 3B according to IEC 60825-1.



| Ref. | Denomination |
|---------------|---------------------------------|
| MOD-37 | 40W Superpulsed IR LASER Module |

Size: 90 x 20 mm

Weight: 28 g, Supplied by one in carton box

Superpulsed IR LASER Module with optic fiber:

Especial designed for treatment in mouth (dental, stomato...) and

For any application where the sterilization of the tip is necessary.

Invisible LASER emission. Classification: class 3B according to IEC 60825-1.

0,7 Joule/30 seconds. Supplied with 2 removable Optic Fibers for sterilization.

NEW



| Réf. | Denomination |
|-------------------|---|
| MOD-37F | 30W Superpulsed IR LASER module with optic fiber (supplied with 2 Fibres Ø 6mm) |
| MOD-37F-01 | Optic fiber for MOD-37F (fiber Ø 6mm) |

Size: 160 x 20 mm

Weight: 31 g, Supplied by one in carton box

Optical fiber alone, Size: 30 x 7 mm, Weight: 8 g

Modularity, evolutivity of Modulo[®] 100 and Modulo[®] 200

You can create and customize your Modulo[®] by selecting modules, accessories, a base and battery charger:

Pulsed MicroLaser Module :

A genuine LASER, micro power. (equivalence Girlase)

Invisible emission laser. Classification: class 3R according IEC 60825-1.



| Réf. | Denomination |
|---------------|----------------------------|
| MOD-21 | Infrared MicroLaser Module |

Size: 90 x 20 mm

Weight: 28 g

Supplied by one in carton box

LASER safety glasses:

LASER safety glasses are indispensable for use with LASER modules.



| Réf. | Denomination |
|---------------|----------------------|
| MOD-52 | LASER safety glasses |

Size: 600 x 650 x 20 mm

Weight: 44 g

Supplied by one in carton box

Light Modules:

Wavelength choice, depending on application.



| Réf. | Denomination |
|---------------|-----------------------|
| MOD-22 | Infrared light Module |
| MOD-23 | Red light Module |
| MOD-24 | White light Module |

Size: 90 x 20 mm

Weight: 28 g

Supplied by one in carton box

Stand and battery charger:

Essential to recharge the batteries of Modulo[®] 100 and Modulo[®] 200.

It permits also to hold 4 modules.



| Réf. | Denomination |
|---------------|-------------------------------------|
| MOD-01 | Battery charger Modulo [®] |

Size: 240 x 185 x 70 mm

Weight: 423 g

Supplied by one in carton box

Modularity and evolutivity of Modulo® 100 and Modulo® 200

Examples of Modulo®'s personalisation:

| Example of kits Modulo® | | | | |
|---|--|--|---|---|
| | Modulo® Auriculotherapy Electric point finder & Stimulator and LASER | Modulo® Acupuncture Reflexotherapy | Modulo® Auriculomedicine according NOGIER | Modulo® Frequential Therapies Micro currents |
| Réf.: | Modulo100 AuriculoLASER | Modulo100Reflexo Acuneedles | Modulo200 Auriculomedecine | Modulo200 Acutrodes |
| Bases | | | | |
| Base Modulo® 100 | ● | ● | | |
| Base Modulo® 200 | | | ● | ● |
| Accessories | | | | |
| Battery charger | ● | ● | ● | ● |
| External Grounding Electrode | ○ | ○ | | ○ |
| LASER safety glasses | ● | | | |
| Modules | | | | |
| Direct Electric Acupuncture | | ● | | ● |
| Differential Electric Auriculotherapy | ● | | | |
| Auto GS Differential Electric Auriculotherapy | ○ | | | |
| Acutrodes micro-currents | | ○ | | ● |
| Acuneedles micro-currents | | ● | | ○ |
| 15 W Pulsed LASER IR | ● | ○ | | |
| 30 W Pulsed LASER IR | | | | |
| 75 W Pulsed LASER IR | | | | |
| 40 W Super Pulsed LASER IR | | | | |
| 30 W Super Pulsed LASER IR with 2 optic fibres | | | | |
| Micro-Laser IR | | | ● | |
| Infrared Light | | | ○ | |
| Red Light | | | ● | |
| White Light | ○ | | ● | |

● Included in the price

○ Recommended option

To equip and enrich your Modulo® 100 or your Modulo® 200, other options and interchangeable accessories are available. **See the spare parts of Modulo® range page 39.**

Modularity, evolutivity of Modulo[®] 100 and Modulo[®] 200 LASER

Examples of Modulo[®] LASER's personalisation:



| Example of kits Modulo [®] | | | |
|--|---|--|---|
| | Modulo[®] 100 Acupuncture Multi IR LASER | Modulo[®] 200 Physiotherapy IR LASER SP | Modulo[®] 200 Stomato Dental Applications |
| Réf.: | Modulo100 Multi LASER | Modulo200 Physio SP IR LASER | Modulo200 Stomato SP IR LASER |
| Bases | | | |
| Base Modulo [®] 100 | ● | | |
| Base Modulo [®] 200 | | ● | ● |
| Accessories | | | |
| Battery charger | ● | ● | ● |
| External Grounding Electrode | | | |
| LASER safety glasses | ● | ● | ● |
| Modules | | | |
| Direct Electric Acupuncture | | | |
| Differential Electric Auriculotherapy | | | |
| Auto GS Differential Electric Auriculotherapy | | | |
| Acutrodes micro-currents | | | |
| Acuneedles micro-currents | | | |
| 15 W Pulsed LASER IR | ● | ○ | ● |
| 40 W Pulsed LASER IR | ○ | | |
| 75 W Pulsed LASER IR | ● | | |
| 40 W Super Pulsed LASER IR | ● | ● | |
| 30 W Super Pulsed LASER IR with 2 optic fibres | | ○ | ● |
| Micro-Laser IR | | | |
| Infrared Light | | | |
| Red Light | | | |
| White Light | | | |

● Included in the price

○ Recommended option

To equip and enrich your Modulo[®] 100 or your Modulo[®] 200, many options and interchangeable accessories are available. **See the spare parts of Modulo[®] range page 39.**

Instruments for Iridology

Iriscope:

Used in Iridology, it permits to do a very precise diagnosis thanks to its two interchangeable magnifying glass.

Two enlargements available: **5x** or **8x**

Iriscope is delivered in a leather case with a spare bulb.

| Réf. | Denomination |
|---------------|------------------|
| IRIH5X | Iriscope H 5 x * |
| IRIH8X | Iriscope H 8 x * |

(*) Limited to our stocks

| Réf. | Spare parts Denomination |
|----------------|------------------------------|
| IRI-011 | Magnifying glass 5 x * |
| IRI-012 | Magnifying glass 8 x * |
| IRI-01 | Iriscope head + bulb (038) * |
| IRI-013 | Bulb (038) |
| IRI-02 | Leather case |

(*) Limited to our stocks

Instruments for sympathicotherapy (nasotherapy)

Endo-nasal stylets and Tipped Applicator:


They are used in Sympathicotherapy and allow the reflex points in the nasal mucosa to be stimulated gently.

The endo-nasal stylets are available in flexible or semi-rigid versions.

The sterile applicators (cotton tip and wooden handle) are packaged in boxes of 200 units.

| Réf. | Denomination |
|----------------|----------------------------------|
| STYSOU | Flexible endo-nasal stylet x 1 |
| STYSER | Semi-rigid endo-nasal stylet x 1 |
| ECOU200 | Cotton tipped applicator x 200 |

Auriculotherapy, Acupuncture Charts and others



| Réf. | Authors - Titres - Editors |
|----------------|--|
| PLBEIJ01-FR | BEIJING INSTITUTE Official localisation of Acupuncture's points 3 charts 100x80cm + 1 booklet 40 pages |
| *PLBFNS01-AN | B.L. FRANK / N.E.SOLIMAN Auricular Therapy Reflex Chart - English - 60x90 cm * |
| PLBFNS02-AL | B.L. FRANK / N.E. SOLIMAN Reflexzonen für Ohrakupunktur - German - 21x27 cm |
| PLBFNS02-AN | B.L. FRANK / N.E.SOLIMAN Auriculartherapy Reflex Chart - English - 21x27 cm |
| *PLBFNS02-ES | B.L. FRANK / N.E. SOLIMAN Areas refleja para terapia auricular - Spanish - 21x27 cm* |
| PLBFNS02-FR | B.L. FRANK / N.E.SOLIMAN Planche d'Auriculothérapie - French - 21x27 cm |
| *PLBFNS02-PO | B.L. FRANK / N.E. SOLIMAN - Auricular Therapy Reflex Chart - Portuguese - 21x27 cm* |
| PLBFNS03-AN | B.L. FRANK / N.E.SOLIMAN Hand Acupuncture Reflex Chart - English - 28x43 cm |
| PLCERS01-FR | CERS TA Nasosympathicothérapie - Somatotopie interne |
| PLDALI01-FR | D. ALIMI Nomenclature et Cartographie Internationales d'Auriculothérapie (set of 4 charts) Edit.Sedatelec - French - 35x50cm |
| *PLDBART01-TL | ACUPUNCTURE ARTISTIQUE (Jeu de 4) * |
| *PLJVBA0113-AN | JAN VAN BAARLE Nose Reflex Zone - English - plasticised 15x21 cm* |
| *PLJVBA0113-FR | JAN VAN BAARLE Zones Réflexes du Nez 15x21 cm* |
| *PLJVBA0114-AN | JAN VAN BAARLE Teeth and body - English - plasticised 15x21 cm* |
| *PLJVBA0114-FR | JAN VAN BAARLE Dentition et Corps 15x21 cm* |
| *PLJVBA0118-FR | JAN VAN BAARLE Aromatherapy 16x24 cm* |
| *PLJVBA0210-AN | JAN VAN BAARLE Hand reflexology - English- plasticised 21x30 cm* |
| *PLJVBA0210-FR | JAN VAN BAARLE Zones Réflexes Main 21x30 cm* |
| *PLJVBA0211-AN | JAN VAN BAARLE Foot Reflexology - English - plasticised 21x30 cm* |
| *PLJVBA0212-AN | JAN VAN BAARLE Hand and Foot Reflexology - English - plasticised 21x30 cm* |
| *PLJVBA0212-FR | JAN VAN BAARLE Zones Réflexes Main et Pied 21x30 cm* |
| *PLJVBA0215-AN | JAN VAN BAARLE Reiki - English - plasticised 21x30cm* |
| *PLJVBA0215-FR | JAN VAN BAARLE Reiki 21x30 cm* |
| *PLJVBA0216-AN | JAN VAN BAARLE Colors Therapy - English - plasticised 21x30 cm* |
| *PLJVBA0216-FR | JAN VAN BAARLE Therapies by Colours 21x30 cm* |
| *PLJVBA0217-AN | JAN VAN BAARLE Iridology - English - plasticised 21x30 cm* |
| PLJVBA0217-FR | JAN VAN BAARLE Iridology 21x30 cm |
| *PLJVBA021-AN | JAN VAN BAARLE Meridians - English - plasticised 21x30 cm* |
| PLJVBA021-FR | JAN VAN BAARLE Meridians 21x30 cm |
| *PLJVBA022-AN | JAN VAN BAARLE Acupuncture - Alarme - MU -JU Points - English - plasticised 21x30 cm* |
| PLJVBA022-FR | JAN VAN BAARLE Points Acupuncture 21x30 cm |
| *PLJVBA023-AN | JAN VAN BAARLE Special Points - English - plasticised 21x30 cm* |
| *PLJVBA023-FR | JAN VAN BAARLE Points Spéciaux 21x30 cm* |
| *PLJVBA025-AN | JAN VAN BAARLE Sympto Acu Points - English - plasticised 21x30 cm* |
| *PLJVBA026-AN | JAN VAN BAARLE Acu Points on Face - English - plasticised 21x30 cm* |

Charts (suite)

| Réf. | Authors - Titres - Editors |
|----------------|--|
| * PLJVBA026-FR | JAN VAN BAARLE Points Acu Faciaux 21x30 cm* |
| * PLJVBA027-AN | JAN VAN BAARLE Five Elements - English - plasticised 21x30 cm* |
| * PLJVBA027-FR | JAN VAN BAARLE Cinq Eléments 21x30 cm* |
| * PLJVBA028-AN | JAN VAN BAARLE Chakras - English - plasticised 21x30 cm* |
| PLJVBA028-FR | JAN VAN BAARLE Chakras 21x30 cm |
| * PLJVBA029-AN | JAN VAN BAARLE Ear Reflexology - English - plasticised 21x30 cm* |
| PLJVBA029-FR | JAN VAN BAARLE Zones Réflexes Oreille 21x30 cm |
| * PLJVBA031-AN | JAN VAN BAARLE The Meridians - English - plasticised * 42x60 cm |
| PLJVBA031-FR | JAN VAN BAARLE Méridiens 40x60 cm |
| * PLJVBA032-LL | JAN VAN BAARLE Chronobioenergetica (FR/EN/DE/HO) 40x60 cm* |
| PLPNRN01-AL | P.F.M NOGIER presented par R. NOGIER Cartographie de l'oreille (Darstellung des Ohres), 1 planche Edit. SEDATELEC – German - 40 x 60 cm |
| PLPNRN01-AN | P.F.M NOGIER presented par R. NOGIER Cartographie de l'oreille (Ear chart), 1 planche Edit. SEDATELEC – English - 40 x 60 cm |
| PLPNRN01-ES | P.F.M NOGIER presented par R. NOGIER Cartographie de l'oreille (Cartografia de la oreja), 1 planche Edit. SEDATELEC – Spanish- 40 x 60 cm |
| PLPNRN01-FR | P.F.M.NOGIER presented par R.NOGIER Cartographie de l'oreille Edit.SEDATELEC – French - 40x60cm |
| PLPNRN01-HO | P.F.M NOGIER presented par R. NOGIER Cartographie de l'oreille (Oorkaart), 1 planche Edit. SEDATELEC – Dutch - 40 x 60 cm |
| PLPNRN01-IT | P.F.M NOGIER presented par R. NOGIER Cartographie de l'oreille (Mappa auricolare), 1 planche Edit. SEDATELEC – Italian- 40 x 60 cm |
| PLTOLE01-AN | T. OLESON Ear reflex points chart Edit. HCA-USA 75x120 cm |

(*) Limited to our stocks

Supplied by one in box roll

Notes

| |
|--|
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |
| |

Auriculotherapy, Acupuncture Books and others



In acupuncture

| Réf. | Authors - Titres - Editors |
|-----------------|--|
| * LIAFAU01-FR | A. FAUBERT Traité didactique d'acupuncture Edit. TREDANIEL* |
| * LIALAD01-FR | A. LADE Images et fonction points d'acupuncture 374p. Edit. SATAS* |
| * LIBAPN01-FR | B. AUTROCHE, P. NAVAILH Acupuncture en gynécologie et obstétrique Edit. MALOINE* |
| * LIBDES01-FR | B. DESOUTTER Merveilleux vaisseaux 2° Edit. SAURAMPS MEDICAL* |
| * LIBMAIAL01-FR | B. MAIRE C. MOUGLALIS A. HUCHET Ce qui marche...en Acupuncture 126p. Edit. JOSETTE LYON* |
| * LIBSAL01-FR | B. SALAGNAC Naissance et acupuncture Edit. SATAS* |
| * LIGMAC01-FR | G. MACIOCIA Examen de la langue en médecine chinoise Edit. SATAS* |
| * LIHSOL01-FR | H. SOLINAS Atlas d'acupuncture Chinoise - Topographie Méridiens Edit. MALOINE* |
| * LIJBFG01-FR | J. BOSSY, F. GUEVIN, TAI THU Acupuncture et Médecine traditionnelle Edit. SATAS* |
| * LIJBOSAL01-FR | J. BOSSY J-L. LAFONT J-CI. MAUREL Sémiologie en acupuncture (230p.) Edit. DOIN* |
| * LIJLAF01-FR | JL LAFONT Emergence - Origine et évolution de l'acu... Edit. SATAS* |
| * LIJROS02-FR | J. ROSS Zang fu : Organes et entrailles en médecine chinoise Edit. SATAS* |
| * LILVSU01-FR | L. VANGERMEERSCH, SUN PEI LIN syndromes de bi: rhumatologie en med. chinoise Edit. SATAS* |
| * LIMPIQ01-AN | M. PIQUEMAL Advance in Acupuncture with bioelectricity 122p. – anglais-Edit. Y ARTE SA* |
| * LINTAI01-FR | N. THAI THU Sémiologie et analgésie en acupuncture Edit. RESURGENCE* |
| * LIPDEA01-FR | P. DEADMAN Manuel d'acupuncture Edit. SATAS* |
| * LIREAM01-FR | REVUE ACUPUNCTURE ET MOXIBUSTION 2003 (3 n°) ou 2004 (4 n°)* |
| * LITKAP01-FR | T.J. KAPTCHUK Comprendre la médecine chinoise Edit. SATAS* |
| LITYAM01-AN | T. YAMAMOTO Yamamoto New Scalp Acupuncture 200p. English -Edit. YAMAMOTO REHAB.CLINIC |
| * LIUHAL01-AL | U. HECKER und Al. Lehrbuch und Repetitorium Akupunktur (mit Cdom)(500p) - Edit.HIPPOKRATES* |
| LIYROU05-FR | Y. ROUXEVILLE Abrégé de physiologie à l'usage des acupuncteurs et des réflexothérapeutes (232p.) Edit. Springer |

(*) Limited to our stocks

Auriculotherapy, Acupuncture Books and others

In Auriculotherapy and Auriculomedicine

| Réf. | Authors - Titres - Editors |
|---------------|--|
| * LIALEN02-LL | A. LENTZ Bibliographie d'Acupuncture Auriculaire 350p. Edit. GUENIOT* |
| * LIAPRI01-FR | A. PRINTEMPS Pourquoi se soigner par Auriculomédecine 153p. Edit. GUENIOT* |
| * LIBFNS01-AN | B. L. FRANK N.E. SOLIMAN Auricular Therapy - A comprehensive texte 370p. – English - Edit. AAP* |
| * LIBFNS02-AN | B. L. FRANK N.E. SOLIMAN Atlas of Auricular Therapy / Medicine 100p. 2°– English - Edit. AAP* |
| LIBFNS03-AN | B. L. FRANK N.E. SOLIMAN Pocket Atlas of Auricular Therapy / Medicine – English - Edit. AAP |
| * LIBLEC01-FR | B. LECLERC Précis d'Auriculothérapie Tome I (170 p.) Edit. LECLERC* |
| * LIBLEC04-FR | B. LECLERC Auriculomedicine théorie et pratique 1998 276p Edit. LECLERC* |
| * LIBSTR01-AN | B. STRITTMATTER Ear Acupuncture - Pocket atlas (382 p.) Edit. THIEME* |
| * LIBSTR02-AN | B. STRITTMATTER Identifying and treating blockages to healing (257 p.) – English - Edit. THIEME* |
| LICTER01-FR | C. TERRAL Douleur et Acupuncture Edit. SAURAMP MED. |
| LICVUL01-FR | Ch. VULLIEZ Auriculothérapie appliquée à l'art dentaire Edit. SAURAMPS MEDICAL |
| LICVUL02-FR | Ch. VULLIEZ Intérêt de l'étude neurophy. de la langue en Auriculo 50p.Edit.SEDATELEC |
| * LIFSCO01-IT | F. SCOPPA Lineamenti di Auricoloterapia (175 p.) – Italian- Edit. MARTINA* |
| LIGLEM01-LL | ANNALES DU GLEM 1995 à 2007 115p. |
| * LIGLEM02-LL | RESUME CONFERENCES SYMPOSIUM LYON Mai 2000 textes en version originale 280p. * |
| * LIGLEM03-LL | RESUME CONFERENCES SYMPOSIUM PORTO-RICO Nov. 2002 114p. * |
| * LIJVIL01-ES | J. R. VILLAVERDE Tratado dictatico de auriculomedicina 288p. – Spanish - Edit. MANDALA* |
| * LIKLAN01-AN | K. LANDGREN Ear Acupuncture - Practical Guide 230 p. – English - Edit. ELSEVIER* |
| LIMROM01-AN | M. ROMOLI Auricular Acupuncture Diagnosis – English - Edit. CHURCHILL LIVINSTONE |
| LIPNOG01-AL | P. NOGIER Praktische Einführung in die Aurikulotherapie - German - (163 p.) Edit. SATAS |
| LIPNOG01-AN | P. NOGIER Handbook to Auriculotherapy (165 p.) – English - Edit. SATAS |
| LIPNOG01-FR | P. NOGIER Introduction pratique à l'Auriculothérapie 169p. Edit.SATAS |
| * LIPNOG01-IT | P. NOGIER Introduzione pratica all'auricoloterapia (168 p.)– Italian- Edit. LIBRERIA CORTINA - Torino* |
| * LIPNOG02-FR | P. NOGIER Traité d'Auriculothérapie 326p .Edit. MAISONNEUVE* |
| * LIPNOG04-FR | P. NOGIER L'homme dans l'oreille 256p. Edit. MAISONNEUVE* |
| LIPNOG08-LL | P. NOGIER frequential Correspondences (4 tongues) (30 p.) Edit. SAURANPS MEDICAL |
| LIPNOGAL07-FR | P. NOGIER ET AL Revue auriculomédecine 1975 à 1984 Edit. MAISONNEUVE |
| * LIRBOU02-FR | R. BOURDIOL L'Auriculo-somatologie 264p. Edit. MAISONNEUVE* |
| * LIRHMF01-AL | R. HELLING M. FELDMEIER Aurikulomedizin nach Nogier (97 p.) - German - Edit. HIPPOKRATES* |
| * LIRNOG01-ES | R. NOGIER Introduccion practica a la medicina auricular. La fotopercepcion cutanea (145 p.) Edit. Pax México-México* – Spanish - |
| LIRNOG01-FR | R. NOGIER Intro. Prat. à l'Auriculomédecine/Photoperception cutanée 152p. Edit. HAUG |
| * LIRNOG01-IT | R. NOGIER Introduzione pratica all'auricolomedicina La fotopercezione cutanea (151 p.) – Italian- Edit. LIBRERIA CORTINA - Torino* |
| * LIRNOG04-AL | R. NOGIER Aurikulotherapie 1. Grades (97 p.) - German - Edit. DR NOGIER/DR BAUDET* |
| LIRNOG04-FR | R. NOGIER Auriculothérapie 1er degré Edit. SAURAMPS MEDICAL |
| LIRNOG06-FR | R. NOGIER Auriculothérapie 2ème degré 80p. Edit. SAURAMPS MEDICAL |
| LIRNOG07-AN | R. NOGIER Auriculotherapy – English - Edit. THIEME |
| LIRSAN03-FR | R. SANANES Guide de l'Arthrose - Homéo Edit. RESURGENCE |
| * LITOLE03-AN | T. OLESON Auriculotherapy Manual 359 p. – English - 3° Edit. * |
| * LIVMAR01-ES | V. MARTINEZ El desafío de la auriculoterapia (196 p.) – Spanish - Edit. MANDALA-LIBR. VERDE* |
| LIYMYR01-FR | Y. MEAS/ Y. ROUXEVILLE Ce qui marche...en AURICULOthérapie 118p Edit. JOSETTE LYON |
| LIYMYR02-FR | Y. MEAS/ Y. ROUXEVILLE Panorama de l'auriculothérapie et de l'auriculomédecine 385p Edit. SPRINGER |
| LIYROU01-FR | Y. ROUXEVILLE Acupuncture auriculaire personnalisée Edit. SAURAMPS MEDICAL |
| * LIYROU02-FR | Y. ROUXEVILLE ET AL. Actes d'Alba 2003 216p. Edit.SAURAMPS MEDICAL * |
| LIYRMB01-EN | Y. ROUXEVILLE / M. LEBEL A course in Auriculomedicine 2012 105p. Edit. Dux Lucis Books |
| LICAMAR01-FR | ICAMAR French Review 2010 – 2011 237p. |
| LIYROU04-FR | Y. ROUXEVILLE ET AL. Auriculothérapie - Acupuncture Auriculaire 323p. Edit. SPRINGER |

Auriculotherapy, Acupuncture Books and others

In Laser

| Réf. | Authors - Titres - Editors |
|---------------|---|
| * LICVRN01-FR | Ch. VULLIEZ-R. NOGIER Précis de laserthérapie en odonto-stomatologie Edit. JOLLOIS* |
| * LIDBAX01-AN | G. D. BAXTER Therapeutic lasers, theory and practice (260 p.) –English- Edit. CHURCHILL LIVINGSTONE* |
| * LILHJT01-AN | L. HODE, J. TUNER The Laser Therapy Handbook (590 p.) – English - Edit. PRIMA BOOKS* |

In Homeopathy

| Réf. | Authors - Titres - Editors |
|---------------|--|
| * LICVUL03-FR | Ch. VULLIEZ Homéopathie pour le soin de la dent Edit. MARCO PIETTER |
| * LILVAN01-FR | L. VANNIER L'homéopathie sans peine Edit. RESURGENCE |
| * LIMDEN01-FR | M. DENIS Homéo Dermatologie Edit. MAISONNEUVE |
| * LIRSAN01-FR | R. SANANES Manifestations allergiques Edit. RESURGENCE |
| * LIRSAN02-FR | R. SANANES Penser et écrire en Homéopathie Edit. RESURGENCE |

In Osteopathy

| Réf. | Authors - Titres - Editors |
|---------------|---|
| * LICAGE01-FR | C. AGERON MARQUE guide pratique d'osteopathie en gynecologie Edit. SATAS* |
| * LIETOR01-ES | E. TORTI Terapeutica humana - Manipulacion vertebral (155 p.) Edit. TORTI* |
| * LIJGES01-FR | J. GESRET Acupuncture et Ostéopathie vérité neurophysiologique* |
| * LIJUPL01-FR | J.E. UPLEDGER-VR Thérapie cranio-sacree tome 1 Edit. SATAS* |
| * LIJUPL02-FR | J.E. UPLEDGER-VR Thérapie cranio-sacree tome 2 Edit. SATAS* |
| * LIPTRI01-FR | P. TRICOT Ce qui marche...en Osthéopathie 126p. Edit. JOSETTE LYON* |
| * LITLTD01-FR | T. LIEM, T.K. DOBLER Guide pratique d'ostéopathie - Techniques pariétales Edit. MALOINE* |
| * LIWSUT01-FR | W.G. SUTHERLAND Enseignement dans la science de l'ostéopathie Edit. SATAS* |

Others

| Réf. | Authors - Titres - Editors |
|----------------|---|
| * LIAVDB01-FR | A. VAN DER BURG Magnétothérapie 346p. Edit. RESURGENCE* |
| * LICKRE01-FR | C. KREPS La kinésiologie 560 p. Edit. RESURGENCE* |
| LIDHEN01-FR | D. HENRY La médecine des ventouses TOME 1 - 165p Edit. GUY TREDANIEL |
| LIDHEN02-FR | D. HENRY La médecine des ventouses TOME 2 - 190p Edit. GUY TREDANIEL |
| * LIELUB01-FR | E. LUBINIC Manuel pratique d'Aromathérapie Edit. VIGOT* |
| * LIGLIP01-FR | G. LIPPERT Les corps cathédrales RESURGENCE* |
| * LIJC DP01-FR | J-M. CLAY, D.M. POUNDS Massotherapie clinique Edit. MALOINE* |
| * LIJC JL01-FR | J.R. CHENARD, J. LI Maux de dos - Med. chinoise - Med. moderne Edit. MALOINE* |
| * LIJTRA01-FR | J.G. TRAVELL Douleur trouble fonctionnels myofaciaux vol I Edit. HAUG* |
| * LIJTRA02-FR | J.G. TRAVELL Douleur trouble fonctionnels myofaciaux vol II Edit. HAUG* |
| * LIMDOS01-FR | M. DOSCH Atlas des techniques de neuralthérapie Edit. HAUG* |
| * LIMDOS02-FR | M. DOSCH Introduction aux techniques de neuralthérapie Edit. HAUG* |
| * LIMGJD01-FR | M. GARNIER, J. DELAMARE Dictionnaire illustre des termes de médecine Edit. MALOINE* |
| LINVOL01-FR | N. VOLF Etre jeune à tout âge Edit. XO |
| LIPLEU01-FR | P.F. de LEUZE De l'aiguille au coeur de mes patients Edit. LES 2 ENCRE |
| * LIRGAUT01-FR | R. GAUTIER Votre GSM - Votre santé Edit. RESURGENCE* |
| * LIRNOG02-FR | R. NOGIER Ce lait qui menace les femmes 230p. Edit. du ROCHER* |
| * LIRNOG03-FR | R. NOGIER Des Allergies cachées aux maladies du sein 152p. Edit. SAURAMPS MEDICAL* |
| * LIRNOG05-FR | R. NOGIER Le bonheur est dans l'assiette 260p. Edit. AMBRE* |
| LIRNOG08-FR | R. NOGIER Alimentation personnalisée Allergie et intolérances alimentaires 320p. Edit. Sully |
| * LIRSAN04-FR | R. SANANES Insomnie - Douleur de vivre Edit. TREDANIEL* |
| * LISDEB02-FR | S. DEBARD L'électro galvanisme buccal 103p. Edit. DHTM* |

(*) Limited to our stocks

Index

| | |
|---|------|
| Adhesive plasts ASP® Original | p.12 |
| AGISCOP D® and AGISCOP D® Gold Silber | p.14 |
| AGISCOP DT® and AGISCOP DT® Gold Silber | p.14 |
| AGISTIM duo | p.31 |
| ASE, Electrolipolysis sterile needle for AGISTIM Duo | p.33 |
| Anotest – Test-Rings for filters, ATEST 10 | p.25 |
| ASP® Original CLASSIC | p.11 |
| ASP® Original Gold | p.11 |
| ASP® Original Titanium | p.12 |
| ASP® Original Test box | p.12 |
| ASP® Perma Regular (Ø 0,8 mm) | p.13 |
| ASP® Perma Small (Ø 0,6 mm) | p.13 |
| Base Modulo® 100 without stand and battery charge | p.38 |
| Base Modulo® 200 without stand and battery charge | p.38 |
| Books for Acupuncture | p.47 |
| Books for Auriculomedecine | p.48 |
| Books for Auriculotherapy | p.48 |
| Books (others) | p.49 |
| Charts (Auriculotherapy, Acupuncture and others) | p.45 |
| CD Audio Dr Paul NOGIER | p.9 |
| Concentrator for Premio® 10 moxa | p.27 |
| DB 165, detector with beating | p.19 |
| Detector (Black – White) | p.17 |
| Detector (Gold – Aluminium) | p.17 |
| Detector (Gold – Silver) | p.17 |
| Detector (North – South) | p.17 |
| Dim bright lamp (Heine Lamp ou LEP) | p.17 |
| Ear stamps (wood or self – inking) | p.7 |
| Electrical or Micro – currents modules for Modulo® 100 or Modulo 200® | p.39 |
| Endo – nasal stylets | p.44 |
| External Grounding handle for Modulo® 100 or Modulo 200® | p.40 |
| GIR 30 + | p.20 |
| Glass rods | p.8 |
| Heine Lamp or LEP (Dim bright lamp) | p.17 |
| Iriscope | p.44 |
| Kit 4 needles module for AGISTIM Duo | p.33 |
| L Filters | p.23 |
| Large magnet | p.18 |
| LASER Safety glasses | p.41 |
| Light Modules for Modulo® 100 or Modulo 200® | p.41 |
| Modulo® 100 Acu 15 W Pulsed IR LASER | p.29 |
| Modulo® 100 Acu 40 W Pulsed IR LASER | p.30 |
| Modulo® 100 Acu 40 Superpulsed IR LASER | p.30 |
| Modulo® 100 Acu 75 W Pulsed IR LASER | p.30 |
| Modulo® 100 Acupuncture | p.28 |
| Modulo® 100 Acupuncture Multi IR LASER | p.43 |

Index (following)

| | |
|--|------|
| Modulo [®] 100 Auriculo GS | p.16 |
| Modulo [®] 100 Auriculo GS Automatic | p.16 |
| Modulo [®] 100 Micro-currents Reflexo Acutrodes | p.29 |
| Modulo [®] 100 Micro-currents Reflexo Needles | p.29 |
| Modulo [®] 200 Auriculomedecine | p.21 |
| Modulo [®] 200 Dental Fiber Superpulsed IR LASER | p.36 |
| Modulo [®] 200 Micro-currents Acutrodes | p.37 |
| Modulo [®] 200 Micro-currents Needles | p.37 |
| Modulo [®] 200 Physio Superpulsed IR LASER | p.36 |
| Modulo [®] 200 Pulsed IR LASER 15 W | p.34 |
| Modulo [®] 200 Pulsed IR LASER 40 W | p.35 |
| Modulo [®] 200 Pulsed IR LASER 75 W | p.35 |
| Modulo [®] 200 Stomato Superpulsed IR LASER | p.43 |
| Modulo [®] 200 Superpulsed IR LASER 40 W | p.35 |
| Micro-currents Acuneedles Module for Modulo [®] 100 or Modulo 200 [®] | p.39 |
| Micro-currents Acutrodes Module for Modulo [®] 100 or Modulo 200 [®] | p.39 |
| Needles connectors for AGISTIM Duo | p.33 |
| NOGIER colour filters | p.22 |
| Optical fiber for 30W Superpulsed IR LASER module | p.36 |
| Patient grounding clip | p.9 |
| PCFPN, Frequencial Colour Programme Paul Nogier | p.23 |
| Polaroid plates | p.24 |
| Polartrons North, South, Unidirectional North, Unidirectional South | p.18 |
| Premio [®] 10 moxa | p.26 |
| Pressure feelers, black 130 g, blue 250 g, orange 400 g, beige 650 g | p.8 |
| Protractors (30°, anti – tobacco) | p.8 |
| Pulsed LASER Modules for Modulo [®] 100 or Modulo 200 [®] | p.40 |
| Pulsed MicroLaser Module for Modulo [®] 100 or Modulo 200 [®] | p.41 |
| Reusable stimulation electrode for AGISTIM Duo | p.33 |
| Rubber ears (left and right) | p.7 |
| Safety LASER goggles for Modulo [®] 100 or Modulo 200 [®] | p.41 |
| Software and CD-ROM NATOM | p.9 |
| Specific filters | p.24 |
| Stand and battery charger for Modulo [®] 100 or Modulo 200 [®] | p.41 |
| Storage box to store 20 filters | p.25 |
| Styl Modul for AGISTIM Duo (Acupuncture, Auriculartherapy) | p.32 |
| Superpulsed LASER Modules for Modulo [®] 100 or Modulo 200 [®] | p.40 |
| Superpulsed LASER Module with optic fiber for Modulo [®] 100 or Modulo 200 [®] | p.40 |
| Tipped applicator | p.44 |
| VAS clip with cable | p.18 |
| VAS clip without cable | p.18 |
| Charts (Auriculotherapy, Acupuncture and others) | p.45 |

Sales Conditions

Our list "rates and references products" **free of taxes** is valid for all countries **recognizing the CE marking**, and the Canada. Only the Original ASP needles, auriculo instruments (pressure feelers, rubber ears ...) and the PREMIO 10 moxa are available in the United States. *This last intended use is "PREMIO 10 moxa, a non-invasive infrared lamp intended to provide topical heating. It is indicated for temporary relief of minor muscle and joint pain, arthritis and muscle spasm, relieving stiffness, promoting relaxation of muscle tissue, and to temporarily increase local blood circulation where heat is indicated."*

GUARANTEE

Our machines are guaranteed for 2 years parts and labour against any defect of manufacture, except the parts wear (wires, sensors, batteries...) and for use in conformity with the operating instructions. Especially are not covered, breakage due to falls or crashes, dumping or reversal of liquids, as well as spreading coulures of batteries in the box. Small instruments are guaranteed 1 year.

RESERVATION OF OWNERSHIP

The merchandises delivered remain Sedatelec's property until full payment of the price in principal and incidental costs.

CLAUSE ATTRIBUTING JURISDICTION

In case of dispute, the Lyon Commercial Court is the sole competent.

INTERNATIONAL CONDITIONS

Shipment conditions are stipulated in «price list 2012 - 2013»

PAYMENT

Payment is due at the order / proforma are established at request.

↪ Bank transfer:

| Address | | | |
|---|--|--------------------|-----------|
| BNP PARIBAS | | | |
| 164 Grande Rue – F-69600 OULLINS | | | |
| Banque | Guichet | N° de compte | Clé |
| 30004 | 00694 | 00024094070 | 22 |
| IBAN | FR76 3000 4006 9400 0240 9407 079 | | |
| BIC | BNPAFRPPTAS | | |

↪ Postal bank transfer:

| Address | | | |
|---|--|--------------------|-----------|
| LA BANQUE POSTALE | | | |
| Centre de Lyon – F-69900 LYON Cedex 20 | | | |
| Banque | Guichet | N° de compte | Clé |
| 20041 | 01007 | 0621038P038 | 22 |
| IBAN | FR78 2004 1010 0706 2103 8P03 822 | | |
| BIC | PSSTFRPLYO | | |

↪ Credit cards: VISA, Mastercard, Amex...

↪ Bank cheque with order

↪ Cash against commercial documents (please supply full bank details)

CONDITIONS in FRANCE

Routing conditions are mentioned in our valid price list 2012-2013.

PAYMENT MODES

↪ Payment can be done by cheque or Visa Card with order.

↪ Any invoice not paid on its due date shall lead, as of right and without any prior official notification, to the payment of a legal penalty of one and a half times the legal interest rate, the collection costs being payable by the debtor (Sapin law 1994).



Purchase Order

Date: .. / .. / Place:

Name: First Name :

Profession: Speciality :

Address:

Postal Code: Town: Country:

Telephone: Mobile:

Email Address:

| Product Ref. | PRODUCT DENOMINATION | Quantity | Unit Price (without tax) €* | Total Price (without tax) € * |
|--------------|----------------------|----------|--------------------------------|----------------------------------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

* Ask for last active price list excluding tax

SHIPMENT + €
TAXES (V.A.T.) + €

| | |
|----------------------|--|
| TOTAL AMOUNT: | |
|----------------------|--|

PAYMENT MODE:

- Cash: , €
 Cheque n°: Bank:
 Credit card: N°
 Exp : Cryptograms :

COMMENTS:

SIGNATURE :

Thank you for having ordered Sedatelec products



Complementary Medecine

Chemin des Mûriers • F-69540 Irigny • (Lyon) France
Tél : +33 (0) 472 663 322 / Fax : +33 (0) 478 508 903
contact@sedatelec.com

www.sedatelec.com