



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY



DX TWIN

The latest innovation in mechanical palpate & roll technology



AESTHETIC

Stretch marks
Facial treatment
Scarring
Cellulite
Neckline
Treatments for men



RELAXATION

Drainage
Heavy legs
Relaxation



SPORT

Tendinitis
Sprain
Muscle spasm
Muscle tear
Stiffness
Performance



THERAPEUTIC

Rheumatology
Sprain
Vascular problems
Scarring
Oedema
Haematomas
Burns

Revolutionary design

The DXtwin has been granted two new patents producing a revolutionary device designed for treating cellulite and redefining body contours.

Triple action

1. **Lipolysis**: that is the melting of fat cells through the elimination of fibrous tissue and the reduction of adipose masses
2. **Circulatory stimulation**: draining and micro-stimulation help eliminate toxins and unblock subcutaneous tissue
3. **Redefinition**: the stimulation of elastin fibres leads to a smoother-looking and better toned skin body contours

Connective tissue mastery

The DXtwin is an indispensable tool in the treatment of cellulite resistant to calorie-control and weight loss

Perfect control


The control panel provides complete and instantaneous regulation of pulse frequency, pressure and treatment cycle.

Multi-purpose

Four modes of treatment cover a wide range of indications

Customisation

Personalise your DXtwin by creating your own design



www.alfaspa.co.uk



DX TWIN - What's in the box?

Features and technical specifications:
 Dimensions (h,w,d): 177.5 x 55.5 x 50.8cm
 Weight: 60kg
 115/230V-50/60Hz-200W
 2-year warranty parts & labour
 International patent
 Registered Trademark & Design

DXtwin

EVOsphere

Double suction cups-double mobilisation

Hands-free application
 Deep elimination of fibrosis
 Toning and firming effect on tissues
 Lipolysis through dipocyte reduction
 Ideal for thighs and buttocks



EVO3D

Toning-Lipolysis-Detox

Palpate and roll action
 3 rollers at different heights, creating 2 skin folds
 Thin rollers, so low suction pressure needed
 Removable rollers for easy maintenance



NEOcup

Drainage-Stimulation-Mobilisation

New, ergonomic shape for easy manipulation
 3 sizes for flexibility of use
 Transparent to allow monitoring of skin reaction



NEOlift

Toning-Drainage-Regeneration

3 suction heads dedicated to facial treatments
 Precision targetting of scars, wrinkles and stretch marks
 Pulsed mode optimised for anti-aging treatment



Professional treatment cream 1L
 2 foam filter disks
 Power cable
 user manual
 certificate of warranty
 treatment protocols book
 Banner



www.alfaspa.co.uk



SPO

The original professional "palpate & roll" device

Indications:

Cellulite
Lymphatic drainage
Heavy legs
Relaxing body massage
Facial massage
Post-exercise/sport

Mobile and easy to use:

The SPO comes with a convenient mobile trolley and shelves that house the treatment heads and accessories.

How it works:

The treatment heads generate a vacuum on the surface of the skin, which is sucked in between the 2 inner rollers of the treatment head. These mobile rollers automatically and permanently adapt to the thickness of the skin fold created.

Palpate & roll:

While sliding the treatment head across, the skin fold is rolled and pressed between the mobile cylinders, which reproduce the palpate & roll massage method.

Benefits:

Slimming and contouring
Reduced appearance of cellulite
Improved circulation
Soothing and relaxing


Lymphatic drainage:

The two spatulas - situated front and back of the mobile rollers - make the resulting lymphatic drainage all the more effective. Finally, the two outer mobile rollers make it easier to slide the treatment head across the treated area.

Targeted treatment:

The basic suction cups mobilise the tissues in the targeted zones, which can be worked either in static or sliding mode.

The 2 patented double suction cups perform a static circular palpate & roll massage, particularly efficient for localised beauty treatments.



www.alfaspa.co.uk



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY



SPO - What's in the box?

Features and technical specifications:

Dimensions (h,w,d): 104.5 x 47.5 x 40.5 cm

Weight: 40kg

115V/230V 50Hz 190W / 115V/230V 60Hz 210W

Continuous or pulsed mode

0-900mbar

Parameters: timer, frequency, cycle, maximum vacuum

2-year warranty parts & labour

International patent

Registered Trademark & Design

2 facial suction cups

These have been specifically designed to treat the face and scars. They promote skin tissue regeneration by stimulating synthesis of collagen and elastin fibres.



2 treatment heads

These rollers constantly adapt to the skin fold created, working on every kind of cellulite and morphology. The 2 spatulas, on the front and back of the rollers, make lymphatic drainage more effective.




4 basic suction cups

The 4 difference sizes are precision designed for each treatment stage: mobilising, defibrolising, relief, drainage



4 double suction cups

These perform a static circular palpate & roll massage on targeted areas, used simultaneously to improve results. Effective for localised beauty treatments, such as old fibrous cellulite.



Main connecting tube

Double suction tube connection

Power cable

Professional massage cream 1l

SPO user manual

Warranty certificate



www.alfaspa.co.uk



S'ONE

The professional "palpate & roll" mobile device

Indications:

Cellulite
Lymphatic drainage
Heavy legs
Relaxing body massage
Facial massage
Post-exercise/sport

Benefits:

Slimming and contouring
Reduced appearance of cellulite
Improved circulation
Soothing and relaxing

How it works:

The S'ONE generates a vacuum on the surface of the skin, which is sucked in between the 2 inner rollers of the treatment head. These mobile rollers automatically and permanently adapt to the thickness of the skin fold created.

Palpate & roll:

While sliding the treatment head across, the skin fold is rolled and pressed between the mobile cylinders, which reproduce the palpate & roll massage method.


Lymphatic drainage:

The two spatulas - situated front and back of the mobile rollers - make the resulting lymphatic drainage all the more effective. Finally, the two outer mobile rollers make it easier to slide the treatment head across the treated area.

Targeted treatment:

The basic suction cups mobilise the tissues in the targeted zones, which can be worked either in static or sliding mode.

The 2 patented double suction cups perform a static circular palpate & roll massage, particularly efficient for localised beauty treatments.





ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY

S'ONE - What's in the box?

The S'ONE - features and technical specifications:

Device dimensions: 285x240x130mm

Case dimensions : 525x156x375mm

Device weight: 3.6kg

115V 60Hz 165W / 230V 50Hz 140W

Continuous or pulsed vacuum

Maximum vacuum: 600mbar

Parameters: timer, frequency, cycle, maximum vacuum

2-year warranty parts & labour

International patent

Registered Trademark & Design



4 basic suction cups

Suitable for all indications;
varying sizes for different
body parts



2 facial suction cups

These have been
specifically designed
for treating the face
and scars.


Power cable

Connection tube

Double suction cups
connection tube

Spare filter

S'ONE® User Manual



Body contouring massage cream

- Combines an active complex of lecithin-caffeine & guarana
- Lubricative and effective for prolonged massages
- Water soluble for ease of maintenance of treatment heads



2 double suction cups

Effective for localised beauty
treatments, such as old fibrous
cellulite.

2 treatment heads

The small head is for
arms and calves. The
larger one is for
extended areas such
as the thighs,
buttocks and back.




www.alfaspa.co.uk



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY



PULSTAR PSE & PULSTAR PSX

Pressure Therapy – Activate / Reshape / Recover

Indications:

Cellulite
Varicose veins
Heavy legs
Post-exercise/sport
Muscle and joint pain
Peripheral oedema
Post-lymph node surgery
Constipation

How it works:


By inflating sleeves that wrap around the body and limbs, PULSTAR creates continuous graduated pressure waves which gently massage interstitial fluid and lymph back into the circulation and help return blood to the heart.

Effects on cellulite:

In stimulating lymphatic circulation, PULSTAR increases the elimination of toxins and fatty deposits, reducing the appearance of unsightly cellulite.

Automatic control:

Unlike other systems, an automatic double pressure gradient is regulated by the PULSTAR unit and inflatable sleeve accessories. The volume of each inflation chamber varies according to the form and structure of the area being treated.



Benefits:

Slimming and contouring
Reduced appearance of cellulite
Personalised treatment
Improved circulation
Soothing and relaxing

Personalised treatment:

A further pressure control automatically regulates the flow of air into each compartment, creating a pressure gradient from top to bottom. This means the device automatically adapts to each client's body morphology, delivering personalised treatment at the touch of a button.

Lymphatic stimulation:

Furthermore, PULSTAR's patented system of connections allows it to inflate and deflate accessories significantly more rapidly than similar products, creating increased circulatory stimulation and drainage.

PULSTAR Technology®:

PULSTAR Technology® was developed by STARVAC from 40 years of experience, knowledge, expertise and patented research in medicine, beauty and well-being, based on techniques pioneered by Dr Emil Vodder. The PULSTAR PSE uses the very best technology in the field of veino-lymphatic drainage.

Active Wave Technology - PULSTAR PSX

The active wave generates waves that ripple from bottom to top of the accessory and back again. The ripple effect generated by the active wave provides a rhythmic pumping motion along the limb being treated, which acts directly on the micro-circulation (the capillaries and small blood vessels), reinforcing the action of the pressure therapy treatment. These effects are enhanced by regular pressure surges.




www.alfaspa.co.uk



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY



What's in the box?



PULSTAR PSE


Inflation & piloting connection tubes, power cable
User manual & warranty certificate

Accessories supplied with the PSE:

5-chamber boots for treatment of legs and thighs
(length 95cm, circumference 75cm)

Optional extras:

2-chamber belt for treatment of the abdomen, hips
and buttocks (length 110cm, height 36cm)
4 chamber sleeve for treatment of arms and forearms
(length 78cm, circumference 63cm)



Features and technical specifications of PSE:

Pressure : 20-80 mmHg
Dimensions (h,l,d) : 12 x 24 x 28cm
Weight : 3.7kg
230V - 50/60Hz - 65W
2-year warranty on parts and labour

PULSTAR PSX


Inflation & piloting connection tubes, power cable
User manual & warranty certificate

Accessories supplied with the PSX:

9-chamber (92cm) or 10-chamber (103cm) boots and
one of the following optional extras

Optional extras:

8-chamber sleeves for arms and forearms (78cm)
3-chamber belt which connects to the boots to treat the
abdomen, hips and buttocks (110cm)
5-chamber belt used alone to treat thighs (110cm)
9-chamber abdominal cushion to treat constipation
(adjustable)
25 pairs of hygienic plastic covers (for boots and sleeves)
Cryoactive drainage cream 400ml bottle (professional size)



Features and technical specifications of PSX:

Pressure : 20-120 mmHg
Dimensions (h,l,d) : 12 x 24 x 28cm
Weight : 4kg
230V - 50/60Hz - 65W
2-year warranty on parts and labour



www.alfaspa.co.uk



ALFA SPA
HEALTH • TECHNOLOGY • BEAUTY



L300

The face wizard



1. Detoxification

Deep-seated impurities are loaded with ions from the cleansing lotion and purifying gel. The current emitted by the VCE heads drives the impurities out of the skin.



2. Ultrasonic peeling

Microvibrations, accompanied by the soothing lotion, gently exfoliate the superficial layers of the skin and eliminate the impurities released, while refining skin pores.



3. Toning

The toning current, emitted by the Vibratory Cellular Electrotherapy (VCE) heads, ensures maximal penetration of the active ingredients of the restructuring gel by making cell membranes more permeable, allowing deep regeneration and visible rejuvenation.



4. Laser smoothing and lifting

Powerful light rays stimulate the dermis, improving fibroblast function and hence elastin and collagen production, reducing wrinkles and producing a lifting effect.

Benefits:

Visible results after the first session

Brows lifted and arches widened

Wrinkles smoothed

Cheekbones lifted

Skin toned

Complexion cleared

Redness diminished



www.alfaspa.co.uk



L300 - What's in the box?

Features and technical specifications:

L300 dimensions (h,w,d): 151 x 380 x 269mm

Stand dimensions (h,w): 915 x 470mm, painted aluminium

Weight: 15.7kg (L300 3.7kg, stand 12kg)

96-240V 50/60Hz 40W

2-year warranty parts & labour

L300
l'expert visage



4 treatment heads

2 VCE heads
& 2 o-rings

Ultrasonic peeling head

Laser head



Creams and gels

IONS+ toning cleansing lotion 400ML

DETOXIFICATION cell purifying gel 400ML

TONIFICATION cell restructuring gel 400ML

ULTRA ACTIVE protective Rejuvenating cream 200ML

REGENERANTE anti-wrinkle serum 200ML

RAFFERMISSANTE firming soothing lotion 400ML

For retail (home kits)

4 x IONS+ toning cleansing lotion 250ML bottle

4 x ULTRA ACTIVE Protective Rejuvenating cream 50ML tube

4 x REGENERANTE anti-wrinkle serum 50ML tube

4 x RAFFERMISSANTE firming soothing lotion 250ML bottle

2 pairs of laser protection goggles

1 return electrode, case & connection cable

2 packs of 100 cotton patches for VCE treatment heads

8 tablets of 6 compact face masks

Trolley, power cable, user manual, treatment protocols



www.alfaspa.co.uk



UX100

Slimming and toning in one device

Indications:

Cellulite
Lymphatic drainage

The UX100 is a completely automatic, hands-free slimming and remodelling device with rapid results. It redefines curves and reduces cellulite through lipolysis, defibrosing fatty deposits and enhancing toxins and fat by improving lymphatic drainage.


The UX100 brings together three treatment techniques in one device:

1. **Low-frequency ultrasound**, which has a defibrosing effect on fatty deposits and stimulates adipocyte lysis.
2. **Extraction currents**, which provoke sustained localised muscle contraction, toning the region through gentle passive exercise and improving the effect of the low-frequency ultrasound.
3. **Muscle fibre stimulation currents**, which ensure that fat and toxins released by lipolysis are effectively cleared and drained by the lymphatic and venous systems.


Benefits:

Slimming, toning and contouring
Reduced appearance of cellulite
Soothing and relaxing

The user-friendly, ergonomic interface allows selection of 31 prerecorded protocols specifically designed for different types of cellulite, e.g.



Soft cellulite, found around the buttocks and thighs



Tonic cellulite, located over the abdomen and waist



Oedemic cellulite, found above the knees and waistline

The ultrasound plates and electrodes are applied with our UX100 Ultrasonic Treatment Gel, specially formulated to enhance the slimming and toning effects of this treatment.

Each plate has an individual on/off switch, allowing exquisite control over treatment duration and location. Warning LEDs light up when the plates are in use.



www.alfaspa.co.uk



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY

UX100 - What's in the box?

Features and technical specifications:

Dimensions (h,w,d): 110 x 57 x 48cm

Weight: 32kg

115V/230V 50Hz 190W / 115V/230V 60Hz 210W

31 pre-recorded protocols

2-year warranty parts & labour

International patent

Registered Trademark & Design

**UX
100**

Un plus en esthétique



Low-frequency ultrasound:

10 low-frequency ultrasound plates with 34 cells

- 4 plates with 4 cells each

(112 x 60mm)

- 6 plates with 3 cells each

(170 x 60mm)

Frequency: 30kHz

Cell power: 0-2W/cm²

Pulsed mode: 12Hz - 3Hz - 1.4Hz



Extraction and stimulation:

8 stimulation electrodes

4 stimulation canals, 8 cells

Output current per canal 100mA

UX100 Ultrasound Treatment Gel

Power cable

UX100 user manual

Warranty certificate



www.alfaspa.co.uk



CREAMS & GELS

Enhanced slimming and toning



la slim
SLIMMING

DXTWIN
SPO
S'ONE
PSE/PSX



Guarana



Horse
chestnut



Ivy

The combination of guarana and horse chestnut promotes deconsolidation of fat tissue and reduction of lipid.

Ivy has powerful agents that firm and smooth skin reducing the peau d'orange effect.

la cryo
DRAINAGE

DXTWIN
SPO
S'ONE
PSE/PSX



Menthol



Camphor



Horse
chestnut

Menthol freshness provides a cooling effect.

Combined with camphor, a toning agent, it relieves heavy legs and swollen feet.

The use of horse chestnut increases capillary resistance, stimulates venous tone, and facilitates absorption of oedema.

la neutre
MASSAGE

DXTWIN
SPO
S'ONE



Sunflower
oil



Vitamin E

Sunflower oil combined with vitamin E helps to rebuild the lipid film of the skin and has softening and soothing properties.

It lubricates and provides a pleasant sensation to the lipomassage.



www.alfaspa.co.uk



CREAMS & GELS

Enhanced slimming and toning



legel
CLARIFYING

L300
DXTWIN
SPO
S'ONE



Orange

Lemon



Green tea

leserum
RE-DEFENSIFYING

L300
DXTWIN
SPO
S'ONE



Ginseng



Horsetail

Hyaluronic acid

lacrème
REGENERATING

L300
DXTWIN
SPO
S'ONE
PSE/PSX



Argan oil



Avocado



Royal jelly

Marine collagen



lotion
SOOTHING

L300
DXTWIN
SPO
S'ONE
PSE/PSX



Chamomile



Aloe vera



Orange

Evens skin texture, stimulates cells and restores radiance.

Orange extract restores brilliance and clarity.

Green tea purifies and activates microcirculation.

Lemon extract evens and illuminates.

Toning and firming, anti-aging, it restores youth to treated skin.

Hyaluronic acid has smoothing and anti-aging effects.

Ginseng stimulates collagen production and repairs.

Horsetail promotes tissue elasticity, firming and remineralization.

Argan oil regenerates skin cells.

Avocado oil has a moisturizing and restructuring power.

Royal jelly stimulates the production of collagen and offers deep nourishment.

Marine collagen improves the flexibility and elasticity of the skin.

Calming and soothing, reduces redness and refreshes the face.


Chamomile soothes and relaxes the muscles.

Aloe vera provides moisture and freshness

Orange extract gently soothes and refreshes the skin.



www.alfaspa.co.uk



ALFA SPA

HEALTH • TECHNOLOGY • BEAUTY

www.alfaspa.co.uk