

HOME

- [Our Principles](#)
- [Sales Orders & Enquiries](#)
- [Customer Support](#)
- [mybiomerieux.com & QC Certificates](#)
- [History bioMérieux](#)
- [Quality & Environment](#)
- [bioMérieux Goes Green](#)
- [Terms & Conditions](#)
- [Legal](#)
- [How to use this Product List](#)



● BACTERIOLOGY ● IMMUNOASSAYS ● MOLECULAR DIAGNOSTICS

*Welcome to the bioMérieux
Clinical PDF Product List.*

*Use the Head Tabs at the top of each page to
skip to each section, then click the titles in the
Tab Menu to find each product group.*



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Our Principles



Customers come First

We must satisfy our customers
- clinicians and laboratory professionals -
with high quality, innovative products and
services to protect and improve the health
of patients and consumers worldwide.

Employees are our Greatest Asset

We must respect and recognise the work and
the ideas of all employees, no matter what
their job. We work together towards a
common goal: satisfy our customers. We must
provide employees with the training and
resources they need to achieve this.

We Belong to a Community

We have social, environmental and economic
responsibilities towards the local communities
where our sites are located. We also have
responsibilities towards the global
community as a partner in public health.

Our Shareholders are Partners

We must communicate company information in
a transparent, thorough and timely manner.
We must grow bioMérieux's sales and
ensure that shareholder investments
provide a proper return.

HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List

Sales Orders & Enquiries

Our telephone lines are open Monday to Friday (except Bank Holidays) according to the info of office addresses below:

 Denmark

Order Team
Telephone: 70 10 84 02
Fax: 70 10 84 01
order.dk@biomerieux.com

General Enquiries
Telephone: 70 10 84 00
Fax: 70 10 84 01
contact.dk@biomerieux.com

Opening Hours
Monday to Thursday
from 8.30am to 4.00pm
Fridays 8.30am to 3.00pm

bioMérieux Denmark ApS
Smedeholm 13C
2730 Herlev

 Finland

Order Team
Telephone: 09 8545 6040
Fax: 09 8545 6045
order.fi@biomerieux.com

General Enquiries
Telephone: 09 8545 600
Fax: 09 8545 6045
contact.fi@biomerieux.com

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Suomi Oy
Konalantie 47C
00390 Helsinki

 Norway

Order Team
Telephone: 2337 5552
Fax: 2337 5551
order.no@biomerieux.com

General Enquiries
Telephone: 2337 5550
Fax: 2337 5551
contact.no@biomerieux.com

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Norge AS
Ulvenveien 75A, 3 etg.
0581 Oslo

 Sweden

Order Team
Telephone: 031 68 84 96
Fax: 031 68 48 48
order.se@biomerieux.com


General Enquiries
Telephone: 031 68 84 90
Fax: 031 68 48 48
contact.se@biomerieux.com

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Sweden AB
Hantverkarsvägen 15
43633 Askim

Quotations

If you have received a quotation from your representative please have the quote reference to hand when placing an order to ensure the correct pricing is set up.



Customer Support

Support enquiry e-mails are logged by our Customer Support team. Our telephone hotlines are open Monday to Friday (except Bank Holidays) between the times listed below.

Application Support

- To help with all of your enquiries we have a team of in-house or field Application Specialists
- If our team can deal with your enquiry over the telephone they will endeavour to do so
- If this is not possible or you require more hands-on help with an evaluation, validation or training then your enquiry will be escalated to our field Application Specialists who may then arrange to visit your laboratory

Technical Support

- Our field Service Engineers install equipment, perform preventative maintenance (servicing) and provide emergency breakdown cover in your laboratory
- The service team perform instrument pre-delivery inspections, repairs and airIDEAL® air flow calibrations

**Denmark**

Hotline: 44 508 408

Opening Hours
Monday to Thursday
from 8.30am to 4.00pm
Fridays 8.30am to 3.00pm

bioMérieux Denmark ApS
Smedeholm 13C
2730 Herlev

**Finland**

Hotline: 09 8545 6020

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Suomi Oy
Konalantie 47C
00390 Helsinki

**Norway**

Hotline: 2337 5555

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Norge AS
Ulvenveien 75A, 3 etg.
0581 Oslo

**Sweden**


Hotline: 031 68 58 58

Opening Hours
Monday to Friday
from 8.30am to 4.00pm

bioMérieux Sweden AB
Hantverkarsvägen 15
43633 Askim

E-mail address for all Nordic Countries

customersupport.nordic@biomerieux.com



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List

mybiomerieux.com


bioMérieux customer portal www.mybiomerieux.com is designed to be part of the future e-platform.

www.mybiomerieux.com offers an easy account creation and hosts Technical Library as well as coming applications like validation packages, distance learning, and others.

Once you have registered, you can access:

- Package Inserts
- Quality Control Certificates
- Material Safety Data Sheets in English (more about MSDS in Nordic languages in "Quality and Environment")
- Thermal Shock Validations
- Instrument Manuals
- Validation guides and addenda

To find Instrument user manuals or validation addenda, choose BROWSE BY SYSTEMS from the BROWSE BY PRODUCT RANGE drop down menu and this will bring you to a list of our instruments.



Technical Library page for www.mybiomerieux.com

Quality Control Certificates


Our pre-poured media products must be passed by our Quality Control Departments and a certificate produced before they are released for sale. The QC Certificates for these and other products can be accessed directly via www.mybiomerieux.com. You can search for these using the product name, reference number or the lot number.

Expiry Dates

For pre-poured media and kits, the expiry date for a specific batch is on the top of the QC Certificate. Product can be used on the day of expiry.

LyfoCults® Plus

To obtain Certificates of Quality for our LyfoCults Plus products, visit www.pmlmicro.com



QC Certificate



History bioMérieux

Over a Century in the Service of Global Public Health

PAGE 1 OF 2

PAGE 2 OF 2

The commitment of the Mérieux family to work in the field of biology goes back to 1897, with the establishment by Marcel Mérieux, a pupil of Louis Pasteur and Emile Roux, of Institut Mérieux. This marks the start of an extraordinary adventure in both biology and industry, which was to have a profound impact on vaccinology and on in vitro diagnostics at the global level.

Initially under the direction of Doctor Charles Mérieux and later of Alain Mérieux, Institut Mérieux was to become a world leader in the field of human and veterinary vaccines. It gave birth to many bodies which are still major players today in the public health arena:

- In human medicine, Pasteur Mérieux Connaught which later became Aventis Pasteur and then Sanofi Pasteur
- In veterinary medicine, IFFA (Institut Français de Fièvre Aphteuse) which became Rhône Mérieux, and later Merial

In parallel to this work on vaccines, in 1963 Alain Mérieux set up bioMérieux to work on in vitro diagnostics.

The family's corporate holding was later enlarged through the acquisition of Transgene, working in immunotherapy, and then the American industrial control company Silliker, followed by ABL (a North American company researching into virology).



Marcel Mérieux, a student of Louis Pasteur, creates the Institut Mérieux (Year 1897).



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List

History bioMérieux

Over a Century in the Service of Global Public Health

The name of Mérieux Alliance was chosen to reflect the full scope of these activities. Under the impetus of Alain Mérieux and his sons, Christophe, Rodolphe and Alexandre, the group developed further, fostering a policy of scientific innovation and international openness by developing a network of sites all over the world.

In 2009, Mérieux Alliance is changing its name and again taking up the standard of “**Institut Mérieux**” to communicate its vision of biology and to take on the new challenges in global public health.

Institut Mérieux has thus set up Mérieux Développement, a company investing in the area of health. 2010 has seen the Institute increase its commitment to Nutrition/Health with the structuring of Mérieux NutriSciences grouping together Silliker, Biofortis and Bioagri. It has also created IMAccess, a company dedicated to biology without borders between the North and South.

Human Health

Vaccines - Biotherapies

2000 > Sanofi Pasteur*

2000 > Aventis Pasteur*


1994 > Pasteur Mérieux Connaught*

Nutrition/Health

Pasteur Mérieux
Sérums & Vaccins

1985 > Pasteur Vaccins

Diagnosis



Immunotherapy

1947 > Creation of IFFA 1968 > IFFA Mérieux 1983 > 1994 > Rhône-Mérieux

Animal Health

Vaccines - Médicaments

1997 >

2004 > Merial*

* Companies founded by the Mérieux family but no longer part of Institut Mériuex since 1994.



Quality & Environment

QMS

QA

QC


MSDS

Quality Policy

At bioMérieux, we never forget that every day, somewhere in the world, a clinician makes a decision for a treatment, or a food product is released for consumption, based on the quality of a test we are providing.

Our tests must provide an optimum level of performance because our customers are themselves responsible for the quality of the test results they provide to those who depend on them. This is why we consider quality to be everybody's responsibility.

bioMérieux's Quality Policy guides our internal processes and drives good business outcomes and efficiencies. In line with regulatory requirements, it also draws upon input from the company's sites and subsidiaries worldwide to ensure that the policy is understood and embraced by all employees.



Quality by Design

We have our minds set on doing better every day and when it comes to our Quality Policy, we adopt a holistic approach. Our goal is quality by design in everything we do: our aim is not to create quality by checking, but quality right from the start by integrating it into every aspect of our work.

We have a [Global Quality Management System](#) to control our processes and manage risk worldwide, certified against ISO 9001. This approach helps assure fulfilment of regulatory requirements. It leads to decreased production cycle times, reduced product and service nonconformities and an optimised supply chain. Products can thus be launched rapidly on the market and processes can be adapted to satisfy the needs of our customers.

Quality is Everybody's Responsibility

Quality is not the job of the Quality Organisation only: it is a shared responsibility by all employees, affecting everybody from R&D and manufacturing to sales and marketing, as well as support functions.

Performance indicators engage all employees in the achievement of quality objectives. Continuous improvement initiatives are in place, and best practices are exchanged among the company's teams worldwide.

It is our responsibility to help customers do their job by understanding their needs and making sure that we meet them. Our Quality Policy is key to achieving this.



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Quality & Environment

QMS

QA

QC

MSDS



To meet the high expectations of our customers, our Quality Assurance system is designed to ensure the reproducibility of product performance throughout a product's lifecycle. Each new product is validated following stringent performance criteria required by our customers and by regulatory bodies such as the FDA, AFSSAPS and other organisations such as LNE/G-MED, DNV and specified by pharmacopoeias, food accreditation bodies and ISO standards.

Quality Regulations

Our Quality Assurance and Regulatory Affairs Departments must approve all products before release. All development and validation documentation is archived and can be accessed during an audit. The following are requirements applicable to all bioMérieux entities:

- ISO 9001:2008, Quality management systems – Requirements
- ISO 13485:2003, Medical-Devices – Quality Management Systems – Requirements for regulatory purposes
- ISO 14971:2007, Application of Risk Management to Medical Devices
- FDA cGMP 21 CFR Part 820 (current Good manufacturing Practises) Quality System Regulation
- Culture media ranges manufactured for use in the pharmaceutical industry are compliant with the US/EU/J pharmacopoeia requirements
- All our products meet, as applicable, the requirements of:
 - European Medical Device Directive 93/42/EEC
 - European In Vitro Diagnostic Directive 98/79/EC
 - US/EP/JP Pharmacopoeia
- Our subsidiaries are compliant with:
 - European Packaging and Packaging Waste Directive 94/62/EC
 - Waste Electrical and Electronic Equipment Directive 2002/96/EC

Copies of our certificates can be obtained upon request.

Inspections at Our Production Sites & Subsidiaries

- 25–30 customer audits a year
- International health authorities such as the FDA and AFSSAPS
- ISO 9000 family and ISO 13485 audits

Raw Material Management

- Strict selection and control of raw materials to maintain the reproducibility of performance between batches
- Control of the supply chain for raw materials
- Raw materials are blended for each product
- Full traceability from raw material to finished product
- TSE/BSE Status Certification

Culture Media Production

- Plate production in clean rooms with environmental monitoring and controls
- Filling areas are operated to ISO 5 classification
- Special manufacture for customer specified formulations, filling volumes and presentations for plates and bottles



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Quality & Environment

QMS

QA

QC

MSDS

Performance & Shelf-life

- Evaluated for each product according to:
 - EN 13640:2002 Stability Testing of in-vitro diagnostic reagents
 - EN 1322:1999 In-vitro Diagnostic Devices – Culture Media for Microbiology – Performance Criteria for Culture Media
- Using validated methods
- Inter-lot and Intra-lot variability analysed to prevent any drift in performance (USP 1117)
- Stability testing:
 - In real-time on 3 batches
 - Includes thermal shocks to guarantee shipment and temperature variations will not impact on performance
 - Pharmaceutical culture media – wild and pharmacopoeia strains tested in development
 - Routinely re-evaluated after launch
 - Extends to beyond the expiry given so that **media can be used on the last day of shelf-life** with no impact on performance

Quality Control

- Microbiological and physical criteria for raw materials and final product
- Includes: growth promotion (recovery rates, colony morphology), pH, colour, agar strength
- Pharmaceutical culture media are tested using microbial strains defined by the USP/EP/JP
- >50% recovery rate compared with a gold standard media
- Unique ≥ 75% recover rate compared with a previously passed batch of the same media for all 3P™ media
- Non-sterile pharmaceutical testing culture media complies with the harmonised chapters of the pharmacopoeias:
 - Microbiological examination of non-sterile products: microbial enumeration tests (EP 2.6.12, USP 61, JP 35.1)
 - Microbiological examination of non sterile products: tests for specified microorganisms (EP 2.6.13, USP 62, JP 35.2)

HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Quality & Environment

- QMS
- QA
- QC
- MSDS

Material Safety Data Sheets

An MSDS provides lab personnel with procedures for handling or working with dangerous products in a safe manner, and include information such as physical data, toxicity, health effects, storage, disposal, protective equipment, and waste-handling procedures.

bioMérieux reagents conform to all European and national regulations applicable to chemical substances that might have an impact on health or environment.

When required MSDS exist for our products. They can easily be retrieved from www.mybiomerieux.com under Technical Library in English. MSDS in Danish, Finnish, Norwegian and Swedish can be found at www.ecoonline.com where customer own report creation tools also exist.

To find our Nordic MSDS easily they can be searched by article number and language wanted.



bioMérieux Goes Green – Environmental Commitment

Five Key Areas

Throughout the world, bioMérieux is going green. On each site, in every subsidiary, a green champion is in charge of ensuring we achieve our environmental goals in the following five areas of focus.

Energy

Reduce our Consumption of Non-Renewable Energy Sources.

Conscious that continued use of non-renewable energy sources will cause greater environmental problems, we want to use energy as responsibly and efficiently as possible. We have implemented an eco-efficient program focused on three major areas: reduce consumption, increase equipment efficiency, experiment with alternative energy sources.

For examples of how we are achieving our goal, [click here](#).

Water

Conserve Water Through Good Management of Resources.

Water is the resource we use the most in our manufacturing processes. To reduce consumption, we favour closed loop circuits and are actively replacing once-through cooling systems.

For examples of how we are achieving our goal, [click here](#).

Paper

Reduce Paper Usage in Processes, Packages and Products.

As we are constantly growing, we are faced with the issue of increased paper consumption. We are meeting this challenge by changing our product packaging and recycling paper.

For examples of how we are achieving our goal, [click here](#).

Waste

Minimise Adverse Impacts on Health and the Environment.

We have been optimising our waste management for a number of years, reducing the weight of landfill waste, recycling and sorting non-hazardous waste. Specialised companies handle hazardous waste.

For examples of how we are achieving our goal, [click here](#).

Emissions

Minimise Adverse Impacts on Health and the Environment.

Global warming due to greenhouse gas emissions is a major problem. Because of the close link between CO2 emissions and global warming, we are focusing our efforts on high CO2 emission sources that have a significant impact on the atmosphere.

For examples of how we are achieving our goal, [click here](#).



HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Terms & Conditions

1. ORDERS

- To avoid errors, each order should specify:
 - product name and reference number
 - the addresses for delivery and invoicing
- bioMérieux's acceptance of any order shall be conditional upon creditworthiness of the Buyer. In particular, bioMérieux reserves the right to suspend the order if the situation of the Buyer visibly presents a risk with regards to the recovery of the outstanding debts of bioMérieux, or in case of unavailability of the Product. Concerning the participation in handling costs and delivery prices:
 - bioMérieux charges an administration fee on each order
 - bioMérieux charges a dangerous goods fee for shipments that include products with UN numbers
 - Express delivery costs, when such deliveries are asked by the customer, are charged to the customer

2. SUBSCRIPTION CONDITIONS

- Adjustable subscriptions are available for usual and regular consumer goods.
- Notwithstanding the above, subscription shall be required for special products, cells, blood and red cells of animals.
- Given the special treatment of these Products, an additional cost per shipment shall be applied.

3. DELIVERIES

- Sales of Products shall be made DDP (Delivered Duty Paid, Incoterms 2010 ICC) to the delivery address specified by the customer in the order. The Products will not be delivered to specific departments or rooms inside the building, unless a special agreement has been made with bioMérieux Sales department.
- bioMérieux undertakes to provide the Buyer with the Products (instruments, reagents and disposables) indicated on the order accepted by bioMérieux. Products will be delivered and installed as specified in Buyer's purchase order accepted by bioMérieux.
- bioMérieux undertakes to deliver goods from stock without delay. In case of temporary unavailability, the customer will receive this information from bioMérieux. The shipment of the outstanding balance of the order shall be done within the frame.

4. PRICES

- Prices of Products are quoted net ex-works (Incoterms 2010 ICC).
- VAT shall be applicable to all sales for re-export unless evidence of exemption is supplied at the date of the order. Reagents, disposables and accessories prices are invoiced at the price in force at the date of delivery and in accordance to agreed prices with bioMérieux Sales department.
- Unless otherwise expressly agreed by the Parties, all shipment costs and expenses are prepaid by bioMérieux and re-invoiced to the Buyer, even for partial shipments.

5. INVOICING

- Products are invoiced on the shipping date, at the price in force at the date of order. Unless otherwise informed by the Buyer, bioMérieux reserves its right to establish a deferred invoicing of the Products.
- Delay by the Buyer in challenging the invoiced amount within a period of 30 days from the date of the invoice shall be deemed acceptance thereof by the Buyer.

6. PAYMENT

- Payment shall be made in cash in advance or, in case of bioMérieux's express prior consent, 30 days net following the invoicing date, by cheque or by bank transfer.
- A penalty is due by Buyer in case of late payment and is to be calculated from the day following the date of payment stated on the invoice. The amount of such penalty is calculated at the interest rate applied by the European Central Bank at its most recent refinancing operation with a ten point increase.
- Failure by the Buyer to fulfil its payment obligations, bioMérieux reserves the right:
 - to remove immediately all the payment facilities and special commercial conditions previously granted
 - to suspend or cancel, without notice and compensation, any current order
 - to require, in order to perform any subsequent delivery, cash payment before each shipment, or any other means of payment chosen by bioMérieux
 - to require immediate payment of the entire outstanding balance
- Payment should be made by bank transfer to the bank account indicated on the invoice.

HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Terms & Conditions

PAGE 1 OF 5

PAGE 2 OF 5

PAGE 3 OF 5

PAGE 4 OF 5

PAGE 5 OF 5

7. LEASE PURCHASE

In case of lease purchase, the funding agency shall assume the rights and obligations hereunder as Buyer, and shall guarantee the rights and obligations of the end-user.

8. OWNERSHIP RIGHTS AND PASSING OF RISK

- 8.1 All products remain bioMérieux's property until full payment thereof by Buyer. In the event of non-payment by Buyer, bioMérieux reserves the right to recover or to make recover possession of the unpaid goods. bioMérieux shall retain the part of the price already paid as a compensation for the use of the goods. The promise to pay shall not be considered, for the purposes of this provision, as payment.
- 8.2 Notwithstanding the above, the Buyer is liable for all damage occurring to the Products or damage caused by the Products. In accordance with the Incoterms as per Article 3, risk of loss or damage are transferred from bioMérieux to the Buyer after delivery of such Products to the Buyer. The Buyer shall therefore ensure that its insurance covers the Products delivered from this time.

9. TRACEABILITY

The technical specificities of products sold by bioMérieux require compliance with traceability rules. Consequently, bioMérieux reserves the right to sell its Products solely to professionals authorized to provide biological analysis results. The Buyer undertakes to comply with all traceability rules which apply to the Products ordered. bioMérieux will not be liable for non-compliance with traceability rules after delivery of its products to the site specified by Buyer.

10. CLAIMS

- 10.1 The Buyer shall have checked the incoming parcel and within three (3) days from the delivery report to the bioMérieux order department any missing or damaged products.
- 10.2 The Buyer shall check the incoming parcel directly on arrival and if any damage is visible from outside it must be written down on the delivery note at arrival.
- 10.3 Claims must be sent to bioMérieux Customer service as soon as possible upon notice. No product returns shall be accepted without bioMérieux's prior written consent. bioMérieux shall decide, at its sole judgement, if the Products shall be replaced.

11. WARRANTIES

- 11.1 Reagents:
bioMérieux's reagents are intended for in vitro use. The conformance of the reagents to the specifications indicated in the package insert is guaranteed until their expiry date.
- 11.2 Instrument:
bioMérieux warrants the Instrument to the Buyer against defects in material and workmanship and defects arising from failure to conform to the user's manual applicable on the date of installation, and as provided to the Buyer, for the period of one (1) year from the date of installation, when the certificate of installation or equivalent is signed.
bioMérieux agrees to correct or have corrected by a third party designated by bioMérieux, either by repair or at its election, by replacement, any such defect found on examination to have occurred, when the Instrument has been used under normal operating conditions and maintained in accordance with the instructions provided by bioMérieux such as, but not limited to, written instructions, package insert and user's manuals, during such one (1) year warranty period, provided bioMérieux is promptly notified in writing upon discovery of such defect by the Buyer.
Disposables and replacement items with a normal life expectancy of less than one (1) year such as, but not limited to, batteries, lamps and tubing are excluded from this warranty.
bioMérieux warranty is valid to the Buyer's original site of installation only. Any transhipment VOIDS this warranty unless bioMérieux facilitates the relocation of the Instrument.
- 11.3 Software:
bioMérieux warrants that the Software (other than the one embedded in the Instrument) will be in accordance with the descriptions and specifications in the documentation accompanying such Software for a period of ninety (90) days from the date of installation under normal use and conditions. During this warranty period, bioMérieux agrees to provide technical support of the Software by correcting any deficiency of the Software and, in no case, by redesigning the Software. Software developed by and the property of third parties is warranted to be in accordance with the descriptions and specifications in the documentation accompanying such Software by bioMérieux and by its software publishers. Any other warranty, such as covering performance or results linked with their use, or such as warranty for compatibility with the Buyer's hardware, is expressly excluded.

HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Terms & Conditions

PAGE 1 OF 5

PAGE 2 OF 5

PAGE 3 OF 5

PAGE 4 OF 5

PAGE 5 OF 5

11. WARRANTIES (CONTINUED)

11.3 Software (continued):

bioMérieux will evaluate warranty claims upon receipt of written notification. Firstly, bioMérieux will make a remote diagnosis and if possible, will correct the defect. Secondly, if needed, onsite evaluations will be initiated. Thirdly, if repair cannot be made on the Buyer's site, the Instrument shall be repaired in bioMérieux premises. Transportation costs will be at bioMérieux's expense. The return of the Instrument to bioMérieux can only be done after bioMérieux's written approval.

The Buyer will provide access to bioMérieux, in a timely manner, to any technical support, facilities, hardware, software or information in the Buyer's possession necessary for bioMérieux to complete such work. bioMérieux will correct such stated defects.

In any event, the Buyer is under an obligation to insure that the Instrument and/or Software is maintained in accordance with the instructions provided by bioMérieux such as, but not limited to, written instructions, package insert and user's manuals and in accordance with local standards during the warranty period of the concerned Products.

Preventive maintenance is not included in the warranty but rather a part of the imposed maintenance associated with the obligations of the Buyer.

Notwithstanding the foregoing, this warranty will not apply to the extent the Instrument and/or Software does not conform to and/or function in accordance with the applicable written technical documentation as a result of a defect arising from:

- any act or omission of the Buyer
- any corporation, joint stock company or any other entity or organization of any kind (other than bioMérieux or any person under the express direction of bioMérieux) making any revisions or modifications to the Instrument and/or Software after its provision to the Buyer
- the malfunction of any Buyer-supplied or Buyer-authorized third-party software or equipment
- Buyer operation of the Instrument and/or Software other than in accordance with applicable documentation or design, or on hardware not recommended, supplied or approved by bioMérieux
- the occurrence of any Force Majeure event such as, but not limited to, fire, electrical current variation, flood and defect in the installation environment of the Instrument and/or Software

THE WARRANTY OF bioMérieux SET FORTH ABOVE AND THE OBLIGATIONS AND LIABILITIES OF bioMérieux THEREUNDER ARE EXCLUSIVE AND IN LIEU OF ALL OTHER REMEDIES, WARRANTIES, GUARANTEES OR LIABILITIES, EXPRESSED OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO ANY PRODUCTS DELIVERED HEREUNDER bioMérieux WARRANTY SHALL ONLY APPLY TO PRODUCTS PURCHASED BY THE BUYER FROM bioMérieux.

This warranty shall not be extended or altered except by a written agreement signed by bioMérieux.

12. LIABILITY

12.1 bioMérieux accepts no liability for any damage incurred through misuse or incorrect storage of Products.

As such, Buyer shall ensure the respect of:

- the storage conditions written on the Products packaging or set forth in the package inserts

12.2 The Buyer or its representative shall take all measures to ensure that instruments, reagents and software are used in compliance with legal requirements and standards, including the installation of instruments and performance analysis. As an expert/a man skilled in the art, the Buyer shall use the instruments in strict compliance with the recommendations of bioMérieux and shall remain liable for any interpretation and any use of results provided. The Buyer guarantees bioMérieux against any third party action on the consequences of using non-compliant Products.

12.3 bioMérieux's liability for any and all losses or damages resulting from any cause whatsoever, including alleged negligence, shall in no event exceed the purchase price of the Product in respect of which the claim is made, or at the election of bioMérieux the repair or replacement of such Product plus associated freight charges. For the avoidance of doubt, bioMérieux's warranty expressly excludes any possible labour and transport costs involved in the repair of an Instrument. IN NO EVENT SHALL bioMérieux BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES WHETHER BUYER'S CLAIM IS IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE.

HOME

- Our Principles
- Sales Orders & Enquiries
- Customer Support
- mybiomerieux.com & QC Certificates
- History bioMérieux
- Quality & Environment
- bioMérieux Goes Green
- Terms & Conditions
- Legal
- How to use this Product List



Terms & Conditions

PAGE 1 OF 5

PAGE 2 OF 5

PAGE 3 OF 5

PAGE 4 OF 5

PAGE 5 OF 5

13. WASTE MATERIAL

Unless anything to the contrary whether expressly agreed by bioMérieux or due to a mandatory local rule, the Buyer is responsible for notifying bioMérieux so that the collection, treatment, recovery and environmental sound disposal for the Waste Electrical and Electronic Equipment can be organised, as provided for WEEE directive (2002/96/EC) and any local implementations thereof in the country in which the bioMérieux selling company is located.

The obligations of the Buyer (in particular in connection with decontamination) to which the Buyer will necessarily have to conform, are defined in the user's manual of the concerned system. The conditions of disposal of these systems, as well as the price rates related to these services will be communicated to the Buyer at his request. The Buyer as keeper of the equipment will be responsible for the integrity of the system and for its disposal at bioMérieux benefit. Every expense incurred or damage undergone caused by the Buyer's negligence and/or any breach of its abovementioned obligations will entitle bioMérieux to a compensation.

The Buyer expressly agrees that in case of resale of any Instrument supplied by bioMérieux and/or transfer of the Instrument in any way, within the framework of a loan or a rent for free or for a fee, or a donation or any other terms, on the national territory or any other territory, then bioMérieux will be released de facto, without any other condition, of its obligations of financing and organization of the elimination of the aforementioned waste material. bioMérieux undertakes to communicate any necessary information about the concerned Instruments. The Buyer will be completely responsible for the good documentary, technical and administrative management. The Buyer will guarantee it towards bioMérieux and will hold bioMérieux harmless for any damaging consequence, expenses, legal actions from third party including involved authorities in case of wrongful execution of this obligation.

14. SOFTWARE

- 14.1 Software included in the price list is subject to a non exclusive license agreement: software programs can neither be reproduced, modified nor transferred without bioMérieux's agreement.
- 14.2 Software sold by bioMérieux is free from any known viruses at the date of manufacture. The Buyer shall take all necessary precautions thereafter.
- 14.3 Moreover, bioMérieux assumes no responsibility for the software compatibility with the hardware of the Buyer.
- 14.4 On the delivery date, Software supplied by bioMérieux are virus free. bioMérieux interventions are done through tools protected by updated antivirus software. It is the Buyer's responsibility to set up a protection plan for all the data of its information system (computing resources of the Buyer) against future risks.
- 14.5 The Buyer has to make sure of the compatibility of his computer hardware with those proposed by bioMérieux. bioMérieux cannot be held responsible for any contamination by virus. In addition, bioMérieux reserves the right to invoice complementary services due to contamination existing in the Buyer's computer hardware.

15. PERSONAL DATA PROTECTION

Personal data of the Buyer (in particular its address and telephone details) collected are subject to processing's which have been declared to bioMérieux. The Buyer accepts that its personal data are used exclusively by bioMérieux and other entities of the bioMérieux group, in conformity with the provisions of The European Data Protection Directive (95/46/EC). The Buyer has the right to access, modify and rectify such data. For this purpose, the Buyer can contact its account manager.

Within the framework of warranty operations of the systems sold by bioMérieux or the quality control of these systems, the Buyer grants access to bioMérieux to the patients personal data for a required period of time in order to perform the abovementioned operation. According to the patient confidentiality regulations applicable in the country in which the bioMérieux selling company is located and The European Data Protection Directive above and after reception of the written instruction of the Buyer, bioMérieux undertakes to set up, from the beginning of the operation, appropriate devices and safety procedures, so as to guarantee strictly the safety; the integrity and the confidentiality of the data pertaining to the Buyer.



Terms & Conditions

16. ASSIGNMENT

- 16.1 The Buyer shall not assign or transfer the agreement or any of its rights and obligations arising hereunder, or any order covered hereunder, to any other party, including to an Affiliate, without prior written consent of bioMérieux. Any change in the ownership or control of the Buyer shall be deemed an assignment, which requires written consent from bioMérieux.
- 16.2 The Buyer acknowledges that bioMérieux can transfer the order or be subject to a change of control for the benefit of a third party.

17. APPLICABLE LAW

- 17.1 These conditions of sales are governed by and construed in accordance with the laws of Sweden, Norway, Denmark and Finland.
- 17.2 Any dispute shall be settled by the Courts of the country in which the bioMérieux selling company is located, which shall have exclusive jurisdiction, even in case of multiple defendants or claim for contribution from a third party.

How to use this Product List

HOW TO USE THIS PRODUCT LIST

- Overview
- Layout & Navigation
- Finding Products

Overview

In support of the bioMérieux Goes Green initiative, we have produced this interactive PDF Product List.

The following notes have been produced with the intention of assisting you to get the best results out of this PDF Product List.

As a progressive company, we consider all feedback valuable. If you have any comments or are having problems, please contact us. Our telephone lines are open Monday to Friday (except Bank Holidays) according to the info of office addresses below:



Denmark

Order Team

Telephone: 70 10 84 02

Fax: 70 10 84 01

order.dk@biomerieux.com

General Enquiries

Telephone: 70 10 84 00

Fax: 70 10 84 01

contact.dk@biomerieux.com

Opening Hours

Monday to Thursday

from 8.30am to 4.00pm

Fridays 8.30am to 3.00pm

bioMérieux Denmark ApS

Smedeholm 13C

2730 Herlev



Finland

Order Team

Telephone: 09 8545 6040

Fax: 09 8545 6045

order.fi@biomerieux.com

General Enquiries

Telephone: 09 8545 600

Fax: 09 8545 6045

contact.fi@biomerieux.com

Opening Hours

Monday to Friday

from 8.30am to 4.00pm

bioMérieux Suomi Oy

Konalantie 47C

00390 Helsinki



Norway

Order Team

Telephone: 2337 5552

Fax: 2337 5551

order.no@biomerieux.com

General Enquiries

Telephone: 2337 5550

Fax: 2337 5551

contact.no@biomerieux.com

Opening Hours

Monday to Friday

from 8.30am to 4.00pm

bioMérieux Norge AS

Ulvenveien 75A, 3 etg.

0581 Oslo



Sweden

Order Team

Telephone: 031 68 84 96

Fax: 031 68 48 48

order.se@biomerieux.com

General Enquiries

Telephone: 031 68 84 90

Fax: 031 68 48 48

contact.se@biomerieux.com

Opening Hours

Monday to Friday

from 8.30am to 4.00pm

bioMérieux Sweden AB

Hantverkarsvägen 15

43633 Askim



HOW TO USE THIS PRODUCT LIST

- Overview
- Layout & Navigation
- Finding Products

Layout & Navigation

To help you navigate this Product List, every page is separated into these three areas:



1. Head Tabs

There are ten Head Tabs across the top of every page.

When the cursor is placed over these Head Tabs, the arrow changes to a hand . By clicking the mouse when the hand icon is over the required Head Tab, you navigate to this section, and the selected tab changes to White to confirm your choice. All of the unselected tabs are shown in Orange.

2. Tab Menus

For every Head Tab, there is a specific Tab Menu shown in the left hand pane.

This Tab Menu presents the complete list of contents within the selected Head Tab.

If you select a title within the Tab Menu, you advance to the specified page within the Product List. Further Tab Menu options may then be shown if additional pages are available.

As with Head Tabs, your selection is highlighted to show your current position.

3. Content

The content for the required Head Tab and Tab Menu is shown in the largest portion of the page.

The footer bar contains contact details, and will automatically launch e-mail and our web site.



HOW TO USE THIS PRODUCT LIST

Overview

Layout & Navigation

Finding Products

Finding Products

There are two ways to locate products within this PDF Product List:

1. Using the Head Tabs & Tab Menus

You can locate a product by selecting the relevant Head Tab and Tab Menu and also by browsing forwards or backwards using navigation controls.

For instance to find Etest® products, as with previous years these are located in the 'Identification & Susceptibility Testing' Tab. Click the 'Identification & Susceptibility Testing' Head Tab to select and then click on 'Etest' in the Tab Menu.

2. Using the Adobe Reader Search Tool

You can use the Adobe Reader built in search tool to find products by name or reference.

The search box is usually located in the toolbar above. To open the search box if it is not currently available, choose the 'Edit' menu and then select 'Find' or press:

Ctrl + **F** if you're using a Windows PC

  +  if you're using an Apple Mac

Type the word, words or part of the word you want to search for in to the search box.

The way in which you cycle through each result depends on the software being used. In most software the first result will be highlighted on the page, to advance to the next result click the arrow next to the search box.

This search facility is good for finding product references, as well as every incidence of say, Etest. This makes it the most flexible search tool available within this PDF Product List.

Note, Adobe Reader searches from the page currently viewed through to the last page, then continues the search from the front cover until it reaches the page you began on.



BLOOD CULTURE / MYCOBACTERIAL CULTURE

Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options

Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles

Feature – Myla™

Feature – Next Generation of BacT/ALERT

Blood Culture

- BacT/ALERT Computer Systems
- BacT/ALERT Culture Media
- BacT/ALERT Consumables & Accessories

Mycobacterial Culture

- BacT/ALERT Computer Systems
- BacT/ALERT Media
- BacT/ALERT Consumables & Accessories

PREVI™ Fluo TB

- PREVI Fluo TB Reagents



BacT/ALERT® 3D Culture System

The BacT/ALERT 3D provides proven recovery of a wide range of clinically significant organisms, including bacteria, yeasts and mycobacteria. This single system can be used to culture from blood and all normally sterile body fluids using our extensive range of culture media. The same system can be used to culture mycobacteria from decontaminated clinical samples using our Mycobacterial Process bottles or blood using the Mycobacteria Blood bottles.

BacT/ALERT 3D System Options

The BacT/ALERT 3D system is available in a number of different modules that can be configured to provide a perfect fit no matter what your workload. Larger systems comprise a Controller module with a touch screen for text free control over the basic functions of the system. This unit can be linked to a maximum of six 240 Incubator modules, which would provide sufficient capacity for an annual workload of 84,000 bottles from a single PC.

For the laboratory requiring extra space saving the BacT/ALERT 3D combination module combines the computer function and touch screen of a controller module with two incubator drawers. A further three Incubator modules can be linked to this module. Finally the smallest BacT/ALERT 3D 60 can handle an annual workload of up to 3,500 bottles per annum.

BacT/ALERT 3D™



Figure 1:
String Configuration



Figure 2:
Stacked Configuration

Configuration

The Incubator modules can be configured in either a Stacked or String formation. The String configurations are supplied on purpose built bench carts that have cupboards underneath which are left empty for the department to utilise. This configuration is the most user friendly as all the cells are situated conveniently for the loading and unloading of bottles. The Stacked configuration has the smallest footprint of any blood culture system currently available and is ideal where floor space is at a premium.

Simplicity

Routine use of the BacT/ALERT 3D system is accessed through the touch screen on the Controller/Combination Module. The software is intuitive, all tasks are easy to perform and the system is designed to minimise the possibility of users making any mistakes. This ensures that staff can rapidly be trained in day-to-day use.

BLOOD CULTURE / MYCOBACTERIAL CULTURE

Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options

Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles

Feature – Myla™

Feature – Next Generation of BacT/ALERT

Blood Culture

- BacT/ALERT Computer Systems
- BacT/ALERT Culture Media
- BacT/ALERT Consumables & Accessories

Mycobacterial Culture

- BacT/ALERT Computer Systems
- BacT/ALERT Media
- BacT/ALERT Consumables & Accessories

PREVI™ Fluo TB

- PREVI Fluo TB Reagents



Polycarbonate Blood Culture Bottles

Our users also benefit from bioMérieux’s revolutionary, virtually unbreakable, polycarbonate blood culture bottles.

These offer enhanced safety as glass is removed from the workplace. The bottles are also significantly lighter than glass ensuring that waste disposal costs are minimal. The bottles are also validated for transport through your pneumatic air system without any need for additional padding or transport containers. This greatly reduces the time taken for these important samples to reach the laboratory.

Colorimetric Technology

The BacT/ALERT® 3D blood culture system uses the most up to date colorimetric technology. This offers improved time to detection over alternative systems in two ways.

First, they provide visual notification of positive bottles as they arrive in the department. This means that bottles may be sub-cultured before they are loaded onto the system thereby improving the turn around time to your clinicians.

Also delayed entry bottles will flag positive two hours after loading. There is no need to wait until these positive bottles enter the death phase of the bacterial life cycle before flagging.

FAN Plus Bottles

During 2012 the BacT/ALERT Bottles will also see the world wide launch of the new range of FAN media FA Plus, FN Plus and PF Plus. This media will incorporate Adsorbent Polymer Beads (APB) in place of the activated charcoal to remove antibiotics and naturally occurring bacterial inhibitors, which may be present in the patient’s blood sample.

The key benefits are:

- An increase in overall organism recovery
- Improved ability to neutralise a broader spectrum of antimicrobials
- Faster Time to Detection compared to the current bottle types
- Smaller particles prevent any issues that may be seen with particles blocking sub-culture devices
- Clear, easy to read Gram stains

BacT/ALERT 3D™



BLOOD CULTURE / MYCOBACTERIAL CULTURE

- Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options
- Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles
- Feature – Myla™
- Feature – Next Generation of BacT/ALERT Blood Culture

BacT/ALERT Computer Systems

BacT/ALERT Culture Media

BacT/ALERT Consumables & Accessories
- Mycobacterial Culture

BacT/ALERT Computer Systems

BacT/ALERT Media

BacT/ALERT Consumables & Accessories
- PREVI™ Fluo TB

PREVI Fluo TB Reagents



Myla™

bioMérieux are proud to introduce you to Myla, an intelligent microbiology middleware solution that is designed by microbiologists for microbiologists.

Myla is designed to be accessible from any device linked to the network. This means that it is possible to have blood culture systems sited in satellite hospitals, removing the need to transport bottles between sites and reducing the time taken between blood culture collection and incubation. Staff would be notified to the presence of a positive bottle via a variety of alerting mechanisms. These include standard on screen alerts, optional notification lights, and the ability to send text or e-mail messages.

Myla would also alert users to any disruptions in the interface or when error codes are generated on an instrument (e.g. a temperature alert due to a drawer not completely closed on a BacT/ALERT® instrument).

With the launch of a new version of Myla software, the functionality with the BacT/ALERT systems was improved. Amongst other features these include a range of workflow indicators:

- Average time between blood culture order and when the bottle is loaded
- Average time to positivity
- Average time between positive detection and bottle unloading for follow up

See [page 29](#) for more details of Myla.



BLOOD CULTURE / MYCOBACTERIAL CULTURE

- Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options
- Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles
- Feature – Myla™
- Feature – Next Generation of BacT/ALERT

- Blood Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Culture Media
 - BacT/ALERT Consumables & Accessories
- Mycobacterial Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Media
 - BacT/ALERT Consumables & Accessories
- PREVI™ Fluo TB
 - PREVI Fluo TB Reagents




Next Generation of BacT/ALERT®



As part of bioMérieux’s ongoing Research and Development we continue to lead the way with blood culture technology. The next generation of BacT/ALERT will be launched in 2013.

The new system will provide many benefits over the existing blood culture systems including automated loading and unloading of blood culture bottles and an even more compact footprint.



BLOOD CULTURE / MYCOBACTERIAL CULTURE

Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options

Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles

Feature – Myla™

Feature – Next Generation of BacT/ALERT

Blood Culture

BacT/ALERT Computer Systems
BacT/ALERT Culture Media
BacT/ALERT Consumables & Accessories

Mycobacterial Culture

BacT/ALERT Computer Systems
BacT/ALERT Media
BacT/ALERT Consumables & Accessories

PREVI™ Fluo TB

PREVI Fluo TB Reagents



Reference	Description	Presentation
-----------	-------------	--------------

Blood Culture



BacT/ALERT® Computer Systems

413641	OBSERVA® Data Management System	One
411401 NEW	Myla™ Middleware	One

BacT/ALERT Culture Media

259789	BacT/ALERT SA Standard Aerobic Culture bottles	100 Bottles
259790	BacT/ALERT SN Standard Anaerobic Culture bottles	100 Bottles
259791	BacT/ALERT FA FAN Aerobic Culture bottles	100 Bottles
410851 NEW	BacT/ALERT FA Plus Aerobic Culture bottles	100 Bottles
259793	BacT/ALERT FN FAN Anaerobic Culture bottles	100 Bottles
410852 NEW	BacT/ALERT FN Plus Anaerobic Culture bottles	100 Bottles
259794	BacT/ALERT PF Paediatric FAN Culture bottles	100 Bottles
410853 NEW	BacT/ALERT PF Plus Paediatric Culture bottles	100 Bottles

BacT/ALERT Consumables & Accessories

233766	BacT/ALERT Airway Sub vent units	100
279012	BacT/ALERT Adaptor caps (blue)	120
279013	BacT/ALERT Adaptor inserts (blue)	60
413251 NEW	Blood collector set (21 gauge needle, 30cm tubing)	240 sets
413201 NEW	Blood collector set (23 gauge needle, 30cm tubing)	240 sets
413200 NEW	Blood collector holder	180 pieces
940-0046-01 NEW	Bottle Transport Rack	48 Wells
259662	Generic Barcode Labels	2000
206543	Reflectance Standards	1 x 4 Standard



BLOOD CULTURE / MYCOBACTERIAL CULTURE

- Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options
- Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles
- Feature – Myla™
- Feature – Next Generation of BacT/ALERT

- Blood Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Culture Media
 - BacT/ALERT Consumables & Accessories

- Mycobacterial Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Media
 - BacT/ALERT Consumables & Accessories

- PREVI™ Fluo TB
 - PREVI Fluo TB Reagents



Reference	Description	Presentation
-----------	-------------	--------------

Mycobacterial Culture



BacT/ALERT® Computer Systems

413641	OBSERVA® Data Management System	One
411401 NEW	Myla™ Middleware	One

BacT/ALERT Media

259797	BacT/ALERT MP Process Bottle (Polycarbonate)	100 Bottles
259760	MAS with Vancomycin. Antibiotic supplement	5 x 20
251011	BacT/ALERT MB Blood Bottles	25 Bottles
259877	BacT/ALERT Blood Enrichment Fluid	5 x 5

BacT/ALERT Consumables & Accessories

210361	BacT/ALERT Adaptor caps (white)	120
210362	BacT/ALERT Adaptor inserts (white)	60
259787	BacT/ALERT Bottle Reseals	100
413251 NEW	Blood collector set (21 gauge needle, 30cm tubing)	240 sets
413201 NEW	Blood collector set (23 gauge needle, 30cm tubing)	240 sets
413200 NEW	Blood collector holder	180 pieces

BLOOD CULTURE / MYCOBACTERIAL CULTURE


- Feature – BacT/ALERT® 3D Culture System & BacT/ALERT 3D System Options
- Feature – Polycarbonate Blood Culture Bottles, Colorimetric Technology & FAN Plus Bottles
- Feature – Myla™
- Feature – Next Generation of BacT/ALERT

- Blood Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Culture Media
 - BacT/ALERT Consumables & Accessories
- Mycobacterial Culture
 - BacT/ALERT Computer Systems
 - BacT/ALERT Media
 - BacT/ALERT Consumables & Accessories

PREVI™ Fluo TB

PREVI Fluo TB Reagents

Reference	Description	Presentation
PREVI™ Fluo TB PREVI Fluo TB Reagents 411010 NEW PREVI Fluo TB kit Up to 200 Tests 411011 NEW Fixative 1000ml 411015 NEW Auramine 500ml 411014 NEW Carbolic Solution 500ml 411012 NEW Degommier Discolouring 1000ml 411009 NEW Thiazine Red 500ml 411013 NEW Carbolic Solution for Thiazine Red 500ml		



FULL MICROBIOLOGY LAB AUTOMATION

Feature – Full Microbiology Lab Automation & Intelligent Automation

Software
Myla™ Middleware

PREANALYTICAL

Urinanalysis
Chromogenic Media
Direct Identification of E. coli, Proteus & Enterococcus
Screening of Multidrug-Resistant Organisms
Urinary Bacterial Count

PREVI™ Isola
PREVI Isola Plate System
PREVI Isola Consumables & Accessories

PREVI Color
PREVI Color Gram System
PREVI Color Gram Reagents
PREVI Color Gram Consumables & Accessories



FMLA Full Microbiology Lab Automation


Full Microbiology Lab Automation focuses on a faster time to result. FMLA brings together bioMérieux’s advances in the automation of operations in Microbiology to deliver greater efficiency, productivity and quality, while retaining the ability to tailor solutions to individual laboratory requirements and enabling the laboratory to effectively respond to increasing workloads.

For more details of FMLA, [click here](#).

Intelligent Automation

Our automated solutions offer maximum modularity without compromising result quality and are designed to link to our bespoke intelligent middleware solution, Myla™ to further improve connectivity, data and real-time workflow management. Myla has been designed by microbiologists for microbiologists and plays a central role in the FMLA approach.

To find out how this new lab partner helps you simplify your life and recaptures your time, [click here](#).



FULL MICROBIOLOGY LAB AUTOMATION

Feature – Full Microbiology Lab Automation & Intelligent Automation

Software
Myla™ Middleware

PREANALYTICAL

Urinanalysis
Chromogenic Media
Direct Identification of E. coli, Proteus & Enterococcus
Screening of Multidrug-Resistant Organisms
Urinary Bacterial Count

PREVI™ Isola
PREVI Isola Plate System
PREVI Isola Consumables & Accessories

PREVI Color
PREVI Color Gram System
PREVI Color Gram Reagents
PREVI Color Gram Consumables & Accessories

Reference	Description	Presentation																														
Software Myla™ Middleware <table> <tr> <td>411401 NEW</td><td>Myla Middleware</td><td>One</td></tr> <tr> <td>411613 NEW</td><td>Myla VITEK® 2 Connection Kit</td><td>1 Kit</td></tr> <tr> <td>411615 NEW</td><td>Myla VITEK MS Connection Kit</td><td>1 Kit</td></tr> <tr> <td>411616 NEW</td><td>Myla BacT/ALERT® Connection Kit</td><td>Per System</td></tr> <tr> <td>411405 NEW</td><td>Myla Monitor LCD 17LN HP L1710</td><td></td></tr> <tr> <td>410494 NEW</td><td>Signal Tower Kit</td><td>1 Kit</td></tr> <tr> <td>410493 NEW</td><td>Serial 485 to Ethernet Adapter</td><td>1 Adapter</td></tr> <tr> <td>410491 NEW</td><td>Cellular Modem Kit</td><td>1 Kit</td></tr> <tr> <td>410151 NEW</td><td>Myla Alert / Notification Option</td><td></td></tr> <tr> <td>410152 NEW</td><td>Myla Multiuser Licence</td><td>One Licence</td></tr> </table>			411401 NEW	Myla Middleware	One	411613 NEW	Myla VITEK® 2 Connection Kit	1 Kit	411615 NEW	Myla VITEK MS Connection Kit	1 Kit	411616 NEW	Myla BacT/ALERT® Connection Kit	Per System	411405 NEW	Myla Monitor LCD 17LN HP L1710		410494 NEW	Signal Tower Kit	1 Kit	410493 NEW	Serial 485 to Ethernet Adapter	1 Adapter	410491 NEW	Cellular Modem Kit	1 Kit	410151 NEW	Myla Alert / Notification Option		410152 NEW	Myla Multiuser Licence	One Licence
411401 NEW	Myla Middleware	One																														
411613 NEW	Myla VITEK® 2 Connection Kit	1 Kit																														
411615 NEW	Myla VITEK MS Connection Kit	1 Kit																														
411616 NEW	Myla BacT/ALERT® Connection Kit	Per System																														
411405 NEW	Myla Monitor LCD 17LN HP L1710																															
410494 NEW	Signal Tower Kit	1 Kit																														
410493 NEW	Serial 485 to Ethernet Adapter	1 Adapter																														
410491 NEW	Cellular Modem Kit	1 Kit																														
410151 NEW	Myla Alert / Notification Option																															
410152 NEW	Myla Multiuser Licence	One Licence																														



FULL MICROBIOLOGY LAB AUTOMATION

Feature – Full Microbiology Lab Automation & Intelligent Automation

Software
Myla™ Middleware

- PREANALYTICAL**

Urinalysis
Chromogenic Media
Direct Identification of E. coli, Proteus & Enterococcus
Screening of Multidrug-Resistant Organisms
Urinary Bacterial Count

PREVI™ Isola
PREVI Isola Plate System
PREVI Isola Consumables & Accessories

PREVI Color
PREVI Color Gram System
PREVI Color Gram Reagents
PREVI Color Gram Consumables & Accessories

Reference	Description	Presentation
-----------	-------------	--------------

Urinalysis



Chromogenic Media

Direct identification of E. coli, Proteus & Enterococcus

43821 NEW	chromID™ CPS	20 Plates
43829 NEW	chromID CPS	100 Plates
43541	chromID CPS	20 Plates
43549	chromID CPS	100 Plates
43473	chromID CPS/Columbia CNA + 5% sheep blood	20 Plates
411617 NEW	chromID CPS/Columbia CNA + 5% sheep blood	20 Plates
56541	ID indole TDA reagents	2 x 200 Tests

Screening of Multidrug-Resistant Organisms

43481	chromID ESBL	20 Plates
-------	--------------	-----------

Urinary Bacterial Count

56507	URILINE™	10 Dipslides
56508	URILINE	100 Dipslides



FULL MICROBIOLOGY LAB AUTOMATION

Feature – Full Microbiology Lab Automation & Intelligent Automation

Software
Myla™ Middleware

- PREANALYTICAL**
- Urinalysis**
Chromogenic Media
Direct Identification of E. coli, Proteus & Enterococcus
Screening of Multidrug-Resistant Organisms
Urinary Bacterial Count
- PREVI™ Isola**
PREVI Isola Plate System
PREVI Isola Consumables & Accessories
- PREVI Color**
PREVI Color Gram System
PREVI Color Gram Reagents
PREVI Color Gram Consumables & Accessories



Reference	Description	Presentation
-----------	-------------	--------------

PREVI™ Isola



PREVI Isola Plate System

29500	PREVI Isola	One
-------	-------------	-----

PREVI Isola Consumables & Accessories

29508	PREVI Isola tips	50 x 96 tips
29509	PREVI Isola Applicators	8 Cartridges
29800	PREVI Isola Rack 1	One
29705	PREVI Isola Rack 2	One
410012	PREVI Isola Rack 3 UF	One
29706	PREVI Isola Rack 4	One
29707	PREVI Isola Rack 5	One
29704	PREVI Isola Rack 6 Swab	One
413958	PREVI Isola Rack 10 Enrichment broth	One
29718	PREVI Isola Waste Bin	10 Waste Bins
411143	PREVI Isola Labels	1500 Labels
29712 NEW	PREVI Isola Input Cassette	One
29713 NEW	PREVI Isola Output Cassette	One



FULL MICROBIOLOGY LAB AUTOMATION

Feature – Full Microbiology Lab Automation & Intelligent Automation

Software
Myla™ Middleware

PREANALYTICAL

Urinanalysis
Chromogenic Media
Direct Identification of E. coli, Proteus & Enterococcus
Screening of Multidrug-Resistant Organisms
Urinary Bacterial Count

PREVI™ Isola
PREVI Isola Plate System
PREVI Isola Consumables & Accessories

PREVI Color
PREVI Color Gram System
PREVI Color Gram Reagents
PREVI Color Gram Consumables & Accessories



Reference	Description	Presentation
-----------	-------------	--------------

PREVI™ Color



PREVI Color Gram System

29502	PREVI Color Gram 12	One
29503	PREVI Color Gram 30	One
29506	PREVI Color Cytocentrifuge Rotor	One

PREVI Color Gram Reagents

29521	Acetone Fuchsin-A	500ml
29519	Acetone Safranin-A	500ml
29524	Crystal Violet-C	500ml
29522	Fuchsin-A	500ml
29523	Iodine-B	500ml
29520	Safranin-A	500ml
29525	Nozzle Cleaning Solution	5 Litres

PREVI Color Gram Consumables & Accessories

29559	D&E 500ml Empty Bottle	One
29558	D&E 5L Empty Bottle	One



QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains



Qualitative QC Organisms for the
Quality Control of API® Identification Strips
& VITEK® 2 Cards



LyfoCults® Plus are high quality, patented, disposable devices containing lyophilized microorganisms derived from ATCC® cultures. The unique, glass-free, all-in-one device stores and hydrates the freeze-dried microorganisms and has a built-in inoculation wand to streak onto culture media.

LyfoCults Plus QC Range

ATCC microorganisms packaged in convenient sets for the Quality Control of bioMérieux test kits:

- VITEK 2 Identification Cards
- VITEK 2 Susceptibility Cards
- API Identification Strips
- ATB™ Susceptibility Strips
- Individual QC Strains

Requirements vary between sites, countries and regulatory authorities but the Clinical and Laboratory Standards Institute (CLSI) recommends that identification products be QC tested once per lot and once per shipment.

 Click play to view LyfoCults Plus video on YouTube

Streamlined QC Sets

- Recommended for the routine QC testing of API Identification Strips and VITEK 2 Identification Cards
- Maximum of two organisms chosen to test the least robust and therefore most vulnerable to the effects of shipping, substrate in the Card or Strip
- One product reference per set for ease of ordering



Comprehensive QC Sets

- Recommended for the validation of a new VITEK 2 system – Performance Qualification
- Validation of new API Identification Strips or new VITEK 2 Identification Cards
- Organisms are chosen to reliably test each well or cupule of the Identification Card or Strip for a positive and negative reaction
- One product reference per set for ease of ordering

End User Agreement (EUA)

Please Read before ordering from the LyfoCults Plus range.

- All users of ATCC-derived products must agree to the terms of use by signing the EUA prior to receiving products from bioMérieux
- The ATCC EUA is a contract that governs the transfer of ATCC-licensed products between two organisations, when the recipient intends to use the licensed products for QC testing
- EUA protects the rights of all parties to maintain a chain-of-custody from ATCC, to bioMérieux and to the End User, and to assure that the QC microorganisms are traceable to their source
- It also assures that the product will be used by qualified personnel and handled according to governing laboratory protocols

For the LyfoCults Plus End User Agreement, [click here](#).



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains

Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for VITEK® 2 Identification Cards (page 1 of 5)

VITEK 2 Gram Negative QC Sets

410586	LyfoCults Plus VITEK 2 GN Streamlined QC set Kit contains: <i>Enterobacter hormaechei</i> <i>Stenotrophomonas maltophilia</i>	ATCC® 700323™ ATCC 17666™	4 Vials 2 Vials 2 Vials
410601	LyfoCults Plus VITEK 2 GN Comprehensive QC set Kit contains: <i>Acinetobacter baumannii</i> <i>Enterobacter hormaechei</i> <i>Elizabethkingia meningoseptica</i> <i>Klebsiella oxytoca</i> <i>Ochrobactrum anthropi</i> <i>Proteus vulgaris</i> <i>Pseudomonas aeruginosa</i> <i>Pseudomonas aeruginosa</i> <i>Shigella sonnei</i> <i>Stenotrophomonas maltophilia</i>	ATCC BAA-747™ ATCC 700323™ ATCC 13253™ ATCC 700324™ ATCC BAA-749™ ATCC 6380™ ATCC 9721™ ATCC BAA-1744™ ATCC 25931™ ATCC 17666™	20 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.



QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for VITEK® 2 Identification Cards (page 2 of 5)

VITEK 2 Gram Positive QC Sets

410616	LyfoCults Plus VITEK 2 GP Streamlined QC set Kit contains: <i>Enterococcus casseliflavus</i> <i>Staphylococcus saprophyticus</i>	ATCC® 700327™ ATCC BAA-750™	4 Vials 2 Vials 2 Vials
410584	LyfoCults Plus VITEK 2 GP Comprehensive QC set Kit contains: <i>Enterococcus casseliflavus</i> <i>Enterococcus saccharolyticus</i> <i>Kocuria kristinae</i> <i>Listeria monocytogenes</i> <i>Staphylococcus saprophyticus</i> <i>Staphylococcus sciuri ssp. sciuri</i> <i>Streptococcus equi ssp. zooepidemicus</i> <i>Streptococcus pneumoniae</i> <i>Streptococcus thermophilus</i>	ATCC 700327™ ATCC 43076™ ATCC BAA-752™ ATCC BAA-751™ ATCC BAA-750™ ATCC 29061™ ATCC 43079™ ATCC 49619™ ATCC 19258™	18 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
- VITEK 2 Gram Negative QC Sets
- VITEK 2 Gram Positive QC Sets
- VITEK 2 YST QC Sets
- VITEK 2 NH QC Sets
- VITEK 2 ANC QC Sets
- QC for VITEK 2 Susceptibility Cards
- VITEK 2 AST QC Sets
- QC for API® Identification Strips
- API 20 Strep QC Sets
- API 20 A QC Set
- API 20 C AUX QC Sets
- API 20 E QC Sets
- API ID 32 GN QC Set
- API Coryne QC Sets
- API 20 NE QC Sets
- API NH QC Sets
- RapiD 20 E™ QC Set
- rapid ID 32 E QC Set
- API Staph QC Sets
- QC for ATB™ Susceptibility Strips
- ATB Gram Negative QC Set
- ATB Gram Positive QC Set
- ATB ANA QC Set
- ATB HAEMO QC Set
- ATB FUNGUS QC Set
- Individual QC Strains

Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus QC for VITEK® 2 Identification Cards (page 4 of 5) VITEK 2 NH QC Sets			
410628	LyfoCults Plus VITEK 2 NH Streamlined QC set Kit contains: <i>Eikenella corrodens</i>	ATCC® BAA-1152™	2 Vials
410611	LyfoCults Plus VITEK 2 NH Comprehensive QC set Kit contains: <i>Aggregatibacter aphrophilus</i> <i>Eikenella corrodens</i> <i>Enterobacter aerogenes</i> <i>Haemophilus influenzae</i> <i>Neisseria gonorrhoeae</i> <i>Neisseria lactamica</i> <i>Oligella urethralis</i> <i>Paenibacillus polymyxa</i> <i>Staphylococcus epidermidis</i>	ATCC 33389™ ATCC BAA-1152™ ATCC 13048™ ATCC 9007™ ATCC 19424™ ATCC 23970™ ATCC 17960™ ATCC 7070™ ATCC 12228™	2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.



QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
 - VITEK 2 Gram Negative QC Sets
 - VITEK 2 Gram Positive QC Sets
 - VITEK 2 YST QC Sets
 - VITEK 2 NH QC Sets
 - VITEK 2 ANC QC Sets

- QC for VITEK 2 Susceptibility Cards
 - VITEK 2 AST QC Sets

- QC for API® Identification Strips
 - API 20 Strep QC Sets
 - API 20 A QC Set
 - API 20 C AUX QC Sets
 - API 20 E QC Sets
 - API ID 32 GN QC Set
 - API Coryne QC Sets
 - API 20 NE QC Sets
 - API NH QC Sets
 - RapiD 20 E™ QC Set
 - rapid ID 32 E QC Set
 - API Staph QC Sets

- QC for ATB™ Susceptibility Strips
 - ATB Gram Negative QC Set
 - ATB Gram Positive QC Set
 - ATB ANA QC Set
 - ATB HAEMO QC Set
 - ATB FUNGUS QC Set

Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus QC for VITEK® 2 Identification Cards (page 5 of 5) VITEK 2 ANC QC Sets LyfoCults™ +++ PLUS			
410604	LyfoCults Plus VITEK 2 ANC Streamlined QC set Kit contains: <i>Bacteroides ovatus</i> <i>Clostridium septicum</i>	ATCC® BAA-1296™ ATCC 12464™	4 Vials 2 Vials 2 Vials
410605	LyfoCults Plus VITEK 2 ANC Comprehensive QC set Kit contains: <i>Bacteroides ovatus</i> <i>Bacteroides vulgatus</i> <i>Clostridium perfringens</i> <i>Clostridium septicum</i> <i>Clostridium sordellii</i> <i>Corynebacterium striatum</i> <i>Parabacteroides distasonis</i>	ATCC BAA-1296™ ATCC 8482™ ATCC 13124™ ATCC 12464™ ATCC 9714™ ATCC BAA-1293™ ATCC BAA-1295™	14 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

- VITEK 2 Gram Negative QC Sets
- VITEK 2 Gram Positive QC Sets
- VITEK 2 YST QC Sets
- VITEK 2 NH QC Sets
- VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

- API 20 Strep QC Sets
- API 20 A QC Set
- API 20 C AUX QC Sets
- API 20 E QC Sets
- API ID 32 GN QC Set
- API Coryne QC Sets
- API 20 NE QC Sets
- API NH QC Sets
- RapiD 20 E™ QC Set
- rapid ID 32 E QC Set
- API Staph QC Sets

QC for ATB™ Susceptibility Strips

- ATB Gram Negative QC Set
- ATB Gram Positive QC Set
- ATB ANA QC Set
- ATB HAEMO QC Set
- ATB FUNGUS QC Set

Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for VITEK® 2 Susceptibility Cards

VITEK 2 AST QC Sets

410609	LyfoCults Plus VITEK 2 AST-GN QC set Kit contains: <i>Escherichia coli</i> <i>Escherichia coli</i> <i>Klebsiella pneumoniae</i> <i>Pseudomonas aeruginosa</i>	ATCC® 25922™ ATCC 35218™ ATCC 700603™ ATCC 27853™	8 Vials 2 Vials 2 Vials 2 Vials 2 Vials
410610	LyfoCults Plus VITEK 2 AST-GP QC set Kit contains: <i>Enterococcus faecalis</i> <i>Enterococcus faecalis</i> <i>Escherichia coli</i> <i>Staphylococcus aureus</i> <i>Staphylococcus aureus</i> <i>Staphylococcus aureus</i> <i>Staphylococcus aureus</i> <i>Streptococcus pneumoniae</i>	ATCC 29212™ ATCC 51299™ ATCC 35218™ ATCC 29213™ ATCC BAA-976™ ATCC BAA-977™ ATCC BAA-1026™ ATCC 49619™	16 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials 2 Vials
410613	LyfoCults Plus VITEK 2 AST-YS QC set Kit contains: <i>Candida krusei</i> <i>Candida parapsilosis</i>	ATCC 6258™ ATCC 22019™	4 Vials 2 Vials 2 Vials



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.



QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for API® Identification Strips (page 1 of 4)

API 20 Strep QC Sets

411002	LyfoCults Plus API 20 Strep Streamlined QC set Kit contains: <i>Streptococcus equi subsp. zooepidemicus</i>	ATCC® 700400™	2 Vials
410590	LyfoCults Plus API 20 Strep Comprehensive QC set Kit contains: <i>Streptococcus equi subsp. zooepidemicus</i> <i>Streptococcus uberis</i>	ATCC 700400™ ATCC 700407™	4 Vials 2 Vials 2 Vials

API 20 A QC Set

410614	LyfoCults Plus API 20 A QC set Kit contains: <i>Bacteroides ovatus</i> <i>Clostridium perfringens</i> <i>Clostridium sordellii</i>	ATCC 8483™ ATCC 13124™ ATCC 9714™	6 Vials 2 Vials 2 Vials 2 Vials
--------	--	---	--

API 20 C AUX QC Sets

411004	LyfoCults Plus API 20 C AUX Streamlined QC set Kit contains: <i>Candida glabrata</i> <i>Cryptococcus laurentii</i>	ATCC 15126™ ATCC 18803™	4 Vials 2 Vials 2 Vials
410589	LyfoCults Plus API 20 C AUX Comprehensive QC set Kit contains: <i>Candida glabrata</i> <i>Candida guilliermondii</i> <i>Cryptococcus laurentii</i>	ATCC 15126™ ATCC 6260™ ATCC 18803™	6 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
 - VITEK 2 Gram Negative QC Sets
 - VITEK 2 Gram Positive QC Sets
 - VITEK 2 YST QC Sets
 - VITEK 2 NH QC Sets
 - VITEK 2 ANC QC Sets
- QC for VITEK 2 Susceptibility Cards
 - VITEK 2 AST QC Sets
- QC for API® Identification Strips
 - API 20 Strep QC Sets
 - API 20 A QC Set
 - API 20 C AUX QC Sets
 - API 20 E QC Sets
 - API ID 32 GN QC Set
 - API Coryne QC Sets
 - API 20 NE QC Sets
 - API NH QC Sets
 - RapiD 20 E™ QC Set
 - rapid ID 32 E QC Set
 - API Staph QC Sets
- QC for ATB™ Susceptibility Strips
 - ATB Gram Negative QC Set
 - ATB Gram Positive QC Set
 - ATB ANA QC Set
 - ATB HAEMO QC Set
 - ATB FUNGUS QC Set
- Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for API® Identification Strips (page 2 of 4)

API 20 E QC Sets

411003	LyfoCults Plus API 20 E Streamlined QC set Kit contains: <i>Proteus mirabilis</i>	ATCC® 35659™	2 Vials 2 Vials
410624	LyfoCults Plus API 20 E Comprehensive QC set Kit contains: <i>Enterobacter cloacae</i> <i>Escherichia coli</i> <i>Klebsiella pneumoniae</i> <i>Proteus mirabilis</i> <i>Stenotrophomonas maltophilia</i>	ATCC 13047™ ATCC 25922™ ATCC 35657™ ATCC 35659™ ATCC 51331™	10 Vials 2 Vials 2 Vials 2 Vials 2 Vials

API ID 32 GN QC Set

410588	LyfoCults Plus API ID 32 GN QC set Kit contains: <i>Delftia acidovorans</i> <i>Sphingobacterium multivorum</i>	ATCC 43868™ ATCC 35656™	4 Vials 2 Vials 2 Vials
--------	---	----------------------------	-------------------------------

API Coryne QC Sets

411001	LyfoCults Plus API Coryne Streamlined QC set Kit contains: <i>Corynebacterium renale</i>	ATCC 19412™	2 Vials 2 Vials
410625	LyfoCults Plus API Coryne Comprehensive QC set Kit contains: <i>Cellulosimicrobium cellulans</i> <i>Corynebacterium renale</i> <i>Microbacterium testaceum</i>	ATCC 27402™ ATCC 19412™ ATCC 15829™	6 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
 - VITEK 2 Gram Negative QC Sets
 - VITEK 2 Gram Positive QC Sets
 - VITEK 2 YST QC Sets
 - VITEK 2 NH QC Sets
 - VITEK 2 ANC QC Sets
- QC for VITEK 2 Susceptibility Cards
 - VITEK 2 AST QC Sets
- QC for API® Identification Strips
 - API 20 Strep QC Sets
 - API 20 A QC Set
 - API 20 C AUX QC Sets
 - API 20 E QC Sets
 - API ID 32 GN QC Set
 - API Coryne QC Sets
 - API 20 NE QC Sets
 - API NH QC Sets
 - RapiD 20 E™ QC Set
 - rapid ID 32 E QC Set
 - API Staph QC Sets
- QC for ATB™ Susceptibility Strips
 - ATB Gram Negative QC Set
 - ATB Gram Positive QC Set
 - ATB ANA QC Set
 - ATB HAEMO QC Set
 - ATB FUNGUS QC Set
- Individual QC Strains

Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for API® Identification Strips (page 3 of 4)

API 20 NE QC Sets

411006	LyfoCults Plus API 20 NE Streamlined QC set Kit contains: <i>Aeromonas hydrophila</i> <i>Alcaligenes faecalis</i>	ATCC® 35654™ ATCC 35655™	4 Vials 2 Vials 2 Vials
410626	LyfoCults Plus API 20 NE Comprehensive QC set Kit contains: <i>Aeromonas hydrophila</i> <i>Alcaligenes faecalis</i> <i>Pseudomonas aeruginosa</i> <i>Sphingobacterium multivorum</i>	ATCC 35654™ ATCC 35655™ ATCC 27853™ ATCC 35656™	8 Vials 2 Vials 2 Vials 2 Vials 2 Vials

API NH QC Sets

411005	LyfoCults Plus API NH Streamlined QC set Kit contains: <i>Neisseria gonorrhoeae</i>	ATCC 31426™	2 Vials 2 Vials
410622	LyfoCults Plus API NH Comprehensive QC set Kit contains: <i>Haemophilus influenzae</i> <i>Aggregatibacter aphrophilus</i> <i>Neisseria gonorrhoeae</i>	ATCC 10211™ ATCC 49917™ ATCC 31426™	6 Vials 2 Vials 2 Vials 2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.



QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
- VITEK 2 Gram Negative QC Sets
- VITEK 2 Gram Positive QC Sets
- VITEK 2 YST QC Sets
- VITEK 2 NH QC Sets
- VITEK 2 ANC QC Sets
- QC for VITEK 2 Susceptibility Cards
- VITEK 2 AST QC Sets
- QC for API® Identification Strips
- API 20 Strep QC Sets
- API 20 A QC Set
- API 20 C AUX QC Sets
- API 20 E QC Sets
- API ID 32 GN QC Set
- API Coryne QC Sets
- API 20 NE QC Sets
- API NH QC Sets
- RapiD 20 E™ QC Set
- rapid ID 32 E QC Set
- API Staph QC Sets
- QC for ATB™ Susceptibility Strips
- ATB Gram Negative QC Set
- ATB Gram Positive QC Set
- ATB ANA QC Set
- ATB HAEMO QC Set
- ATB FUNGUS QC Set
- Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



QC for API® Identification Strips (page 4 of 4)

RapiD 20 E™ QC Set

410621	LyfoCults Plus RapiD 20 E QC set		6 Vials
	Kit contains:		
	<i>Escherichia coli</i>	ATCC® 11775™	2 Vials
	<i>Klebsiella pneumoniae</i>	ATCC 35657™	2 Vials
	<i>Proteus hauseri</i>	ATCC 13315™	2 Vials

rapid ID 32 E QC Set

410619	LyfoCults Plus rapid ID 32 E QC set		6 Vials
	Kit contains:		
	<i>Klebsiella oxytoca</i>	ATCC 43863™	2 Vials
	<i>Morganella morganii</i>	ATCC 25829™	2 Vials
	<i>Proteus mirabilis</i>	ATCC 35659™	2 Vials

API Staph QC Sets

411000	LyfoCults Plus API Staph Streamlined QC set		2 Vials
	Kit contains:		
	<i>Staphylococcus capitis</i>	ATCC 35661™	2 Vials
410620	LyfoCults Plus API Staph Comprehensive QC set		6 Vials
	Kit contains:		
	<i>Staphylococcus capitis</i>	ATCC 35661™	2 Vials
	<i>Staphylococcus lentus</i>	ATCC 700403™	2 Vials
	<i>Staphylococcus xylosus</i>	ATCC 700404™	2 Vials



The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus
LyfoCults Plus

- QC for VITEK® 2 Identification Cards
- VITEK 2 Gram Negative QC Sets
- VITEK 2 Gram Positive QC Sets
- VITEK 2 YST QC Sets
- VITEK 2 NH QC Sets
- VITEK 2 ANC QC Sets
- QC for VITEK 2 Susceptibility Cards
- VITEK 2 AST QC Sets
- QC for API® Identification Strips
- API 20 Strep QC Sets
- API 20 A QC Set
- API 20 C AUX QC Sets
- API 20 E QC Sets
- API ID 32 GN QC Set
- API Coryne QC Sets
- API 20 NE QC Sets
- API NH QC Sets
- RapiD 20 E™ QC Set
- rapid ID 32 E QC Set
- API Staph QC Sets
- QC for ATB™ Susceptibility Strips
- ATB Gram Negative QC Set
- ATB Gram Positive QC Set
- ATB ANA QC Set
- ATB HAEMO QC Set
- ATB FUNGUS QC Set
- Individual QC Strains



Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus QC for ATB™ Susceptibility Strips ATB Gram Negative QC Set LyfoCults™ +++ PLUS			
410997	LyfoCults Plus ATB Gram Negative QC set Kit contains: <i>Escherichia coli</i> <i>Escherichia coli</i> <i>Pseudomonas aeruginosa</i>	ATCC® 25922™ ATCC 35218™ ATCC 27853™	6 Vials 2 Vials 2 Vials 2 Vials
ATB Gram Positive QC Set			
410995	LyfoCults Plus ATB Gram Positive QC set Kit contains: <i>Enterococcus faecalis</i> <i>Streptococcus pneumoniae</i> <i>Staphylococcus aureus</i> <i>Staphylococcus aureus</i>	ATCC 29212™ ATCC 49619™ ATCC 29213™ ATCC 43300™	8 Vials 2 Vials 2 Vials 2 Vials 2 Vials
ATB ANA QC Set			
410999	LyfoCults Plus ATB ANA QC set Kit contains: <i>Bacteroides fragilis</i>	ATCC 25285™	2 Vials 2 Vials
ATB HAEMO QC Set			
410996	LyfoCults Plus ATB HAEMO QC set Kit contains: <i>Haemophilus influenzae</i>	ATCC 49247™	2 Vials 2 Vials
ATB FUNGUS QC Set			
410998	LyfoCults Plus ATB FUNGUS QC set Kit contains: <i>Candida krusei</i> <i>Candida parapsilosis</i>	ATCC 6258™ ATCC 22019™	4 Vials 2 Vials 2 Vials

ATCC Licensed
Derivative

®

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4



Reference	Description	Equivalent to Strain	Presentation
-----------	-------------	----------------------	--------------

LyfoCults® Plus



Individual QC Strains (page 1 of 4)

300772	LyfoCults Plus <i>Acinetobacter baumannii</i>	ATCC® BAA-747™	2 Vials
300750	LyfoCults Plus <i>Aeromonas hydrophila</i>	ATCC 35654™	2 Vials
301215	LyfoCults Plus <i>Aggregatibacter aphrophilus</i>	ATCC 49917™	2 Vials
301192	LyfoCults Plus <i>Aggregatibacter aphrophilus</i>	ATCC 33389™	2 Vials
300765	LyfoCults Plus <i>Alcaligenes faecalis</i>	ATCC 35655™	2 Vials
300776	LyfoCults Plus <i>Aneurinibacillus aneurinolyticus</i>	ATCC 11376™	2 Vials
410270	LyfoCults Plus <i>Arcanobacterium haemolyticum</i>	ATCC BAA-1784™	2 Vials
303485	LyfoCults Plus <i>Bacillus badius</i>	ATCC 14574™	2 Vials
300803	LyfoCults Plus <i>Bacillus circulans</i>	ATCC 61™	2 Vials
300836	LyfoCults Plus <i>Bacillus megabacterium</i>	ATCC 14581™	2 Vials
303489	LyfoCults Plus <i>Bacillus pumilus</i>	ATCC BAA-1434™	2 Vials
300826	LyfoCults Plus <i>Bacillus subtilis</i>	ATCC 6633™	2 Vials
300850	LyfoCults Plus <i>Bacteroides fragilis</i>	ATCC 25285™	2 Vials
303524	LyfoCults Plus <i>Bacteroides ovatus</i>	ATCC BAA-1296™	2 Vials
300862	LyfoCults Plus <i>Bacteroides ovatus</i>	ATCC 8483™	2 Vials
303592	LyfoCults Plus <i>Bacteroides ovatus</i>	ATCC BAA-1304™	2 Vials
300868	LyfoCults Plus <i>Bacteroides vulgatus</i>	ATCC 8482™	2 Vials
300784	LyfoCults Plus <i>Brevibacillus agri</i>	ATCC 51663™	2 Vials
300805	LyfoCults Plus <i>Brevibacillus laterosporus</i>	ATCC 64™	2 Vials
301497	LyfoCults Plus <i>Burkholderia cepacia</i>	ATCC 25416™	2 Vials
300905	LyfoCults Plus <i>Campylobacter jejuni ssp. jejuni</i>	ATCC 33291™	2 Vials
300908	LyfoCults Plus <i>Campylobacter jejuni ssp. jejuni</i>	ATCC BAA-1153™	2 Vials
300927	LyfoCults Plus <i>Candida albicans</i>	ATCC 14053™	2 Vials
300921	LyfoCults Plus <i>Candida albicans</i>	ATCC 10231™	2 Vials
300914	LyfoCults Plus <i>Candida albicans</i>	ATCC 90028™	2 Vials
300953	LyfoCults Plus <i>Candida glabrata</i>	ATCC MYA-2950™	2 Vials
300947	LyfoCults Plus <i>Candida glabrata</i>	ATCC 15126™	2 Vials
300945	LyfoCults Plus <i>Candida guilliermondii</i>	ATCC 6260™	2 Vials
300959	LyfoCults Plus <i>Candida krusei</i>	ATCC 6258™	2 Vials
300962	LyfoCults Plus <i>Candida lusitanae</i>	ATCC 34449	2 Vials
300965	LyfoCults Plus <i>Candida parapsilosis</i>	ATCC 22019™	2 Vials
300983	LyfoCults Plus <i>Candida utilis</i>	ATCC 9950	2 Vials

ATCC Licensed Derivative

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4



Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus Individual QC Strains (page 2 of 4)			
410279	LyfoCults Plus <i>Cellulosimicrobium cellulans</i>	ATCC® BAA-1816™	2 Vials
410282	LyfoCults Plus <i>Cellulosimicrobium cellulans</i>	ATCC BAA-1817™	2 Vials
300996	LyfoCults Plus <i>Cellulosimicrobium cellulans</i>	ATCC 27402	2 Vials
300986	LyfoCults Plus <i>Clostridium difficile</i>	ATCC 9689™	2 Vials
300993	LyfoCults Plus <i>Clostridium perfringens</i>	ATCC 13124™	2 Vials
300999	LyfoCults Plus <i>Clostridium septicum</i>	ATCC 12464	2 Vials
301000	LyfoCults Plus <i>Clostridium sordellii</i>	ATCC 9714™	2 Vials
410276	LyfoCults Plus <i>Corynebacterium renale</i>	ATCC BAA-1785™	2 Vials
301016	LyfoCults Plus <i>Corynebacterium renale</i>	ATCC 19412™	2 Vials
303528	LyfoCults Plus <i>Corynebacterium striatum</i>	ATCC BAA-1293™	2 Vials
410273	LyfoCults Plus <i>Corynebacterium urealyticum</i>	ATCC 43044™	2 Vials
301039	LyfoCults Plus <i>Cryptococcus laurentii</i>	ATCC 18803™	2 Vials
410291	LyfoCults Plus <i>Curtobacterium pusillum</i>	ATCC 19096™	2 Vials
303600	LyfoCults Plus <i>Delftia acidovorans</i>	ATCC 43868™	2 Vials
301064	LyfoCults Plus <i>Eikenella corrodens</i>	ATCC BAA-1152™	2 Vials
301153	LyfoCults Plus <i>Elizabethkingia meningoseptica</i>	ATCC 13253™	2 Vials
301053	LyfoCults Plus <i>Enterobacter aerogenes</i>	ATCC 13048™	2 Vials
301072	LyfoCults Plus <i>Enterobacter hormaechei</i>	ATCC 700323™	2 Vials
301067	LyfoCults Plus <i>Enterobacter cloacae</i>	ATCC 13047™	2 Vials
301105	LyfoCults Plus <i>Enterococcus casseliflavus</i>	ATCC 700327™	2 Vials
301090	LyfoCults Plus <i>Enterococcus faecalis</i>	ATCC 19433™	2 Vials
301095	LyfoCults Plus <i>Enterococcus faecalis</i>	ATCC 29212™	2 Vials
301107	LyfoCults Plus <i>Enterococcus faecalis</i> (VanB)	ATCC 51299™	2 Vials
303502	LyfoCults Plus <i>Enterococcus saccharolyticus</i>	ATCC 43076™	2 Vials
301145	LyfoCults Plus <i>Escherichia coli</i>	ATCC 8739™	2 Vials
301131	LyfoCults Plus <i>Escherichia coli</i>	ATCC 25922™	2 Vials
301140	LyfoCults Plus <i>Escherichia coli</i>	ATCC 35218™	2 Vials
301122	LyfoCults Plus <i>Escherichia coli</i>	ATCC 11775™	2 Vials
300873	LyfoCults Plus <i>Geotrichum capitatum</i>	ATCC 28576™	2 Vials
301182	LyfoCults Plus <i>Haemophilus influenzae</i>	ATCC 49247™	2 Vials
301201	LyfoCults Plus <i>Haemophilus influenzae</i>	ATCC 9007™	2 Vials
301198	LyfoCults Plus <i>Haemophilus influenzae</i>	ATCC 10211™	2 Vials

ATCC Licensed Derivative

®

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4



Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus Individual QC Strains (page 3 of 4) 301188 LyfoCults Plus <i>Haemophilus influenzae</i> ATCC® 49766™ 2 Vials 301221 LyfoCults Plus <i>Klebsiella oxytoca</i> ATCC 700324™ 2 Vials 301236 LyfoCults Plus <i>Klebsiella oxytoca</i> ATCC 43863™ 2 Vials 301250 LyfoCults Plus <i>Klebsiella pneumoniae</i> ATCC 35657™ 2 Vials 301247 LyfoCults Plus <i>Klebsiella pneumoniae</i> ATCC 700603™ 2 Vials 301281 LyfoCults Plus <i>Kloeckera apiculata</i> var. <i>apis</i> ATCC 32857™ 2 Vials 303509 LyfoCults Plus <i>Kloeckera japonica</i> ATCC 58370™ 2 Vials 301277 LyfoCults Plus <i>Kocuria kristinae</i> ATCC BAA-752™ 2 Vials 301283 LyfoCults Plus <i>Listeria monocytogenes</i> ATCC BAA-751™ 2 Vials 410288 LyfoCults Plus <i>Microbacterium liquefaciens</i> ATCC BAA -1819™ 2 Vials 410285 LyfoCults Plus <i>Microbacterium paraoxydans</i> ATCC BAA-1818™ 2 Vials 301407 LyfoCults Plus <i>Microbacterium testaceum</i> ATCC 15829™ 2 Vials 301324 LyfoCults Plus <i>Morganella morganii</i> ATCC 25829™ 2 Vials 301316 LyfoCults Plus <i>Neisseria gonorrhoeae</i> ATCC 19424™ 2 Vials 301328 LyfoCults Plus <i>Neisseria gonorrhoeae</i> ATCC 43069™ 2 Vials 301341 LyfoCults Plus <i>Neisseria gonorrhoeae</i> ATCC 49226™ 2 Vials 301321 LyfoCults Plus <i>Neisseria gonorrhoeae</i> ATCC 31426™ 2 Vials 301351 LyfoCults Plus <i>Neisseria lactamica</i> ATCC 23970™ 2 Vials 301375 LyfoCults Plus <i>Ochrobactrum anthropi</i> ATCC BAA-749™ 2 Vials 301382 LyfoCults Plus <i>Oligella ureolytica</i> ATCC 43534™ 2 Vials 301385 LyfoCults Plus <i>Oligella urethralis</i> ATCC 17960™ 2 Vials 301399 LyfoCults Plus <i>Paenibacillus macerans</i> ATCC 8509™ 2 Vials 301445 LyfoCults Plus <i>Paenibacillus polymyxa</i> ATCC 7070™ 2 Vials 303492 LyfoCults Plus <i>Paenibacillus validus</i> ATCC 29948™ 2 Vials 303531 LyfoCults Plus <i>Parabacteroides distasonis</i> ATCC BAA-1295™ 2 Vials 301448 LyfoCults Plus <i>Proteus hauseri</i> ATCC 13315™ 2 Vials 301436 LyfoCults Plus <i>Proteus mirabilis</i> ATCC 35659™ 2 Vials 301455 LyfoCults Plus <i>Proteus vulgaris</i> ATCC 6380™ 2 Vials 301473 LyfoCults Plus <i>Prototheca wickerhamii</i> ATCC 16529™ 2 Vials 301480 LyfoCults Plus <i>Pseudomonas aeruginosa</i> ATCC 27853™ 2 Vials 303498 LyfoCults Plus <i>Pseudomonas aeruginosa</i> ATCC BAA-1744™ 2 Vials 301490 LyfoCults Plus <i>Pseudomonas aeruginosa</i> ATCC 9721™ 2 Vials			

ATCC Licensed Derivative

®

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC.

bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

Denmark: order.dk@biomerieux.com Finland: order.fi@biomerieux.com Norway: order.no@biomerieux.com Sweden: order.se@biomerieux.com www.biomerieux.com

48

QC ORGANISMS

Feature – LyfoCults® Plus

LyfoCults Plus

QC for VITEK® 2 Identification Cards

VITEK 2 Gram Negative QC Sets

VITEK 2 Gram Positive QC Sets

VITEK 2 YST QC Sets

VITEK 2 NH QC Sets

VITEK 2 ANC QC Sets

QC for VITEK 2 Susceptibility Cards

VITEK 2 AST QC Sets

QC for API® Identification Strips

API 20 Strep QC Sets

API 20 A QC Set

API 20 C AUX QC Sets

API 20 E QC Sets

API ID 32 GN QC Set

API Coryne QC Sets

API 20 NE QC Sets

API NH QC Sets

RapiD 20 E™ QC Set

rapid ID 32 E QC Set

API Staph QC Sets

QC for ATB™ Susceptibility Strips

ATB Gram Negative QC Set

ATB Gram Positive QC Set

ATB ANA QC Set

ATB HAEMO QC Set

ATB FUNGUS QC Set

Individual QC Strains

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4



Reference	Description	Equivalent to Strain	Presentation
LyfoCults® Plus Individual QC Strains (page 4 of 4) 301526 LyfoCults Plus <i>Salmonella typhimurium</i> ATCC® 14028™ 2 Vials 301579 LyfoCults Plus <i>Shigella sonnei</i> ATCC 25931™ 2 Vials 301583 LyfoCults Plus <i>Sphingobacterium multivorum</i> ATCC 35656™ 2 Vials 303514 LyfoCults Plus <i>Sporobolomyces salmonicolor</i> ATCC MYA-4550™ 2 Vials 301589 LyfoCults Plus <i>Staphylococcus aureus</i> ATCC 29213™ 2 Vials 301616 LyfoCults Plus <i>Staphylococcus aureus</i> ATCC BAA-1026™ 2 Vials 301586 LyfoCults Plus <i>Staphylococcus aureus</i> ATCC 25923™ 2 Vials 301627 LyfoCults Plus <i>Staphylococcus aureus</i> ATCC BAA-977™ 2 Vials 301624 LyfoCults Plus <i>Staphylococcus aureus</i> ATCC BAA-976™ 2 Vials 301621 LyfoCults Plus <i>Staphylococcus aureus</i> (MRSA) ATCC 43300™ 2 Vials 301630 LyfoCults Plus <i>Staphylococcus capitis</i> ATCC 35661™ 2 Vials 301634 LyfoCults Plus <i>Staphylococcus epidermidis</i> ATCC 12228™ 2 Vials 301640 LyfoCults Plus <i>Staphylococcus epidermidis</i> ATCC 14990™ 2 Vials 301650 LyfoCults Plus <i>Staphylococcus lentus</i> ATCC 700403™ 2 Vials 301658 LyfoCults Plus <i>Staphylococcus saprophyticus</i> ATCC BAA-750™ 2 Vials 301662 LyfoCults Plus <i>Staphylococcus sciuri ssp. sciuri</i> ATCC 29061™ 2 Vials 301664 LyfoCults Plus <i>Staphylococcus xylosus</i> ATCC 700404™ 2 Vials 301697 LyfoCults Plus <i>Stenotrophomonas maltophilia</i> ATCC 17666™ 2 Vials 301668 LyfoCults Plus <i>Stenotrophomonas maltophilia</i> ATCC 51331™ 2 Vials 301701 LyfoCults Plus <i>Streptococcus equi subsp. Zooepidemicus</i> ATCC 700400™ 2 Vials 301707 LyfoCults Plus <i>Streptococcus equi subsp. Zooepidemicus</i> ATCC 43079™ 2 Vials 301728 LyfoCults Plus <i>Streptococcus pneumoniae</i> ATCC 49619™ 2 Vials 301731 LyfoCults Plus <i>Streptococcus pneumoniae</i> ATCC 6303™ 2 Vials 301744 LyfoCults Plus <i>Streptococcus pyogenes</i> ATCC 19615™ 2 Vials 301762 LyfoCults Plus <i>Streptococcus sp</i> ATCC 12392™ 2 Vials 301757 LyfoCults Plus <i>Streptococcus thermophilus</i> ATCC 19258™ 2 Vials 301767 LyfoCults Plus <i>Streptococcus uberis</i> ATCC 700407™ 2 Vials 301773 LyfoCults Plus <i>Trichophyton mentagrophytes</i> ATCC 9533™ 2 Vials 301780 LyfoCults Plus <i>Trichosporon mucoides</i> ATCC 204094™ 2 Vials 301783 LyfoCults Plus <i>Virgibacillus pantothenicus</i> ATCC 14576™ 2 Vials 303518 LyfoCults Plus <i>Zygosaccharomyces bailii</i> ATCC MYA-4549™ 2 Vials			

ATCC Licensed Derivative

®

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. bioMérieux is licensed to use these trademarks and to sell products derived from ATCC cultures.

CULTURE MEDIA

- Transport Media & Dipslides
 - All Purpose Transport Medium
 - ESwab Collection & Transport Medium
 - Specific Transport Medium
 - Urinary Bacterial Count
 - Transport & Culture of Gonococci
 - Transport & Culture of Fungi
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Transport Media & Dipslides

All Purpose Transport Medium

42105	PORTAGERM™ agar (tube without swab)	20 Tubes
41998	PORTAGERM Amies agar & swab - sterile zone	Pack of 10 Units
41999	PORTAGERM Amies agar + swab	Pack of 50 Units
41995	PORTAGERM vial for liquid samples	Pack of 10 Vials

ESwab Collection & Transport Medium

41992	ESwab - nose, throat, vagina, rectum, wound	Pack of 50 Units
41990	ESwab - eye, ear, nasal passages, throat, urogenital tract	Pack of 50 Units
41991	ESwab - naso pharynx and pediatric samples	Pack of 50 Units

Specific Transport Medium

56528	Gonoline® DUO 2 + CO2 generators	10 Dipslides
56525	Mycoline®	10 Dipslides
42041	PORTAGERM pylori	8 Bottles x 2.5ml

Urinary Bacterial Count

56507	URILINE™	10 Dipslides
56508	URILINE	100 Dipslides

Transport & Culture of Gonococci

56528	Gonoline DUO 2 + CO2 generators	10 Dipslides
-------	---------------------------------	--------------

Transport & Culture of Fungi

56525	Mycoline	10 Dipslides
-------	----------	--------------



CULTURE MEDIA

Transport Media & Dipslides

Pre-Poured Media – Plates

Chromogenic Media

- Direct identification of *E. coli*, *Proteus* & *Enterococcus*
- Direct identification of *Candida albicans*
- Detection of *Salmonellas*
- Isolation & identification of *Staphylococci*
- Screening of *S. agalactiae*
- Screening of Multidrug-Resistant Organisms
- Direct identification of *P. aeruginosa*
- Direct identification of *Vibrio* sp

General Purpose Media

Specific Media

Yeasts & Moulds

Media for Antibiotic Susceptibility:

Disc Diffusion Methods

Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles

Pre-Poured Media – Tubes & Bottles

Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

Air Sampling

Dehydrated Media

Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates



Chromogenic Media (page 1 of 2)

Direct identification of *E. coli*, *Proteus* & *Enterococcus*

43821 NEW	chromID™ CPS	20 Plates
43829 NEW	chromID CPS	100 Plates
43541	chromID CPS	20 Plates
43549	chromID CPS	100 Plates
43473	chromID CPS/Columbia CNA + 5% sheep blood	20 Plates
411617 NEW	chromID CPS/Columbia CNA + 5% sheep blood	20 Plates
56541	ID indole TDA reagents	2 x 200 Tests

Direct identification of *Candida albicans*


43631	chromID Candida	20 Plates
43639	chromID Candida	100 Plates
43464	chromID Candida/Sabouraud Gentamicin Chloramphenicol 2	20 Plates

Detection of *Salmonellas*

43621	chromID Salmonella	20 Plates
43465	chromID Salmonella/Hektoen	20 Plates

Isolation & identification of *Staphylococci*

43371	chromID <i>S. aureus</i>	20 Plates
43466 NEW	chromID MRSA/chromID <i>S. aureus</i>	20 Plates



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
 - Chromogenic Media
 - Direct identification of E. coli, Proteus & Enterococcus
 - Direct identification of Candida albicans
 - Detection of Salmonellas
 - Isolation & identification of Staphylococci
 - Screening of S. agalactiae
 - Screening of Multidrug-Resistant Organisms
 - Direct identification of P. aeruginosa
 - Direct identification of Vibrio sp
 - General Purpose Media
 - Specific Media
 - Yeasts & Moulds
 - Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - Divided Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates



Chromogenic Media (page 2 of 2)

Screening of S. agalactiae

43461	chromID™ Strepto B	20 Plates
-------	--------------------	-----------

Screening of Multidrug-Resistant Organisms

43871	NEW chromID C. difficile	20 Plates
43861	NEW chromID CARBA agar	20 Plates
43481	chromID ESBL	20 Plates
43451	chromID MRSA	20 Plates
43459	chromID MRSA	100 Plates
43466	NEW chromID MRSA/chromID S. aureus	20 Plates
43002	chromID VRE	20 Plates
43004	NEW chromID VRE	20 Plates

Direct identification of P. aeruginosa

43462	chromID P. aeruginosa	20 Plates
-------	-----------------------	-----------

Direct identification of Vibrio sp

43761	chromID Vibrio	20 Plates
43762	NEW chromID Vibrio CE	20 Plates



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
 - Chromogenic Media
 - General Purpose Media
 - Non Selective agar
 - Non Selective blood agar
 - Selective blood agar
 - Lactose agar
 - Specific Media
 - Yeasts & Moulds
 - Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - Divided Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates

General Purpose Media (page 1 of 2)

Non Selective agar

33606	Nutrient agar	20 Plates
43011	Trypcase Soy agar	20 Plates

Non Selective blood agar

33633	Blood agar (horse blood)	20 Plates
44282	Brain Heart Infusion agar	20 Plates
44117	Columbia agar chocolated	20 Plates
43101	Chocolate agar + PolyViteX™	20 Plates
43109	Chocolate agar + PolyViteX	100 Plates
43050	Columbia agar + 5% horse blood	20 Plates
43059	Columbia agar + 5% horse blood	100 Plates
44813	Columbia agar + 5% horse blood + NAD	20 Plates
44724	Columbia agar + 5% horse blood + NAD	100 Plates
43041	Columbia agar + 5% sheep blood	20 Plates
43049	Columbia agar + 5% sheep blood	100 Plates
43401	Schaedler agar + 5% sheep blood	20 Plates
43061	Trypcase Soy agar + 5% horse blood *	20 Plates
43001	Trypcase Soy agar + 5% sheep blood	20 Plates
44279	Yeast Extract agar	20 Plates

Selective blood agar

44196	Columbia agar + Nalidixic Acid + 5% horse blood	20 Plates
43071	Columbia CNA agar + 5% sheep blood	20 Plates
43079	Columbia CNA agar + 5% sheep blood	100 Plates
43223	Schaedler Neomycin Vancomycin agar + 5% sheep blood	10 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
 - Chromogenic Media
 - General Purpose Media
 - Non Selective agar
 - Non Selective blood agar
 - Selective blood agar
 - Lactose agar
 - Specific Media
 - Yeasts & Moulds
 - Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - Divided Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------


Pre-Poured Media – Plates

General Purpose Media (page 2 of 2)

Lactose agar

43331	CLED medium	20 Plates
43339	CLED medium	100 Plates
44113	CLED medium (with Andrade Indicator)	20 Plates
43341	Drigalski agar *	20 Plates
43231	Endo agar *	20 Plates
43081	Eosin Methylene Blue agar	20 Plates
44153	MacConkey agar (with salt)	20 Plates
43141	MacConkey agar (with salt and crystal violet)	20 Plates
44309	MacConkey agar No.2 *	20 Plates
33648	MacConkey agar without salt	20 Plates
43021	Purple Lactose agar *	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA

Transport Media & Dipslides

Pre-Poured Media – Plates

Chromogenic Media

General Purpose Media

Specific Media

- Actinomyces
- Anaerobes
- Bacillus cereus
- Bordetella
- Campylobacter
- Clostridium difficile
- Corynebacterium
- Dermatophytes
- E. coli O157
- Enterobacteriaceae
- Enterococcus
- Gardnerella
- Haemophilus
- Helicobacter pylori
- Legionella
- Listeria
- Urogenital Mycoplasmas
- Neisseria
- Pseudomonas
- Salmonella, Shigella
- Staphylococcus
- Streptococcus
- Vibrio
- Yersinia

Yeasts & Moulds

Media for Antibiotic Susceptibility:

Disc Diffusion Methods

Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles

Pre-Poured Media – Tubes & Bottles

Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

Air Sampling

Dehydrated Media

Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates

Specific Media (page 1 of 5)

Actinomyces

44814	Actinomyces agar	10 Plates
-------	------------------	-----------

Anaerobes

44100	Anaerobe agar + 5% horse blood	20 Plates
44102	Anaerobe agar + Nalidixic acid + Tween + 5% horse blood	20 Plates
33624	Anaerobe agar + Neomycin + 5% horse blood	20 Plates
44196	Columbia agar + Nalidixic Acid + 5% horse blood	20 Plates
43401	Schaedler agar + 5% sheep blood	20 Plates
43223	Schaedler Neomycin Vancomycin agar + 5% sheep blood *	10 Plates

Bacillus cereus

44105	Bacillus cereus selective agar	20 Plates
-------	--------------------------------	-----------

Bordetella

44110	Bordet Gengou agar	10 Plates
44358	Charcoal agar + Cephalixin	10 Plates

Campylobacter

33627	Campylobacter Selective agar (CCDA)	20 Plates
43361	Campylosel agar	20 Plates
43369	Campylosel agar	100 Plates

Clostridium difficile

43431	Clostridium difficile agar	20 Plates
-------	----------------------------	-----------

Corynebacterium

44136	Hoyle medium	10 Plates
-------	--------------	-----------

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

Transport Media & Dipslides

Pre-Poured Media – Plates

Chromogenic Media
General Purpose Media
Specific Media

- Actinomyces
- Anaerobes
- Bacillus cereus
- Bordetella
- Campylobacter
- Clostridium difficile
- Corynebacterium
- Dermatophytes
- E. coli O157
- Enterobacteriaceae
- Enterococcus
- Gardnerella
- Haemophilus
- Helicobacter pylori
- Legionella
- Listeria
- Urogenital Mycoplasmas
- Neisseria
- Pseudomonas
- Salmonella, Shigella
- Staphylococcus
- Streptococcus
- Vibrio
- Yersinia

Yeasts & Moulds
Media for Antibiotic Susceptibility:
Disc Diffusion Methods
Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles

Pre-Poured Media – Tubes & Bottles

Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

Air Sampling

Dehydrated Media

Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates



Specific Media (page 2 of 5)

Dermatophytes

43464	chromID™ Candida/Sabouraud Gentamicin Chloramphenicol 2	20 Plates
44815	Dermatophyte Test medium	10 Plates
44125	Dermatophyte Test medium + Phenol Red	10 Plates
56525	Mycoline®	10 Dipslides
44176	Sabouraud Dextrose + Chloramphenicol + Actidione	20 Plates
44175	Sabouraud Dextrose agar (with Chloramphenicol)	20 Plates
43651	Sabouraud Gentamicin Chloramphenicol 2 agar *	20 Plates
43659	Sabouraud Gentamicin Chloramphenicol 2 agar	100 Plates

E. coli O157

43391	CT-SMAC agar	20 Plates
44154	MacConkey agar with Sorbitol (SMAC)	10 Plates

Enterobacteriaceae

43821 NEW	chromID CPS	20 Plates
43829 NEW	chromID CPS	100 Plates
43541	chromID CPS	20 Plates
43549	chromID CPS	100 Plates
43481	chromID ESBL	20 Plates
43331	CLED medium	20 Plates
43339	CLED medium	100 Plates
44113	CLED medium (with Andrade Indicator)	20 Plates
43341	Drigalski agar *	20 Plates
43231	Endo agar *	20 Plates
43239	Endo agar *	100 Plates
43081	Eosin Methylene Blue agar *	20 Plates
44811	Membrane Lauryl Sulphate agar	20 Plates
44153	MacConkey agar (with salt)	20 Plates
43141	MacConkey agar (with salt and crystal violet)	20 Plates
43149	MacConkey agar (with salt and crystal violet)	100 Plates
33648	MacConkey agar without salt	20 Plates
43021	Purple Lactose agar *	20 Plates
44827	Tryptone Bile X	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

Transport Media & Dipslides

Pre-Poured Media – Plates

Chromogenic Media
General Purpose Media
Specific Media

- Actinomyces
- Anaerobes
- Bacillus cereus
- Bordetella
- Campylobacter
- Clostridium difficile
- Corynebacterium
- Dermatophytes
- E. coli O157
- Enterobacteriaceae
- Enterococcus
- Gardnerella
- Haemophilus
- Helicobacter pylori
- Legionella
- Listeria
- Urogenital Mycoplasmas
- Neisseria
- Pseudomonas
- Salmonella, Shigella
- Staphylococcus
- Streptococcus
- Vibrio
- Yersinia

Yeasts & Moulds
Media for Antibiotic Susceptibility:
Disc Diffusion Methods
Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles

Pre-Poured Media – Tubes & Bottles

Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

Air Sampling

Dehydrated Media

Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates



Specific Media (page 3 of 5)

Enterococcus

43002	chromID™ VRE	20 Plates
43004 NEW	chromID VRE	20 Plates
43151	D-Coccosel (Bile Aesculin) agar	20 Plates

Gardnerella

43411	Gardnerella agar (human blood)	20 Plates
-------	--------------------------------	-----------

Haemophilus

43101	Chocolate agar + PolyViteX™	20 Plates
43681	Chocolate Haemophilus agar 2	20 Plates
43689	Chocolate Haemophilus agar 2	100 Plates
44813	Columbia agar + 5% horse blood + NAD	20 Plates
44117	Columbia agar chocolate	20 Plates

Helicobacter pylori

43263	Pylori agar	10 Plates
-------	-------------	-----------

Legionella

44143	Legionella agar (BCYE) *	20 Plates
43013	Legionella BCYE with L-cysteine	10 Plates
43023	Legionella BCYE without L-cysteine	10 Plates
44146	Legionella Selective agar (MWY) *	20 Plates
43031	Legionella GVPC medium	20 Plates
43032	Legionella GVPC medium	100 Plates

Listeria

43560	Listeria Selective agar (Oxford)	20 Plates
43559	Palcam agar *	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

Transport Media & Dipslides

Pre-Poured Media – Plates

Chromogenic Media
General Purpose Media
Specific Media

Actinomyces
Anaerobes
Bacillus cereus
Bordetella
Campylobacter
Clostridium difficile
Corynebacterium
Dermatophytes
E. coli O157
Enterobacteriaceae
Enterococcus
Gardnerella
Haemophilus
Helicobacter pylori
Legionella
Listeria

Urogenital Mycoplasmas
Neisseria
Pseudomonas
Salmonella, Shigella
Staphylococcus
Streptococcus
Vibrio
Yersinia

Yeasts & Moulds
Media for Antibiotic Susceptibility:
Disc Diffusion Methods
Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles

Pre-Poured Media – Tubes & Bottles

Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

Air Sampling

Dehydrated Media

Supplements



ReferenceDescriptionPresentation

Pre-Poured Media – Plates



Specific Media (page 4 of 5)

Urogenital Mycoplasmas

43003	A7 agar *	10 Plates
-------	-----------	-----------

Neisseria

43101	Chocolate agar + PolyViteX™	20 Plates
43611	Chocolate agar + PolyViteX + VCAT 3	20 Plates
33625	New York City medium + LCAT	20 Plates
44163	New York City medium + VCAT	20 Plates
44824	New York City medium + VCT	20 Plates

Pseudomonas

44347	Cepacia Selective agar	10 Plates
43565	Cetrimide agar *	20 Plates
43462	chromID™ P. aeruginosa	20 Plates
44170	Pseudomonas Selective agar (C-F-C)	20 Plates
44171	Pseudomonas Selective agar (C-N)	20 Plates

Salmonella, Shigella

43588	Brilliant Green agar modified	20 Plates
43621	chromID Salmonella	20 Plates
43629	chromID Salmonella *	100 Plates
43465	chromID Salmonella/Hektoen	20 Plates
44127	DCA agar	20 Plates
44128	DCA Hynes agar	20 Plates
44124	DCLS	20 Plates
43231	Endo agar *	20 Plates
43111	Hektoen agar	20 Plates
43119	Hektoen agar	100 Plates
43091	Salmonella Shigella agar *	20 Plates
43099	Salmonella Shigella agar *	100 Plates
43563	XLD agar	20 Plates
43564	XLD agar *	100 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

Transport Media & Dipslides
Pre-Poured Media – Plates

- Chromogenic Media
- General Purpose Media
- Specific Media
 - Actinomyces
 - Anaerobes
 - Bacillus cereus
 - Bordetella
 - Campylobacter
 - Clostridium difficile
 - Corynebacterium
 - Dermatophytes
 - E. coli O157
 - Enterobacteriaceae
 - Enterococcus
 - Gardnerella
 - Haemophilus
 - Helicobacter pylori
 - Legionella
 - Listeria
 - Urogenital Mycoplasmas
 - Neisseria
 - Pseudomonas
 - Salmonella, Shigella
 - Staphylococcus
 - Streptococcus
 - Vibrio
 - Yersinia

- Yeasts & Moulds
- Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - Divided Plates

Pre-Poured Media – Plates, Tubes & Bottles
Pre-Poured Media – Tubes & Bottles
Urogenital Mycoplasmas
Incubation System: GENbag-GENbox
Air Sampling
Dehydrated Media
Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates



Specific Media (page 5 of 5)

Staphylococcus

43521	Baird-Parker agar	20 Plates
43451	chromID™ MRSA	20 Plates
43459	chromID MRSA	100 Plates
43466 NEW	chromID MRSA/chromID S. aureus	20 Plates
43371	chromID S. aureus	20 Plates
43379	chromID S. aureus	100 Plates
44122	DNase agar	10 Plates
44123	DNase agar (with Methyl Green) *	10 Plates
43671	Mannitol Salt 2 agar	20 Plates
43679	Mannitol Salt 2 agar	100 Plates
44165	Nutrient agar with 5% salt *	20 Plates

Streptococcus

43461	chromID Strepto B	20 Plates
43071	Columbia CNA agar + 5% sheep blood	20 Plates
44135	G.B.S (Islam) medium *	10 Plates
43712	Granada agar	20 Plates
43467	Granada/Columbia CNA agar + 5% sheep blood	20 Plates
44177	Streptococcus Selective agar (COBA)	10 Plates

Vibrio

43761	chromID Vibrio	20 Plates
44178	T.C.B.S. medium	10 Plates

Yersinia

43421	Yersinia CIN	20 Plates
43209	Yersinia CIN	100 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Yeasts & Moulds
 - Non Selective agar
 - Selective agar
- Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - Divided Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates

Yeasts & Moulds

Non Selective agar

44156	Malt Extract agar *	20 Plates
44715	Malt Extract irradiated agar	20 Contact Plates
44174	Sabouraud Dextrose agar	20 Plates

Selective agar

43631	chromID™ Candida	20 Plates
43464 NEW	chromID Candida/Sabouraud Gentamicin Chloramphenicol 2	20 Plates
44815	Dermatophyte Test medium	10 Plates
44125	Dermatophyte Test medium + Phenol Red	10 Plates
56525	Mycoline®	10 Dipslides
43555 NEW	Sabouraud Dextrose agar	20 Plates
44717	Sabouraud Dextrose irradiated agar	20 Contact Plates
44175	Sabouraud Dextrose agar (with Chloramphenicol)	20 Plates
43596 NEW	Sabouraud Dextrose Chloramphenicol agar	20 Plates
44176	Sabouraud Dextrose + Chloramphenicol + Actidione agar	20 Plates
43651	Sabouraud Gentamicin Chloramphenicol 2 agar *	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Yeasts & Moulds
- Media for Antibiotic Susceptibility:
 - Disc Diffusion Methods
 - BSAC compliant media
 - NCCLS compliant media
 - Divided Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates

Media for Antibiotic Susceptibility: Disc Diffusion Methods

BSAC compliant media

44321	Columbia agar + 2% NaCl	10 Plates
33620	Sensitivity Test agar	20 Plates
44138	Sensitivity Test agar + 5% horse blood	20 Plates
44198	Sensitivity Test agar + 5% horse blood + NAD	10 Plates
44721	Sensitivity Test agar + 5% horse blood + NAD	100 Plates

NCCLS compliant media

33612	Haemophilus Test medium (HTM) *	20 Plates
43302	Mueller Hinton 2 agar (150mm) *	20 Plates
43511	Mueller Hinton 2 agar (120 x 120mm) *	20 Square Plates
43301	Mueller Hinton 2 agar *	20 Plates
43309	Mueller Hinton 2 agar *	100 Plates
43321	Mueller Hinton 2 agar + 5% sheep blood *	20 Plates
43901	Mueller Hinton agar + 5% horse blood agar + NAD	20 Plates

Divided Plates

44807	Anaerobe agar + Neomycin + 5% horse blood/ Columbia + 5% horse blood	20 Plates
43464	chromID™ Candida/Sabouraud Gentamicin Chloramphenicol 2	20 Plates
43473	chromID CPS/Columbia CNA + 5% sheep blood	20 Plates
43466 NEW	chromID MRSA/chromID S. aureus	20 Plates
43465	chromID Salmonella/Hektoen	20 Plates
41167	chromID4 CPS/Columbia CNA + 5% sheep blood	20 Plates
43467	Granada/Columbia CNA agar + 5% sheep blood	20 Plates
44808	New York City medium + LCAT/Sabouraud agar + Chloramphenicol	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA


- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
 - Pharmaceutical Applications
 - Irradiated media
 - Page 1 of 2
 - Page 2 of 2
 - Sterility Testing media
 - Pre-Poured Media – Tubes & Bottles
 - Urogenital Mycoplasmas
 - Incubation System: GENbag-GENbox
 - Air Sampling
 - Dehydrated Media
 - Supplements

Reference

Description

Presentation

Pre-Poured Media – Plates, Tubes & Bottles



Pharmaceutical Applications (page 1 of 3)

Irradiated media (page 1 of 2)

96301	Bi-Box (double wrapped and irradiated) †	10 x 10 Contact Plate Boxes
43501 NEW	Count-Tact™ agar	20 Contact Plates
43491	Count-Tact irradiated agar †	20 Contact Plates
43492	Count-Tact irradiated agar †	100 Contact Plates
43593	Count-Tact irradiated agar (isolator) †	20 Contact Plates
43691 NEW	Count-Tact 3P™ irradiated agar †	20 Contact Plates
43699 NEW	Count-Tact 3P irradiated agar †	100 Contact Plates
96300	Count-Tact Applicator †	One
43580 NEW	Count-Tact Sabouraud Dextrose agar with Chloramphenicol and neutralisers †	20 Contact Plates
43581	Count-Tact Sabouraud Dextrose irradiated agar with Chloramphenicol and neutralisers †	20 Contact Plates
43812 NEW	Count-Tact Sabouraud Dextrose 3P irradiated agar with neutralizers †	20 Contact Plates
43582 NEW	Count-Tact Trypcase Soy agar †	20 Contact Plates
44229	Environmental Monitoring irradiated agar *	20 Plates
44300	Environmental Monitoring irradiated agar isolator *	20 Plates
44715	Malt Extract irradiated agar	20 Contact Plates
44717	Sabouraud Dextrose irradiated agar	20 Contact Plates
44343	Sabouraud Dextrose irradiated agar (140mm) *	4 x 5 Plates
44301	Sabouraud Dextrose irradiated agar isolator *	20 Plates
43554	Sabouraud Dextrose 3P irradiated agar	20 Plates
43814 NEW	Sabouraud Dextrose 3P irradiated agar with neutralizers	20 Plates
43595	Sabouraud Dextrose Chloramphenicol irradiated agar	20 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

† FIX system™



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
 - Pharmaceutical Applications
 - Irradiated media
 - Page 1 of 2
 - Page 2 of 2
 - Sterility Testing media
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates, Tubes & Bottles


Pharmaceutical Applications (page 2 of 3)

Irradiated media (page 2 of 2)

43557	Trypcase Soy irradiated agar	100 Plates
D1113	Trypcase Soy irradiated agar †	100 Contact Plates
43131	Trypcase Soy irradiated agar †	20 Plates
43556	Trypcase Soy irradiated agar (140mm)	20 Plates
43594	Trypcase Soy irradiated agar (isolator)	20 Plates
43562	Trypcase Soy irradiated agar with neutralisers	100 Plates
D1106	Trypcase Soy irradiated agar with Penicillinase * †	20 Contact Plates
43711 NEW	Trypcase Soy 3P™ irradiated agar	20 Plates
43169	Trypcase Soy 3P irradiated agar	100 Plates
43811 NEW	Trypcase Soy 3P irradiated agar with neutralizers	20 Plates
43819 NEW	Trypcase Soy 3P irradiated agar with neutralizers	100 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

† FIX system™



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
 - Pharmaceutical Applications
 - Irradiated media
 - Sterility Testing media
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Plates, Tubes & Bottles

Pharmaceutical Applications (page 3 of 3)

Sterility Testing media

43558 NEW	Plate Count agar	20 Plates
42616	Rinsing Solution irradiated	4 x 300ml
42044	Rinsing Solution (EP, USP fluid A) irradiated	12 x 100ml
44502	Thioglycollate broth	10 x 110ml
42074	Thioglycollate broth + resazurin	20 Tubes x 9ml
44001	Thioglycollate USP with injectable cap	12 x 100ml
44277	Trypcase Soy broth	10 x 50ml
42100	Trypcase Soy broth	20 Tubes x 9ml
41146	Trypcase Soy broth	6 x 100ml
44288	Trypcase Soy broth	50 x 20ml
44011	Trypcase Soy broth (EP, USP) irradiated	12 x 100ml
44362	Trypcase Soy broth (wnb) *	10 x 90ml

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

wnb – wide necked bottle.



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - E. coli & coliforms
 - General Purpose Media
 - Non Selective broth
 - Non Selective agar
 - Diluent
 - Specific Media
 - Yeasts & Moulds
 - Media for Biochemical Identification
 - Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
 - Media for Mycobacteria
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

Chromogenic Media

E. coli & coliforms

42605	chromID™ O157:H7	6 x 200ml
42017	Coli ID	6 x 200ml

General Purpose Media (page 1 of 2)

Non Selective broth

44260	Brain Heart Infusion broth	50 x 9ml
42081	Brain Heart Infusion broth	20 Tubes x 9ml
44327	Nutrient broth *	50 x 5ml
44271	Nutrient broth *	50 x 9ml
42100	Trypcase Soy broth	20 Tubes x 9ml
44266	Trypcase Soy broth	50 x 9ml
41146	Trypcase Soy broth	6 x 100ml
44288	Trypcase Soy broth	50 x 20ml

Non Selective agar

41536	Chocolate agar	6 x 100ml
44333	Chocolate agar (slopes)	50 x 3.5ml
42079	Chocolate agar + PolyViteX™ (slopes)	20 Tubes
41244	Columbia agar	6 x 200ml
44283	Nutrient agar (slopes)	50 x 5ml
43551 NEW	R2A agar	20 Plates
42101	Trypcase Soy agar (slopes)	20 Tubes
44267	Trypcase Soy agar (slopes)	50 x 15ml
43019 NEW	Trypcase Soy agar	100 Plates
41467	Trypcase Soy agar	6 x 100ml
41466	Trypcase Soy agar	6 x 200ml

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - E. coli & coliforms
 - General Purpose Media
 - Non Selective broth
 - Non Selective agar
 - Diluent
- Specific Media
- Yeasts & Moulds
- Media for Biochemical Identification
- Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
- Media for Mycobacteria
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

General Purpose Media (page 2 of 2)

Diluent

42076 NEW	Maximum Recovery Diluent	20 Tubes x 9ml
42085	Peptone water	20 Tubes x 9ml
44286	Peptone water	50 x 3ml
44287	Purified water	50 x 4.5ml
44332	Saline Solution 0.85%	50 x 3ml
44329	Tryptone water (Report 71) *	50 x 5ml

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Anaerobes
 - Corynebacterium
 - Enterobacteriaceae
 - Enterococcus
 - Legionella
 - Listeria
 - Mycobacteria
 - Mycoplasma hominis, Ureaplasma urealyticum
 - Salmonella, Shigella
 - Staphylococcus
 - Trichomonas
 - Vibrio
 - Streptococci
 - Yeasts & Moulds
 - Media for Biochemical Identification
 - Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
 - Media for Mycobacteria
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

Specific Media

Anaerobes

44272	Cooked Meat medium	50 x 15ml
41936	Schaedler agar (+ Vit. K3)	6 x 100ml
42106	Schaedler broth (+ Vit. K3)	20 Tubes x 13ml
42097	Schaedler broth + 0.02% agar (+ Vit. K3)	20 Tubes x 13ml
42098	Schaedler broth + 0.2% agar (+ Vit. K3)	20 Tubes x 13ml
42074	Thioglycollate broth + resazurin	20 Tubes x 9ml

Corynebacterium

42088	Loeffler medium (slopes)	20 Tubes
-------	--------------------------	----------

Enterobacteriaceae

41594	CLED agar	6 x 200ml
41336	Eosin Methylene Blue agar	6 x 100ml

Enterococcus

44334	Aesculin broth *	50 x 3ml
42083	D-Coccosel (Bile Aesculin) agar	20 Tubes

Legionella

41054	Legionella - BCYE agar base	6 x 200ml
55641	Legionella enrichment supplement	4 x 2ml
55645	Legionella selective mixture	4 x 2ml

Listeria

44302	Buffered Listeria Enrichment broth	10 x 225ml
42072	Fraser broth	20 Tubes x 10ml
42046	Fraser broth	6 x 225ml
42048	Fraser broth (half)	6 x 225ml
44523	Listeria broth (UVM I)	50 x 10ml

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Anaerobes
 - Corynebacterium
 - Enterobacteriaceae
 - Enterococcus
 - Legionella
 - Listeria
 - Mycobacteria
 - Mycoplasma hominis, Ureaplasma urealyticum
 - Salmonella, Shigella
 - Staphylococcus
 - Trichomonas
 - Vibrio
 - Streptococci
- Yeasts & Moulds
- Media for Biochemical Identification
- Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
 - Media for Mycobacteria
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

Specific Media

Mycobacteria

42082	Coletsos medium (slopes)	20 Tubes
42107	Lowenstein-Jensen + TCH (slopes)	20 Tubes
42089	Lowenstein-Jensen medium (slopes)	20 Tubes

Mycoplasma hominis, Ureaplasma urealyticum

42508	Urea-arginine LYO-2 broth	25 Tests
-------	---------------------------	----------

Salmonella, Shigella

42042	Buffered Peptone water	6 x 90ml
42043	Buffered Peptone water	6 x 225ml
41554	Hektoen agar	6 x 200ml
44536	Mannitol Selenite broth	50 x 15ml
42114 NEW	Muller Kauffmann Tetrathionate and novobiocin broth (MKTn)	20 x 10ml
42091	Rappaport broth	20 Tubes x 9ml
42099	Selenite-F broth	20 Tubes x 9ml
44263	Selenite-F broth	50 x 15ml

Staphylococcus

44003	Baird Parker & RPF agar	4 x 100ml
44259	Brain Heart Infusion broth + 7% salt *	50 x 9ml
42080	Mannitol Salt agar (slopes)	20 Tubes

Trichomonas

44526	Diamonds medium	50 x 4.5ml
-------	-----------------	------------

Vibrio

44261	Alkaline Peptone water	50 x 10ml
-------	------------------------	-----------

Streptococci

42116 NEW	Todd Hewitt broth selective enrichment for detection of group B streptococci in pregnant women	20 Tubes x 9ml
------------------	--	----------------

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Yeasts & Moulds
 - Non Selective Media
 - Selective Media
 - Media for Biochemical Identification
 - Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
 - Media for Mycobacteria
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

Yeasts & Moulds

Non Selective Media

44524	Cornmeal agar *	50 x 20ml
42037	Sabouraud 2 agar	20 Tubes
42066	Sabouraud 2 agar	6 x 100ml
44253	Sabouraud Dextrose agar	10 x 250ml
42108	Sabouraud liquid broth	20 Tubes x 9ml

Selective Media

42084	Dermatophyte Test agar (slopes)	20 Tubes
42067	Sabouraud Chloramphenicol 2 agar	6 x 100ml
42038	Sabouraud Chloramphenicol 2 agar	20 Slopes
42094	Sabouraud Chloramphenicol Actidione agar (slopes)	20 Tubes
42031	Sabouraud Gentamicin Chloramphenicol 2 agar	20 Tubes
42056	Sabouraud Gentamicin Chloramphenicol 2 agar	6 x 100ml
42096	Sabouraud Tetrazolium Genta. Chloramphenicol agar (slopes)	20 Tubes
44525	Yeast Extract agar	10 x 100ml

Media for Biochemical Identification

42086	Esculin agar	20 Tubes
42087	Kligler iron agar (stab/slope)	20 Tubes
42085	Peptone water	20 Tubes x 9ml
44284	Urea agar	50 Agar Slopes
55752	Urea Indole medium	10 x 10ml

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
 - Chromogenic Media
 - General Purpose Media
 - Specific Media
 - Yeasts & Moulds
 - Media for Biochemical Identification
 - Media for Antibiotic Susceptibility Testing:
 - Disc Diffusion Methods
 - Media for Mycobacteria
 - Sputum digestion & decontamination
 - Staining
 - Culture
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
-----------	-------------	--------------

Pre-Poured Media – Tubes & Bottles

Media for Antibiotic Susceptibility Testing: Disc Diffusion Methods

41864	Mueller Hinton 2 agar	6 x 200ml
-------	-----------------------	-----------

Media for Mycobacteria

Sputum digestion & decontamination

42104	R2 : benzalkonium chloride	20 Tubes x 10ml
42103	R3 : Phosphate Buffer	20 Tubes x 15ml
55242	R4 : solution of bovine albumin	16 x 10ml

Staining

55531	Gabett solution	1 x 450ml
55521	Kinyoun solution	1 x 450ml

Culture

42082	Coletsos medium (slopes)	20 Tubes
42107	Lowenstein-Jensen + TCH (slopes)	20 Tubes
42089	Lowenstein-Jensen medium (slopes)	20 Tubes



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
 - Culture, Identification Enumeration & Susceptibility Testing
 - Enumeration
 - Detection
 - Collection & Culture
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Urogenital Mycoplasmas

Culture, Identification Enumeration & Susceptibility Testing

Enumeration

42505	Mycoplasma IST 2	25 Tests
-------	------------------	----------


Detection

42508	Urea-arginine LYO-2 broth	25 Tests
-------	---------------------------	----------

Collection & Culture

42507	Mycoplasma Preparation (collection + transport)	8 Tests
43003	A7 agar *	10 Plates

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas

Incubation System: GENbag-GENbox

- GENbag
 - Generators
 - Indicator
 - Accessories
- GENbox
 - Generators
 - Indicator

- Air Sampling
- Dehydrated Media
- Supplements



Reference	Description	Presentation
Incubation System: GENbag-GENbox		
GENbag		
Generators		
45534	GENbag Anaer	20 Sachets
45533	GENbag CO2	20 Sachets
45532	GENbag Microaer	20 Sachets
Indicator		
96118	GENbag/box Anaer indicator	50 Tests
Accessories		
45511	GENbag clip seals	Pack of 10
GENbox		
96127	GENbox 2.5L	12 Plate capacity
96128	GENbox 7L	42 Plate capacity
96129	GENbox Lid 2.5L	1 Lid
96130	GENbox Lid 7L	1 Lid
Generators		
96124	GENbox Anaer	10 Sachets
96126	GENbox CO2	10 Sachets
96125	GENbox Microaer	10 Sachets
Indicator		
96118	GENbag/box Anaer indicator	50 Tests

CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
- Supplements

Implement a Simple Traceability Solution in Your Facility

with *air*IDEAL® 3P™ Traceability & 3P Plates

- air*IDEAL 3P Traceability has all of the benefits as before, but now incorporates features to allow full traceability of all sampling operations
- Implement a simple traceability solution or compliment your existing electronic record keeping
- Complete audit trail that includes technician, sampler, individual agar plate and sample

The full traceability package comprises 3 components:

- air*IDEAL 3P Traceability air sampler unit with indicating LED light
- Remote control with embedded software and built-in 2D-barcode scanner (Datamatrix)
- 21 CFR Part 11 compliant and network-capable software for designing and following air sampling campaigns

ISO 14698 validation for physical and biological efficiency, H₂O₂ isolator compatibility and air circuit decontamination.



Reference	Description	Presentation
-----------	-------------	--------------

Air Sampling

410174	<i>air</i> IDEAL 3P Traceability (air sampler Contact Plates)	Instrument
410175	<i>air</i> IDEAL 3P Traceability (air sampler 90mm Plates)	Instrument
410173	<i>air</i> IDEAL RUID Remote Control	Instrument
410169	<i>air</i> IDEAL Traceability CD Software	Software
96304	Sampling Screen (Contact Plates)	One
96309	Sampling Screen (90mm)	One
96308	Tripod	One



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
 - Page 1 of 2
 - Page 2 of 2
- Supplements

Reference	Description	Presentation
Dehydrated Media		
51024	Azide dextrose broth (Rothe)	500g
51074	Baird Parker agar base	500g
51065	Bismuth Sulphite agar (Wilson Blair)	500g
51039	Blood agar base	500g
51009	Brain Heart Infusion	500g
51045	Brilliant Green agar	500g
51003	Brilliant Green bile broth	500g
51052	CLED agar	500g
51026	Columbia agar	500g
51025	D-Coccosel agar - Aesculia Bile agar	500g
51029	Desoxycholate citrate agar	500g
51027	Desoxycholate lactose agar	500g
51032	Endo agar	500g
51033	Eosin Methylene Blue agar	500g
51022	Ethyl violet azide broth (Litsky)	500g
51034	Eugon agar	500g
51013	Eugon broth	500g
51050	Hektoen agar	500g
51059	Kligler iron agar	500g
51014	Lactose broth	500g
51036	MacConkey agar without crystal violet	500g
51015	MacConkey broth	500g
51053	Mannitol Salt agar	500g
51075	Mueller Hinton 2 agar	500g
51035	Purple Lactose agar	500g
51077	Sabouraud agar	500g
51078	Sabouraud Chloramphenicol agar	500g
51063	Sabouraud liquid broth	500g
51043	Salmonella Shigella agar	500g
51076	Schaedler broth	500g
51017	Selenite-F broth	500g
51072	Standard methods agar (PCA)	500g



CULTURE MEDIA

- Transport Media & Dipslides
- Pre-Poured Media – Plates
- Pre-Poured Media – Plates, Tubes & Bottles
- Pre-Poured Media – Tubes & Bottles
- Urogenital Mycoplasmas
- Incubation System: GENbag-GENbox
- Air Sampling
- Dehydrated Media
 - Page 1 of 2
 - Page 2 of 2
- Supplements

Reference	Description	Presentation
-----------	-------------	--------------

Dehydrated Media

51018	Thioglycollate broth + resazurin	500g
51047	Triple sugar iron agar	500g
51044	Trypcase Soy agar	500g
51019	Trypcase Soy broth	500g
51048	TSN agar	500g
51005	Violet Red Bile Lactose agar	500g
51049	XLD agar	500g

Supplements

55691	Clostridium difficile mixture	4 x 2ml
55673	CNA mixture	8 x 2ml
55682	Kanamycin - Vancomycin mixture	8 x 1ml (Lyophilised)
55641	Legionella enrichment supplement	4 x 2ml
55645	Legionella selective mixture	4 x 2ml
55651	PolyViteX™ mixture	4 x 1ml + Diluent
55652	PolyViteX mixture	4 x 10ml + Diluent
55501	Potassium tellurite solution	6 x 5ml
55663	VCN mixture	8 x 1ml



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
 - Staining
 - Catalase, Oxidase
 - Other Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

Orientation Tests

Staining

55542	Color Gram 2, kit	4 x 240ml
55545	Color Gram 2, R1 Oxalate Crystal Violet	2 Litre Bottle
55546	Color Gram 2, R2 Lugol PVP	2 Litre Bottle
55547	Color Gram 2, R3 Decolorizer	2 Litre Bottle
55548	Color Gram 2, R4 Safranin	2 Litre Bottle
55531	Gabett solution	1 x 450ml
55521	Kinyoun solution	1 x 450ml

Catalase, Oxidase

55561	ID color Catalase	2 x 5ml
55635	Oxidase Reagent	50 Ampoules x 0.75ml
54991	Paper discs (for use with Oxidase Reagent)	50 discs

Other Tests

50110	API® OF medium	10 Ampoules x 5ml
55631	Kovacs reagent	1 x 25ml
55601	ONPG	30 discs
55912	Optochin test	2 x 30 discs
55501	Potassium tellurite solution	6 x 5ml
55181	Rabbit plasma	8 x 0.5ml
55182	Rabbit plasma	1 x 3.5ml



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- SLIDEX®: Latex Agglutination
- Staphylococcus
- Streptococcus
- Streptococcus pneumoniae
- Meningitis
- RAPIDEC®: Identification in 2 Hours
- Direct IF: Immunofluorescence Antibody
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

Rapid Antigen Detection & Identification

SLIDEX®: Latex Agglutination

Staphylococcus

73115	SLIDEX Staph-Plus	50 Tests
73116	SLIDEX Staph-Plus	250 Tests (Stirring Sticks not included)

Streptococcus


58811	SLIDEX Strepto-Plus	40-50 Tests
58818	SLIDEX Strepto-Plus A	40 Tests
58819	SLIDEX Strepto-Plus B	40 Tests

Streptococcus pneumoniae

58821	SLIDEX pneumo-Kit	50 Tests
-------	-------------------	----------

Meningitis

58831	SLIDEX meningite Strepto B	30 Tests
58803	SLIDEX meningite-Kit 5	25 Tests



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- SLIDEX®: Latex Agglutination
- RAPIDEC®: Identification in 2 Hours
- RAPIDEC staph
- RAPIDEC ur
- RAPIDEC Listeria monocytogenes
- Direct IF: Immunofluorescence Antibody
- Chlamydia
- RSV (Respiratory Syncytial Virus)
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

Rapid Antigen Detection & Identification

RAPIDEC®: Identification in 2 Hours

RAPIDEC staph

03300	RAPIDEC staph	30 Tests
70500	API® FVB	5ml Bottle
70700	Suspension medium (2ml)	100 Ampoules

RAPIDEC ur

03200	RAPIDEC ur	30 Tests
70700	Suspension medium (2ml)	100 Ampoules
70500	API FVB	5ml Bottle
70542	JAMES Reagent	2 x 5ml Ampoules

RAPIDEC Listeria monocytogenes

03700	RAPIDEC Lmono	20 Tests
-------	---------------	----------

Direct IF: Immunofluorescence Antibody

Chlamydia

55331	Chlamydia Direct IF prelevement (collection kit)	20 Kits
55321	Chlamydia Direct IF identification	50 Tests
75741	Direct IF slides	100 Slides
55311	Chlamydia Direct IF control +/-	10 Slides

RSV (Respiratory Syncytial Virus)

55341	RSV Direct IF identification	50 Tests
55351	RSV Direct IF control +/-	10 Slides

Rapid Antigen Detection of Clostridium difficile toxin with VIDAS® and mini VIDAS®, see [page 110](#).



IDENTIFICATION & SUSCEPTIBILITY TESTING

Orientation Tests
 Rapid Antigen Detection & Identification

API® & ID 32 Identification

Identification Software

API & RapiD Strips

Enterobacteriaceae

Page 1 of 2

Page 2 of 2

Non Enteric Gram Negative rods

Staphylococcus

Streptococcus

Neisseria & Haemophilus

Corynebacterium

Campylobacter

Anaerobes

Listeria

Bacillus

Lactobacillus

Yeasts

Enzymes

ID 32 Strips

Consumables & Accessories

Etest®

ATB™ Susceptibility Testing

VITEK® 2: Automated Method

VITEK 2 Compact: Automated Method

VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



Identification Software

40011	APIWEB®	Licence
40012	APIWEB Stand Alone	CD

API & RapiD Strips (page 1 of 7)

Enterobacteriaceae (page 1 of 2)



API 20 E

20100	API 20 E	25 Strips
20120	API 20 E reagent kit *	6 x 5ml Ampoules (w/o dropper-bottles)
	Kit contains:	
	JAMES Reagent	
	API NIT 1 & NIT 2	
	API VP 1 & VP 2	
	API TDA	
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
50110	API OF medium	10 Ampoules x 5ml
70402	API TDA	2 x 5ml Ampoules
70422	API VP 1 & VP 2	2 x 5ml VP 1; 2 x 5ml VP 2
70380	API Zn	2 x 10g
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
20230	NaCl 0.85% (5ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml
20150	Suspension medium (5ml)	100 Ampoules
411003	LyfoCults® Plus API 20 E Streamlined QC set ‡	2 Vials
410624	LyfoCults Plus API 20 E Comprehensive QC set ‡	10 Vials

* Standing orders are recommended. This product may have a 6 week lead time, please contact us to check for availability.

‡ See LyfoCults Plus [page 42](#) for strains contained within each set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
 - Identification Software
 - API & RapiD Strips
 - Enterobacteriaceae
 - Page 1 of 2
 - Page 2 of 2
 - Non Enteric Gram Negative rods
 - Staphylococcus
 - Streptococcus
 - Neisseria & Haemophilus
 - Corynebacterium
 - Campylobacter
 - Anaerobes
 - Listeria
 - Bacillus
 - Lactobacillus
 - Yeasts
 - Enzymes
 - ID 32 Strips
 - Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 2 of 7)

Enterobacteriaceae (page 2 of 2)

API 10 S

10100	API 10 S	50 Strips
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70402	API TDA	2 x 5ml Ampoules
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
20230	NaCl 0.85% (5ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml
20150	Suspension medium (5ml)	100 Ampoules

RapiD 20 E™ (Identification in 4 hours)

20701	RapiD 20 E	25 Strips
70422	API VP 1 & VP 2	2 x 5ml VP 1; 2 x 5ml VP 2
70542	JAMES Reagent	2 x 5ml Ampoules
20070	NaCl 0.85% medium (2ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml
20710	RapiD 20 E Suspension medium	25 x 1.25ml
410621	LyfoCults® Plus RapiD 20 E QC set ‡	6 Vials

API 50 CH

50300	API 50 CH	10 Strips
50430	API 50 CHB/E medium	10 Ampoules
70902	Ampoule guards (50 CH range)	Pack of 3
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules
20150	Suspension medium (5ml)	100 Ampoules

‡ See LyfoCults Plus [page 44](#) for strains contained within the set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Identification Software
- API & RapiD Strips
- Enterobacteriaceae
- Non Enteric Gram Negative rods
- Staphylococcus
- Streptococcus
- Neisseria & Haemophilus
- Corynebacterium
- Campylobacter
- Anaerobes
- Listeria
- Bacillus
- Lactobacillus
- Yeasts
- Enzymes
- ID 32 Strips
- Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 3 of 7)

Non Enteric Gram Negative rods

API 20 NE

20050	API 20 NE	25 Strips & media
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70380	API Zn	2 x 10g
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
20070	NaCl 0.85% medium (2ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml
411006	LyfoCults® Plus API 20 NE Streamlined QC set ‡	4 Vials
410626	LyfoCults Plus API 20 NE Comprehensive QC set ‡	8 Vials

Staphylococcus

API Staph

20500	API Staph	25 Strips & media
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70422	API VP 1 & VP 2	2 x 5ml VP 1; 2 x 5ml VP 2
70494	API ZYM A	2 x 8ml Ampoules
70493	API ZYM B	2 x 8ml Ampoules
70100	Mineral oil	125ml Bottle
411000	LyfoCults Plus API Staph Streamlined QC set ‡	2 Vials
410620	LyfoCults Plus API Staph Comprehensive QC set ‡	6 Vials

‡ See LyfoCults Plus [pages 43 and 44](#) for strains contained within each set.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Identification Software
- API & RapiD Strips
- Enterobacteriaceae
- Non Enteric Gram Negative rods
- Staphylococcus
- Streptococcus
- Neisseria & Haemophilus
- Corynebacterium
- Campylobacter
- Anaerobes
- Listeria
- Bacillus
- Lactobacillus
- Yeasts
- Enzymes
- ID 32 Strips
- Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 4 of 7)



Streptococcus

API 20 Strep

20600	API 20 Strep	25 Strips & media
70491	API NIN	2 x 5ml Bottle
70422	API VP 1 & VP 2	2 x 5ml VP 1; 2 x 5ml VP 2
70494	API ZYM A	2 x 8ml Ampoules
70493	API ZYM B	2 x 8ml Ampoules
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules
411002	LyfoCults® Plus API 20 Strep Streamlined QC set ‡	2 Vials
410590	LyfoCults Plus API 20 Strep Comprehensive QC set ‡	4 Vials

Neisseria & Haemophilus

API NH

10400	API NH	10 Strips & media
70100	Mineral oil	125ml Bottle
411005	LyfoCults Plus API NH Streamlined QC set ‡	2 Vials
410622	LyfoCults Plus API NH Comprehensive QC set ‡	6 Vials

Corynebacterium

API Coryne

20900	API Coryne	12 Strips & media
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70492	API PYZ	2 x 5ml Bottle
70494	API ZYM A	2 x 8ml Ampoules
70493	API ZYM B	2 x 8ml Ampoules
70100	Mineral oil	125ml Bottle
411001	LyfoCults Plus API Coryne Streamlined QC set ‡	2 Vials
410625	LyfoCults Plus API Coryne Comprehensive QC set ‡	6 Vials

‡ See LyfoCults Plus [pages 41 to 43](#) for strains contained within each set.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Identification Software
- API & RapiD Strips
- Enterobacteriaceae
- Non Enteric Gram Negative rods
- Staphylococcus
- Streptococcus
- Neisseria & Haemophilus
- Corynebacterium
- Campylobacter
- Anaerobes
- Listeria
- Bacillus
- Lactobacillus
- Yeasts
- Enzymes
- ID 32 Strips
- Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 5 of 7)

Campylobacter

API Campy

20800	API Campy	12 Strips & media
70562	API FB	2 x 5ml Ampoules
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70100	Mineral oil	125ml Bottle
55635	Oxidase Reagent	50 Ampoules x 0.75ml

Anaerobes

API 20 A

20300	API 20 A	25 Strips & media
70510	API BCP	5ml Bottle
70520	API EHR	5ml Bottle
70530	API XYL	5ml Bottle
70100	Mineral oil	125ml Bottle
410614	LyfoCults® Plus API 20 A QC set ‡	6 Vials

Listeria

API Listeria

10300	API Listeria	10 Strips & media
-------	--------------	-------------------

‡ See LyfoCults Plus [page 41](#) for strains contained within the set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Identification Software
- API & RapiD Strips
- Enterobacteriaceae
- Non Enteric Gram Negative rods
- Staphylococcus
- Streptococcus
- Neisseria & Haemophilus
- Corynebacterium
- Campylobacter
- Anaerobes
- Listeria
- Bacillus
- Lactobacillus
- Yeasts
- Enzymes
- ID 32 Strips
- Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 6 of 7)

Bacillus

API 50 CH

50300	API 50 CH	10 Strips
50430	API 50 CHB/E medium	10 Ampoules
70902	Ampoule guards (50 CH range)	Pack of 3
70100	Mineral oil	125ml Bottle

Lactobacillus

API 50 CH

50300	API 50 CH	10 Strips
50410	API 50 CHL medium	10 Ampoules
70902	Ampoule guards (50 CH range)	Pack of 3
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules
20150	Suspension medium (5ml)	100 Ampoules



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Identification Software
- API & RapiD Strips
- Enterobacteriaceae
- Non Enteric Gram Negative rods
- Staphylococcus
- Streptococcus
- Neisseria & Haemophilus
- Corynebacterium
- Campylobacter
- Anaerobes
- Listeria
- Bacillus
- Lactobacillus
- Yeasts
- Enzymes
- ID 32 Strips
- Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



API & RapiD Strips (page 7 of 7)

Yeasts

API 20 C AUX

20210	API 20 C AUX	25 Strips & media
70700	Suspension medium (2ml)	100 Ampoules
411004	LyfoCults® Plus API 20 C AUX Streamlined QC set ‡	4 Vials
410589	LyfoCults Plus API 20 C AUX Comprehensive QC set ‡	6 Vials

API Candida


10500	API Candida	10 Strips & media
-------	-------------	-------------------

Enzymes

API ZYM

25200	API ZYM	25 Strips
70494	API ZYM A	2 x 8ml Ampoules
70493	API ZYM B	2 x 8ml Ampoules

‡ See LyfoCults Plus [page 41](#) for strains contained within each set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
 - Identification Software
 - API & RapiD Strips
 - ID 32 Strips
 - Enterobacteriaceae
 - Staphylococcus
 - Streptococcus
 - Anaerobes
 - Yeasts
 - Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



ID 32 Strips (page 1 of 2)

Enterobacteriaceae

ID 32 E

32400	ID 32 E	25 Strips
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
20070	NaCl 0.85% medium (2ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml

rapid ID 32 E (Identification in 4 hours)

32700	rapid ID 32 E	25 Strips
70510	API BCP (for automatic reading only)	5ml Bottle
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
20070	NaCl 0.85% medium (2ml)	100 Ampoules
55635	Oxidase Reagent	50 Ampoules x 0.75ml
410619	LyfoCults® Plus rapid ID 32 E QC set ‡	6 Vials

Staphylococcus

ID 32 Staph

32500	ID 32 STAPH	25 Strips
70562	API FB	2 x 5ml Ampoules
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70572	API VP A & VP B	1 x 5ml VP A; 1 x 5ml VP B
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules

‡ See LyfoCults Plus [page 44](#) for strains contained within the set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
 - Identification Software
 - API & RapiD Strips
 - ID 32 Strips
 - Enterobacteriaceae
 - Staphylococcus
 - Streptococcus
 - Anaerobes
 - Yeasts
 - Consumables & Accessories
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

API® & ID 32 Identification



ID 32 Strips (page 2 of 2)

Streptococcus

rapid ID 32 STREP (Identification in 4 hours)

32600	rapid ID 32 STREP	25 Strips
70562	API FB	2 x 5ml Ampoules
70491	API NIN	2 x 5ml Bottle
70572	API VP A & VP B	1 x 5ml VP A; 1 x 5ml VP B
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules

Anaerobes

rapid ID 32 A (Identification in 4 hours)

32300	rapid ID 32 A	25 Strips
70562	API FB	2 x 5ml Ampoules
70442	API NIT 1 & NIT 2	2 x 5ml NIT 1; 2 x 5ml NIT 2
70542	JAMES Reagent	2 x 5ml Ampoules
70100	Mineral oil	125ml Bottle
70700	Suspension medium (2ml)	100 Ampoules

Yeasts

32200	ID 32 C	25 Strips & media
70700	Suspension medium (2ml)	100 Ampoules

Consumables & Accessories

70901	Ampoule guards (API/ID 32 range)	Pack of 3
70902	Ampoule guards (50 CH range)	Pack of 3
99234	Densimat	One
15512	Densimat Control Kit (for use with Densimat 99234)	5 Ampoules
99597	Lithium batteries (for use with Densimat 99234)	Two
70900	McFarland Standard Set	6 Ampoules
70250	PSIpettes (5ml) sterile	Pack of 400
70200	Stand for 12 ampoules	12 Ampoule Stand
70610	Sterile cotton swabs	Pack of 100
70760	Sterile loops	50 x 10

IDENTIFICATION & SUSCEPTIBILITY TESTING

Orientation Tests

Rapid Antigen Detection & Identification

API® & ID 32 Identification

- Etest®
- Antimicrobial Susceptibility Testing

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4

Antifungal Susceptibility Testing

Antimycobacterial Susceptibility Testing

Confirmation of Extended Spectrum

Beta-Lactamase (ESBL)

Detection of Metallo Beta-Lactamase (MBL)

Glycopeptide Resistance Detection (GRD)

AmpC Detection

Cefinase™ for Detection of Beta-Lactamase

Accessories

Biotoools™

Biotoools Accessories

ATB™ Susceptibility Testing

VITEK® 2: Automated Method

VITEK 2 Compact: Automated Method

VITEK MS: MALDI-TOF Mass Spectrometer



Description

Etest®



Etest®

Antimicrobial Susceptibility Testing (page 1 of 4)

	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack	30 Strips Foil Pack	100 Strips Foil Pack
Amikacin (AK 0.016-256) NEW	—	—	412219	—	501318
Amoxicillin (ACL 0.002-32 (low))	—	532608	—	—	—
Amoxicillin (AC 0.016-256 (high)) NEW	—	—	412243	—	500918
Amoxicillin-clavulanic acid (2/1) (XL 0.016-256) NEW	—	—	412241	—	501018
Ampicillin (AM 0.016-256) NEW	—	—	412253	—	501518
Ampicillin-sulbactam (2/1) (AB 0.016-256) AVAILABLE 2013	—	—	412251	—	501818
Azithromycin (AZ 0.016-256) AVAILABLE 2013	—	—	412257	—	501618
Aztreonam (AT 0.016-256) AVAILABLE 2013	—	—	412259	—	501718
Bacitracin (BA 0.016-256) AVAILABLE 2013	—	—	412261	—	—
Benzylpenicillin (PGL 0.002-32 (low))	—	—	412265	—	502618
Benzylpenicillin (PG 0.016-256 (high))	—	—	412263	—	502518
Cefaclor (CF 0.016-256)	504500	504508	—	504510	504518
Cefepime (PM 0.016-256) AVAILABLE 2013	—	—	412273	—	505018
Cefixime (IX 0.016-256) AVAILABLE 2013	—	—	412275	—	529918
Cefoperazone/sulbactam (2/1) (CPS 0.016-256) AVAILABLE 2013	—	—	412277	—	529318
Cefotaxime (CTL 0.002-32 (low)) NEW	—	—	412281	—	505618
Cefotaxime (CT 0.016-256 (high)) NEW	—	—	412279	—	505518
Cefotetan (CN 0.016-256) AVAILABLE 2013	—	—	412283	—	—
Cefoxitin (FX 0.016-256) AVAILABLE 2013	—	—	412285	—	506518
Cefpirome (CR 0.016-256)	506400	506408	—	—	—
Cefpodoxime (PX 0.016-256) AVAILABLE 2013	—	—	412289	—	505818

Blister packaging (storage at max. -20 °C).

Single Push Blister packaging (storage at max. -20 °C to +8 °C).

Foam packaging (storage at max. +8 °C and +22 °C).

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

Antimicrobial Susceptibility Testing

Page 1 of 4

Page 2 of 4

Page 3 of 4

Page 4 of 4

Antifungal Susceptibility Testing

Antimycobacterial Susceptibility Testing

Confirmation of Extended Spectrum

Beta-Lactamase (ESBL)

Detection of Metallo Beta-Lactamase (MBL)

Glycopeptide Resistance Detection (GRD)

AmpC Detection

Cefinase™ for Detection of Beta-Lactamase

Accessories

Biotoools™

Biotoools Accessories

ATB™ Susceptibility Testing

VITEK® 2: Automated Method

VITEK 2 Compact: Automated Method

VITEK MS: MALDI-TOF Mass Spectrometer

Description



Etest®

Etest®

Antimicrobial Susceptibility Testing (page 2 of 4)

	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack	30 Strips Foil Pack	100 Strips Foil Pack
Ceftaroline (CPT 0.002-32) AVAILABLE 2013	537500	537508	412291	—	—
Ceftazidime (TZ 0.016-256) NEW	—	—	412293	—	506718
Ceftizoxime (CZ 0.002-32)	527300	527308	—	—	—
Ceftrobiprole (BPRL 0.002-32 (low))	536000	536008	—	—	—
Ceftrobiprole (BPR 0.016-256 (high))	535300	535308	—	—	—
Ceftriaxone (TXL 0.002-32 (low)) NEW	—	—	412303	—	507018
Ceftriaxone (TX 0.016-256 (high)) AVAILABLE 2013	—	—	412301	—	506618
Cefuroxime (XM 0.016-256) AVAILABLE 2013	—	—	412305	—	506918
Cephalothin (CE 0.016-256) AVAILABLE 2013	—	—	412307	—	503518
Chloramphenicol (CL 0.016-256) AVAILABLE 2013	—	—	412309	—	507518
Ciprofloxacin (CI 0.002-32) NEW	—	—	412311	—	508618
Clarithromycin (CH 0.016-256)	508700	508708	—	508710	508718
Clindamycin (CM 0.016-256) NEW	—	—	412315	—	509518
Colistin (CO 0.016-256) AVAILABLE 2013	—	—	412317	—	—
Dalbavancin (DLC 0.002-32)	—	535208	—	—	—
Daptomycin (DPC 0.016-256) AVAILABLE 2013	—	—	412324	—	535018
Doripenem (DOR 0.002-32) AVAILABLE 2013	—	—	412326	—	535918
Doxycycline (DC 0.016-256) AVAILABLE 2013	—	—	412328	—	509718
Enrofloxacin (EF 0.002-32)	528900	528908	—	—	—
Ertapenem (ETP 0.002-32) AVAILABLE 2013	—	—	412332	—	531618
Erythromycin (EM 0.016-256) NEW	—	—	412334	—	510518

Blister packaging (storage at max. -20 °C).

Single Push Blister packaging (storage at max. -20 °C to +8 °C).

Foam packaging (storage at max. +8 °C and +22 °C).



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

Antimicrobial Susceptibility Testing

- Page 1 of 4
- Page 2 of 4
- Page 3 of 4
- Page 4 of 4

- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotoools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Description

Etest®



Etest®

Antimicrobial Susceptibility Testing (page 3 of 4)

Description	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack	30 Strips Foil Pack	100 Strips Foil Pack
Fosfomycin (FM 0.064-1024)	529100	529108	—	—	—
Fusidic Acid (FU 0.016-256) AVAILABLE 2013	—	—	412358	—	511518
Gatifloxacin (GA 0.002-32)	530200	530208	—	530210	530218
Gemifloxacin (GEM 0.002-32) AVAILABLE 2013	—	—	412364	—	530318
Gentamicin (GML 0.016-256 (low)) NEW	—	—	412368	—	512518
Gentamicin (GM 0.064-1024 (high)) AVAILABLE 2013	—	—	412366	—	—
Imipenem (IP 0.002-32)	—	—	412374	—	513618
Kanamycin (KM 0.016-256)	527800	527808	—	527810	527818
Levofloxacin (LE 0.002-32) NEW	—	—	412393	—	527418
Linezolid (LZ 0.016-256) NEW	—	—	412396	—	531318
Mecillinam (MM 0.016-256)	513700	513708	—	—	—
Meropenem (MP 0.002-32)	—	—	412402	—	513818
Metronidazole (MZ 0.016-256) NEW	—	—	412404	—	530018
Minocycline (MC 0.016-256) AVAILABLE 2013	—	—	412409	—	516018
Moxifloxacin (MX 0.002-32) AVAILABLE 2013	—	—	412411	—	529018
Mupirocin (MU 0.064-1024) AVAILABLE 2013	—	—	412417	—	—
Nalidixic Acid (NA 0.016-256)	516500	516508	—	—	—
Netilmicin (NC 0.016-256)	517500	517508	—	517510	517518
Nitrofurantoin (NI 0.032-512) AVAILABLE 2013	—	—	412426	—	—
Norfloxacin (NX 0.016-256) AVAILABLE 2013	—	—	412428	—	—
Ofloxacin (OF 0.002-32)	519600	519608	—	519610	519618

- Blister packaging (storage at max. -20 °C).
- Single Push Blister packaging (storage at max. -20 °C to +8 °C).
- Foam packaging (storage at max. +8 °C and +22 °C).

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

Antimicrobial Susceptibility Testing

- Page 1 of 4
- Page 2 of 4
- Page 3 of 4
- Page 4 of 4

- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Description

Etest®



Etest®

Antimicrobial Susceptibility Testing (page 4 of 4)

Description	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack	30 Strips Foil Pack	100 Strips Foil Pack
Oxacillin (OX 0.016-256) NEW	—	—	412432	—	520518
Piperacillin (PP 0.016-256) AVAILABLE 2013	—	—	412436	—	521518
Piperacillin-tazobactam (4 mg/L) (PTc 0.016-256) NEW	—	—	412434	—	521418
Polymyxin B (PO 0.064-1024) AVAILABLE 2013	—	—	412438	—	—
Quinupristin-dalfopristin (QDA 0.002-32)	528700	528708	—	528710	528718
Rifampicin (RI 0.002-32) AVAILABLE 2013	—	—	412450	—	526018
Spectinomycin (SC 0.064-1024) AVAILABLE 2013	—	—	412452	—	529218
Streptomycin (SM 0.064-1024) AVAILABLE 2013	—	—	412454	—	—
Sulbactam (SUL 0.016-256)	535100	535108	—	—	—
Sulfamethoxazole (SX 0.064-1024) AVAILABLE 2013	—	—	412458	—	534118
Teicoplanin (TP 0.016-256) NEW	—	—	412461	—	522018
Telavancin (TLV 0.002-32)	534500	534508	—	—	—
Temocillin (TMO 0.064-1024) *	535600	535608	—	—	—
Tetracycline (TC 0.016-256) NEW	—	—	412471	—	522518
Ticarcillin-clavulanic acid (2 mg/L) (TLc 0.016-256) AVAILABLE 2013	—	—	412473	—	522618
Tigecycline (TGC 0.016-256) NEW	—	—	412475	—	533518
Tobramycin (TML 0.016-256 (low)) AVAILABLE 2013	—	—	412479	—	522718
Tobramycin (TM 0.064-1024 (high))	533100	533108	—	—	—
Trimethoprim (TR 0.002-32) AVAILABLE 2013	—	—	412483	—	523618
Trimethoprim-sulphamethoxazole (1/19) (TS 0.002-32) AVAILABLE 2013	—	—	412481	—	524418
Vancomycin (VA 0.016-256)	—	—	412488	—	525518

- Blister packaging (storage at max. -20 °C).
- Single Push Blister packaging (storage at max. -20 °C to +8 °C).
- Foam packaging (storage at max. +8 °C and +22 °C).
- * For investigational use only, not for clinical use.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

- Antimicrobial Susceptibility Testing
- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotoools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Description	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack	30 Strips Foil Pack	100 Strips Foil Pack
-------------	---------------------------	----------------------------	--	------------------------	-------------------------

Etest®



Antifungal Susceptibility Testing

Amphotericin B (AP 0.002-32) AVAILABLE 2013	—	—	412249	—	526318
Anidulafungin (AND 0.002-32)	532000	532008	—	—	—
Caspofungin (CS 0.002-32) AVAILABLE 2013	—	—	412269	—	532418
Fluconazole (FL 0.016-256) AVAILABLE 2013	—	—	412350	—	510818
Flucytosine (FC 0.002-32) AVAILABLE 2013	—	—	412352	—	510918
Itraconazole (IT 0.002-32) AVAILABLE 2013	—	—	412380	—	525818
Ketoconazole (KE 0.002-32)	525900	525908	—	525910	525918
Micafungin (MYC 0.002-32)	535700	535708	—	—	—
Posaconazole (POS 0.002-32) AVAILABLE 2013	—	—	412440	—	532118
Voriconazole (VO 0.002-32) AVAILABLE 2013	—	—	412490	—	532818

Antimycobacterial Susceptibility Testing

Ethambutol (EB 0.016-256)	527700	—	—	—	—
Ethionamide (ET 0.016-256)	527500	—	—	—	—
Isoniazide (IZ 0.016-256)	527900	—	—	—	—

- Blister packaging (storage at max. -20 °C).
- Single Push Blister packaging (storage at max. -20 °C to +8 °C).
- Foam packaging (storage at max. +8 °C and +22 °C).

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

- Antimicrobial Susceptibility Testing
- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotoools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer





Etest®

Description	30 Strips Blister Pack	100 Strips Blister Pack	30 Strips Single Push Blister Pack
Confirmation of Extended Spectrum Beta-Lactamase (ESBL)			
ESBL - Cefepime/Cefepime + Clavulanic acid (4 mg/L) AVAILABLE 2013 (PM/PML 0.25-16/0.064-4)	—	—	412338
ESBL - Cefotaxime/Cefotaxime + Clavulanic acid (4 mg/L) AVAILABLE 2013 (CT/CTL 0.25-16/0.016-1)	—	—	412336
ESBL - Ceftazidime/Ceftazidime + Clavulanic acid (4 mg/L) AVAILABLE 2013 (TZ/TZL 0.5-32/0.064-4)	—	—	412340
Detection of Metallo Beta-Lactamase (MBL)			
MBL - Imipenem/Imipenem + EDTA (IP/IPI 256/1-64)	534200	534208	—
MBL - Meropenem/Meropenem + EDTA (MP/MPI 8/0.032-2)	411361	411362	—
Glycopeptide Resistance Detection (GRD)			
GRD - Vancomycin/Teicoplanin (VA/TP 0.5-32/0.5-32)	537200	537208	—
AmpC Detection			
AmpC - Cefotetan/Cefotetan + Cloxacillin (CN/CNI 0.5-32/0.5-32) *	537100	537108	—

Reference	Description	Presentation
55622	Cefinase	Pack of 50

Blister packaging (storage at max. -20 °C).

Single Push Blister packaging (storage at max. -20 °C to +8 °C).

* For investigational use only, not for clinical use.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

- Antimicrobial Susceptibility Testing
- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotoools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
Etest®		Etest®
Accessories		
559807	Applicator (1 applicator and 5 strips of tape)	One
559806	Applicator tape	Pack of 5
559809	Clamps - For Resealing Etest Package	Pack of 5
559709	Desiccant Capsules	Pack of 8
559904	Measuring Sticks	Pack of 5
411200	MiniGrip Vacuum Pen	One
559801	Storage Box, Maxi (12 trays, 12 desiccant strips, 1 applicator, 1 strip of tape and an instruction manual)	One
559902	Storage Box, Maxi (12 desiccant strips, 1 applicator, 1 strip of tape and an instruction manual)	One
559901	Storage Box, Midi (6 desiccant strips, 1 applicator, 1 strip of tape and an instruction manual)	One
559900	Storage Box, Mini (2 desiccant strips, 1 applicator, 1 strip of tape and an instruction manual)	One
501603	Storage Container	One

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

- Antimicrobial Susceptibility Testing
- Antifungal Susceptibility Testing
- Antimycobacterial Susceptibility Testing
- Confirmation of Extended Spectrum Beta-Lactamase (ESBL)
- Detection of Metallo Beta-Lactamase (MBL)
- Glycopeptide Resistance Detection (GRD)
- AmpC Detection
- Cefinase™ for Detection of Beta-Lactamase
- Accessories
- Biotoools™
- Biotoools Accessories

- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

Etest®

Etest®

Biotoools™

559802	Simplex C76 - For Automatic Application of Etest Strips	One
559804	Nema C88 - Vacuum Pen For Application of Etest Strips	One
559803	Retro C80 - For Inoculation of Agar Plates	One

Biotoools Accessories

5100700700	Biotoools AC Adapter	One
5100700000	Biotoools Cord Europe	One
5100700200	Biotoools Petridish Holder 90mm	One
5100700300	Biotoools Petridish Holder 150mm	One
5100700500	Nema C88 Vacuum Pen Set	One
5100599300	Nema/Simplex Suction Cups	Pack of 2
5100700400	Simplex C76 Tray Holder	One



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®

ATB™ Susceptibility Testing

- ATB Galleries
 - Gram Negative rods
 - Pseudomonas & non-Fermenters
 - Urines
 - Staphylococcus
 - Streptococcus
 - Enterococcus
 - Haemophilus
 - Anaerobes
 - Yeasts
 - Veterinary Isolates

Consumables & Accessories

- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

ATB™ Susceptibility Testing

ATB Galleries (page 1 of 2)

Gram Negative rods

14315	ATB G- 5 ◇	25 Strips
14960	ATB medium	100 Ampoules
70700	Suspension medium (2ml)	100 Ampoules
410997	LyfoCults® Plus ATB Gram Negative QC set ‡	6 Vials

Pseudomonas & non-Fermenters

14345	ATB PSE 5 ◇	25 Strips
14960	ATB medium	100 Ampoules
70700	Suspension medium (2ml)	100 Ampoules

Urines

14335	ATB UR 5 ◇	25 Strips
-------	------------	-----------

Staphylococcus

14325	ATB STAPH 5 ◇	25 Strips & ATB Na medium
14960	ATB medium	100 Ampoules
70700	Suspension medium (2ml)	100 Ampoules

Streptococcus

14355	ATB STREP 5 (Streps & Pneumos ONLY) ◇	10 Strips & ATB S medium
70700	Suspension medium (2ml)	100 Ampoules

Enterococcus

14365	ATB ENTEROC 5 (Enterococci ONLY) ◇	25 Strips & ATB S medium
70700	Suspension medium (2ml)	100 Ampoules

‡ See LyfoCults Plus [page 45](#) for strains contained within the set.

◇ For full details of strip compositions please call your local Product Specialist.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
 - ATB Galleries
 - Gram Negative rods
 - Pseudomonas & non-Fermenters
 - Urines
 - Staphylococcus
 - Streptococcus
 - Enterococcus
 - Haemophilus
 - Anaerobes
 - Yeasts
 - Veterinary Isolates
 - Consumables & Accessories

- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

ATB™ Susceptibility Testing

ATB Galleries (page 2 of 2)

Haemophilus

14299	ATB HAEMO (Haemophilus + Branhamella ONLY)	10 Strips + ATB S medium
20070	NaCl 0.85% medium (2ml)	100 Ampoules
410996	LyfoCults® Plus ATB HAEMO QC set ‡	2 Vials

Anaerobes

14269	ATB ANA ◊	10 Strips and ATB S medium
20070	NaCl 0.85% medium (2ml)	100 Ampoules
410999	LyfoCults Plus ATB ANA QC set ‡	2 Vials

Yeasts

14204	ATB FUNGUS 3	25 Strips and ATB F2 medium
20070	NaCl 0.85% medium (2ml)	100 Ampoules
410998	LyfoCults Plus ATB FUNGUS QC set ‡	4 Vials

Veterinary Isolates

14289	ATB VET (G+ and G- veterinary isolates) ◊	25 Strips
14960	ATB medium	100 Ampoules
14931	ATB S medium	10 Ampoules
70700	Suspension medium (2ml)	100 Ampoules

‡ See LyfoCults Plus [page 45](#) for strains contained within each set.

◊ For full details of strip compositions please call your local Product Specialist.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
 - ATB Galleries
 - Gram Negative rods
 - Pseudomonas & non-Fermenters
 - Urines
 - Staphylococcus
 - Streptococcus
 - Enterococcus
 - Haemophilus
 - Anaerobes
 - Yeasts
 - Veterinary Isolates
- Consumables & Accessories
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
ATB™ Susceptibility Testing		
Consumables & Accessories		
70901	Ampoule guards (API/ID 32 range)	Pack of 3
70902	Ampoule guards (50 CH range)	Pack of 3
14960	ATB medium	100 Ampoules
14931	ATB S medium	10 Ampoules
70700	Suspension medium (2ml)	100 Ampoules
70640	ATB Suspension medium (3ml)	100 Ampoules
99234	Densimat	One
15512	Densimat Control Kit (for use with Densimat 99234)	5 Ampoules
99597	Lithium batteries (for use with Densimat 99234)	Two
70900	McFarland Standard Set	6 Ampoules
70250	PSIpettes (5ml) sterile	Pack of 400
70200	Stand for 12 ampoules	12 Ampoule Stand
70610	Sterile cotton swabs	Pack of 100
70760	Sterile loops	50 x 10



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
 - VITEK 2 Analysers
 - Myla™ Middleware Software
 - VITEK 2 Identification Cards
 - VITEK 2 Antibiotic Susceptibility Cards
 - VITEK 2 Other Consumables & Accessories
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

VITEK® 2: Automated Method



VITEK 2 Analysers

27226	VITEK 2 Configuration	Instrument
27228	VITEK 2 XL Configuration	Instrument
27202	VITEK 2 additional module	Instrument
413641	OBSERVA® Data Management System	One

Myla™ Middleware Software

411401 NEW	Myla Middleware	One
411613 NEW	Myla VITEK 2 Connection Kit	1 Kit
410152 NEW	Myla Multiuser Licence	One Licence



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
 - VITEK 2 Analysers
 - Myla™ Middleware Software
 - VITEK 2 Identification Cards
 - Gram Negatives (Enterobacteriaceae / non-fermenters)
 - Gram Positive
 - Neisseria, Haemophilus & Fastidious Bacteria
 - Anaerobes
 - Yeasts
 - VITEK 2 Antibiotic Susceptibility Cards
 - VITEK 2 Other Consumables & Accessories
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

VITEK® 2: Automated Method



VITEK 2 Identification Cards

Gram Negatives (Enterobacteriaceae / non-fermenters)

21341	VITEK 2 GN	20 Cards
410586	LyfoCults® Plus VITEK 2 GN Streamlined QC set ‡	4 Vials
410601	LyfoCults Plus VITEK 2 GN Comprehensive QC set ‡	20 Vials

Gram Positive

21342	VITEK 2 GP	20 Cards
410616	LyfoCults Plus VITEK 2 GP Streamlined QC set ‡	4 Vials
410584	LyfoCults Plus VITEK 2 GP Comprehensive QC set ‡	18 Vials

Neisseria, Haemophilus & Fastidious Bacteria

21346	VITEK 2 NH	20 Cards
410628	LyfoCults Plus VITEK 2 NH Streamlined QC set ‡	2 Vials
410611	LyfoCults Plus VITEK 2 NH Comprehensive QC set ‡	18 Vials


Anaerobes

21347	VITEK 2 ANC	20 Cards
410604	LyfoCults Plus VITEK 2 ANC Streamlined QC set ‡	4 Vials
410605	LyfoCults Plus VITEK 2 ANC Comprehensive QC set ‡	14 Vials

Yeasts

21343	VITEK 2 YST	20 Cards
410612	LyfoCults Plus VITEK 2 YST Streamlined QC set ‡	2 Vials
410608	LyfoCults Plus VITEK 2 YST Comprehensive QC set ‡	22 Vials

‡ See LyfoCults Plus [pages 35 to 39](#) for strains contained within each set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
 - VITEK 2 Analysers
 - Myla™ Middleware Software
 - VITEK 2 Identification Cards
 - VITEK 2 Antibiotic Susceptibility Cards
 - Gram Negative rods
 - Gram Positive cocci
 - Yeasts
 - VITEK 2 Other Consumables & Accessories
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
VITEK® 2: Automated Method VITEK 2 Antibiotic Susceptibility Cards Gram Negative rods  		
Nordic		
412938	VITEK 2 AST-N209	20 Cards
413080	VITEK 2 AST-N217	20 Cards
413081	VITEK 2 AST-N218	20 Cards
413083	VITEK 2 AST-N222	20 Cards
413147	VITEK 2 AST-N230	20 Cards
410609	LyfoCults® Plus VITEK 2 AST-GN QC set ‡	8 Vials
Gram Positive cocci		
Nordic		
22294	VITEK 2 AST-P594 (Replaced by 414576 by end of 2012)	20 Cards
414576 NEW	VITEK 2 AST-P629	20 Cards
410610	LyfoCults Plus VITEK 2 AST-GP QC set ‡	16 Vials
Global		
22216	VITEK 2 AST-P576	20 Cards
22233	VITEK 2 AST-P580	20 Cards
22276	VITEK 2 AST-P586	20 Cards
22287	VITEK 2 AST-P592	20 Cards
XXXXXX	VITEK 2 AST-P5XXX	20 Cards
410028 NEW	VITEK 2 AST-ST01 (Standard Streptococcus Card)	20 Cards
410610	LyfoCults Plus VITEK 2 AST-GP QC set ‡	16 Vials
Yeasts		
412610 NEW	VITEK 2 AST-YS06	20 Cards
410613	LyfoCults Plus VITEK 2 AST-YS QC set ‡	4 Vials

‡ See LyfoCults Plus [page 40](#) for strains contained within each set.

IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
 - VITEK 2 Analysers
 - Myla™ Middleware Software
 - VITEK 2 Identification Cards
 - VITEK 2 Antibiotic Susceptibility Cards
 - VITEK 2 Other Consumables & Accessories
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
VITEK® 2: Automated Method VITEK 2 Other Consumables & Accessories VITEK 2 Other Consumables & Accessories		
V1204	Saline 0.45%	3 x 500ml
V1200	Saline Dispenser	One
70610	Sterile cotton swabs	Pack of 100
69285	Test tubes (uncapped)	2000
21218	Saline reservoir bags	12 x 1000ml
21219	Pipetter / Diluter accessory kit	330 Pipette Straws & 1 Tube
27209	Smart Carrier Station (SCS)	One
27700	Smart Carrier Cassette	One cassette with memory chip
21250	DensiCHEK™ Plus	One
21255	DensiCHEK Plus - Standard Reference Kit	4 Vials
93059	DensiCHEK - Standard Reference	1 Vial
96128	VITEK 2 Smart Box 7L	One



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
 - VITEK 2 Compact Analyser
 - VITEK 2 Compact Identification Cards
 - VITEK 2 Compact Antibiotic Susceptibility Cards
 - VITEK 2 Compact Other Consumables & Accessories
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
VITEK® 2 Compact: Automated Method VITEK 2 Compact Analyser 27560 VITEK 2 Compact 60 Instrument 413641 OBSERVA® Data Management System One		



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
 - VITEK 2 Compact Analyser
 - VITEK 2 Compact Identification Cards
 - Gram Negatives (Enterobacteriaceae / non-fermenters)
 - Gram Positive
 - Neisseria, Haemophilus & Fastidious Bacteria
 - Anaerobes
 - Yeasts
 - VITEK 2 Compact Antibiotic Susceptibility Cards
 - VITEK 2 Compact Other Consumables & Accessories
- VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
-----------	-------------	--------------

VITEK® 2 Compact: Automated Method



VITEK 2 Compact Identification Cards

Gram Negatives (Enterobacteriaceae / non-fermenters)

21341	VITEK 2 GN	20 Cards
410586	LyfoCults® Plus VITEK 2 GN Streamlined QC set ‡	4 Vials
410601	LyfoCults Plus VITEK 2 GN Comprehensive QC set ‡	20 Vials

Gram Positive

21342	VITEK 2 GP	20 Cards
410616	LyfoCults Plus VITEK 2 GP Streamlined QC set ‡	4 Vials
410584	LyfoCults Plus VITEK 2 GP Comprehensive QC set ‡	18 Vials

Neisseria, Haemophilus & Fastidious Bacteria

21346	VITEK 2 NH	20 Cards
410628	LyfoCults Plus VITEK 2 NH Streamlined QC set ‡	2 Vials
410611	LyfoCults Plus VITEK 2 NH Comprehensive QC set ‡	18 Vials


Anaerobes

21347	VITEK 2 ANC	20 Cards
410604	LyfoCults Plus VITEK 2 ANC Streamlined QC set ‡	4 Vials
410605	LyfoCults Plus VITEK 2 ANC Comprehensive QC set ‡	14 Vials

Yeasts

21343	VITEK 2 YST	20 Cards
410612	LyfoCults Plus VITEK 2 YST Streamlined QC set ‡	2 Vials
410608	LyfoCults Plus VITEK 2 YST Comprehensive QC set ‡	22 Vials

‡ See LyfoCults Plus [pages 35 to 39](#) for strains contained within each set.



IDENTIFICATION & SUSCEPTIBILITY TESTING

Orientation Tests

Rapid Antigen Detection & Identification

API® & ID 32 Identification

Etest®

ATB™ Susceptibility Testing

VITEK® 2: Automated Method

VITEK 2 Compact: Automated Method

VITEK 2 Compact Analyser

VITEK 2 Compact Identification Cards

VITEK 2 Compact Antibiotic Susceptibility Cards

Gram Negative rods

Gram Positive cocci

Yeasts

VITEK 2 Compact Other Consumables & Accessories

VITEK MS: MALDI-TOF Mass Spectrometer



Reference	Description	Presentation
-----------	-------------	--------------

VITEK® 2 Compact: Automated Method



VITEK 2 Compact Antibiotic Susceptibility Cards

Gram Negative rods

Nordic

412938	VITEK 2 AST-N209	20 Cards
413080	VITEK 2 AST-N217	20 Cards
413081	VITEK 2 AST-N218	20 Cards
413083	VITEK 2 AST-N222	20 Cards
413147	VITEK 2 AST-N230	20 Cards
410609	LyfoCults® Plus VITEK 2 AST-GN QC set ‡	8 Vials

Gram Positive cocci

Nordic

22294	VITEK 2 AST-P594 (Replaced by 414576 by end of 2012)	20 Cards
414576 NEW	VITEK 2 AST-P629	20 Cards
410610	LyfoCults Plus VITEK 2 AST-GP QC set ‡	16 Vials

Global

22216	VITEK 2 AST-P576	20 Cards
22233	VITEK 2 AST-P580	20 Cards
22276	VITEK 2 AST-P586	20 Cards
22287	VITEK 2 AST-P592	20 Cards
XXXXXX	VITEK 2 AST-P5XXX	20 Cards
410028 NEW	VITEK 2 AST-ST01 (Standard Streptococcus Card)	20 Cards
410610	LyfoCults Plus VITEK 2 AST-GP QC set ‡	16 Vials

Yeasts

412610 NEW	VITEK 2 AST-YS06	20 Cards
410613	LyfoCults Plus VITEK 2 AST-YS QC set ‡	4 Vials

‡ See LyfoCults Plus [page 40](#) for strains contained within each set.

- IDENTIFICATION & SUSCEPTIBILITY TESTING
- Orientation Tests

Rapid Antigen Detection & Identification

API® & ID 32 Identification

Etest®

ATB™ Susceptibility Testing

VITEK® 2: Automated Method

VITEK 2 Compact: Automated Method

VITEK 2 Compact Analyser

VITEK 2 Compact Identification Cards

VITEK 2 Compact Antibiotic Susceptibility Cards

VITEK 2 Compact Other Consumables & Accessories

VITEK MS: MALDI-TOF Mass Spectrometer

Reference	Description	Presentation
VITEK® 2 Compact: Automated Method VITEK 2 Compact Other Consumables & Accessories 		
V1204	Saline 0.45%	3 x 500ml
V1200	Saline Dispenser	One
70610	Sterile cotton swabs	Pack of 100
69285	Test tubes (uncapped)	2000
V1221	Pipette 145 µl	One
30507	Pipette tips (0.5-250 µl) Non-sterile	10 Racks x 96 tips
V1222	Pipette 280 µl	One
30501	Pipette tips (100-1000 µl) Non-sterile	10 Racks x 60 tips
21250	DensiCHEK™ Plus	One
21255	DensiCHEK Plus - Standard Reference Kit	4 Vials
93059	DensiCHEK - Standard Reference	1 Vial
96128	VITEK 2 Smart Box 7L	One



IDENTIFICATION & SUSCEPTIBILITY TESTING

- Orientation Tests
- Rapid Antigen Detection & Identification
- API® & ID 32 Identification
- Etest®
- ATB™ Susceptibility Testing
- VITEK® 2: Automated Method
- VITEK 2 Compact: Automated Method
- VITEK MS: MALDI-TOF Mass Spectrometer

VITEK MS Analyser & Accessories

VITEK MS Consumables

Reference	Description	Presentation
-----------	-------------	--------------

VITEK® MS: MALDI-TOF Mass Spectrometer

VITEK MS Analyser & Accessories

410895	VITEK MS Mass Spectrometer with Myla™	Instrument
411655	VITEK MS Prep Station	Accessory
411402	VITEK MS-ID database	Accessory

VITEK MS Consumables

410893	Target Slides	32 Disposable Slides (8 x 4)
411071	CHCA Matrix Solution	5 x 0.5ml
411072	Extraction fluid for yeasts	5 x 0.5ml
411721 NEW	VITEK MS Silica Gel	500g



IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests

- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



VIDAS® – Simplicity, Flexibility & Accuracy



As a health care provider, you face new challenges every day. VIDAS is constantly adapting to help you overcome them by providing simplicity, flexibility and accuracy.

The VIDAS kits include all reagents in a single-dose and ready-to-use format together with separate standards and controls. Making it easy and cost-effective to run one or several patient samples at a time.



Reference	Description	Presentation
-----------	-------------	--------------

VIDAS

VIDAS Instruments

410416	mini VIDAS® Blue System	
410417	VIDAS Blue System	

VIDAS Consumables & Accessories

99091	mini VIDAS Printer paper	5 Rolls per Pack
99087	VIDAS decontamination sponges	Pack of 150
99187	VIDAS MLE card holder	One
99098	VIDAS Loading tray	One
99099	mini VIDAS Loading tray	One
93152	mini VIDAS bar code reader	One
99746	VIDAS Blue bar code reader	One
99747	mini VIDAS Blue bar code reader	One
30531	Dacron swabs	100 x 2

IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests

- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



Reference	Description	Presentation
VIDAS® VIDAS Reagents for Clinical Specimens (page 1 of 4) Hepatitis A 30307 VIDAS HAV IgM 30 Tests 30312 VIDAS HAV total 30 Tests Hepatitis B 30314 VIDAS Anti-HBc Total II 60 Tests 30238 VIDAS Anti-HBs Quick 60 Tests 30439 VIDAS HBc IgM II 30 Tests 30305 VIDAS HBe/Anti-HBe 30 Tests 30317 VIDAS HBs Ag Ultra confirmation 1 x 0.5ml 30315 VIDAS HBs Ag Ultra 60 Tests Hepatitis C 30308 VIDAS HCV 60 Tests HIV 30443 VIDAS HIV DUO Ultra 60 Tests 30447 VIDAS HIV DUO Quick 60 Tests 30117 VIDAS p24 II 30 Tests 30444 VIDAS p24 Confirmation 60 Tests TORC 30204 VIDAS CMV IgG 60 Tests 30205 VIDAS CMV IgM 30 Tests 30203 VIDAS CMV IgG Avidity 30 Tests 30221 VIDAS RBG II 60 Tests 30214 VIDAS RUB IgM 30 Tests 30211 VIDAS TOXO Competition 60 Tests 30210 VIDAS TOXO IgG II 60 Tests 30202 VIDAS TOXO IgM 60 Tests 30222 VIDAS TOXO Avidity 30 Tests		



IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests

- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



Reference	Description	Presentation																																																															
VIDAS® VIDAS Reagents for Clinical Specimens (page 2 of 4) Other Serology <tr> <td>30235</td><td>VIDAS EBNA IgG</td><td>30 Tests</td></tr> <tr> <td>30237</td><td>VIDAS EBV VCA IgM</td><td>30 Tests</td></tr> <tr> <td>30236</td><td>VIDAS EBV VCA/EA IgG</td><td>30 Tests</td></tr> <tr> <td>30320</td><td>VIDAS Lyme IgG</td><td>60 Tests</td></tr> <tr> <td>30319</td><td>VIDAS Lyme IgM</td><td>60 Tests</td></tr> <tr> <td>30219</td><td>VIDAS Measles IgG</td><td>60 Tests</td></tr> <tr> <td>30218</td><td>VIDAS Mumps IgG</td><td>60 Tests</td></tr> <tr> <td>30217</td><td>VIDAS Varicella Zoster IgG</td><td>60 Tests</td></tr> <tr> <td>30192</td><td>VIDAS H. pylori IgG</td><td>30 Tests</td></tr> Sepsis <tr> <td>30450</td><td>VIDAS B.R.A.H.M.S. Procalcitonin</td><td>60 Tests</td></tr> Antigen detection <tr> <td>30118</td><td>VIDAS CDA/B</td><td>60 Tests</td></tr> <tr> <td>30107</td><td>VIDAS Rotavirus</td><td>60 Tests</td></tr> Thyroid <tr> <td>30461</td><td>VIDAS Anti-TPO</td><td>30 Tests</td></tr> <tr> <td>30462</td><td>VIDAS Anti-Tg</td><td>30 Tests</td></tr> <tr> <td>30402</td><td>VIDAS FT3</td><td>60 Tests</td></tr> <tr> <td>30459</td><td>VIDAS FT4</td><td>60 Tests</td></tr> <tr> <td>30403</td><td>VIDAS T3</td><td>60 Tests</td></tr> <tr> <td>30404</td><td>VIDAS T4</td><td>60 Tests</td></tr> <tr> <td>30400</td><td>VIDAS TSH</td><td>60 Tests</td></tr> <tr> <td>30441</td><td>VIDAS TSH3 (Ultra sensitive)</td><td>60 Tests</td></tr> Anaemia <tr> <td>30411</td><td>VIDAS Ferritin</td><td>60 Tests</td></tr>			30235	VIDAS EBNA IgG	30 Tests	30237	VIDAS EBV VCA IgM	30 Tests	30236	VIDAS EBV VCA/EA IgG	30 Tests	30320	VIDAS Lyme IgG	60 Tests	30319	VIDAS Lyme IgM	60 Tests	30219	VIDAS Measles IgG	60 Tests	30218	VIDAS Mumps IgG	60 Tests	30217	VIDAS Varicella Zoster IgG	60 Tests	30192	VIDAS H. pylori IgG	30 Tests	30450	VIDAS B.R.A.H.M.S. Procalcitonin	60 Tests	30118	VIDAS CDA/B	60 Tests	30107	VIDAS Rotavirus	60 Tests	30461	VIDAS Anti-TPO	30 Tests	30462	VIDAS Anti-Tg	30 Tests	30402	VIDAS FT3	60 Tests	30459	VIDAS FT4	60 Tests	30403	VIDAS T3	60 Tests	30404	VIDAS T4	60 Tests	30400	VIDAS TSH	60 Tests	30441	VIDAS TSH3 (Ultra sensitive)	60 Tests	30411	VIDAS Ferritin	60 Tests
30235	VIDAS EBNA IgG	30 Tests																																																															
30237	VIDAS EBV VCA IgM	30 Tests																																																															
30236	VIDAS EBV VCA/EA IgG	30 Tests																																																															
30320	VIDAS Lyme IgG	60 Tests																																																															
30319	VIDAS Lyme IgM	60 Tests																																																															
30219	VIDAS Measles IgG	60 Tests																																																															
30218	VIDAS Mumps IgG	60 Tests																																																															
30217	VIDAS Varicella Zoster IgG	60 Tests																																																															
30192	VIDAS H. pylori IgG	30 Tests																																																															
30450	VIDAS B.R.A.H.M.S. Procalcitonin	60 Tests																																																															
30118	VIDAS CDA/B	60 Tests																																																															
30107	VIDAS Rotavirus	60 Tests																																																															
30461	VIDAS Anti-TPO	30 Tests																																																															
30462	VIDAS Anti-Tg	30 Tests																																																															
30402	VIDAS FT3	60 Tests																																																															
30459	VIDAS FT4	60 Tests																																																															
30403	VIDAS T3	60 Tests																																																															
30404	VIDAS T4	60 Tests																																																															
30400	VIDAS TSH	60 Tests																																																															
30441	VIDAS TSH3 (Ultra sensitive)	60 Tests																																																															
30411	VIDAS Ferritin	60 Tests																																																															



IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests
- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



Reference	Description	Presentation																																																									
VIDAS® VIDAS Reagents for Clinical Specimens (page 3 of 4) Reproduction / Fertility <tr> <td>30431</td><td>VIDAS Estradiol II</td><td>60 Tests</td></tr> <tr> <td>30407</td><td>VIDAS FSH</td><td>60 Tests</td></tr> <tr> <td>30405</td><td>VIDAS HCG</td><td>60 Tests</td></tr> <tr> <td>30406</td><td>VIDAS LH</td><td>60 Tests</td></tr> <tr> <td>30409</td><td>VIDAS Progesterone</td><td>60 Tests</td></tr> <tr> <td>30410</td><td>VIDAS Prolactin</td><td>60 Tests</td></tr> <tr> <td>30418</td><td>VIDAS Testosterone</td><td>30 Tests</td></tr> Tumour markers <tr> <td>30413</td><td>VIDAS AFP</td><td>60 Tests</td></tr> <tr> <td>30426</td><td>VIDAS CA 125 II™</td><td>30 Tests</td></tr> <tr> <td>30429</td><td>VIDAS CA 15-3™</td><td>30 Tests</td></tr> <tr> <td>30427</td><td>VIDAS CA 19-9™</td><td>30 Tests</td></tr> <tr> <td>30440</td><td>VIDAS FPSA</td><td>60 Tests</td></tr> <tr> <td>30428</td><td>VIDAS TPSA</td><td>30 Tests</td></tr> Drugs <tr> <td>30603</td><td>VIDAS Digoxin</td><td>60 Tests</td></tr> Cardiovascular <tr> <td>30421</td><td>VIDAS CK MB</td><td>30 Tests</td></tr> <tr> <td>30446</td><td>VIDAS Myoglobin</td><td>30 Tests</td></tr> <tr> <td>30449</td><td>VIDAS NT-proBNP</td><td>60 Tests</td></tr> <tr> <td>30448</td><td>VIDAS Troponin I Ultra</td><td>60 Tests</td></tr> Quality Control <tr> <td>30706</td><td>VIDAS QCV</td><td>60 Tests</td></tr>			30431	VIDAS Estradiol II	60 Tests	30407	VIDAS FSH	60 Tests	30405	VIDAS HCG	60 Tests	30406	VIDAS LH	60 Tests	30409	VIDAS Progesterone	60 Tests	30410	VIDAS Prolactin	60 Tests	30418	VIDAS Testosterone	30 Tests	30413	VIDAS AFP	60 Tests	30426	VIDAS CA 125 II™	30 Tests	30429	VIDAS CA 15-3™	30 Tests	30427	VIDAS CA 19-9™	30 Tests	30440	VIDAS FPSA	60 Tests	30428	VIDAS TPSA	30 Tests	30603	VIDAS Digoxin	60 Tests	30421	VIDAS CK MB	30 Tests	30446	VIDAS Myoglobin	30 Tests	30449	VIDAS NT-proBNP	60 Tests	30448	VIDAS Troponin I Ultra	60 Tests	30706	VIDAS QCV	60 Tests
30431	VIDAS Estradiol II	60 Tests																																																									
30407	VIDAS FSH	60 Tests																																																									
30405	VIDAS HCG	60 Tests																																																									
30406	VIDAS LH	60 Tests																																																									
30409	VIDAS Progesterone	60 Tests																																																									
30410	VIDAS Prolactin	60 Tests																																																									
30418	VIDAS Testosterone	30 Tests																																																									
30413	VIDAS AFP	60 Tests																																																									
30426	VIDAS CA 125 II™	30 Tests																																																									
30429	VIDAS CA 15-3™	30 Tests																																																									
30427	VIDAS CA 19-9™	30 Tests																																																									
30440	VIDAS FPSA	60 Tests																																																									
30428	VIDAS TPSA	30 Tests																																																									
30603	VIDAS Digoxin	60 Tests																																																									
30421	VIDAS CK MB	30 Tests																																																									
30446	VIDAS Myoglobin	30 Tests																																																									
30449	VIDAS NT-proBNP	60 Tests																																																									
30448	VIDAS Troponin I Ultra	60 Tests																																																									
30706	VIDAS QCV	60 Tests																																																									



IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests

- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



Reference	Description	Presentation
-----------	-------------	--------------

VIDAS®



VIDAS Reagents for Clinical Specimens (page 4 of 4)

Immuno-Haemostatis

30455	VIDAS D-Dimer Exclusion II	60 Tests
30115	VIDAS Protein C	30 Tests
30436	VIDAS von Willebrand Factor (vWFAg)	30 Tests

Other tests

66581	Serum free (for urinary assay of VIDAS Cortisol)	1 x 10ml
30420	VIDAS Beta-2 Microglobulin	30 Tests
30451	VIDAS Cortisol S	60 Tests

IMMUNOASSAYS

- VIDAS®
- VIDAS Instruments
- VIDAS Consumables & Accessories
- VIDAS Reagents for Clinical Specimens
- Hepatitis A
- Hepatitis B
- Hepatitis C
- HIV
- TORC
- Other Serology
- Sepsis
- Antigen detection
- Thyroid
- Anaemia
- Reproduction / Fertility
- Tumour markers
- Drugs
- Cardiovascular
- Quality Control
- Immuno-Haemostatis
- Other tests

- Microelisa
- Microbiology Reagent Kits
- Hepatitis
- HIV
- Syphilis



Reference	Description	Presentation
-----------	-------------	--------------

Microelisa

Microbiology Reagent Kits

Hepatitis


284132	HEPANOSTIKA™ – HBsAg Ultra	24 x 8 Tests
284133	HEPANOSTIKA – HBsAg Ultra	72 x 8 Tests
280253	HEPANOSTIKA – HBsAg Ultra Confirmation	25 Tests
284144	HEPANOSTIKA – anti-HBc Uni-Form	24 x 8 Tests

HIV

284017	VIRONOSTIKA™ HIV Uni-Form II plus 0	24 x 8 Tests
284018	VIRONOSTIKA HIV Uni-Form II plus 0	72 x 8 Tests
259851	VIRONOSTIKA HIV Ag/Ab	24 x 8 Tests
259852	VIRONOSTIKA HIV Ag/Ab	72 x 8 Tests

Syphilis

285035	TREPANOSTIKA™-TP recombinant	72 x 8 Tests
285034	TREPANOSTIKA-TP recombinant (includes 1ml positive control, 2ml negative control)	24 x 8 Tests



MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence
- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINAKit
- Control Slides
- Other
- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus

Reference	Description	Presentation
Serodiagnosis		
Leishmaniasis		
75931	Leishmania-Spot IF	10 x 10 Circles
Malaria		
72751	Falciparum-Spot IF	10 x 10 Circles
Rheumatic Fever (Streptococcal Infection)		
Antistreptodornase B		
72331	ASD-Kit	20 Tests
Antistreptolysin		
72351	ASL-Kit	20 Tests
Rheumatoid Factors		
72013	Arthri-Slidex	100 Tests



MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence

- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINAKit
- Control Slides
- Other

- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
-----------	-------------	--------------

Serodiagnosis

Syphilis

RPR inhibition

280447	RPR - nosticon II	100 Tests
--------	-------------------	-----------

Immunofluorescence

75661	Sorbent	5ml
75681	Trepo-Spot IF	10 x 10 Circles

Positive control serum

72191	Serotrol	1ml
-------	----------	-----

Toxoplasmosis

Agglutination

75361	Toxo ISAGA	Kit for 192 Tests
79322	Toxo-ISAGA IgA	192 Wells : 48 Screening Tests
75481	Toxo-Screen DA	Complete Kit with 4 Microtitration plates

Immunofluorescence

75471	Toxo-spot IF™	10 x 10 Circles
-------	---------------	-----------------

MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence

- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINakit
- Control Slides
- Other

- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
-----------	-------------	--------------

Immunology

Liquid Globulins for Immunofluorescence

75692	Fluoline-G - anti-human IgG	1ml Bottle
75603	Fluoline-H - anti-human Immunoglobulins	1ml Bottle
75672	Fluoline-M - anti-human IgM	1ml Bottle
75521	Fuoprep - mounting medium	15ml Bottle

Labware

75751	Slides of 10 x 6 mm Circles	72 Slides
-------	-----------------------------	-----------

Argene™ Immunology Reagents (page 1 of 4)



Herpes Virus Antibodies

14-002	Anti CMV CINApool + F(ab')2 FITC	0.5ml (20x)
11-002	Anti CMV CINApool purified	0.5ml (20x)
11-004	Anti HCMV I.E.A + E.A purified	0.5ml (50x)
11-003	Anti HCMV I.E.A purified	0.5ml (50x)
17-090	anti HSV-1 + 2 FITC	80 Tests
12-090	Anti HSV-1 + 2 FITC concentrated	7.5ml
11-088	Anti HSV-1 purified	7.5ml
11-089	Anti HSV-2 purified	7.5ml
11-090	Anti HSV-1 + 2 purified	7.5ml
17-017	anti VARICELLA ZOSTER VIRUS FITC	80 Tests
11-018	Anti VZV 1U1 purified	0.5ml (40x)
11-017	Anti VZV 2013 purified	0.5ml (40x)
11-002G	Dilution buffer for CINApool	4.5ml
11-002I	Evans blue 1% for CINApool	250µl

MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence
- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINAKit
- Control Slides
- Other
- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
-----------	-------------	--------------

Immunology



Argene™ Immunology Reagents (page 2 of 4)

Respiratory Virus Antibodies

11-091	Anti ADENO-INFLU-PARA purified	0.5ml (40x)
17-020	anti ADENOVIRUS GROUP FITC	80 Tests
12-020	Anti ADENOVIRUS GROUP FITC concentrated	0.5ml (40x)
11-020	Anti ADENOVIRUS GROUP purified	0.5ml (40x)
17-091	anti ADENOVIRUS, INFLUENZA, PARAINFLUENZA FITC	80 Tests
17-092	anti ADENOVIRUS, INFLUENZA, PARAINFLUENZA, RSV FITC	80 Tests
11-030	Anti FLU A GROUP purified	0.5ml (40x)
11-035	Anti FLU B GROUP purified	0.5ml (40x)
17-270	anti human METAPNEUMOVIRUS FITC	80 Tests
17-030	anti INFLUENZA A GROUP FITC	80 Tests
12-030	Anti INFLUENZA A GROUP FITC concentrated	0.5ml (40x)
17-035	anti INFLUENZA B GROUP FITC	80 Tests
12-035	Anti INFLUENZA B GROUP FITC concentrated	0.5ml (40x)
17-036	anti PARAINFLUENZA 1 FITC	80 Tests
11-036	Anti PARAINFLUENZA 1 purified	0.5ml (40x)
17-037	anti PARAINFLUENZA 2 FITC	80 Tests
11-037	Anti PARAINFLUENZA 2 purified	0.5ml (40x)
17-038	anti PARAINFLUENZA 3 FITC	80 Tests
12-038	Anti PARAINFLUENZA 3 FITC concentrated	0.5ml (40x)
11-038	Anti PARAINFLUENZA 3 purified	0.5ml (40x)
17-040	anti PARAINFLUENZA GROUP FITC	80 Tests
12-040	Anti PARAINFLUENZA GROUP FITC concentrated	1ml (20x)
11-040	Anti PARAINFLUENZA purified	0.5ml (40x)
17-042	anti RSV GROUP FITC	80 Tests
12-042	Anti RSV GROUP FITC concentrated	0.5ml (40x)
11-042	Anti RSV GROUP purified	0.5ml (40x)

MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence
- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINAKit
- Control Slides
- Other
- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
-----------	-------------	--------------

Immunology



Argene™ Immunology Reagents (page 3 of 4)

CINAKit

19-002	CINAKit complete kit HCMV ppUL83 (pp65) I.F. Antigenemia	100 Tests
19-0028	CINAKit complete kit HCMV ppUL83 (pp65) I.F. Antigenemia	200 Tests
19-0028US	CINAKit complete kit HCMV ppUL83 (pp65) I.F. Antigenemia	200 Tests
19-002K	Dextran 6% for CINAKit 19-002	60ml
19-002C	Fixative solution for CINAKit 19-002	110ml
33-041	Mounting medium for immunofluorescence	15ml
19-002E	Permeabilization solution for CINAKit 19-002	110ml
19-002B	Rapid lysis solution for CINAKit 19-002	45ml

Control Slides

40-061	Control slide ADENOVIRUS	5 Slides
40-022	Control slide CYTOMEGALOVIRUS IEA-EA	5 Slides
40-017	Control slide HERPES SIMPLEX TYPE 1 & 2	5 Slides
40-071	Control slide INFLUENZA A	5 Slides
40-081	Control slide INFLUENZA B	5 Slides
40-161	Control slide PANEL ADV, IA, IB, VRS, P1, P2, P3	5 Slides
40-091	Control slide PARAINFLUENZA 1	5 Slides
40-101-1	Control slide PARAINFLUENZA 2	5 Slides
40-111	Control slide PARAINFLUENZA 3	5 Slides
40-121	Control slide RSV	5 Slides
40-041	Control slide VARICELLA - ZOSTER	5 Slides

MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence
- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINakit
- Control Slides
- Other
- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
-----------	-------------	--------------

Immunology



Argene™ Immunology Reagents (page 4 of 4)

Other

14-002	Anti CMV CINApool + F(ab')2 FITC	
11-033	Anti HAV 1009 purified	250ug 0.5ml
11-034	Anti HAV 813 purified	250ug 0.5ml
12-090	Anti HSV-1 + 2 FITC concentrated	
11-045	Anti MEASLES purified	0.5ml (40x)
50-010	anti MOUSE IgG + IgM "human ads." FITC	0.5ml
50-011	anti MOUSE IgG + IgM "human ads." FITC	1ml
50-012	anti MOUSE IgG + IgM "human ads." FITC	2ml
50-021	anti MOUSE IgG + IgM "human ads." PEROXIDASE	1ml
50-022	anti MOUSE IgG + IgM "human ads." PEROXIDASE	2ml
11-046	Anti MUMPS purified	0.5ml (40x)
33-030	Evans blue 1%	2ml
51-010	F(ab')2 anti MOUSE IgG + IgM "human ads." FITC	0.5ml
51-011	F(ab')2 anti MOUSE IgG + IgM "human ads." FITC	1ml
31-010	IF slides - Teflon coated - 6mm wells	100 Slides
31-030	"Kova glasstic slides" for antigenemia	100 Slides
33-041	Mounting medium for immunofluorescence	15ml
33-011	PBS pH 7.2 without Ca ²⁺ and Mg ²⁺	4 Litres

MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis
- Leishmaniasis
- Malaria
- Rheumatic Fever (Streptococcal Infection)
- Antistreptodornase B
- Antistreptolysin
- Rheumatoid Factors
- Syphilis
- RPR inhibition
- Immunofluorescence
- Positive control serum
- Toxoplasmosis
- Agglutination
- Immunofluorescence
- Immunology
- Liquid Globulins for Immunofluorescence
- Labware
- Argene™ Immunology Reagents
- Herpes Virus Antibodies
- Respiratory Virus Antibodies
- CINAKit
- Control Slides
- Other

- Rapid Tests
- Faecal Occult Blood
- Bladder Tumour Antigen
- Troponin
- C-Reactive Protein
- Influenza
- HIV
- Hepatitis B
- Rotavirus / Adenovirus



Reference	Description	Presentation
Rapid Tests Faecal Occult Blood 410594 bioNexia® FOBplus Test Rapid (5 minute) immunological point of care test for detection of human haemoglobin in stools 25 Tests 414275 bioNexia FOBplus Lab Kit Rapid (5 minute) immunological point of care test for detection of human haemoglobin in stools, kit specially developed for laboratory personnel 25 Tests 414276 bioNexia FOBplus Patient Set Patient collection set complementary kit to ref 410594 bioNexia FOBplus Test and 414275 bioNexia FOBplus Lab Kit 25 Tests		
Bladder Tumour Antigen 410591 bioNexia BTA Test Rapid (5 minute) point of care test for bladder tumour antigen in urine 10 Tests		
Troponin 410598 bioNexia Troponin I Test Rapid (10 minute) point of care test for Troponin I in whole blood or serum/plasma 10 Tests		
C-Reactive Protein 410592 bioNexia CRPplus Test Rapid (10 minute), semi-quantitative point of care test for C-Reactive Protein in whole blood 10 Tests		
Influenza 410918 bioNexia Influenza A+B Test Rapid (10 minute), test for detection of Influenza A and B in nasopharyngeal aspirates 10 Tests		

All bioNexia kits contain everything needed for sample collection and test processing, e.g. lancets, swabs etc.



MANUAL SEROLOGY / RAPID TEST / IMMUNOLOGY

- Serodiagnosis

Leishmaniasis

Malaria

Rheumatic Fever (Streptococcal Infection)

Antistreptodornase B

Antistreptolysin

Rheumatoid Factors

Syphilis

RPR inhibition

Immunofluorescence

Positive control serum

Toxoplasmosis

Agglutination

Immunofluorescence
- Immunology

Liquid Globulins for Immunofluorescence

Labware

Argene™ Immunology Reagents

Herpes Virus Antibodies

Respiratory Virus Antibodies

CINAKit

Control Slides

Other
- Rapid Tests

Faecal Occult Blood

Bladder Tumour Antigen

Troponin

C-Reactive Protein

Influenza

HIV

Hepatitis B

Rotavirus / Adenovirus

Reference	Description	Presentation
Rapid Tests HIV 31112 VIKIA® HIV 1/2 25 Test Devices Hepatitis B 31124 VIKIA HBsAg 25 Test Devices Rotavirus / Adenovirus 31111 VIKIA Rota/Adeno 20 Test Devices		



MOLECULAR BIOLOGY

- NucliSENS®
Feature – NucliSENS easyMAG®
NucliSENS Magnetic Extraction Systems
 easyMAG Automated Extraction System
 miniMAG® Semi Automated Extraction System
NucliSENS EasyQ® System
 Real Time Amplification & Detection System
NucliSENS EasyQ Reagents, Consumables & Accessories
 Enterovirus
 HSV
 MRSA
 Consumables
- Argene™
Real-Time PCR
 Complete Kits
 Respiratory Multi Well System
 Optimised Primers & Probes
 Hybridowell™
 Ancillary Reagents
- DiversiLab®
Feature – DiversiLab
DiversiLab Equipment
DiversiLab Fingerprinting Kits
DiversiLab Test Accessories
DiversiLab Extraction Accessories



NucliSENS® easyMAG®
Automated Extraction System

High Extraction Quality

- Magnetic extraction based on the gold standard BOOM technology (silica particles)
- High concentration factor: up to 100
- High quality RNA and DNA extraction: the easyMAG has been validated with a wide range of clinical commercial kits for real-time amplification and detection

Flexibility/Versatility

Optimise your Workflow and Productivity, having the possibility to perform several applications in the same run:

- 40 minutes for 24 extractions
- Generic protocol approach for DNA and RNA extraction
- One set of reagents whatever the application
- High flexibility in term of sample input and output volumes



MOLECULAR BIOLOGY

- NucliSENS®
- Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables
- Argene™
- Real-Time PCR

Complete Kits

Repiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents
- DiversiLab®
- Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories

Reference	Description	Presentation
NucliSENS® NucliSENS Magnetic Extraction Systems easyMAG® Automated Extraction System		
280140	easyMAG Automated Magnetic Extraction System	One
278303	Strip Plate (Greiner)	100 Plates
280141	BioHit Pipettor	One
280151	BioHit adaptor	One
280146	BioHit Pipette tips	960
280130	NucliSENS easyMAG Extraction Buffer 1	4 x 1 Litre
280131	NucliSENS easyMAG Extraction Buffer 2	4 x 1 Litre
280132	NucliSENS easyMAG Extraction Buffer 3	4 x 1 Litre
280133	NucliSENS easyMAG Magnetic Silica	384 Extractions
280134	NucliSENS easyMAG Lysis Buffer	4 x 1 Litre
280135	NucliSENS easyMAG Disposables	384 Positions
miniMAG® Semi Automated Extraction System		
200305	NucliSENS miniMAG Extractor	One
200299	Magnetic Rack 12 hole	One
200292	NucliSENS Lysis Buffer 2ml	48 Tubes
200293	NucliSENS Magnetic Extraction Reagents	48 Tests
200294	Micro Tubes 1.5ml with cap	500



MOLECULAR BIOLOGY

- NucliSENS®
- Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables

- Argene™
- Real-Time PCR

Complete Kits

Repiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents

- DiversiLab®
- Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories



Reference	Description	Presentation
-----------	-------------	--------------

NucliSENS®



NucliSENS EasyQ® System

Real Time Amplification & Detection System

285052	EasyQ System for Real Time Amplification and Detection using Molecular Beacons. Includes EasyQ Reader, PC Workstation, Director Software, Incubator, Printer & Tube-strip Centrifuge	One
285056	Tube-strip Centrifuge	One
270677	Mo Bio Adaptor FOR Genie Vortex	One

NucliSENS EasyQ® Reagents, Consumables & Accessories

Enterovirus

410311	NucliSENS EasyQ Enterovirus v1.1	48 Tests
410312	NucliSENS EasyQ Enterovirus Controls	6 x Positive and 6 x Negative

HSV

280110	NucliSENS EasyQ HSV 1 + 2	48 Tests
280104	NucliSENS EasyQ Basic Kit v2	48 Tests

MRSA

280115	NucliSENS EasyQ MRSA Kit	48 Tests
280101	COPAN Sterile dry flocked swabs	Pack of 100

Consumables

285048	0.2ml 8-tube strips	125
285051	0.2ml 8-tube caps	125



MOLECULAR BIOLOGY

- NucliSENS®
- Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables

- Argene™
- Real-Time PCR

Complete Kits

Respiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents

- DiversiLab®
- Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories



Reference	Description	Presentation
-----------	-------------	--------------

Argene™



Real-Time PCR (page 1 of 3)

Complete Kits

69-010B	ADENOVIRUS R-gene™ - Quantification kit	90 Tests
69-013B	BK Virus R-gene Quantification Kit	90 Tests
69-011B	BORDETELLA R-gene - Real-Time PCR Kit	60 Tests
69-100B	CMV, HHV6, 7, 8 R-gene Quantification Kit	140 Tests
69-003B	CMV R-gene Quantification Kit	90 Tests
69-002B	EBV R-gene Quantification Kit	90 Tests
69-005B	ENTEROVIRUS R-gene - Real-Time PCR Kit	90 Tests
69-004B	HSV1 HSV2 VZV R-gene Quantification Kit	180 Tests

Respiratory Multi Well System

71-043	Adenovirus/human Bocavirus r-gene™	60 Tests
71-044	Chla/Myco Pneumo r-gene	60 Tests
71-045	Coronavirus/Parainfluenza virus r-gene	60 Tests
71-040	Influenza A/B r-gene	60 Tests
71-042	Rhinovirus, Enterovirus/Cell control r-gene	60 Tests
71-041	RSV/hMPV r-gene	60 Tests

Real-Time PCR Assays for Transplant & Respiratory Medicine

➤ Validated for

- All major Real-Time PCR Instruments
- Multiple automatic extraction systems
- Various types of samples
- Uniform procedure

➤ Included

- All controls to monitor the entire process
- Ready-to-use reagents
- Quantification standards
- Turnkey solution



MOLECULAR BIOLOGY

- NucliSENS®**
 Feature – NucliSENS easyMAG®
 NucliSENS Magnetic Extraction Systems
 easyMAG Automated Extraction System
 miniMAG® Semi Automated Extraction System
 NucliSENS EasyQ® System
 Real Time Amplification & Detection System
 NucliSENS EasyQ Reagents, Consumables & Accessories
 Enterovirus
 HSV
 MRSA
 Consumables
- Argene™**
 Real-Time PCR
 Complete Kits
 Respiratory Multi Well System
 Optimised Primers & Probes
 Hybridowell™
 Ancillary Reagents
- DiversiLab®**
 Feature – DiversiLab
 DiversiLab Equipment
 DiversiLab Fingerprinting Kits
 DiversiLab Test Accessories
 DiversiLab Extraction Accessories



Reference	Description	Presentation
Argene™ Real-Time PCR (page 2 of 3) Optimised Primers & Probes		
71-010	Adenovirus r-gene™ Primers/Probes	60 Reactions
71-005	BK Virus r-gene Primers/Probe	60 Reactions
71-009	Bocavirus r-gene Primers/Probe	60 Reactions
71-012	Bordetella parapertussis r-gene	60 Reactions
71-001	Chlamydomphila pneumoniae r-gene Primers/Probe	60 Reactions
69-100R6	HHV6 amplification premix	60 Reactions
69-100R7	HHV7 amplification premix	20 Reactions
69-100R8	HHV8 amplification premix	20 Reactions
71-008	hMPV A/B r-gene Primers/Probes	60 Reactions
71-015	HSV1 r-gene	60 Reactions
71-016	HSV2 r-gene	60 Reactions
71-300	Influenza A(M) Group & H1N1 2009 r-gene	60 Tests
71-004	JC Virus r-gene Primers/Probe	60 Reactions
71-003	Legionella species r-gene Primers/Probe	60 Reactions
71-002	Mycoplasma pneumoniae r-gene Primers/Probe	60 Reactions
71-007	RSV A/B r-gene Primers/Probes	60 Reactions
71-017	VZV R-gene™	60 Reactions



MOLECULAR BIOLOGY

- NucliSENS®
- Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables
- Argene™
- Real-Time PCR

Complete Kits

Repiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents
- DiversiLab®
- Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories

Reference	Description	Presentation
-----------	-------------	--------------

Argene™



Real-Time PCR (page 3 of 3)

Hybridowell™

67-065BC	ADENOVIRUS CONSENSUS amplification/detection	50 Samples
67-025BC	CHLAMYLEGE amplification/detection	50 Samples
67-080BC	ENTEROVIRUS CONSENSUS amplification/detection	50 Samples
67-090BC	HERPES CONSENSUS GENERIC amplification/detection	50 Samples
67-050	Hybridowell HERPES Identification (for 67-090BC)	96 Tests
67-012	Hybridowell Universal, complete kit	96 Tests
67-130BC	JC/BK CONSENSUS amplification/detection	50 Samples

Ancillary Reagents

68-009	Bocavirus r-gene™ positive control	10 Reactions
71-106	CELL Control r-gene	100 Reactions
68-020	CMV Clear QC Panel	3 x 1ml
71-103	Colour Compensation r-gene	3 x 60µl
71-100	DICO Ampli r-gene	100 Reactions
71-101	DICO Extra r-gene	100 Reactions
68-006	Quanti FluA QS r-gene	30 Reactions
68-008	Quanti HHV8 QS r-gene	30 Reactions
71-105	RICO Extra r-gene	100 Reactions
65-001	UNG r-gene	500 Tests



MOLECULAR BIOLOGY

NucliSENS®

Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables

Argene™

Real-Time PCR

Complete Kits

Respiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents

DiversiLab®

Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories



DiversiLab® – Automated Strain Typing

DiversiLab is an automated platform using rep-PCR technology that provides standardised, reproducible DNA fingerprinting of bacterial and fungal samples for complete isolate characterisation. It rapidly and accurately distinguishes your isolates to strain level.

With DiversiLab technology, strain typing becomes much easier and possible in routine use, even with numerous specimens. Intra-laboratory and Inter-laboratory reproducibility are key to accurately track HAI. The DiversiLab platform has been shown to have:

- <1% variability for inter-laboratory reproducibility
- <2% variability for intra-laboratory reproducibility

The synergy between this innovative genotyping system and bioMérieux's other systems supports clinical decisions that can ultimately reduce the number of HAI and ensure the appropriate use of antibiotic therapy.



MOLECULAR BIOLOGY

- NucliSENS®
- Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables
- Argene™
- Real-Time PCR

Complete Kits

Repiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents
- DiversiLab®
- Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories

Reference	Description	Presentation
DiversiLab® DiversiLab Equipment 		
	DiversiLab System	Instrument
DiversiLab Fingerprinting Kits		
410946	DiversiLab <i>Acinetobacter</i> Kit	24 Tests
410948	DiversiLab <i>Archaeobacteria</i> Kit	24 Tests
410949	DiversiLab <i>Aspergillus</i> Kit	24 Tests
410962	DiversiLab <i>Bacillus</i> Kit	24 Tests
411007	DiversiLab Bacteria (Gram Positive & Gram Negative) Kit	24 Tests
410963	DiversiLab <i>Bifidobacterium</i> (& <i>Rhodococcus</i> & <i>Nocardi spp</i>) Kit	24 Tests
410964	DiversiLab <i>Campylobacter</i> Kit	24 Tests
410965	DiversiLab <i>Candida</i> Kit	24 Tests
410966	DiversiLab <i>Clostridium</i> Kit	24 Tests
410967	DiversiLab <i>Clostridium perfringens</i> Kit	24 Tests
410968	DiversiLab <i>Enterobacter</i> Kit	24 Tests
410969	DiversiLab <i>Enterococcus</i> Kit	24 Tests
410980	DiversiLab <i>Escherichia</i> (& <i>Shigella spp</i>) Kit	24 Tests
410993	DiversiLab <i>Francisella</i> Kit	24 Tests
411008	DiversiLab Fungal (Yeasts & Molds) Kit	24 Tests
410981	DiversiLab <i>Klebsiella</i> Kit	24 Tests
410982	DiversiLab <i>Lactobacillus</i> Kit	24 Tests
410983	DiversiLab <i>Listeria</i> Kit	24 Tests
410984	DiversiLab <i>Mycobacterium</i> Kit	24 Tests
410986	DiversiLab <i>Mycobacterium tuberculosis</i> Kit	24 Tests
410985	DiversiLab <i>Mycoplasma</i> Kit	24 Tests
410987	DiversiLab <i>Propionibacterium</i> (& <i>Legionella spp</i>) Kit	24 Tests
410988	DiversiLab <i>Pseudomonas</i> Kit	24 Tests
410994	DiversiLab <i>Saccharomyces</i> Kit	24 Tests
410989	DiversiLab <i>Salmonella</i> Kit	24 Tests
410990	DiversiLab <i>Serratia</i> Kit	24 Tests
410991	DiversiLab <i>Staphylococcus</i> Kit	24 Tests
410992	DiversiLab <i>Streptococcus</i> Kit	24 Tests



MOLECULAR BIOLOGY

- NucliSENS®

Feature – NucliSENS easyMAG®

NucliSENS Magnetic Extraction Systems

easyMAG Automated Extraction System

miniMAG® Semi Automated Extraction System

NucliSENS EasyQ® System

Real Time Amplification & Detection System

NucliSENS EasyQ Reagents, Consumables & Accessories

Enterovirus

HSV

MRSA

Consumables
- Argene™

Real-Time PCR

Complete Kits

Repiratory Multi Well System

Optimised Primers & Probes

Hybridowell™

Ancillary Reagents
- DiversiLab®

Feature – DiversiLab

DiversiLab Equipment

DiversiLab Fingerprinting Kits

DiversiLab Test Accessories

DiversiLab Extraction Accessories

Reference	Description	Presentation
-----------	-------------	--------------

DiversiLab®



DiversiLab Test Accessories

270670	DiversiLab LabChip Kit 25 Chips	325 Tests
	Kit contains:	
	DiversiLab LabChip Reagents	325 Tests
	DiversiLab LabChip	25 LabChips x 13 Tests
	DiversiLab Syringe	Two
270715	DiversiLab LabChip Priming Station	One
270673	DiversiLab Gel Standards	1 Set
270656	DiversiLab Thermal Cycler QC Kit	48 Tests
411164	DiversiLab Training Kit	24 Tests

DiversiLab Extraction Accessories

270676	DiversiLab Mo Bio UltraClean™	250 Tests
270675	DiversiLab Mo Bio UltraClean	50 Tests
270677	Mo Bio Adaptor FOR Genie Vortex	One

