

USER MANUAL



ELECTRIC HEATING PAD

SY•AL30

WELCOME

Dear Customer,

Welcome to Sytech! Thank you for purchasing the Electric Heating Pad SY-AL30. Please read this manual carefully before operating your product and pay close attention to the precautions which are mentioned below. Store this user manual in a safe place for future reference. We are sure that this new Electric Heating Pad will satisfy your needs.

SAFETY GUIDELINE

1. Do not fold heating pad
2. Do not puncture, cut or otherwise damage any part of the heated electric pad
3. Keep away from water, moisture and high humidity
4. Do not attempt to modify, repair, disassemble or reassemble any component(s) of the heated electric pad. Repairs should be carried out by the manufacturer or other qualified persons.
5. Not for use by children or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.
6. For domestic use only
7. Warning! Never go to sleep with the heating pad in use. Extended use can cause burns.
8. Do not cut or pierce the heating pad
9. Do not pull, twist or otherwise apply force to the cables.
10. The device should not be connected to the mains via a timer.
11. The switch should not be placed on or under the appliance, nor covered when the appliance is in operation.
12. When in use do not put heavy or large objects on the pad
13. When in use do not place any source of heat over the pad, e.g. hot water bottles, heating pads or similar
14. When not in use, unplug the pad from the socket. The appliance should be stored in a dry place. Take care not to put additional weight on it.
15. Do not allow children being supervised to play with the appliance.
16. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
17. The appliance is not to be used if there are signs of damage. Observe for signs of damage such as rips, stitching coming undone, cords becoming loose or similar.
18. The appliance is not to be used for warming animals.
19. The appliance has to be repaired or replaced if cover is worn.

PARTS OF THE CONTROLLER

- 1) 2 Levels of heat to select from
- 2) LED indicator



OPERATION

- Plug the socket into the mains
- To turn the power on, move the sliding switch upwards so '1' or '2' shows in the window. '1' is low power (lower temperature), '2' is high power (higher temperature). The heated electric pad provides automatic electronic temperature regulation and overheating protection. In the case of failure of the automatic temperature controller, the pad will automatically permanently turn off once the temperature reaches approximately 100 degrees Celsius.
- Automatic temperature control: within about 10 minutes the temperature will rise to 55-60 degrees Celsius. The temperature will then vary between 40-45 and 55-60 degrees Celsius as the pad regulates its own temperature.
- To turn the pad off, move the slider down so '0' shows in the window.

MAINTENANCE

- The heated electric pad should be checked regularly to ensure the control box and the pad is intact. Use of the pad should be stopped immediately if any wiring is exposed or the pad / control box is compromised.
- The heated electric pad should be kept in a cool dry place. DO NOT EXPOSE THE CONTROLLER, CABLES OR HEATING PAD TO WATER, MOISTURE OR HIGH HUMIDITY
- Sponge clean only with luke-warm water and neutral detergent. Rinse with clean water by sponging and allow to air dry before re-use. WARNING: do not allow the electric control to get wet while cleaning. Inner: polyester. Outer: polar fleece.

SPECIFICATIONS

Electric heating Pad – Therapeutic Heat
Long lasting heat – provides great relief for muscular pain
2 Temperature settings: High, Low
LED indicator on switch
Over heating protection system
Fixed connector + VDE plug
Removable and washable cover
Material: Cotton
Size: 30*40cm
Voltage: 220V-240V~ 50-60HZ
Power: 45W

DECLARATION OF CONFORMITY

Manufacturer Name: Satyatrade S. L.
Address: Pol.Ind.La Raya. C/ Guadalquivir, 2.
Camarma de Esteruelas, 28816 Madrid
Tel: 902 430 967
Fax: 91 8864285
NIF: B83254763
Manufacturing country: China
Sytech® SY-AL30 Complies with the following standards:
EMC Standards: Directive 2004/108/CE
LVD Standards: 2006/95/EC
ROHS standard: 2011/65/EC Restriction of Hazardous Substance

Signed: Ajeet Nebhwani Utamchandani
Sole Administrator

