



SURPLUS MEDICAL EQUIPMENT

FROM

THE SWEDISH NATIONAL BOARD OF HEALTH AND WELFARE



SEMS AB
Henrikdalsvägen 2
SE-386 92 FÄRJESTADEN

Telefon/Phone:
0485-48480
Int. + 46 485 48480

Telefax:
0485-48490
Int- + 46 485 48490

Org.nummer
556576-3660
VAT:SE556576366001

Bankgiro
5547-6501
E-mail: sems@sems.se

Postgiro:
2 49 99 -5

<i>Name/Function</i>	Autoclave
<i>Manufacture typ</i>	Getinge GE-406
<i>Manufacture</i>	Getinge
<i>Made in</i>	Sweden
<i>Part number.</i>	81136000
<i>Production year</i>	1982
<i>Comment´s</i>	Needs a connection to water pipe and outflow.



Contents list: Autoclave 81136000

<u>Pairt number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
Edition 1 1982	Contents list	1	
	Instruction manual for Autoclave GE-406	1	
	Test certificate from the government inspection unit	1	
	Autoclave Getinge GE-406	1	
	Upper grating	1	
	Lower grating	1	
460 40 21-48	Graph paper pack of 100 pcs	1	
470 26 96 -01	Tools for disassembling pumps	1	

Technical data

Width cabinet type	750 mm	Steam consumption	
Width recessed type	1050 mm	/process(135 textiles progr)	4 kg
Height cabinet type	1660 mm	peak value	30kg/h
Height recessed type	1700 mm	Pilot voltage	24 VAC
Depth	930 mm	Fuse pilot circuits 5 x 20 mm	10 A
Loading height	800 mm	Electricity alternative 1	
Weight	360 kg	3 x 380 volts + Beutral, 50 Hz	
Chamber width	455 mm	Electricity alternative 2	
height	460 mm	3 x 220 volts 50 Hz	21 kW
depth	660 mm	Electric rating incl. steamgenerator	0,9
corner radius	75 mm	Electric rating exkl. steamgenerator	kW
volume	135 litre	Wlectric energy consumption/process	
design pressure	300 kpa	Incl. steam generator	3,6 kWh
material SIS 2347		Steam generator	
corresponding to AISI 316		Gross volume	21 litre
Vacuumpump electric rating	0,75 kW	Water volume	10 litre
Water consumption at 110 mbar	3/l/min	Design pressure	350 kpa
Feedwater pump electric rating	0,37 kW	Heat losses	0,5 kW
Feedwater pumpcapacity at 400 kpa		Electric rating	19 kW
Conterpressure	5,6 l/min		
Heat losses,stream generator			
equipped	about 2 kW		
Water consumption/process	200 litre		
Peak value	10 l/min		

<i>Name/Function</i>	Autoclave
<i>Manufacture typ</i>	Webeco
<i>Manufacture</i>	Webeco
<i>Made in</i>	Germany
<i>Part number.</i>	81136200
<i>Production year</i>	1970
<i>Comment's</i>	Does not need water and outflow installation. If power failure should occur, it is possibility to use the enclosed paraffin-stove for heating.



Contents list: Autoclave 81136200

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
M7102-541030-6	Contents list	1	
	Instruction Manual, SoS	1	
	Sterilization chart, SjøvS	1	
	Autoclave Webeco	1	
	Paraffin stove	1	
	Measuring cup	1	
	Holder	1	
	Base plate	1	
	Tray	1	
	Spade	1	
	Cassette, small	1	
	Cassette, large	1	

Technical data

Work pressure and temperature 20kPa(2kg/cm ²); 130C, for instrument, syringe, bandage etc. 100kPa(1kg/cm ²); 120C, for rubber articles, plastics and bottles with solvents.	Outside measurements Breadth: 43 cm Height: 71 cm (including the bottom plate) Deep: 75 cm (total deep which includes the lock device on the door.
Camber Breadth: 16,2 cm Height: 16,2 cm Deep: 49,5 cm	Fittings Plate 45x15x3 cm Cassette 45x14x5 cm Cassette 45x14x12cm Bottom plate Insert rack Spade
Weight (With all fittings) 71 Kg	Glass gauge Paraffin stove
Running-options Electricity; 3kW fuse 16A 220V	
Paraffin: 2,1 kW paraffin stove with three burner.	Spanner three cleaning needles Spare part: Door gasket

<i>Name/Function</i>	Defibrillator
<i>Manufacture typ</i>	LifePac 10
<i>Manufacture</i>	Physicontrol
<i>Made in</i>	USA
<i>Part number.</i>	1168100
<i>Production year</i>	1993
<i>Comment's</i>	With cardioscope display



Contents list: Defibrillator 1168100

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual Lifepac 10	2	
	Instruction Manual AC Auxiliary power module	1	
	Book about defibrillation- What you should know	1	
80 42 00-37	Defibrillator, Lifepac 10	1	
80 42 17-07	AC Auxiliary power module (Charger & Power supply)	1	
09-10 424-09	Battery, fastpac	3	
	Main cabel 220 VAC	1	
80 54 00 01	ECG-cable	1	
	Accessories bag (velcro)	2	
80 04 18	Child size paddles (for attachment to the regular paddles)	2	
8047 00-003	ECG-paper	2	rolls

Technical data

Energy select	Disarm 5,10,20,50,100,200,300,360
Charge time	Charge to 360 J in less than 12 seconds
Paddle area	82 cm ²
Cord Length	2.3 m
Synchronizer	Ensures energy discharge within 20 msec of sync marker on cardioscope (triggers to patient generated QRS complex)
Paddle Controls	Charge(with unticator light) Record (activates strip chart recorder) Energy select (rotating dial; disarm to 360J) Shock (dual energy discharge buttons, one on each paddle)
Battery	3 NiCad batteries, 12 V 1.o amp hours each A new Battery pak registering at least 100% capacity on the Battery support system will minimally provide: -45 minutes of monitoring, or -30 minutes of pacing, or -25 discharges at 360 joules per battery pak
Height	10,4 cm
Width	40,6 cm
Depth	37 cm
Weight	9 kg

<i>Name/Function</i>	Diathermy
<i>Manufacture typ</i>	Elektrotom 200
<i>Manufacture</i>	Berchtold
<i>Made in</i>	Germany
<i>Part number.</i>	81167300
<i>Production year</i>	1993
<i>Comment's</i>	Has got a cutting and spray function. Also bipolar function for use with forceps.



Contents list: Diathermy 81167300

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual	1	
	Elektrotom 200 Diathermy apparatus	1	
	Main cable 220 V	1	
	Ground, yellow/green	1	
	Neutral electrode, multiple use	1	
	Cable to neutral electrode	1	
	Handle with cable	1	
	Blade electrode, straight	1	
	Spark electrode, ø 4mm	1	
	Tweezers	1	
	Cable to tweezers	1	

Technical data

Power requirement	110 / 220 V AC +/- reconnectable
Current consumption	max 440 VA
Nominal frequency	500 kc/s
Impulse frequency of modulation	33 kc/s
CREST-FACTOR	
Current type spray-coagulation	7.6
Current type contact- coagulation	3.0
Current type pure cut	1.6
Current type blend I	2.2
Current type blend II	2.6
Current type blend III	3.0
Current type bipolar coagulation	1.8
HF-OUTPUT	
Current type spray-coagulation	80 Watts +/- 20 % on 300 Ω
Current type contact- coagulation	175 Watts +/- 20 % on 200 Ω
Current type pure cut	175 Watts +/- 20 % on 300 Ω
Current type blend I	175 Watts +/- 20 % on 300 Ω
Current type blend II	175 Watts +/- 20 % on 300 Ω
Current type blend III	175 Watts +/- 20 % on 300 Ω
Current type bipolar coagulation	50 Watts +/- 20 % on 100 Ω
Weight	8,8 kg
Dimensions	130 mm height
	400 mm width
	380 mm depth

<i>Name/Function</i>	Development machine
<i>Manufacture typ</i>	Curix 60
<i>Manufacture</i>	Agfa-Gevaert
<i>Made in</i>	Germany
<i>Part number.</i>	1169600
<i>Production year</i>	1994
<i>Comment´s</i>	Film and development equipment are not included



Contents list: Development machine 1169600

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual Agfa, Curix 60	1	
	Developing machine Agfa, Curix 60	1	
	Light-proof lid	1	
7 9462 1081 0	Bottle 2,5 litres	6	
8 9432 1090 4	Valve for bottle	4	
	Cork for bottle	3	

Technical data

Film range/films sizes	All standard brand films of the size 10 x 10 cm – 35 x 43 cm which are suitable for machine processing.
Hourly output	60 films/h of size 24 x 30 cm
Heating time	About 7 min
Processing time	3 min
Processing rate	28 cm/min, 20 cm/min with transformation kit
Chemicals	
Developer solution	G 153
Fixer solution	G 353 and G 354
Contents of replenishment tanks	2.5 l each
Contents of baths	0,91 l each
Bath temperature	
Developer	34° (ex works)
Fixer	34° (ex works)
Water	Not heated
Drying system	Infrared drying with 6 drying levels, Air throughput about 80 m ³ /h
Heat emission	About 2100 kJ/h in continuous Operation with the dryer set at 4
Dimensions	
Width	650 mm
Length	940 mm incl. feed table
Height	370 mm incl. feet
Base area	650 x 800 mm
Weight	54 kg empty, 64 kg full
Operating ambient temperature	15° C to 30° C
Operating humidity	25% to 80%


<i>Name/Function</i>	X-ray and development equipment
<i>Manufacture typ</i>	Curix 60
<i>Manufacture</i>	Agfa-Gevaert
<i>Made in</i>	Germany
<i>Part number.</i>	1169500
<i>Production year</i>	1994
<i>Comment's</i>	Cassette Kodak, X-Omatic LW, Medium & Regular



Contents list: X-ray and development equipment 1169500

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
00 11 009	X-ray cassette 35 x 43 cm	4	
00 11 010	X-ray cassette 24 x 30 cm	6	
00 11 011	X-ray cassette 18 x 24 cm	4	
9196882+6802532	Filter box and filters 35 x 43 cm	1	
9196825+6802474	Filter box and filters 24 x 30 cm	1	
	Label SIN	5	
	Label DX	5	
91 98 979	X-ray apron	2	
66 44 447	X-ray interpretation cupboard	1	

<i>Name/Function</i>	Oxygen supply system
<i>Manufacture typ</i>	
<i>Manufacture</i>	Aspira Medical AB
<i>Made in</i>	Sweden
<i>Part number.</i>	1169700
<i>Production year</i>	1996
<i>Comment's</i>	Support as maximum 16 patients with air from the compressor and oxygen from gas bottle.



Contents list: Oxygen supply system 1169700

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual	1	
	Bottle regulator	10	
	Gas loop consisting of 4 blocks	4	
	Extension tube air pressure	2	
	Extension tube oxygen	2	
	Spray nozzle	10	
	Gas block, air pressure	3	
	Tube, air pressure		
	Tools		
	Spare parts		

<i>Name/Function</i>	Floor lamp, surgery
<i>Manufacture typ</i>	B 250 SS
<i>Manufacture</i>	Martin
<i>Made in</i>	Germany
<i>Part number.</i>	1144900
<i>Production year</i>	1981
<i>Comment´s</i>	2 lamps packed in a box



Contents list: Floor lamp, surgery 1144900

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	2	
	Instruction Manual	2	
	Lamp body with halogen lamp	2	
	Stand with transformer and mains connection	2	
	Handle, washable in autoclave	2	
	Spring	2	
	Foot	2	
	Halogen lamp, extra 24 V/50W	4	
	Glass fuse F 2,5 A	8	
	Allen key 6 mm	2	

Technical data

Bulb	Top sealed 24 V/ 50 W Two lamps in reserve are enclosed
Light intensity	36.000 lux
Transformer	220V/24 V 60 VA S-marked Double isolated (class II) 2,5 A fuse (secondary) Five fuses in reserve are enclosed
Handle	Can be sterilised in chemical disinfection or autoclaved in 120C.
Able to focus	
Weight	27 Kg

<i>Name/Function</i>	Floor lamp, examination
<i>Manufacture typ</i>	
<i>Manufacture</i>	Kifa
<i>Made in</i>	Germany
<i>Part number.</i>	81144800
<i>Production year</i>	1970
<i>Comment´s</i>	



Contents list: Floor lamp, examination 81144800

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Floor stand	1	
	Stand	1	
	Lamp body with arm and electrical cable	1	

<i>Name/Function</i>	Lamp universal
<i>Manufacture typ</i>	3060-50000
<i>Manufacture</i>	Hanau
<i>Made in</i>	Germany
<i>Part number.</i>	1164000
<i>Production year</i>	
<i>Comment's</i>	Can be mounted on the wall or the ceiling



Contents list: Lamp universal 1164000

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Template for drill holes in ceiling	1	
	Ceiling attachment casing	1	
	Attachment with arm	1	
	Adjustable arm	1	
	Lamp housing	1	
	Transformer Tufvasson 240V-24 V AC, protective transformer PVS 75	1	
	Ziploc bag containing key and bolt for adjusting the spring pressure in the arm between the lamp and the ceiling attachment.	1	

<i>Name/Function</i>	Cardioscope
<i>Manufacture typ</i>	Diascope 1
<i>Manufacture</i>	S & W
<i>Made in</i>	Denmark
<i>Part number.</i>	1167400
<i>Production year</i>	1993
<i>Comment's</i>	Singel Channel ECG Monitor



Contents list: Cardioscope 1167400

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual	1	
	Diascope 1	1	
	Main cable 220 V	1	
	Cords for disposable electrodes	1	
	Cords for multiple use electrodes		
	Multiple use electrodes		
	Tape for multiple use electrodes	1	Pack 500 pcs

Technical data

<u>Dimensions</u> H: 125 x W: 200 x D: 270 mm <u>Weight</u> 4.3 kg <u>Power requirements</u> 220-240 VAC +/-10% 50/60 Hz or 100-117 VAC +/-10% 50/60 Hz <u>Consumption</u> 24 VA <u>Operating temperature</u> +5 C to + 40 C <u>Display</u> CRT screen dim. H 80 mm, W: 100 mm. Display speed: 25 or 50 mm/s. Max. display gain: 1 mV input equals 2 cm on the screen <u>Digital display</u> 3 digit LED Digit size: H: 11 mm, W 6 mm. Updating cycle: 2.1 s simultaneously <u>Cathode ray tube</u> Size: 5.5" (80 x 100 mm). Mode: Non-fade, rolling mode. <u>Light indications</u> Alarm off. <u>Controls</u> Stop/run for freezing and starting the ECG. Alarm on/off for enabling or disabling the alarm. Main function selector for selecting the input of the digital display, ECG LOW/HIGH alarm. Up/down button for increasing/decreasing the parameters selected by the main functions selector. QRS-volume/ECG trace for adjusting the QRS-beep and the ECG trace speed. Mains on/off. Soft functions test selector. <u>Connectors</u> ECG input: Colour coded S&W standard connectors Interface output: 15 way subminiature D-type	<u>ECG Amplifier</u> <u>Input</u> Mode: Differential Impedance: 10 mOhm DC, 4 mOhm at 35 Hz Difference input max: 360 mV CMRR: 80 db Defi shock protection: 5 KV Filter 50/60 Hz: 40 db Defibrillator marker, input typ: 1 V/100 ms <u>Amplification</u> Electrosurgery noise filter Frequency response: 0.3 – 35 Hz (-3 db) Overload recovery time 6 s. <u>Pace discriminator</u> Trigger level: +/-2.3 mV, slew rate sensitive <u>Output</u> Heart rate test generator: 120 beats/min. Analogue heart rate: 10mV/beat. ECG direct (1 mV inp.): 1 Vpp ECG delayed (1 mV inp.): 1 Vpp Alarm: =V when active, -11 V when off Alarm block: 0 V when active, -11 V when off 0 V Ds 1 on: +5 V when on, -11 V when off. <u>Memory back-up</u> The values set by the operator are retained During short power-off periods, i.e ECG-gain, alarm on/off, alarm limits The memory time is minimum 15 s. <u>Patient safety system</u> Leakage current between chassis and any patient terminal at 50 Hz <0.04 uA/V <u>QRS detector</u> Trigger level 0.3 mV (ref. input). <u>Tachometer</u> Mode: Averaging. Range: 20-300 beats/min. <u>Classification</u> IEC 6001 class 1 Type CF with defibrillator protection Vibration IEC 68
---	--

<i>Name/Function</i>	Air compressor, artificial respiration
<i>Manufacture typ</i>	
<i>Manufacture</i>	Biab
<i>Made in</i>	Sweden
<i>Part number.</i>	81150300
<i>Production year</i>	1982
<i>Comment's</i>	The compressor produces clean air to be used with the Ventilator UV 705 or other medical air machine. At least ten Ventilators can be used at the same time.



Contents list: Air compressor, artificial respiration 81150300

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual National Board of Health and Welfare	1	
	SA-certificate for first inspection	1	
	SA-certificate pressure vessel	1	
	Compressor Biab/Haug 75/100	1	
	Surveyance apparatus	1	
	Intake air filter	2	
	Piston rings	2	set
	Piston seal rings	2	
	Valve seals	2	
	Hook wrench 110mm – 115mm	1	

Technical data

Storage mark	M2350-815511-4	
Storage name	Kompaggr 1 vent MT	
Manufacturer	BIAB Tryckluft AB	
Manufacture mark	BIAB/HAUG 75/100	
Voltage demands	3x 220V 50 Hz	
Fuse demands	16 A	
Dimensions:	Without the box	With the box
Length	970 mm	120mm
Breadth	710 mm	800 mm
Height	830 mm	1420 mm
Weight	165 kg	280 kg
Compressor, type	HAUG BOZ 75-D4	
Cylinder, number	2	
Free stated air quantity at four bar	4,1 l/s, 245l/min	
Maximum work pressure	0,8 MPa, 8 bar	
Engine:		
Type	Elmo RA 100LB/4-2769	
Effect	2,2 kW	
Number of revolution	1450 rpm	
Voltage	380/220 V +- 25%	
Frequency	50 HZ +- 6%	
Cooling method	IP 54	
Air container		
Volume	100l	
Calculations-and maximum pressure	0,8 MPa, 8 bar	
Testing pressure	1,1 MPa, 11bar	
Surface treatment	Silicone inside and outside	
Muffs and plugs	Stainless	

<i>Name/Function</i>	Medicine technical service kit
<i>Manufacture type</i>	
<i>Manufacture</i>	
<i>Made in</i>	
<i>Part number.</i>	1169900
<i>Production year</i>	1995
<i>Comment's</i>	The service kit contents test instruments. About sixty different tools and general spare parts for the medical technician.



Contents list: Medicine technical service kit 1169900

<u>Number</u>	<u>Description</u>	<u>Qty</u>	<u>Number</u>	<u>Description</u>	<u>Qty</u>
11 69 940	Test pin, red, for 4 mm B	1	11 69 901	CSC 5-56 400 ml	1
11 69 941	Trimmer	1	11 69 902	Work top soft	1
11 69 942	Cutting pliers	1	11 69 903	Hammer drill MAK	1
11 69 943	Vice universal 2B	1	11 69 904	Drill kit Projahn	1
11 69 944	Water level Typ 446 mm	1	11 69 905	Installation drill kit	1
11 69 945	Electrical tape PVC	5	11 69 906	Thermometer DIDI-Therm	1
11 69 946	Thread tape PTFE	2	11 69 907	ESD Wrist strap	1
11 69 947	Seam tape type 1160	3	11 69 908	Flash light Mini- Maglite	1
11 69 948	Multiway socket with 6 sockets	1	11 69 909	Snipe nose pliers 7890	1
11 69 949	Lubricating grease silicone 7220	1	11 69 910	Threading kit Projahn	1
11 69 950	Tin soldering Multicore 2,3 mm 500G	1	11 69 911	Hole saw ca 40 – 83 mm	1
11 69 951	Pincers 200 mm	1	11 69 912	Instrument skrewdriver kit 1 – 3	1
11 69 952	Rubber seal kit Loctite 1622	1	11 69 913	Cable clip kit with pliers	1
11 69 953	Light machine oil 100 ml ME	1	11 69 914	Lab cable red 100 cm	4
11 69 954	Compass saw blade	3	11 69 915	Lab cable black 100 cm	4
11 69 955	Steel wire 1.0 mm 600 G	1	11 69 916	Lab cable red 200 cm	2
11 69 956	Protective goggles	1	11 69 917	Lab cable black 200 cm	2
11 69 957	Block gas seal bell	1	11 69 918	Soldering iron 300 W	1
11 69 958	Cleaning valve	1	11 69 919	Soldering station Temtronic WE	1
20 00 730	Tool equipment RePairatör	1	11 69 920	Tin solder 1.0 mm 500G	1
90 00 479	Service Manual Servoventilator	1	11 69 921	Mini pliers CHY 932 C 600	1
90 00 480	Service Manual FX-3101	1	11 69 922	Fluke Scope Master 99B	1
90 00 481	Service Manual Diascope 1	1	11 69 923	Multimeter CHY 14 INKL	1
90 00 482	Service Manual Lifepak 10	1	11 69 924	Needle file kit 470 Sandvik	1
90 00 483	Service Manual Curix 60	1	11 69 925	Pliers 426G 0 – 22	1
90 00 484	Service Manual Pneupak ventilators	1	11 69 926	Sidecutting pliers 7190	1
90 00 485	Service Manual Elektrotom 200	1	11 69 927	Extension cord grounded 25 m	1
90 00 486	Service Manual Pulsoximeter 515 A	1	11 69 928	Screwdriver Stubby	1
90 00 487	Service Manual Egnell Compact 30	1	11 69 929	Spanner 6 OC for hose clamps	1
90 00 488	Users manual FX 3101	1	11 69 930	Hose clamps, kit 8 – 14,4	5
90 00 489	Users manual Pulsoximeter 515 A	1	11 69 931	Stress detector	1
90 00 490	Spare parts Autoklav	1	11 69 932	Service box with plastic straps	1
90 00 491	Spare parts AGFA Curix	1	11 69 933	Fuse selection Snabb	1
90 00 492	Spare parts Mobilett	1	11 69 934	Tin vacuum Val-U-Vac	1
			11 69 935	Avlödningsfläta 1,5 – 2	2
			11 69 936	Test pin with grabber R	1
			11 69 937	Test pin with grabber S	1
			11 69 938	Test pin with clamp	1
			11 69 939	Test pin with clamp	1

<i>Name/Function</i>	Anaesthetic apparat.
<i>Manufacture typ</i>	AGA military
<i>Manufacture</i>	AGA
<i>Made in</i>	Sweden
<i>Part number.</i>	81153800
<i>Production year</i>	1964
<i>Comment's</i>	It can be handled by hand or together with Ventilator UV 705.



Contents list: Anaesthetic apparat. 81153800

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Documentation		
	Contents list	1	
	Instruction Manual Anaesthesia apparatus	1	
	Instruction Manual Ventilator UV 705	1	
	Maintenance instructions anaesthesia gases		
	Penlon OMV-50 (Triservice)	1	
	Anaesthesia apparatus		
	Anaesthesia apparatus AGA Anestor Militar	1	
	Valve and soda lime container holders	1	
	Soda lime container	2	
	Plastic tube 8/12 x 500 with angle pieces	1	
	Inhalation valve	1	
	Exhalation valve	1	
	Berner valve with ISO cone	1	
	Potential equilibration valve	1	
	Revivator and rubber parts		
	Valve unit revivator	1	
	Plastic distance piece	2	
	Respiration mask no 3	2	
	Y-piece for respiration mask	1	
	Neck strap and metal holder for respiration mask	1	
	Respiration bladder 2, 3, or 4 litres	2	
	Rubber bladder self-filling	1	
	Rubber bladder folded 250 x 22 mm	2	
	Rubber bladder folded 1050 x 22 mm	2	
	Rubber bladder folded 1050 x 26 mm	1	
	Anaesthetic scale		
	Halotane/Trilene		
	Enflurane		
	Isoflurane		
	Penthrane/ Chloroform		

Technical data

AGA Militar	
Weight: Anaesthetic apparatus	40 kg
gas bottles, totally	43 kg
Volume: box	0,2 m3

<i>Name/Function</i>	Operating table
<i>Manufacture typ</i>	OP 1650
<i>Manufacture</i>	Merivaara
<i>Made in</i>	Finland
<i>Part number.</i>	1167600
<i>Production year</i>	1992
<i>Comment's</i>	Positioning functions are: elevation, trendelenburg position, lateral-tilt position, head section position, back section position and foot/leg section



Contents list: Operating table 1167600

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Users Manual	1	
	Service Manual	1	English
	Operating instructions	1	English
	List of accessories	1	
	Basic equipment		
	Operating-table	1	
	Headrest and mattress	1	
	Dividable footrest and mattress	1	
	Mattress for seat and back	1	
	Narcosis stand with vice	1	
	Armrest with vice and mattress	1	
	Wristrest with vice	1	Pair
	Legrest with vice	1	Pair
	Accessories		
	Legrest with accessories	1	Pair
	Rinser with accessories	1	
	Locating screw for armrest	1	
	Locating screw for footrest	2	
	Locating screw for rinsing bowl	1	
	Locating screw for headrest	1	Pair
	Locating screw for leg stirrups	2	
	Locating screw for arm straps	2	
	Transport lock	1	

Technical data

General

Model OP 1650 is manually operated, general purpose surgery table.

The elevation system is actuated by a foot operated, manual/hydraulic pump.

All table posturing with the exception of the foot/leg section is operable from the head end of the table.

Floor Lock/Brake System

The floor lock/brake system consists of two jacks which raise and support the table base off from the swivel castor.

Elevation System

The table top elevation cylinder is actuated by a foot operated, manual- hydraulic pump.

Elevation adjustment range is 730 – 1040 mm.

<i>Name/Function</i>	Operating table, casualties
<i>Manufacture typ</i>	OP 1650
<i>Manufacture</i>	Merivaara
<i>Made in</i>	Finland
<i>Part number.</i>	1167800
<i>Production year</i>	1992
<i>Comment's</i>	Positioning functions are: elevation, trendelenburg position, lateral-tilt position, head section position, back section position and foot/leg section



Contents list: Operating table, casualties 1167800

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Instruction Manual	1	
	Service Manual	1	English
	Operation Instructions	1	English
	List of accessories	1	
	Basic equipment		
	Operating-table	1	
	Headrest and mattress	1	
	Dividable footrest and mattress	1	
	Mattress for seat and back	1	
	Narcosis stand with vice	1	
	Armrest with vice and mattress	1	
	Wristrest with vice	1	Pair
	Legrest with vice	1	Pair
	Accessories		
	Legrest with accessories	1	Pair
	Rinser with accessories	1	
	Locating screw for armrest	1	
	Locating screw for footrest	2	
	Locating screw for rinsing bowl	1	
	Locating screw for headrest	1	Pair
	Locating screw for leg stirrups	2	
	Locating screw for arm straps	2	
	Transport lock	1	

Technical data

General

Model OP 1650 is manually operated, general purpose surgery table.

The elevation system is actuated by a foot operated, manual/hydraulic pump.

All table posturing with the exception of the foot/leg section is operable from the head end of the table.

Floor Lock/Brake System

The floor lock/brake system consists of two jacks which raise and support the table base off from the swivel castor.

Elevation System

The table top elevation cylinder is actuated by a foot operated, manual- hydraulic pump.

Elevation adjustment range is 730 – 1040 mm.

<i>Name/Function</i>	Operating table with extra accessories
<i>Manufacture typ</i>	OP 1650
<i>Manufacture</i>	Merivaara
<i>Made in</i>	Finland
<i>Part number.</i>	1167700
<i>Production year</i>	1992
<i>Comment's</i>	Positioning functions are: elevation, trendelenburg position, lateral-tilt position, head section position, back section position and foot/leg section



Contents list: Operating table with extra accessories 1167700

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Users Manual	1	
	Service Manual	1	English
	Operating instructions	1	English
	List of accessories	1	
	Basic equipment		
	Operating-table	1	
	Headrest and mattress	1	
	Dividable footrest and mattress	1	
	Mattress for seat and back	1	
	Narcosis stand with vice	1	
	Armrest with vice and mattress	1	
	Wristrest with vice	1	Pair
	Legrest with vice	1	Pair
	Accessories		
	Legrest with accessories	1	Pair
	Rinser with accessories	1	
	Locating screw for armrest	1	
	Locating screw for footrest	2	
	Locating screw for rinsing bowl	1	
	Locating screw for headrest	1	Pair
	Locating screw for leg stirrups	2	
	Locating screw for arm straps	2	
	Transport lock	1	
	Extra accessories		
	Shoulder / side protection with vice	2	

Technical data

General

Model OP 1650 is manually operated, general purpose surgery table.

The elevation system is actuated by a foot operated, manual/hydraulic pump.

All table posturing with the exception of the foot/leg section is operable from the head end of the table.

Floor Lock/Brake System

The floor lock/brake system consists of two jacks which raise and support the table base off from the swivel castor.

Elevation System

The table top elevation cylinder is actuated by a foot operated, manual- hydraulic pump.

Elevation adjustment range is 730 – 1040 mm.

<i>Name/Function</i>	X-ray machine mobile
<i>Manufacture typ</i>	Mobilett
<i>Manufacture</i>	Siemens-Elema
<i>Made in</i>	Germany
<i>Part number.</i>	1169400
<i>Production year</i>	1994
<i>Comment's</i>	



Contents list: X-ray machine mobile 1169400

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual Siemens-Elema	1	
	Service Manual Siemens-Elema	1	
	Mobilett Plus	1	
	Transport lock lever	1	
	Exposure trigger	1	
	Angle gauge	1	
	Wrench 16 mm och 18 mm	1	

Technical data

See separate documentation

<i>Name/Function</i>	Suction pump intensive care
<i>Manufacture typ</i>	IVA
<i>Manufacture</i>	Einar Egnell AB
<i>Made in</i>	Sweden
<i>Part number.</i>	81160000
<i>Production year</i>	1981
<i>Comment's</i>	



Contents list: Suction pump intensive care 81160000

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual Pump Egnell	1	
	Assembly instructions Pump Egnell	1	
Model 23	Suction pump, Egnell Minor	1	
	Stand	1	
	Floor stand	1	
	T-piece		
	Holder for 2-litre bottle	1	
17-725	Bottle 2 litres	2	
	Holder for 0,7 litre bottle	1	
17828	Bottle 0,7 litres	2	
17413 S	Lid with overflow protection	1	
17413	Lid	1	
	Filter container	1	
	Connecting tube 7 -10 mm	1	
	Connecting tube 10 -10 mm	1	
	Hose 1,7 m. 8 x 12 mm with 2 female sockets (17 397)	2	
	Hose 30cm. 6 x 8 mm with 2 female sockets (17 397)	4	
421/5	Electrical cord, cast 5 m	1	
FP 50	Filter paper ø 70 mm	1	Pack of 50
	Rubber seal 12,1 x 1,6 mm (Inside suction pump)	4	
	Cotton inlay (Inside suction pump)	1	
	Molykote BR 2 tube. (Inside suction pump)	1	
1198	Rubber non-return valve (Inside suction pump)	4	
	Allen key 6 mm	1	
	U-wrench 17 mm (Inside suction pump)	1	

Technical data

Suction capacity	Approximately 18 litres of air per minute
Suction intensity	Infinitely variable from 0.1 to 0.9 kg/cm ² (10-90 kPa or 75-685 mm Hg)
Power supply	200 – 250 V. 50 c/s
Rating	approx 90 W
Motor speed	Approx 144 r.p.m.
Dimensions	
Length x width x height	400 x 210 x 970 mm
Wheelbase	390 x 390 mm
Weight	approx 18kg

<i>Name/Function</i>	Suction pump surgery
<i>Manufacture typ</i>	Fricare
<i>Manufacture</i>	Siemens-Elema
<i>Made in</i>	Germany
<i>Part number.</i>	81160100
<i>Production year</i>	1981
<i>Comment's</i>	



Contents list: Suction pump surgery 81160100

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual, Fricar suction pump by Siemens-Elema	1	
	Fricar suction pump	1	
	Sealed lid with overflow valve	1	
	Plastic tube (assembled)	1	
	Secretion bottle 4 litres	2	
	Overflow bottle 0,3 litres	1	
	Lid for overflow bottle	1	
	Plastic tube 2 meters	1	
	Timing belt	2	
	Seal for suction bottle	1	
M6141-002010	Phillips skrewdriver Ph no 1	1	
	U-wrench 17 mm	1	
	Allen key 4 mm	1	

Technical data

Motor

Short-circuit one phase motor.


Unit

Rotation unit without valves with automatic ring lubricates.

Net

Alternating current only. Make sure that the voltage agrees with the unit's voltage mark. The current consumption is 500W at 220V. Connect to grounded wall socket.

<i>Name/Function</i>	Ceiling lamp, surgery
<i>Manufacture typ</i>	
<i>Manufacture</i>	Martin
<i>Made in</i>	Germany
<i>Part number.</i>	1160800
<i>Production year</i>	1981
<i>Comment´s</i>	The handle can be disinfected or autoclaved.



Contents list: Ceiling lamp, surgery 1160800

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	
	Instruction Manual	1	
	Lamp body with 2 pcs halogen lamps	1	
	Fittings support	1	
	Handle, washable in autoclave	2	
	Fittings handle	1	
	Ceiling attachment with 6 st nuts and bolts	1	
	Transformer 220/2v V	1	
	Ceiling attachment casing	2	
	Halogen lamp, extra	2	
	Glass fuse F6, 3A	4	
	Allen key 2 mm	1	

Technical data

Lamp type	Two halogen lamps (one in reserve) 24V/100W
Light intensity	Two extra lamps are enclosed 70.000 lux (variable)
Transformer	220V/24V 150 VA S-marked Double isolated (class II) 6,3 A fuse (secondary) Five fuse in reserve are enclosed
Handle	Two sterilised handles are enclosed.
Able to focus	
Weight	28 kg

<i>Name/Function</i>	Ventilator
<i>Manufacture typ</i>	UV 705 std
<i>Manufacture</i>	AGA
<i>Made in</i>	Sweden
<i>Part number.</i>	1162900
<i>Production year</i>	1981
<i>Comment's</i>	Can be used with the anesthetic apparatus AGA 1153800.



Contents list: Ventilator 1162900

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	In lid
	Instruction Manual UV 705, SoS	1	In lid
M7835-292010-4	AGA ventilator UV 705	1	Box A2
F3600-012049-5	Power eliminator for 220 V AC	1	Box A
M7835-292028-6	Patient unit	1	Box A2
F3600-012074-3	Manometer for patient pressure	1	Box A3
	Overflow valve with expansion hose for the circular system	1	Box A6
	Removable writing board	1	Box A8
F3600-012070-1	Compressed air tube 8/14 quick-release, 3 m	1	Box B
	Five-footed stand with feet	1	Box C
	Lower stand with holder	1	Box C
	Upper stand with storage platform	1	Box C
	Attachment for anaesthesia apparatus	2	Box B
M7835-990398-8	Rubber bladder 3 litres	2	Box B
F3600-0120545	Grommet 30 mm (silicon)	2	Box B
M7835-990288-1	Connector piece, plastic ISO	1	Box B
M1536-165308-6	Rubber seal 10. 3 x 2,4	2	Box B
M1536-165405-0	Rubber seal 22. 2 x 2,3	2	Box B
F3600-012073-5	Expansion hose PVC 2/4, 1 meter	1	Box B
F3600-012076-8	Simmering	1	Box B
M2486-840334	Fuse 1A	10	Box B
M2486-840806-9	Fuse 100mA	10	Box B
F3600-012052-9	Membrane bladder 2 x 10 kit	2	Box B
FRP-95-235			
(Wilkerson)	Filter for steam dryer	1	Box B
GRP-96-310			Box B
(Wilkerson)	Holder for steam dryer	1	

Technical data

Compressed air pressure	0,2-0,8 Mpa (2-8 bar)
consumption	
at re-breathing	= minute volume
at non-re-breathing (normal adjustment)	about 2 times minute volume
Electricity supply	220 V 50 Hz
Running time with battery	24 hours
Breathing frequency	8-40/min
Inspirations time, in percent of time per breath	10%-90%
Minute volume	
at re-breathing	2-18 l/min
at non-re-breathing	0-23 l/min
Emptying pressure (non-re-breathing)	70-100cm H ₂ O
Safety valve,	
adjustable	0-70H ₂ O
PEEP (final expiratory pressure at non-re-breathing)	0-20cm H ₂ O
Weight	9,7kg
Size	breadth 410 mm, deep 250 mm, height 155

<i>Name/Function</i>	Ventilator
<i>Manufacture typ</i>	UV 705 intensive care
<i>Manufacture</i>	AGA
<i>Made in</i>	Sweden
<i>Part number.</i>	1162700
<i>Production year</i>	1981
<i>Comment's</i>	Can be used with the anesthetic apparatus AGA 1153800.



Contents list: Ventilator 1162700

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	In lid
	Instruction Manual UV 705, SoS	1	In lid
M7835-292010-4	AGA ventilator UV 705	1	Box A2
F3600-012049-5	Current unit for 220 V AC	1	Box A
F3600-012059-4	Dosage unit for air and O2	1	Box A1
M7835-292028-6	Patient unit	1	Box A2
F3600-012074-3	Manometer for patient pressure	1	Box A3
F3600-012050-3	Valve for exhalation system	1	Box A4
M7835-990358-2	Y-piece, straight	1	Box A5
M7835-990448-1	Y-piece, angled	1	Box A5
	Overflow valve with expansion hose for the circular system	1	Box A6
	Removable writing board	1	Box A8
	Bed support holder	1	Box B
M7863-131010-1	Folded tube 105 cm	4	Box B
69 34 509 (Siemens- Elema)	Servo-humidifier 150	3	Box B
68 08 943 (Siemens-Elema)	Inlay for servo-humidifier	32	Box B
F3600-012070-1	Compressed air hose 8/14 with quick-release, 3 m	1	Box B
	Five-footed stand with feet	1	Box C
	Lower stand with storage platform	1	Box C
	Upper stand with storage platform	1	
	Attachment for anaesthesia apparatus	2	Box B
M7835-292068-2	Support for patient tubes	1	Box C
M7835-990398-8	Rubber bladder 3 litres	2	Box B
F3600-0120545	Grommet 30 mm (silicon)	2	Box B
M7835-990288-1	Connector piece, plastic ISO	1	Box B
M1536-165308-6	Rubber seal 10. 3 x2,4	2	Box B
M1536-165405-0	Rubber seal 22. 2 x 2,3	2	Box B
F3600-012073-5	Expansion hose PVC 2/4, 1 meter	1	Box B
F3600-012076-8	Simmering	1	Box B
M2486-840334	Fuse 1A	10	Box B
M2486-840806-9	Fuse 100mA	10	Box B
F3600-012052-9	Membrane bladder 2 x 10 kit	2	Box B
FRP-95-235 (Wilkerson)	Filter for steam dryer	1	Box B
GRP-96-310 (Wilkerson)	Holder for steam dryer	1	Box B

Technical data

Compressed air pressure	0,2-0,8 Mpa (2-8 bar)
consumption	
at re-breathing	= minute volume
at non-re-breathing (normal adjustment)	about 2 times minute volume
Electricity supply	220 V 50 Hz
Running time with battery	24 hours
Breathing frequency	8-40/min
Inspirations time, in percent of time per breath	10%-90%
Minute volume	
at re-breathing	2-18 l/min
at non-re-breathing	0-23 l/min
Emptying pressure (non-re-breathing)	70-100cm H2O
Safety valve, adjustable	0-70H2O
PEEP (final expiratory pressure at non-re-breathing)	0-20cm H2O
Weight	9,7kg
Size	breadth 410 mm, deep 250 mm, height 155

<i>Name/Function</i>	Ventilator + spiro
<i>Manufacture typ</i>	UV 705 intensive care
<i>Manufacture</i>	AGA
<i>Made in</i>	Sweden
<i>Part number.</i>	1162800
<i>Production year</i>	1981
<i>Comment's</i>	Can be used with the anesthetic apparatus AGA 1153800.



Contents list: Ventilator + spiro 1162800

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Note</u>
	Contents list	1	In lid
	Instruction Manual UV 705, SoS	1	In lid
M7835-292010-4	AGA ventilator UV 705	1	Box A2
F3600-012049-5	Power eliminator for 220 V AC	1	Box A
F3600-012059-4	Dosing equipment for air and O ₂	1	Box A1
M7835-292028-6	Patient unit	1	Box A2
F3600-012074-3	Manometer for patient pressure	1	Box A3
F3600-012050-3	Valve for exhalation system	1	Box A4
M7835-990358-2	Y-piece, straight	1	Box A5
M7835-990448-1	Y-piece, angled	1	Box A5
	Overflow valve with expansion hose for the circular system	1	Box A6
	Removable writing board	1	Box A8
	Bed support holder	1	Box B
M7863-131010-1	Folded tube 105 cm	4	Box B
69 34 509 (Siemens- Elema)	Servo-humidifier 150	3	Box B
68 08 943 (Siemens-Elema)	Inlay for servo-humidifier	32	Box B
F3600-012070-1	Compressed air hose 8/14 with quick-release, 3 m	1	Box B
	Five-footed stand with feet	1	Box C
	Lower stand with storage platform	1	Box C
	Upper stand with storage platform	1	
	Attachment for anaesthesia apparatus	2	Box B
M7835-292068-2	Support for patient tubes	1	Box C
M7835-990398-8	Rubber bladder 3 litres	2	Box B
F3600-0120545	Grommet 30 mm (silicon)	2	Box B
M7835-990288-1	Connector piece, plastic ISO	1	Box B
M1536-165308-6	Rubber seal 10. 3 x2,4	2	Box B
M1536-165405-0	Rubber seal 22. 2 x 2,3	2	Box B
F3600-012073-5	Expansion hose PVC 2/4, 1 meter	1	Box B
F3600-012076-8	Simmering	1	Box B
M2486-840334	Fuse 1A	10	Box B
M2486-840806-9	Fuse 100mA	10	Box B
F3600-012052-9	Membrane bladder 2 x 10 kit	2	Box B
	Spirometer US-1000	1	Box A10
	Folded tube with attachments for spiro	1	Box A
	Fluidistor	1	A7

Technical data

Compressed air pressure consumption	0,2-0,8 Mpa (2-8 bar)
at re-breathing	= minute volume
at non-re-breathing (normal adjustment)	about 2 times minute volume
Electricity supply	220 V 50 Hz
Running time with battery	24 hours
Breathing frequency	8-40/min
Inspirations time, in percent of time per breath	10%-90%
Minute volume	
at re-breathing	2-18 l/min
at non-re-breathing	0-23 l/min
Emptying pressure (non-re-breathing)	70-100cm H ₂ O
Safety valve, adjustable	0-70H ₂ O
PEEP (final expiratory pressure at non-re-breathing)	0-20cmm H ₂ O
Weight	9,7kg
H ₂ O Size	breadth 410 mm, deep 250 mm, height 155

<i>Name/Function</i>	Respirator
<i>Manufacture typ</i>	Engström 300
<i>Manufacture</i>	Gambro Engström AB
<i>Made in</i>	Sweden
<i>Part number.</i>	1156700
<i>Production year</i>	1981
<i>Comment's</i>	Can be used as a “stand alone” ventilator. It does no need air supply.



Contents list: Respirator 1156700

<u>Part number</u>	<u>Description</u>	<u>Qty</u>	<u>Part number</u>	<u>Description</u>	<u>Qty</u>
115 67 00	Respirator Engström 300	1	52-41554-00	Contents of rubber seal set	
			03-37177	Rubber seal	1
52-41553-00	Sterilized filter	1	03-37417	Rubber seal	1x3
52-41552-02	Sterilized filter with holder	1	03-37097	Rubber seal	1x2
	Thermometer Bennet	1	03-41601	Rubber seal	1
	Flex tube	1	03-37327	Rubber seal	1x2
52-46456	Oil	1	03-37147	Rubber seal	1
52-41554-00	Rubber seal set	3	03-37116	Rubber seal	1x8
52-41473-00	Rubber bladder yellow	1	03-37037	Rubber seal	1
52-41474-00	Rubber bladder svart	1	03-37466	Rubber seal	1
52-40551-00	Membrane	1	03-37537	Rubber seal	1
56-10827-00	Y-piece plastic	1	03-37526	Rubber seal	1x10
56-10826-00	Y-piece plastic angled	1	03-37137	Rubber seal	1x6
52-44120-00	T-valve steel	1	03-37437	Rubber seal	1
52-46411-00	Rubber bladder kit, yellow	1	52-40551-00	Membrane	1x4
52-46963-00	Rubber bladder kit, black	1	52-41151-00	Rubber sleeve	1
52-413 80-05	Tube	1	52-40590-00	Rubber sleeve	1
52-44100-00	Tube set	2			
52-17611-00	Stand telescope	1	52-46963-00		
				Contents of rubber sleeve kit, black	
52-17614	Stand frame	1		Rubber bladder with steel support for membrane	1
	Humidifier	1	53-54194-00	Connection ring	1
	Tube, blue with connections	1		Hose with tube nipple	1
	Tube, white with connections	1		Screw	2
	Patient tray	1		Instructions	1
	Instruction Manual	1		Tube PVC	1
	Engström 300 Instruction Manual	1			
	Instruction Manual Cascade	1	52-46411-00		
	Humidifier			Contents of rubber bladder, black	
				Rubber bladder, yellow	1
				Silicone grease MS4	1
			52-40589-50	Valve	1
			52-44100-00	Contents of tube kit	
				Y-piece plastic	1
				Y-piece plastic, angled	1
				V-bladder, plastic	2
				Rubber tube 1100 mm	1
				Tube with screw connection 1100 mm	1
				Tube with screw connection 750 mm	
				Tube with screw connection at both ends 450 mm	1

Technical data

Breathing frequency	12-35
Inspiration- Expiration ratio	1:2
Dosage of air or rebreathing gas	0-30 l/min
Variable end expiratory pressure	0-20 cm H ₂ O
Variable end expiratory resistance	0-∞
Venturi(neg. expiratory pressure)	0-10 cm H ₂ O
Rotameter O ₂	0-15,0-5 l/min
Rotameter N ₂ O	0-10 l/min
Dosage of air or rebreathing gas	0-30 l/min
Emptying pressure	50-110 cm H ₂ O
Weight	130 kg
Power	220-240 V

Contents list: Emergecy kit 81168550

Emergency handbag	1
Stethoscope	1
Laryngoscope	1
Spatula nr 1	1
Spatula nr 2	1
Spatula nr 3	1
Extra lamp to laryngoscope	1
breathing apparatus	1
blood-pressure gauge	1

<i>Name/Function</i>	Instrument hospital ward
<i>Manufacture typ</i>	
<i>Manufacture</i>	
<i>Made in</i>	
<i>Part number.</i>	8116980
<i>Production year</i>	
<i>Comment's</i>	



Contents list: Instrument hospital ward 8116980

<u>Description</u>	<u>Qty</u>
DRAPE CLAMP	3
MOUTH CATCH	1
NAIL RINSE	2
NAIL PLIERS	1
NEEDLE THREAD	3
BROW MIRROR	1
CLAMP STRAIGHT 14 CM	1
CLAMP STRAIGHT 15 CM	1
CLAMP STRAIGHT 18 CM	3
CLAMP CLAW 14 CM	4
CLAMP SEVERAL KINDS	5
TWEEZERS ANATOMICAL 13 CM	3
TWEEZERS ANATOMICAL 17 CM	1
TWEEZERS ANATOMICAL 20 CM	2
TWEEZERS SEVERAL KINDS	10
TWEEZERS CLAW 3x4	3
RAZOR	2
RESECTION HOOK	1
SCISSORS, SURGERY, SEVERAL KINDS 15 CM	7
SCISSORS BANDAGE	3
ACISSORS SEVERAL KINDS	2
THROAT PLIERS	1
BOWL STAINLESS	4
WOUND SYRINGE	2
WOUND LADLE	1
WOUND HOOK	4
TONGUE SPATULA	2