



# ACCESSORIES

for use with the  
**LIFEPAK® 20**  
DEFIBRILLATOR/MONITOR




Genuine accessories from Physio-Control, a division of Medtronic.



# POWER OPTIONS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR



Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know the LIFEPAK reputation for quality and reliability in our defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept any less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 20 Defibrillator/Monitor. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



## NiMH Rechargeable Internal Battery

Rechargeable nickel metal hydride 12-volt battery for use as AC power backup.  
 (One included per unit at no charge)..... **11141-00068**



## AC Power Cord

Powers the LIFEPAK 20 defibrillator/monitor and recharges the internal battery.  
 (One included per unit at no charge)..... **11140-00015**

# ECG MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR



## 5-Wire ECG Cable

Right-angle connector. 4-wire limb, plus 1 chest lead, labeled "C" on the device reports.....11110-000066



## 3-Wire ECG Cable

Right-angle connector. (One included per unit at no charge) .....11110-000029



## LIFE-PATCH® ECG Electrodes

Adult, pregelled (package of 3 electrodes).

11110-000001



## LIFE-PATCH ECG Electrodes

Adult, pregelled (package of 4 electrodes).

11110-000002



## Box of ECG Printer Paper

50mm x 30m (3 rolls per box).

11240-000013

# THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR



## LIFEPAK 20 Defibrillator/Monitor Standard Adult Detachable Hard Paddles

1 pair, for use with the LIFEPAK 20 defibrillator/monitor. Adult paddle plates are removable to reveal pediatric paddles.....**11130-00037**

### Adult Hard Paddle

Replacement attaches to adult hard paddle assembly P/N 11130-00037 .....**21330-001024**

## LIFEPAK® 12 Defibrillator/Monitor Standard Adult Hard Paddles

1 pair, for use with the LIFEPAK 20 or LIFEPAK 12 defibrillator/monitors...**11130-000001**

### External Sterilizable Paddles

1 pair, for use with the LIFEPAK 20 or LIFEPAK 12 defibrillator/monitors...**11134-000003**



### SIGNAGEL® Electrode Gel

One 8.5 oz. tube for use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback, and EMG.

**21300-005847**



### Internal Paddle Handles with Discharge Control

1 pair, for use with the LIFEPAK 20 or LIFEPAK 12 defibrillator/monitors.

**11131-000001**

### QUIK-COMBO® Therapy Cable

For use with the LIFEPAK 20 or LIFEPAK 12 defibrillator/monitors. (One included per unit at no charge)

**11110-000040**

### FAST-PATCH® Adapter Cable

For use with the QUIK-COMBO therapy cable for the LIFEPAK 20 or LIFEPAK 12 defibrillators.

**11110-000052**



### Internal Paddles

1" size (6.25" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000011**

### Internal Paddles

2" size (8.75" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000022**

### Internal Paddles

1.5" size (6.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000011**

### Internal Paddles

2.5" size (5.50" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000013**

### Internal Paddles

1.5" size (9.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000021**

### Internal Paddles

2.5" size (8.50" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000019**

### Internal Paddles

1.5" size (14.25" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000024**

### Internal Paddles

3.5" size (5.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000014**

### Internal Paddles

2" size (5.75" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000012**

### Internal Paddles

3.5" size (8.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000023**

# THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR

## EDGE SYSTEM™ ELECTRODES FOR PACING/DEFIBRILLATION/ECG WITH QUIK-COMBO CONNECTORS

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector ..... 11996-000091



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

11996-000017



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

11996-000090



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only on manual defibrillator/monitors, 12-month minimum shelf life at time of shipment.

11996-000093

## OTHER EDGE SYSTEM ELECTRODES FOR PACING/DEFIBRILLATION/ECG

18-month minimum shelf life at time of shipment

### EDGE System Electrodes with FAST-PATCH Connector

Must be used with QUIK-COMBO to FAST-PATCH cable

(11110-000052)..... 11996-000092

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR

## MASIMO SET® LNCS SENSORS AND CABLES



### Masimo SET LNC-4 LNCS Patient Cable

4-foot reusable connector cable .....11171-00024

### Masimo SET LNC-10 LNCS Patient Cable

10-foot reusable connector cable .....11171-00016

### Masimo SET LNC-14 LNCS Patient Cable

14-foot reusable connector cable .....11171-00025



### Masimo SET LNC-4-Ext LNCS Extension Cable

4-foot extension patient cable (May be used with devices equipped with Nellcor with conventional SpO<sub>2</sub>)

11171-00027



### Masimo SET LNCS DCIP Pediatric Reusable Sensor, for patients 10-50 kg.

Multiuse Sensor for patients 10-50kg.

11171-000018



### Masimo SET LNCS DCI Adult Reusable Sensor

Multiuse sensor for patients >30 kg .....11171-000017



### Masimo SET LNCS Aadx Adult Adhesive Sensor

Single patient use sensor for patients >30 kg (20/box).

11171-000019



### Masimo SET LNCS Pdx Pediatric Adhesive Sensor

Single patient use sensor for patients 10-50 kg (20/box).

11171-000020

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR

## MASIMO SET LNCS SENSORS AND CABLES (continued)



### Masimo SET LNCS Neo-L Neonatal Adhesive Sensor

Single patient use sensor for patients <3kg or >40kg (20/box).

11171-000021



### Masimo SET LNCS NeoPt-L Neonatal Adhesive Sensor

Single patient use sensor for patients <1kg (20/box).

11171-000022



### Masimo SET LNCS Inf-L Infant Adhesive Sensor

Single patient use sensor for patients 3–20 kg (20/box).

11171-000026

## MASIMO SET LNOP SENSORS AND CABLES



### Masimo SET PC-04 LNOP Patient Cable

4-foot reusable connector cable ..... 11171-000006

### Masimo SET PC-08 LNOP Patient Cable

8-foot reusable connector cable ..... 11171-000008

### Masimo SET PC-12 LNOP Patient Cable

12-foot reusable connector cable ..... 11171-000009



### Masimo SET LNOP DCI Adult Reusable Sensor

For patients >30 kg

11171-000007



### Masimo SET LNOP DCIP Reusable Pediatric Sensor

For patients 10-50 kg

11171-000010

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR


## MASIMO SET LNOP SENSORS AND CABLES (continued)



### Masimo SET LNOP Adt Adult Adhesive Sensor

Single patient use sensor for patients >30kg (20/box).

11171-000011



### Masimo SET LNOP Pdt Pediatric Adhesive Sensor

Single patient use sensor for patients 10–50 kg, (20/box).

11171-000012



### Masimo SET LNOP Neo Neonatal Adhesive Sensor

Single patient use sensor for patients <10kg (20/box).

11171-000013



### Masimo SET LNOP Neo-Pt Neonatal Adhesive Sensor

Single patient use sensor for patients <1kg (20/box).

11171-000014

---

## Adapter Cables

### MNC-1 Adapter Cable (10 ft)

Allows LIFEPAK 12 or LIFEPAK® 20 defibrillator/monitors with Masimo® SpO<sub>2</sub> to connect to Nellcor sensors.

11996-000183

### MNC-1 Adapter Cable (4 ft)

Allows LIFEPAK 12 or LIFEPAK 20 defibrillator/monitors with Masimo SpO<sub>2</sub> to connect to Nellcor sensors.

11996-000198



# CARRYING CASES AND MOUNTING OPTIONS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR

## Basic Carrying Case/ Accessories Organizer

Includes set of right and left pouches for additional storage of accessories.

11260-000018



## Top Pouch

For additional storage: insert in place of standard paddles .....11260-000020

## Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware.....11260-000019



## LIFEPAK 20 Defibrillator/Monitor Docking Station

Attaches to crash cart to lock LIFEPAK 20 defibrillator/monitor safely in place.....21330-000996

## Right Side Accessory Pouch for the LIFEPAK 20 Defibrillator/Monitor

Attaches to the LIFEPAK 20 defibrillator/monitor and holds items such as QUIK-COMBO therapy cable, electrodes, and SIGNAGEL electrode gel .....11260-000016

# LITERATURE AND TRAINING MATERIAL

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR

**LIFEPAK 20 Defibrillator/Monitor Operating Instructions .....26500-002441**

## **LIFEPAK 20 Defibrillator/Monitor Inservice Video**

AED defibrillation; VHS format .....26500-001216

## **LIFEPAK 20 Defibrillator/Monitor Inservice Video**

Manual defibrillation; VHS format.....26500-001217

**LIFEPAK 20 Defibrillator/Monitor Service Manual on CD-ROM ..26500-001212**

## **LIFEPAK 20 Defibrillator/Monitor Service Manual**

(Paper Version).....26500-001213

---

## **COMMUNICATION ACCESSORIES**

### **Serial Port Cable—6 foot**

LIFEPAK 20 defibrillator/monitor to PC .....11230-000018

### **Transport Configuration Cable—6 foot**

LIFEPAK 20 defibrillator/monitor to LIFEPAK 20 defibrillator .....11230-000019

# TRAINING TOOLS AND TESTERS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR



## QUIK-COMBO 12-Lead Patient Simulator

Simulates patient heart rhythms. Allows users to practice and demonstrate monitoring with the 12-lead or 5-wire ECG cable and defibrillating with the QUIK-COMBO therapy cable. Simulates: asystole, bradycardia, normal sinus rhythm, ventricular fibrillation (VF), ventricular tachycardia, and acute MI. Also simulates patient motion. Includes QUIK-COMBO cable, batteries, and operating instructions .....**11200-00005**

## QUIK-COMBO 3-Lead Patient Simulator

Simulates patient heart rhythms. Allows users to practice and demonstrate monitoring with the 3-wire ECG cable and defibrillating with the QUIK-COMBO therapy cable. Simulates: asystole, bradycardia, normal sinus rhythm, ventricular fibrillation (VF), and ventricular tachycardia. Also simulates patient motion. Includes QUIK-COMBO cable, batteries, and operating instructions.....**11201-00002**



## Clip-on Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and five (5) pair of reusable training electrodes.

**11250-000052**



## Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

**11998-000060**

## Replacement AED Training Electrodes

Includes five (5) pair of reusable training electrodes. Adhesive electrode portion only.

**11101-000003**

## QUIK-COMBO Therapy Cable Tester

Tests integrity of the QUIK-COMBO therapy cable.

**11998-000160**

## QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

**11113-000002**



© 2007 Physio-Control, Inc., a division of Metronic. LIFEPAK, LIFE•PATCH, FAST-PATCH, and QUIK-COMBO are registered trademarks of Medtronic, Inc. EDGE System and REDI-PAK are trademarks of Medtronic, Inc. MASIMO SET is a registered trademark of MASIMO CORPORATION. SIGNAGEL is a registered trademark of Parker Laboratories. Specifications and pricing are subject to change without notice. Not all products available outside the U.S.

MIN 3207693-000 / CAT 26500-002679