



ACCESSORIES

LIFEPAK[®] 20 DEFIBRILLATOR/MONITOR SERIES



GENUINE ACCESSORIES FROM PHYSIO-CONTROL

**PHYSIO
CONTROL**

POWER OPTIONS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES



Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know our reputation for quality and reliability in LIFEPAK defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you will not accept anything less. Inside this catalog you will find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with both LIFEPAK 20 and 20e defibrillator/monitors unless otherwise noted. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



NiMH Rechargeable Internal Battery for LIFEPAK 20 Defibrillator/Monitor
Rechargeable Nickel Metal Hydride 12-volt battery for use as AC power backup
11141-000068



AC Power Cord
Powers LIFEPAK 20 /20e defibrillator/monitors and recharges the internal battery.
(One included per unit at no charge)
11140-000015



Lithium-ion Rechargeable Internal Battery for LIFEPAK 20e Defibrillator/Monitor
Rechargeable Lithium-ion battery for use as AC power backup
11141-000112

ECG MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES



5-Wire ECG Cable

Right-angle connector. 4-wire limb, plus 1 chest lead, labeled "C" on the device reports

11110-000066



3-Wire ECG Cable

Right-angle connector. (One included per unit at no charge)

11110-000029



LIFE•PATCH® ECG Electrodes

Adult, pregelled (package of 3 electrodes).

11100-000001



LIFE•PATCH ECG Electrodes

Adult, pregelled (package of 4 electrodes).

11100-000002



Box of ECG Printer Paper

50mm x 30m (3 rolls per box).

11240-000013

THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES



Standard Adult Detachable Hard Paddles

1 pair, for use with LIFEPAK 20/20e defibrillator/monitors. Adult paddle plates are removable to reveal pediatric paddles

11130-00037

Adult Hard Paddle Plate Replacement

Replacement attaches to adult hard paddle assembly P/N 11130-00037

21330-001024

External Sterilizable Paddles

1 pair, for use with LIFEPAK 20/20e or LIFEPAK 12 defibrillator/monitors

11134-000003



SIGNAGEL® Electrode Gel

One 8.5 oz. tube for use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback, and EMG.

21300-005847



Internal Paddle Handles with Discharge Control

1 pair, for use with LIFEPAK 20/20e or LIFEPAK 12 defibrillator/monitors.

11131-000001



Internal Paddles

1" size (6.25" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000010

Internal Paddles

1.5" size (6.00" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000011

Internal Paddles

1.5" size (9.00" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000021

Internal Paddles

1.5" size (14.25" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000024

Internal Paddles

2" size (5.75" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000012

Internal Paddles

2" size (8.75" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000022

Internal Paddles

2.5" size (5.50" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000013

Internal Paddles

2.5" size (8.50" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000019

Internal Paddles

3.5" size (5.00" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000014

Internal Paddles

3.5" size (8.00" shaft length), 1 pair, requires Internal Paddle Handles.

11131-000023

THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES

EDGE SYSTEM™ ELECTRODES FOR PACING/DEFIBRILLATION/ECG WITH QUIK-COMBO CONNECTORS

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector
11996-000091



QUIK-COMBO® Therapy Cable
For use with LIFEPAK 20/20e or LIFEPAK 12 defibrillator/monitors. (One included per unit at no charge)
11110-000040

FAST-PATCH® Adapter Cable
For use with the QUIK-COMBO therapy cable for LIFEPAK 20/20e or LIFEPAK 12 defibrillator/monitors.
11110-000052



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System
11996-000017



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector
11996-000090



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector
For use only on manual defibrillator/monitors, 12-month minimum shelf life at time of shipment.
11996-000093

OTHER EDGE SYSTEM ELECTRODES FOR PACING/DEFIBRILLATION/ECG

18-month minimum shelf life at time of shipment

EDGE System Electrodes with FAST-PATCH Connector
Must be used with QUIK-COMBO to FAST-PATCH cable (11110-000052)
11996-000092

SPO₂ MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES

MASIMO® SET® LNCS SENSORS AND CABLES



Masimo SET LNC-4 LNCS Patient Cable
4-foot reusable connector cable
11171-000024

Masimo SET LNC-10 LNCS Patient Cable
10-foot reusable connector cable
11171-000016

Masimo SET LNC-14 LNCS Patient Cable
14-foot reusable connector cable
11171-000025



Masimo SET LNCS DCI Adult Reusable Sensor
Multiuse sensor for patients >30 kg
11171-000017



Masimo SET LNCS Aadx Adult Adhesive Sensor
Single patient use sensor for patients >30 kg (20/box).
11171-000019



Masimo SET LNCS DCIP Pediatric Reusable Sensor
Multiuse Sensor for patients 10 to 50 kg.
11171-000018



Masimo SET LNCS Neo-L Neonatal Adhesive Sensor
Single patient use sensor for patients <3 kg or >40 kg (20/box).
11171-000028



Masimo SET LNCS Pdx Pediatric Adhesive Sensor
Single patient use sensor for patients 10 to 50 kg (20/box).
11171-000020



Masimo SET LNCS NeoPt-L Neonatal Adhesive Sensor
Single patient use sensor for patients <1 kg (20/box).
11171-000029



Masimo SET LNCS Inf-L Infant Adhesive Sensor
Single patient use sensor for patients 3–20 kg (20/box).
11171-000031

SPO₂ MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES

MASIMO SET LNOP SENSORS AND CABLES



Masimo SET PC-04 LNOP Patient Cable

4-foot reusable connector cable
11171-000006

Masimo SET PC-08 LNOP Patient Cable

8-foot reusable connector cable
11171-000008

Masimo SET PC-12 LNOP Patient Cable

12-foot reusable connector cable
11171-000009



Masimo SET LNOP DCI Adult Reusable Sensor

For patients >30 kg
11171-000007



Masimo SET LNOP DCIP Reuseable Pediatric Sensor

For patients 10 to 50 kg
11171-000010



Masimo SET LNOP Adt Adult Adhesive Sensor

Single patient use sensor for patients >30 kg (20/box).
11171-000011



Masimo SET LNOP Pdt Pediatric Adhesive Sensor

Single patient use sensor for patients 10 to 50 kg (20/box).
11171-000012



Masimo SET LNOP Neo Neonatal Adhesive Sensor

Single patient use sensor for patients <10 kg (20/box).
11171-000013



Masimo SET LNOP Neo-Pt Neonatal Adhesive Sensor

Single patient use sensor for patients <1 kg (20/box).
11171-000014

SPO₂ MONITORING ACCESSORIES

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES

NELLCOR® SPO₂ SENSORS AND CABLES



DS100A DURASENSOR®
Adult reusable, includes 1 reusable clip,
1-year warranty.
11996-000060



**Oxisensor II
Disposable Pediatric Sensor**
(box of 24)
11996-000116

DURA-Y® Multisite Sensor
Reusable, includes 1 reusable sensor,
6-month warranty
11996-000106

**Disposable Adhesive
Bandage Wrap**
For OXI-A/N (2 bags of 50)
11996-000048



**Oxisensor II
Disposable Infant Sensor**
(box of 24)
11996-000115

Oxiband® Adult/Neonatal Sensor
Reusable, 3-month warranty (sensor),
includes 50 adhesive bandage wraps
11996-000061

Disposable Adhesive Foam Wrap
For OXI-P/I (2 bags of 50)
(Not available in Canada)
11996-000108

Oxiband Pediatric/Infant Sensor
Reusable, 3-month warranty (sensor),
includes 50 adhesive bandage wraps
11996-000062

**Disposable Adhesive
Bandage Wrap**
For OXI-P/I (2 bags of 50)
11996-000049

Disposable Adhesive Foam Wrap
For OXI-A/N (2 bags of 50)
(Not available in Canada)
11996-000110

Disposable Posey Wrap
For Oxiband (bag of 12)
(Not available in Canada)
11996-000107



**Oxisensor® II Disposable
Adult Sensor**
(box of 24)
11996-000113

Long cable (box of 24)
11996-000114



DEC-4 Cable Extension
Includes 2xN25 and 1xI20
11110-000042

DEC-8 Cable Extension
8-foot extension cable
11110-000176

Oxisensor II Disposable Neonatal Sensor
(box of 24)
11996-000117

Oxisensor II Disposable Adult Sensor Assorted Pack
Includes 2xD25 and 1xD25L
11996-000052

Oxisensor II Disposable Pediatric Sensor Assorted Pack
Includes 2xD20 and 1xI20
11996-000051

Oxisensor II Disposable Neonatal Sensor Assorted Pack
Includes 2xN25 and 1xI20
11996-000050

LIFEPAK 20e DEFIBRILLATOR/MONITOR ADAPTER CABLES

Compatible only with LIFEPAK 20e Defibrillator/Monitor
with the Masimo/Nellcor SpO₂ configuration

MNC-1 Adapter Cable (4 ft)
Allows specially configured LIFEPAK 20e
defibrillator/monitors to connect
to Nellcor sensors.
11996-000198

MNC-1 Adapter Cable (10 ft)
Allows specially configured LIFEPAK 20e
defibrillator/monitors to connect
to Nellcor sensors.
11996-000183

CARRYING CASES AND MOUNTING OPTIONS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES



Basic Carrying Case/ Accessories Organizer

Includes set of right and left pouches for additional storage of accessories.

11260-000018



Top Pouch

For additional storage: insert in place of standard paddles

11260-000020

Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware

11260-000019



Docking Station

Attaches to crash cart to lock LIFEPAK 20/20e defibrillator/monitor safely in place

21330-000996

Right Side Accessory Pouch

Attaches to LIFEPAK 20/20e defibrillator/monitors and holds items such as QUIK-COMBO therapy cable, electrodes, and SIGNAGEL electrode gel

11260-000016

TRAINING TOOLS AND TESTERS

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES



QUIK-COMBO 3-Lead Patient Simulator

Connects directly to your LIFEPAK defibrillator for safe, interactive training. Select from 17 ECG rhythms including fibrillations, tachycardias and bradycardias.

11996-000310

QUIK-COMBO 12-Lead Patient Simulator

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillations, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities. For use with LIFEPAK devices with a 12-lead ECG feature.

11996-000311



Clip-on Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and five (5) pair of reusable training electrodes.

11250-000052



Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

11998-000060

Replacement AED Training Electrodes

Includes five (5) pair of reusable training electrodes. Adhesive electrode portion only.

11101-000003



QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

11113-000002

LITERATURE AND TRAINING MATERIAL

for the LIFEPAK 20 DEFIBRILLATOR/MONITOR SERIES

LIFEPAK 20 Defibrillator/Monitor Operating Instructions

Also available as free online download at www.physio-control.com
26500-002538 (US English)

LIFEPAK 20 Defibrillator/Monitor Operating Instructions

26500-002441 (International English)

LIFEPAK 20e Defibrillator/Monitor Operating Instructions

Also available as free online download at www.physio-control.com
26500-002570 (US English)

LIFEPAK 20e Defibrillator/Monitor Operating Instructions

26500-002571 (International English)

LIFEPAK 20 Defibrillator/Monitor Inservice Video

AED defibrillation; VHS format Supports former AHA Guidelines 2000 Protocols
26500-001216

LIFEPAK 20 Defibrillator/Monitor Inservice Video

Manual defibrillation; VHS format Supports former AHA Guidelines 2000 Protocols
26500-001315

LIFEPAK 20/20e Defibrillator/Monitor Service Manual

on CD-ROM
26500-002703

LIFEPAK 20 Defibrillator/Monitor Service Manual

(Paper Version)
26500-001213

COMMUNICATION ACCESSORIES



Serial Port Cable—

6 foot

LIFEPAK 20/20e defibrillator/monitor
to PC

11230-000018

Transport Configuration Cable—

6 foot

LIFEPAK 20/20e defibrillator/monitor to
LIFEPAK 20/20e defibrillator/monitor

11230-000019

To order accessories for the LIFEPAK 20 defibrillator/monitor visit our webstore at <https://store.physio-control.com> (US only) or contact us at one of the numbers listed below.

**PHYSIO
CONTROL**

Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll Free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, ON
L6Y 0R3
Tel 888 879 0977
Fax 866 430 6115

©2010 Physio-Control, Inc. All rights reserved. All names herein are trademarks or registered trademarks of their respective owners. Specifications subject to change without notice. Not all products available worldwide. Please contact your local sales representative for a list of approved products in your region.
GDR 3300292_D

