



ACCESSORIES

for use with the
LIFEPAK® 500 AED



Genuine accessories from Physio-Control, a division of Medtronic.

**PHYSIO
CONTROL**

POWER OPTIONS / MISCELLANEOUS ACCESSORIES

for the LIFEPAK 500 AED

Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know the LIFEPAK reputation for quality and reliability in our defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept any less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 500 AED. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



Non-rechargeable Lithium Sulfur Dioxide Battery Pak

For use with the LIFEPAK 500 AED.

11141-000013

Non-rechargeable Lithium Sulfur Dioxide Battery Pak

For use with the LIFEPAK 500 DPS AED. Non-rechargeable battery, charcoal color

11141-000069

Rechargeable Sealed Lead-Acid Battery Pak

For use with the LIFEPAK 500 AED **11141-000002**



Battery Charger

For use with LIFEPAK 500 AED rechargeable SLA Battery Pak Power cord included.

11140-000002



AMBU® Res-Cue Mask® First Responder Kit

Includes 2 sets of disposable vinyl gloves, 1 reusable mouth barrier mask with valve and filter, 1 disposable razor, 1 disposable anti-microbial wipe, and 1 pair of trauma scissors. Can be attached to the AED handle.

11998-000321

AMBU Res-Cue Key® First Responder Kit

Includes 2 sets of disposable vinyl gloves, 1 disposable CPR barrier, 1 disposable razor, 1 disposable anti-microbial wipe, and 1 pair of trauma scissors. Can be stored in the AED carrying case

11998-000320

THERAPY DELIVERY / COMMUNICATION ACCESSORIES

for the LIFEPAK 500 AED

EDGE SYSTEM™ ELECTRODES FOR PACING/DEFIBRILLATION/ECG WITH QUIK-COMBO® CONNECTORS

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System..... 11996-000017



Infant/Child Reduced Energy Defibrillation Electrode Starter Kit

For use only with LIFEPAK 500 Biphasic AED with pink connector or LIFEPAK 1000 defibrillator, LIFEPAK EXPRESS® AED, or LIFEPAK CR® Plus AED. For use on children less than 8 years of age or under 55 lbs. 11101-000017

Replacement Infant/Child Reduced Energy Defibrillation Electrodes

11101-000016

COMMUNICATION ACCESSORIES



AED-to-Modem Cable

11150-000001

AED-to-PC Cable

For use with DATA TRANSFER™ 500 or QUIK-VIEW™ 500 programs.

11230-000001

Setup Transfer Cable

To transfer setup option configuration from unit to device for multiple unit systems.

11110-000050

External Modem for Connection to LIFEPAK 500 AED

Requires AED-to-modem cable.

11150-000010

CARRYING CASES / AED LOCATION SIGNS

for the LIFEPAK 500 AED



Complete Soft Shell Carrying Case.....11998-000014

Complete Soft Shell Carrying Case

LIFEPAK 500 DPS AED—with “stealth” non-reflective surface21300-001058

Replacement Top for Carrying Case 21300-004157

Replacement Shoulder Strap for Carrying Case.....11998-000012

AED/Oxygen Carry Bag.....21300-004834

Hard Shell, Water-tight Carrying Case11998-000021

Battery Pouch

Attaches to side of soft shell carrying case.....11220-000025



AED Location Sign

Flat, 8" x 8",
with logo 11210-000036

AED Location Sign

Flat, 10" x 10",
with logo 11210-000037

AED Location Sign

Flat, 8" x 8",
without logo 11210-000038

AED Location Sign

Flat, 10" x 10",
without logo 11210-000039



AED Location Sign

T-mount, 8" x 8",
with logo 11210-000031

AED Location Sign

T-mount, 10" x 10",
with logo 11210-000032

AED Location Sign

T-mount, 8" x 8",
without logo 11210-000040

AED Location Sign

T-mount, 10" x 10",
without logo 11210-000041



CABINETS AND MOUNTING OPTIONS

for the LIFEPAK 500 AED



AED Wall Cabinet with Alarm and Strobe

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel, with stainless steel trim. Surface mounted trim style with 7" return.....

11210-000028



AED Wall Cabinet with Alarm

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 7" return.

11220-000079



AED Wall Cabinet with Alarm

Recessed mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Recessed mounted trim style with 1.5" return.

11998-0000293

AED Wall Cabinet with Alarm, Fire Rated

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim with 3" return.

11210-000026

AED Wall Cabinet with Alarm, Fire Rated

Recessed, square edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Recessed mounted trim style with 1.5" return.

11210-000027



AED Floor Stand Cabinet with Alarm

Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. White epoxy painted steel 52" tall, 14" x 17-3/8" high x 7" deep.....

11210-000028

For additional cabinet options call JL Industries at (952) 835-6850 or visit their Web site at www.jlindustries.com.

CABINETS AND MOUNTING OPTIONS / TRAINING MATERIAL

for the LIFEPAK 500 AED

AED Wall Cabinet with Alarm

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 3" return.

11998-000292

AED Wall Cabinet with Alarm and Strobe

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with white trim. Surface mounted trim style with 7" return.

11220-000083

AED Floor Stand Cabinet with Alarm

Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Gray epoxy painted steel. 52" tall, 14" x 17-3/8" high x 7" deep.

11210-000029



AED Wall Cabinet with Alarm

Surface mount, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel, with stainless steel trim. Surface mounted trim style with 7" return.

11220-000076

AED Wall Cabinet with Alarm

Recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel finish wall cabinet with stainless steel trim. Recessed mounted trim style with 1.5" return.

11220-000078

AED Wall Cabinet with Alarm

Semi-recessed, rolled edges. Works with LIFEPAK 500, LIFEPAK 1000, LIFEPAK CR Plus, or LIFEPAK EXPRESS defibrillators. Steel wall cabinet with stainless steel trim. Semi-recessed mounted trim style with 3" return.

11220-000077

AED Cabinet Window Replacement Kit

21300-006797



Wall Mount Bracket for LIFEPAK 1000 or LIFEPAK 500 defibrillators

11210-000001

LITERATURE AND TRAINING MATERIAL

AED Instruction Card

Laminated easy reference AED and CPR instruction card. Fits inside the LIFEPAK 500 AED soft shell carrying case**26500-000185**

AED/CPR Reference Poster

Poster-sized easy reference AED and CPR instruction guide. Applicable for wall-mounting near an AED**26500-000268**

Operating Instructions: LIFEPAK 500 Automated External Defibrillator26500-002364

AED Challenge, CD-ROM

An interactive refresher training program that simulates cardiac arrest scenarios using LIFEPAK automated external defibrillators. PC training tool provides practice for AED personnel who are just learning to use an AED, and helps experienced AED users keep their skills sharp. Consistent with the American Heart Association's 2005 Guidelines. Not a replacement for initial training. Designed for responders who have completed initial CPR/AED classroom training**26996-000015**

AED Challenge®, Web-Enabled

An interactive refresher training program that simulates cardiac arrest scenarios using LIFEPAK automated external defibrillators. Web-based training tool provides practice for AED personnel who are just learning to use an AED, and helps experienced AED users keep their skills sharp. Consistent with the American Heart Association's 2005 Guidelines. Not a replacement for initial training. Designed for responders who have completed initial CPR/AED classroom training. Includes online tracking, reporting, and administration tools. Annual license for individual users is required. List price is per user/per year. Quantity discounts are available. Web version is currently available to Heart Safe Solutions customers only**26996-000014**

Inservice Video: LIFEPAK 500 Automated External Defibrillator

VHS format, NTSC.....**26500-000037**

Service Manual: LIFEPAK 500 Automated External Defibrillator

CD-ROM, Microsoft® Windows® 3.1, or Windows 95.**26500-000036**

For additional cabinet options call JL Industries at (952) 835-6850 or visit their Web site at www.jlindustries.com.

TRAINING TOOLS AND TESTERS

for the LIFEPAK 500 AED



LIFEPAK 500T AED Training System

Includes one (1) AED trainer, remote control and cable, simulated battery pak, training electrodes, carrying case (soft-sided, gray), AED instruction card, operating instructions, and batteries for AED trainer and remote control.....**11250-000053**



Replacement Remote Control and Cable

(For 11250-000053)

11250-000005



Replacement Simulated Battery Pak

(For 11250-000053)

11250-000006



AED Training Electrode Set

Consists of five (5) pair of reusable AED training electrodes, cable connector assembly, and reusable foil pouch. Reusable training electrodes good for 30 applications.

11101-000004

Replacement AED Training Electrodes

Five (5) pair of reusable training electrodes. Adhesive electrode portion only.

11101-000003

Cable/Connector Assembly and Reusable Foil Pouch for Training Electrodes

Cable, connector, and reusable foil pouch only.

11101-000006



Test Load

Checks the AED charging circuits during a typical charging cycle. During this test, the AED delivers a low energy test shock. Connects directly into the QUIK-COMBO therapy cable on the defibrillator.

11998-000022

Infant/Child AED Training Electrode Set

Consists of five(5) pair of reusable AED training electrodes, cable connector assembly, and reusable foil pouch. Reusable training electrodes good for 30 applications.

11250-000045

Replacement Carrying Case (Softcase)

(For 11250-000053)

11250-000004

Replacement Analyze Key Cover Labels

(For 11250-000053)

21501-000158

Replacement Infant/Child AED Training Electrodes

Five (5) pair of reusable training electrodes. Adhesive electrode portion only. Cable/connector assembly and reusable foil pouch sold separately.

11250-000042

LIFEPAK 500T AED Training System Operating Instructions

(For 11250-000053)

26500-000294



© 2007 Physio-Control, Inc., a division of Medtronic. LIFEPAK, LIFEPAK EXPRESS, LIFEPAK CR, and QUIK-COMBO are registered trademarks of Medtronic, Inc. EDGE System, DATA TRANSFER, QUICK-VIEW and REDI-PAK are trademarks of Medtronic, Inc. AMBU, Res-Cue Key and Res-Cue Mask are registered trademarks of AMBU, Inc. AED Challenge software is a registered trademark of Insight Instructional Media, LLC. Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Specifications and pricing are subject to change without notice. Not all products available outside the U.S.

MIN 3207699-000 / CAT 26500-002682

