

Invacare®
Formula™ PTO
Plus for Pronto®
M71, M91™ & Storm
Series TDX™

Invacare® Formula™ PTO Plus (Power Tilt Only)
for Pronto® M71, M91™ and Storm Series® TDX™



Yes, you can.™



The Invacare® Formula™ PTO Plus (Power Tilt Only) is designed for active consumers who require a basic tilt for positioning and pressure relief. The PTO Plus is available on the Pronto® M71™ and M91™ and the Storm Series® TDX™, while still maintaining sleek style, indoor and outdoor maneuverability and a tight turning radius. With the new and improved features of the PTO Plus on the M71, M91 and TDX consumers are offered the powered tilt they need on the bases they want.



Invacare® Formula™ PTO Plus for Storm Series® TDX™

FEATURES:

- 18" Seat to floor height for easy access to everyday environments
- Weight capacity of 275 lb. or optional 375 lb. when combined with heavy-duty motor options
- Up to 55° of tilt for pressure relief (0° – 55°, 5° – 60°)
- Weight Shifting Technology
- MK5™ Electronics
- Up to 7.5 mph
- SureStep™ for smooth driving over transitions up to 3 inches in height
- Stability Lock for added stability and control



Invacare® Formula™ PTO Plus for Pronto® M91™ with SureStep™

FEATURES:

- 18" Seat to floor height for easy access to everyday environments
- Weight capacity of 275 lb. or optional 375 lb. with heavy-duty motor package
- Up to 55° of tilt for pressure relief (0°–55°, 5°–60°)
- Weight Shifting Technology
- MK5™ NX-80 Electronics
- Up to 6.5 mph
- SureStep for smooth driving to overcome transitions 3" in height



Invacare Formula PTO Plus for Pronto® M71™ with SureStep®

FEATURES:

- 18" Seat to floor height for easy access to everyday environments
- Weight capacity of 275 lb.
- Up to 55° of tilt for pressure relief (0°–55°, 5°–60°)
- Weight Shifting Technology
- MK5 NX-80 Electronics
- Up to 4.25 mph
- SureStep for smooth driving to overcome transitions 2" in height

Battery Access



The PTO Plus integrates into the base of the chairs while still maintaining easy battery access through a simple flip back design



Yes, you can.™

Specifications

	M91™	M71™	TDX™
Tilt Range	0° – 55° 5° – 60°	0° – 55° 5° – 60°	0° – 55° 5° – 60°

Seat-to-floor height*	18"	18"	18"
-----------------------	-----	-----	-----

* At 0° of tilt

Overall dimensions

Length	39"	34"	35"
Width (excluding joystick)	26"	24"	25"
Height*	39"	33.5"	42"

* Height measured without headrests and with a 20° back on the Pronto M71 and M91 and a 24° back on TDX

PTO Plus Seat dimensions

Widths	16" – 22"	16" – 20"	16" – 22"
Depths	16" – 22"	16" – 20"	16" – 22"
Back Height	16" – 24"	16" – 24"	16" – 24"

Weight

With Batteries	296 lb.	231 lb.	361 lb*
----------------	---------	---------	---------

* Product Weight on TDX was measured on a TDX5 with Gearless Brushless GB motors and Group 24 batteries

Wheel and Caster Sizes

Front Casters	6" x 2"	6" x 2"	6" x 2"
Drive Wheels	14" x 3"	10.5" x 3.5"	12" x 2.5"* 14" x 3"
Rear Casters	6" x 2"	6" x 2"	6" x 2"

* On TDX standard wheel size dependent on base pick

Maximum Speed	6.5 mph	4.25 mph	6 - 7.5 mph*
---------------	---------	----------	--------------

* Speed on TDX dependent on base pick

Incline Capability	9°	9°	9°
Weight Capacity	275 lb. 375 lb.	275 lb.	275 lb. 375 lb.

* Weight capacity dependent on motor package

Turning Radius

Center Mount	19"	19"	22"
--------------	-----	-----	-----

Limited Warranties

Base Frame	60 months	60 months	60 months
Seat Frame	36 months	36 months	36 months
Electronics/Actuator	12 months	12 months	12 months

**Due to manufacturing variances, all dimensions, weights, and degrees may be ±.5", ±10%, or ± 3°

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44036-2125
(800) 333-6900

Canada

Invacare Canada Inc.
570 Matheson Blvd E
Unit 8
Mississauga Ontario
L4Z 4G4 Canada

Invacare, the blue medallion design, Pronto and Storm Series are registered trademarks and Yes, you can., Formula, M71, M91, TDX, MK₅, SureStep, NX and GB are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2004

Invacare Corporation

Form No. 04-093



"Invacare says 'Yes, you can', and I couldn't have said it better myself."

Invacare® Formula™ PTO Plus (Power Tilt Only)
for Pronto® M71, M91 and Storm Series TDX™