Troubleshooting:

"the application has failed to start because mcImcrrt76.dll was not found." - this error occurs due to the MCR failing to install properly either from the download or from the CD.

We ask that you try running the SaSAT exe again, and the error still persists then you can try the following:

1.) Obtain and install a separate copy of the MCRInstaller. (Note: a copy of the MCRInstaller should be included with the CD, but may fail to install properly on some systems). Separate copies of the MCR are available separately by request. If required, please contact us for more information or details.

Or,

2.) If you were able to successfully install the MCR but are still encountering this error, your systems environment variables may need to be changed. We provide the link from The Mathworks which provides an explanation of the process (see point 2):

http://www.mathworks.com/support/solutions/data/1-1IW46N.html

If you have any further problems, please do not hesitate to contact us.