## **GPComplete Installation Manual**

- 1. Insert the GPComplete CD.
- 2. If you have autorun enabled then it should start the GPComplete setup program automatically. Otherwise open the CD drive and double-click setup.exe.

		Name +	Size	Туре	Date Modified
File and Folder Tasks	*	🛅 dotnet	an and a starting	File Folder	20/07/2006 4:00 PM
- 11 Jan		irebird		File Folder	20/07/2006 4:00 PM
new folder		🗀 ikey		File Folder	20/07/2006 4:00 PM
Publish this folder to the Web		🛅 jre		File Folder	20/07/2006 4:00 PM
Share this folder		🚞 pki		File Folder	20/07/2006 4:00 PM
Share this folder		💁 _setup.dll	149 KB	Application Extension	17/05/2006 11:21 AM
		📑 autorun.inf	1 KB	Setup Information	7/07/2006 11:31 AM
Other Places	\$	🛃 data 1.cab	4,421 KB	Cabinet File	6/07/2006 4:14 PM
ourer risees		🛅 data 1.hdr	125 KB	HDR. File	6/07/2006 4:14 PM
🛅 Disk Images		🛃 data2.cab	169,369 KB	Cabinet File	6/07/2006 4:14 PM
My Documents		Setup.dll	540 KB	Application Extension	6/07/2006 4:13 PM
Shared Documents		avout bin	1 KB	BIN File	6/07/2006 4:14 PM
😡 My Computer		(Elsetup.exe)	445 KB	Application	24/05/2006 12:10 PM
<b>a</b>		setup.ini	1 KB	Configuration Settings	6/07/2006 4:13 PM
🧐 My Network Places		🖾 setup.inx	244 KB	InstallShield Compil	21/07/2006 11:13 AM
		📅 setup.isn	423 KB	ISN File	30/06/2006 2:11 PM

3. Click Next.

	Welcome to the InstallShield Wizard for GPComplete The InstallShield Wizard will install GPComplete on your computer. To	continue, click Next.
nstalionHeid in (	< Back	Cancel

4. Please read the license agreement and if you accept the terms of the agreement click Next.

THE PRACTICE MANAGEMENT SOFTWARE COMPANY SOFTWARE LICENCE ACREEMENT         I.       The Agreement         I.       The Agreement         This Software Licence Agreement (Agreement) sets the terms and conditions between you and The Practice Management Software Company ("PMSC") in regards to your use of the software you are currently installing on your computer, and any updates, upgrades or maintenance releases, related services, materials, and documentation that PMSC may provide (Software).         1.1       Acceptance of the Agreement         By installing the Software, you confirm that you:       (a)         (a)       accept and agree to be bound by the Agreement;         (b)       acknowledge that, to the maximum extent permitted by applicable law,	
this Agreement sets forth your exclusive remedies as to any claims you may have related to the Software, and  I accept the terms of the license agreement  Etint	<b>-</b> 1

5. Type in the user name and company name and then click Next.

PComplete - InstallShield	Wizard	2
Customer Information Please enter your informatio	n.	
	Please enter your name and the name of the company for which you work. User Name: Company Name:	
n <b>isemanik</b> iated	< Back Next >	Cancel

6. Click the <u>Server and Client</u> button to install the server and client on this machine. (Only install the server on *one* machine in the network.)

Click the <u>Client</u> button to install only the client on a machine- you will need to have already installed the Server on *one* other machine in the network.

GPComplete - InstallShield	Wizard	×
Select Features Select the feature you want	to install.	
	Please select the feature that you want to install.           Image: Client         Client           Install this if you want only the client software on this machine. One server be installed in the network           Image: Client         Server and Client           A server and the client will be installed           Click on Server and Client to install the server and client software on this machine.           Note: You should install the server on only ONE comput if you want to install the software on additional workstat you should install only the client by using the top button	nt
Indratishieldijeld	K Back Next> Ca	ancel

7. Select the folder where you want the software to be installed. Click Next if you want to leave the current location unchanged (Recommended)



8. Click the Install button to begin the installation process.



9. The installation progress bar will appear.



10. Skip this step if you are installing only the Client software (as selected by you in step 5). Choose the folder where you want to install the database. We recommend that the drive where you are installing the database should have at least 1 GB of free space.

Complete - InstallShield	Vizard	
Database Folder		
	Please select a folder to copy database	
	Destination Folder	
	C:\GPCDATA\	BIOWSE
<b>Searce He</b> lderic d	<back next=""></back>	Cancel

11. The installation of GPComplete is finished so click on Finish.



12. You will need to restart the machine for all the changes to take effect.

	InstallShield Wizard Complete Setup finished
	<ul> <li>Yes, I want to restart my computer now.</li> <li>No, I will restart my computer later.</li> <li>Please reboot the computer for GPComplete to complete installation. The program might not work unless you reboot</li> </ul>
Independent de la companya de la company	< Back Finish Cancel