Options and Accessories

- · Crutch and cane holder
- Oxygen holder
- · Safety flag
- Rear basket
- Docking device
- · Elevating legrest kit

Specifications

Overall dimensions	
Length with	21 5" 22 5"
footplate folded	31.5"-33.5"
Length with	4011 4011
footplate down	40"- 42"
Width	23''
Height	39.75"- 45.75"
Seat-to-floor height	19.25", 20", 20.75"
Seat-to-footplate height	14"-15.5"
Seat dimensions	
Width	17''
Depth	17.5"
Back height*	
(with headrest)	21"-27"
Fixed-back angle	
adjustment	80°-115°
Seat cushion thickness	3.5"
Back cushion thickness	3.5"
Arm dimensions	-
Length	14.5"-19.5"
Distance between	
arm pads	15"-21"
Maximum speed	4.5 mph
Incline capability	9° grade
Weight capacity	250 lb.
Overall turning radius	23''
Weight (with batteries)	205 lb.
Warranty	
Base frame	3 years
Seat frame	3 years
Electronics &	
motor/gearbox	1 year

^{*} Back height, 15"; Hendrest, 6"; Hendrest adjustment range, 6".



Yes, you can[™]

Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900


Canada 570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare is a registered trademark and Yes, you can. and Pronto M6 are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2002 Invacare Corporation Form No. 00-270 Rev. 03/02



"Invacare says 'Yes, you can,' and I couldn't have said it better myself."



ere's the independence, comfort and peace of mind you're looking for: Introducing the Invacare® Pronto™ M6™ power wheelchair.

The compact size and stable, full-time, 6-wheel suspension of the Pronto M6 allows you to move around freely and easily in your home. Its true mid-wheel drive system and smooth hand control make the Pronto M6 easier and more natural to drive for first-time users.

Maneuverability, reliability, comfort and solid good looks – the Invacare Pronto M6 has it all. And it all adds up to the independence and respectability you're looking for.

FEATURES

- True mid-wheel drive system and compact size for tightest overall turning radius in its class make it ideal for small homes with narrow hallways and limited space
- Full-time, 6-wheel suspension for outstanding traction/solid feel eliminates the rocking instability of most mid-wheel drive wheelchairs
- Powerful two-motor design and flat-free tires for performance you can depend on
- Invacare's industry-leading electronics; hand control is easy to learn and gives you full control of your wheelchair
- Reinforced carbon steel frame for long-lasting durability
- Van-style padded vinyl seat and numerous seating adjustments to ensure a perfect fit for your all-day comfort
- Adjustable armrests and footplate flip up to make it easier for you to get in and out of wheelchair
- Hand control conveniently retracts out of the way and allows you to drive right up to tables
- On-board charger and maintenance-free batteries for simple overnight charging
- Meets all RESNA standards for product performance





On-board battery charger means $\ensuremath{\mathsf{essy}}$ charging .



Retractable hand control enables you to get closer to tables, while armests easily adjust for confort and ease.



Six wheels on the ground for easy and stable driving over most thresholds.