JIMMETM (UL8WPB, UL8WPBS)



Shoprider[®] introduces another first in powered mobility with its extremely lightweight, portable, userfriendly Jimmie[™] power chair. The Shoprider[®] Jimmie[™] features the latest in wireless take-apart technology, electronics and high efficiency in-line motors. All these advanced features plus many more adjustability features are stylishly packaged into a comfortable power chair designed to move with your active lifestyle today. The Jimmie[™] is Travel-Friendly too with One-Touch effortless disassembly and transportability.

<image>

JIMMIE[™](UL8WPB, UL8WPBS)

- Connectorless technology for easily assembly
- Extra light total weight 93.5 lbs (UL8WPB)
- Easy to dismantle into three major parts: Seat w/ controller: 24.7 lbs (UL8WPB) Front Chassis w/ batts: 35.4 lbs Rear Chassis w/ drive train: 29.9 lbs
- Adjustable seat height
- Padded full function, armrest (width, height & angle adjustable)
- Removable fold up footrest with full function (height, length & angle adjustable)

| JIMMIE [™] Specifications | | UL8WPB | UL8WPBS |
|---|--------|-----------------------|---------|
| Overall Dimensions ($L^{*1} \times W \times H^{*2}$) | in | 33.5 x 21.7 x 19-22 | |
| Number, Size of Tires (front; mid; rear) | ea, in | 2, 6"; 2, 9" | |
| Fish-On/Capt seat w controller assy; | lb | 24.7 | 37.3 |
| Front chassis w/ batteries | lb | 35.4 | |
| Rear chassis w/ drive train | lb | 29.9 | |
| Total Weight (w/Batteries, Seat) | lb | 93.5 | 106.1 |
| Battery pack weight | lb | 19.3 | |
| Battery Capacity | Ah | 12v -12 Ah x 2 | |
| Charger | Amp | 2 amp Off-Board | |
| Output Power of Motor (Reference Only) | hp | 0.45 x2 | |
| Turning Radius | in | 15.5 | |
| Max. Single User Weight (On the Level Road) | lb | 250 | |
| Maximum Speed*3 (On the Level Road) | mph | 3.75 | |
| Suggested Climbing Angle | %, lbs | Shown as illustration | |
| Range *4 (per charge with 12 Ah Battery/ 198 lbs User Weight)(After the Battery and Mechanical Moving Parts fully Break In) | Mile | Up to 10 miles | |

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. eserves the rights of any changes to the unit. Include the anti-tip wheel. 1) Depending on the height of the 21 seat back When proceeding up any *31 incline, please move your seat to the most forward position or if you have the deluxe seat, make sure that it is in the 90 degree (upright) position. *4) The actual driving range varies with the factors shown as below. A. The weight of occupant B. Level ground C. Battery Condition D. Type of Charger E. Ambient Temperature F. The Way of Driving G. ETC.

Length: 33.5" Width: 21.7" Height: 19"-22" (Ground to Seat) Weight: 106lbs (UL8WPBS) 93.5lbs (UL8WPB) Maximum Speed: 3.75 mile Maximum User weight: 250 lbs

Safe Climbing Angle: 10 % (6°) Turning Radius: 15.5" Battery range: Up to 10 miles



_ _ 10° (17), 90 (198) - - - 6° (10), 100 (220)



Toll-Free: 800.743.0772 www.shoprider.com

Single User