

Fast. Easy. Accurate.

TALKING BLOOD PRESSURE MONITOR UA-1030T



3A
TECHNOLOGY
Advanced cuff design
Advanced inflation
Advanced deflation

**WORLD CLASS™
MONITOR**

 Clinically
Validated
for Accuracy

UA-1030T, part of the World Class Monitor family.

The UA-1030T Talking Blood Pressure Monitor improves upon all of the innovations A&D Medical has created over its 30+ year history. This fully featured monitor is equipped with innovations in cuff design, inflation and deflation (3A Technology). To speed up readings users can preset pressure levels or set them to automatic. The UA-1030T offers the option of TriCheck measurement which allows users to take three consecutive readings and receive an average, complying with the American Heart Association's latest recommendation. The UA-1030T has the option to give all readings and instructions in English, French or Spanish.

People resist taking their blood pressure because it can be uncomfortable. As a cuff constricts during a measurement, wrinkles form on the inside of the cuff. The material bunches up against a user's arm and creates discomfort. The new SmoothFit Cuff is designed to seamlessly conform to the user's arm and give the most comfortable experience possible. It comes standard with the UA-1030T Talking Blood Pressure Monitor.

- **TriCheck** – Three consecutive measurements and advanced averaging.
- **Pressure Rating Indicator** – Providing a snapshot of the user's blood pressure classification according to the US JNC7 and WHO guidelines.
- **Irregular Heartbeat Detection** – Alerts the user to the presence of an irregular heartbeat while still providing an accurate reading.
- **Advanced Warning System** – Alerts the user to the most common problems while monitoring, and helps to avoid them.
- **3A Technology** – Advanced cuff design, advanced inflation and advanced deflation, which translates to faster more accurate and more comfortable readings.
- **SmoothFit Cuff** – Seamless conforming cuff to take the pain out of monitoring.


For questions about A&D Medical or LifeSource products, call us toll-free at 1-800-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.



- TriCheck™ Measurement Feature
- 90 Memory Recall
- Readings and Instructions Given Verbally in English, French or Spanish
- Clock, Time & Date Stamp
- Pressure Rating Indicator™
- Irregular Heartbeat Detection
- Adjustable Pressure Setting
- Advanced Warning System
- Carrying Case Included
- AC Adapter Included
- SmoothFit™ Cuff

UA-1030T Talking Blood Pressure Monitor

To obtain an accurate reading, it is important to use the correct cuff size. A cuff that is too large will produce a reading that is lower than the correct blood pressure; a cuff that is too small will produce a reading that is higher than the correct blood pressure. The UA-1030T Talking Blood Pressure Monitor works with many A&D Medical and LifeSource brand cuffs.



PART NUMBER	UPC NUMBER	ACCESSORY
UA-1030T	0-93764-60293-1	SmoothFit Cuff 9.0"-14.6" (23-37cm)

Each unit comes with a trilingual instruction manual with log book, quick start card, storage case, AC Adapter and 4 AA batteries.

Accessories

PART NUMBER	UPC NUMBER	ACCESSORY
UA-289	0-93764-60297-9	SlimFit Small Cuff 6.3"-9.4" (16-24cm)
UA-290	0-93764-60298-6	SlimFit Medium Cuff 9.0"-14.6" (23-37cm)
UA-291	0-93764-60299-3	SlimFit Large Cuff 12.2"-17.6" (31-45cm)
UA-392	0-93764-60300-6	SmoothFit Cuff 9.0"-14.6" (23-37cm)
TB:233	0-93764-60235-1	AC Power Adapter

Technical Specifications

Dimensions	5.5" (w) x 4.1" (l) x 2.4" (h) (140mm x 105mm x 60mm)
Weight	10.6 oz. (300g) without batteries
Display	Digital, 20/10-mm character height. Pressure/pulse displayed simultaneously
Memory	90 readings
Pressurization	Automatic, using micropump
Deflation	Constant air release valve system
Type	Oscillometric
Accuracy	Pressure: ± 3 mmHg or $\pm 2\%$, whichever is greater. Pulse: $\pm 5\%$
Pressure Range	Pressure: 0 mmHg to 299 mmHg. Pulse: 40 pulses to 180 pulses
Power Source	4 type 'AA' (1.5 volt) alkaline batteries or 120V AC Adapter (TB:233)
Recommended Operating Environment	50°F to 100°F (10°C to 40°C). Less than 85% relative humidity
Recommended Storage Environment	-4°F to 140°F (-20°C to 60°C). Less than 95% relative humidity

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the American National Mediums Institute (ANSI/AAMI SP10) for electronic sphygmomanometers. Manufactured by an ISO 9000 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	Dimensions (L x W x H)	Weight	Individual Item Cube
UA-1030T	7.1 x 4.0 x 8.5	2.1	0.14
Unit Case (10 units)	W x D x H (inches)	WT (lb.)	Case Cube (cubic feet)
UA-1030T	21 x 15 x 9.5	22.6	1.73



1756 Automation Parkway
San Jose, CA 95131
1-888-726-9966
www.andmedical.com

LifeSource, Pressure Rating Indicator, TriCheck, SlimFit, SmoothFit, 3A Technology and World Class Monitor are trademarks and registered trademarks of A&D Company and/or A&D Medical. ©2012 A&D Medical. All rights reserved.

Ask your A&D Medical Sales Representative about other A&D Medical & Wellness Connected products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes
- Activity Monitors

For questions about A&D Medical products, call us toll-free at 1-888-726-9966 or visit one of our websites at: www.andmedical.com, www.WellnessConnected.com, www.LifeSourceOnline.com