

Respirator Fit Testing Kits and Accessories



Convenient qualitative fit testing kits and quantitative fit test probes and adapters to make fittesting as easy as possible. All qualitative kits include complimentary instructional video and CD-ROM on how to fit test. The quantitative accessories make for quick conversion of an AOSafety[®] air purifying or airline respirator to fit test mode for those workers who want to be tested in their own masks.

Respirator Fit Testing Kits and Accessories

AOSafety[®] Respirator Fit Testing Kits and Accessories

All kits include training tools for showing you how to do the various fit tests—easy to follow video and CD-ROM! Also, all written instructions, including those for the test subjects, are in English, French and Spanish. Fit test solutions and other fit test kit components can be purchased separately.

Isoamyl Acetate-For use with elastomeric respirators.

Kit includes all components necessary to conduct the test, packaged in a heavy gauge plastic, foampadded carrying case. Just add paper towels! The unique flexible vinyl fit test tent has an easy-toassemble frame that hangs on any door jamb for quick set-up. It fits the largest test subjects, yet disassembles and folds to fit flat in a clear carrying pouch which tucks away neatly inside the carrying case.

Bitrex and Saccharin Fit Test Kits-For use with filtering facepieces or elastomeric respirators.

All items needed for the test are arranged in a foam-padded carrying case with vinyl handle. The head and shoulder fit-test hood is quick to don with a clear window and opening for inserting the nebulizer nozzle for the test. Saccharin and Bitrex kits are identical with the exception of the challenge agent.





AOSafety [®] Qualitative Fit Test Kits						
Product Code	Description					
50007-00000	Isoamyl Acetate Fit Test Kit					
50005-00000	Saccharin Fit Test Kit					
50006-00000	Bitrex Fit Test Kit					

Respirator Quantitative Fit-Testing Accessories.

For those who prefer quantitative fit testing, insertable probes and screw-in adapters are available to quickly convert any AOSafety elastomeric mask to fit-test mode.

AOSafety® El	astomeric Respirator	Accessories for Converting to Quantitative Fit Test Mode		
Product Code	Catalogue Number	Description		
51170-00000	51170	Half Mask Side Port Plug with Probe for Use with PortaCount®		
52023-00000	52023	Full Face Mask Front Fitting Plug with Probe for Use with PortaCount $^{\circ}$		
51171-00000	51171	Quarter, Half, and Full Face Mask Fit Test Adapter Kit for Use with $\text{PortaCount}^{\circledast}$		

For	For technical information:		To order or locate your nearest stocking distributor:				
Pho Fax		Other Countries 1-508-764-5721, 5273 1-508-764-5640	Phone: Fax:	USA 1-800-225-9038 1-800-488-8007	Canada 1-800-387-9681 1-800-488-8007		Countries 1-317-692-6643 1-317-692-6570





Made in the USA

Aearo Company Southbridge Operations WWW.aosafety.com Aearo Company • 8001 Woodland Drive • Indianapolis, IN 46278 A0Safety® is a trademark licensed to Aearo Company. ©2000 Aearo Company. 34051 R3/00 PortaCount® is a registered trademark of TSI.

