Invacare® MNS500 Low Air Loss Mattress Replacement System





Yes, you can:

- Portable, lightweight
  mattress replacement system
- Available in standard, narrow and bariatric sizes

# Invacare® MNS500 Low Air Loss Mattress Replacement System

The Invacare<sup>®</sup> MNS500 is a portable low air loss therapy mattress replacement system that provides the therapeutic benefits of pressure relief and skin moisture management. The MNS500 also provides the benefit of true low air loss for pressure ulcer therapy.



- MNS500 easily attaches to most medical beds
- Full depth mattress with 16 individual air cells
- Built-in lower safety mattress supports the patient in the event of power failure

### **Specifications**

Model no. MNS500-S (Standard width) Model no. MNS500-N (Narrow width) Model no. MNS500-B (Bariatric width)

Flectrical	
Power 115vac, 50/60Hz, 6A	
Leakage current	Less than 100 microamps
Dimensions	
Standard mattress	80" L x 36" W x 8.5" H
Narrow mattress	80" L x 32" W x 8.5" H
Bariatric mattress	80" L x 42" W x 8.5" H
Power unit	11" W x 12.5" H x 5.25" D
Weight	
Standard mattress	22 lb.
Narrow mattress	20 lb.
Bariatric mattress	24 lb.
Power unit	9 lb.
Weight Capacity	
MNS500-S	350lb.
MNS500-N	350lb.
MNS500-B	600lb.
Safety Code approval CE, UL2601, CSA, CTL	
Material	
Mattress cells	200 denier, polyurethane
coated nylon	
CA# 117 Fire Retardency	
Therapy cover Ventex Recovery 5 with a polyurethane,	
breathable membrane layer. Polyester	
	, CA# 117 Fire Retardency
Limited warranty	
Power unit	1 year
Mattress	6 months



#### Invacare Corporation www.invacare.com

#### USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks of Invacare Corporation. Specifications are subject to change without notification.

Manufactured for Invacare Corporation by Genadyne Biotechnologies, Inc.

© 2005 Invacare Corporation

Form No. 04-217 Rev. 10/05



"Invacare says 'Yes, you can."' and I couldn't have said it better myself.

Tal

## Features

- 135 liters per minute of fresh, dry air circulates between the patient and the support surface
- 8" air cell height for deep patient immersion
- Waterproof, vapor permeable therapy pad provides a low friction surface
- Provides pressure relief below capillary closure
- Pulsate option every 45 seconds, the pressure increases and decreases in each cell

### User-Friendly Controls

- · Auto-firm mode for patient transfers or nursing care
- Auto-firm mode automatically cancels after thirty minutes
- · Pressure control setting for easy comfort adjustment
- Fowler button on key pad with LED
- Single connector rapid attach and CPR deflate
- Low pressure and power failure alarm

